Blom-Singer®

voice restoration systems

37-453-01 Rev. D

Laryngectomy Tube	3	ENG
Laryngektomietubus	10	DEU
Tubo de Laringuectomía	17	ESP
Tube de laryngectomie	24	FRA
Tubo per laringectomia	31	ITA
Laryngectomiebujsje	38	NED
Tubo de laringectomia	46	POR

INHEALTH®
T E C H N O L O G I E S

User Information Data Sheet

[Informationsblatt für Benutzer / Hoja de información del Usuario /](#)
[Fiche d'information destinée à l'utilisateur / Scheda informativa per l'utente /](#)
[Informatieblad voor de gebruiker / Folheto de Dados e Informações para o Usuário](#)

MADE IN

USA

LARYNGECTOMY TUBE

Laryngektomietubus
 Tubo de Laringectomía
 Tube de laryngectomie
 Tubo per laringectomia
 Laryngectomiebuísje
 Tubo de laringectomia

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BLOM-SINGER® LARYNGECTOMY TUBE

The Blom-Singer Laryngectomy Tube is a durable, lightweight, silicone tube that is used to keep the tracheostoma open after laryngectomy. The Laryngectomy Tube is available in a variety of sizes and flexibly contours to the postoperative tracheal anatomy. Along the posterior (backside) it features a recessed area with smooth, rounded edges (diagram 1) that can be modified (fenestrated) for use with a Blom-Singer Voice Prosthesis. A pre-fenestrated series (diagram 2) is also available.

This product is intended for single patient use only. Not for use between patients.

INDICATIONS

The Blom-Singer Laryngectomy Tube is designed to maintain tracheostoma patency after laryngectomy to avoid tracheostomal closure.

CONTRAINDICATIONS

The Blom-Singer Laryngectomy Tube is a medical product that should be used only by professionals with experience and training in its use and care. Prolonged infection, especially when resistant to therapy, may contraindicate use of this tube. Radiation therapy and resultant inflammation, bleeding, discharge, or coughing is an additional contraindication.

HOW SUPPLIED

The Blom-Singer Laryngectomy Tube is available in three product offerings: **Sterile** and **Nonsterile** Laryngectomy Tube (non-fenestrated), and a **Nonsterile Fenestrated** Laryngectomy Tube. A laryngectomy tube holder neck strap is available separately.

INSTRUCTIONS FOR USE

The user's hands and the laryngectomy tube should always be thoroughly cleaned to help avoid introducing contaminants into the tracheostoma and respiratory tract. Removal and reinsertion of the device should only be done with the user positioned in front of a mirror with a bright light focused directly on the stoma.

If the package has been opened and the laryngectomy tube has been handled it should be thoroughly cleaned (or resterilized if deemed necessary) adhering to procedures recommended in the Sterilization section of this data sheet.

Please refer to the diagrams at the front of this data sheet.

The following procedure is recommended by Eric D. Blom, Ph.D.

1. Select the largest diameter and length laryngectomy tube that is comfortably tolerated by the patient and permits easy breathing. Insert the smaller tip end of the tube completely into the tracheostoma until its larger funnel shaped end rests against the entrance to the tracheostoma.
2. A very small amount of medical grade, water-soluble lubricant applied to the rim of the tracheostoma may ease insertion of the laryngectomy tube. In some instances it may be necessary to slightly fold or compress the tip portion of the laryngectomy tube, forming a taper, to facilitate introduction into the tracheostoma. The tube is curved and should conform comfortably to the post-laryngectomy contour of the trachea.
3. The laryngectomy tube may be secured around the neck using a laryngectomy tube holder (sold separately). The straps are inserted into the openings in the side of the laryngectomy tube so that the tabs on the straps point away from the tube. The straps are passed around the back of the neck and adjusted so that the straps are secure. The straps should be loose enough to permit one finger breadth between the backside of the straps and the surface of the

neck. After adjusting the straps to a comfortable tension, the user may choose to cut off some of the excess strap length.

FENESTRATING INSTRUCTIONS

A Blom-Singer Laryngectomy Tube used simultaneously with a voice prosthesis usually will require fenestration of the posterior wall to facilitate unrestricted airflow for tracheoesophageal voice. A pre-fenestrated series is available for this purpose. On non-prefenestrated laryngectomy tubes, a recessed area has been molded into the laryngectomy tube in this area (diagram 1) to facilitate secondary modification by a physician or other experienced and trained clinician. The tube should not be in place during the fenestration procedure. Do not make the fenestration with razors or scissors. Tearing or other damage to the device can result if cutting is done improperly. The laryngectomy tube should **never** be structurally modified by the user (patient).

The following steps are recommended by Eric D. Blom, Ph.D., for fenestrating a Blom-Singer Laryngectomy Tube:

1. Standard aseptic techniques should be observed when handling the laryngectomy tube. With the voice prosthesis positioned correctly in the tracheoesophageal puncture, carefully insert the appropriate diameter and length laryngectomy tube following the procedures recommended in the Instructions for Use in this data sheet.
2. With the laryngectomy tube fully inserted, visualize through the transparent posterior wall of the laryngectomy tube the point of airflow entry to the voice prosthesis. In most patients this will be within the recessed area on the tube. Using a pen, carefully mark that point against the interior surface of the laryngectomy tube.
3. Carefully remove the laryngectomy tube from the tracheostoma. Using a Fenestration Punch that has been autoclave sterilized or thoroughly disinfected with an antiseptic prior to use, make a hole corresponding to the ink mark. Only a 1/4" punch should be used. The punch is positioned by sliding it down over the funnel-shaped

end of the laryngectomy tube. After the fenestration is placed, remove the punched disc completely. Be certain that all particles of silicone cut by the punch are removed from the tube.

4. Thoroughly rinse the laryngectomy tube with water and wash with an antiseptic or autoclave sterilize. It is recommended that the laryngectomy tube be reinserted wet. **Never** use facial or toilet tissue to dry a silicone laryngectomy tube as particles of lint or fabric may be deposited onto the laryngectomy tube that can be inhaled into the airway upon reinsertion. If a laryngectomy tube must be dried, use only a lint-free material and rinse the laryngectomy tube again in water before reinsertion. Carefully inspect for loose particles and structural integrity. **Never** use a laryngectomy tube that has tears or other potentially hazardous structural weaknesses.

5. Reinsert the clean, disinfected, or resterilized laryngectomy tube into the tracheostoma according to the Instructions for Use and/or Sterilization Instructions recommended in this data sheet. Observe for adequate fenestration-to-voice prosthesis alignment, and assess the patient's capacity to phonate. If alignment is inadequate, follow the fenestration procedure previously described and enlarge the size of the punch hole as required. Clean and inspect again as described above.

PATIENT CLEANING INSTRUCTIONS

The user's hands should be thoroughly washed with an antiseptic before handling the laryngectomy tube. Remove the laryngectomy tube from neck and rinse it clean with tap water. Wash the laryngectomy tube with an antiseptic and rinse again.

Prior to reinsertion, inspect the laryngectomy tube for structural damage, such as cracks or tears, resulting from prolonged use. Do not use a tube that is damaged.

It is recommended that the laryngectomy tube be reinserted wet. **Never** use facial or toilet tissue to dry a silicone laryngectomy tube as particles of lint or fabric may be deposited onto the laryngectomy

tube that can be inhaled into the airway upon reinsertion. If a laryngectomy tube must be dried, use only a lint-free material and rinse the laryngectomy tube again in water before reinsertion.

Some users experience deposits of bacteria or yeast (*Candida albicans*) on the surface of the laryngectomy tube indicating the need for more effective device disinfection. These deposits appear as patches. In order to prevent these yeast deposits, and to routinely disinfect the laryngectomy tube, alternately wearing and soaking two devices is recommended.

After removing and cleaning the tube as described above, place it in a small container of fresh 3% hydrogen peroxide overnight. After soaking, the tube should be thoroughly rinsed again with tap water, air dried completely, and then placed in a clean resealable plastic bag until its next use. **Note:** Soaking in hydrogen peroxide will not remove yeast that is already present on the laryngectomy tube. It is a preventive measure only. If yeast is present, the laryngectomy tube should be replaced.

Never insert any cleaning device or punch into the laryngectomy tube while the tube is positioned in the neck. Always clean the device *after* it has been removed.

Do not use solvents or petroleum-based products for cleaning or lubricating the device. These materials may damage silicone.

INSTITUTIONAL STERILIZATION INSTRUCTIONS

Clean the laryngectomy tube thoroughly with an antiseptic and water, and wrap in hospital grade muslin or other suitable lint-free autoclave wrapping material. Steam sterilize at 121° C (250° F) for 20 minutes in a standard autoclave unit. Do not sterilize the laryngectomy tube in the original package. Do not autoclave the neck strap. It should be cleaned with a mild detergent and water. **Do Not Use An Ethylene Oxide Sterilization Cycle.**

WARNINGS AND PRECAUTIONS

Never use facial or toilet tissue to dry a silicone laryngectomy tube as particles of lint or fabric may be deposited onto the laryngectomy tube that can be inhaled into the airway upon reinsertion.

Never insert any cleaning device or punch into the laryngectomy tube while the tube is positioned in the neck. Always clean the device after it has been removed.

Never use a laryngectomy tube that has tears or other potentially hazardous structural weaknesses.

The tube should not be in place during the fenestration procedure. Do not make the fenestration with razors or scissors. Tearing or other damage to the device can result if cutting is done improperly.

Tube modifications should not be attempted by the patient.

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician or licensed speech pathologist.

COMPLICATIONS

The following complications have been identified to occur with the use of this device, but are infrequent. They include: Stoma contamination or infection, which may require removal of the laryngectomy tube and/or appropriate antibiotics; Tissue granulation around the stoma; Accidental inhalation of an inappropriately sized laryngectomy tube into the airway, which may require removal by a physician; Obstruction of the tube by mucus necessitating removal and cleaning; Structural failure from excessive usage and/or cleaning.

ORDERING INFORMATION

USA

Blom-Singer products may be ordered directly from International Healthcare Technologies (InHealth) at the following toll free number: (800) 477-5969, Monday — Thursday, 9:30 am—7:00 pm; and Friday, 9:30 am—5:00 pm, Eastern Standard Time. The toll free

number may be used from the continental USA, Alaska, Hawaii, Puerto Rico, and the Virgin Islands. Products also may be ordered at the InHealth website: www.inhealth.com

Canada

In Canada, please contact Auto Control Medical at (800) 461-0991, or Fax (800) 699-5936.

International

For a list of international distributors, please contact: International Healthcare Technologies, 1110 Mark Avenue, Carpinteria, CA 93013-2918, USA. Tel: +01(805) 684-9337, Fax: +01(805) 684-8594.

Consumer Affairs

If you have any questions or dissatisfaction with a product please contact customer service via telephone, fax, post, or e-mail. Tel: (800) 477-5969; Fax: (888) 371-1530; Post: 1110 Mark Avenue, Carpinteria, CA 93013-2918 USA; E-mail: productcomplaints@inhealth.com

RETURNED GOODS POLICY

All returned merchandise must have a Return Merchandise Authorization number (RMA), and must be unopened and undamaged. RMA numbers are obtained from InHealth Customer Service. Returns without an RMA will not be accepted. Products must be returned in unopened packages, with manufacturer's tamper evident seals intact to be accepted for replacement or credit. Products will not be accepted for replacement or credit if they have been in possession of the customer for more than 90 days. Special order products are nonreturnable. There is a 20% restocking fee on all returned merchandise.

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BIBLIOGRAPHY, BIBLIOGRAPHIE, BIBLIOGRAFIA, BIBLIOGRAFIE

1. Blom, E.D.: "Some Comments on the Escalation of Tracheoesophageal Voice Prosthesis Dimensions," **Arch Otolaryngol**, 129: 500-502, 2003.
2. Eerenstein, S.E., Grolman, W., Schouwenburg, P.F.: "Downsizing of Voice Prosthesis Diameter in Patients with Laryngectomy," **Arch Otolaryngol**, 128: 838-841, 2003.
3. Hamaker, R.C., Blom, E.D.: "Botulinum Neurotoxin for Pharyngeal Constrictor Muscle Spasm in Tracheoesophageal Voice Restoration," **Laryngoscope**, 113:1479-1482, 2003.
4. Evans, P.H.R., Blom, E.D.: "Functional Restoration of Speech," **Principles and Practice of Head and Neck Oncology**, Evans, P.H.R., Montgomery, P.W., Gullane, P.J. (Eds.) Taylor and Francis Group, London, 2003.
5. Blom, E.D., Singer M.I., and Hamaker, R.C., **Tracheoesophageal Voice Restoration Following Total Laryngectomy**. San Diego: Singular Publishing Group, Inc.1998.
6. Blom, E.D., Hamaker, R.C.: "Tracheoesophageal Voice Restoration Following Total Laryngectomy," **Cancer of the Head and Neck**, E.N. and Suen, J. (Eds). W.B. Saunders Publishers. 839-852, 1996.
7. Blom, E.D., "Tracheoesophageal Speech," **Seminars in Speech and Language**. New York: Thieme Medical Publishers. 16:3, 191-204, 1995.
8. Pauloski, B.R., Blom, E.D., Logemann, J.A., Hamaker, R.C. "Functional Outcome After Surgery for Prevention of Pharyngospasms in Tracheoesophageal Speakers. Part II: Swallow Characteristics," **Laryngoscope**, 105:10: 1104-1110, 1995.
9. Blom, E.D., Pauloski, B.R., Hamaker, R.C. "Functional Outcome After Surgery for Prevention of Pharyngospasms in Tracheoesophageal Speakers. Part I: Speech Characteristics," **Laryngoscope**, 105:10: 1093-1103, 1995.
10. Blom, E.D., Hamaker, R.C., Freeman, S.B.: "Postlaryngectomy Voice Restoration," Lucente, F.E. (Ed.): **Highlights of the Instructional Courses**, Vol. 7. St. Louis: Mosby-Year Book. pp 3-10, 1994.
11. Blom, E.D., Singer M.I., and Hamaker, R.C. "A Prospective Study of Tracheoesophageal Speech," **Arch Otolaryngol**, 112: 440-447, 1986.
12. Blom, E.D., Singer, M.I. "Disinfection of Silicone Voice Prostheses," **Arch Otolaryngol**, 112: 1986.
13. Blom, E.D., Singer M.I., and Hamaker, R.C. "An Improved Esophageal Insufflation Test," **Arch Otolaryngol**, 111: 211-212, 1985.
14. Hamaker, R.C., Singer M.I., Blom E.D. and Daniels, H.A. "Primary Voice Restoration at Laryngectomy," **Arch Otolaryngol**, 111: 182-186, 1985.
15. Singer, M.I. and Blom, E.D. "A Selective Myotomy for Voice Restoration After Total Laryngectomy," **Arch Otolaryngol**, 107: 670-673, 1981.
16. Singer, M.I. and Blom, E.D. "An Endoscopic Technique for Restoration of Voice After Laryngectomy," **Ann Otol Rhinol Laryngol**, 89:6: 529-533, 1980.

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Manufacturer

Hersteller/Fabricante/Fabricant/
Produttore/Fabrikant

International Healthcare Technologies

A Division of Helix Medical LLC
1110 Mark Avenue
Carpinteria, CA 93013-2918 USA
www.inhealth.com

European Representative

Europäische Vertretung/Representante
para Europa/Représentant européen/
Rappresentante europeo/Europese
vertegenwoordiger/Representante Europeu

EMERGO EUROPE

Molenstraat 15
2513 BH, The Hague, The Netherlands

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Blom-Singer is a registered trademark of Hansa Medical Products, Inc., InHealth is a registered trademark of International Healthcare Technologies, A Division of Helix Medical, LLC, © 2007 International Healthcare Technologies / Blom-Singer ist ein eingetragenes Warenzeichen der Hansa Medical Products Inc., InHealth ist ein eingetragenes Warenzeichen von International Healthcare Technologies, ein Unternehmen der Helix Medical, LLC / Blom-Singer es una marca comercial registrada de Hansa Medical Products, Inc., InHealth es una marca comercial registrada de International Healthcare Technologies, División de Helix Medical, LLC / Blom-Singer est une marque déposée de Hansa Medical Products, Inc., InHealth est une marque déposée de International Healthcare Technologies, filiale de Helix Medical, LLC / Blom-Singer è un marchio depositato di Hansa Medical Products, Inc., InHealth è un marchio depositato di International Healthcare Technologies, una divisione di Helix Medical, LLC / Blom-Singer is een gedeponeerd handelsmerk van Hansa Medical Products, Inc., InHealth is een gedeponeerd handelsmerk van International Healthcare Technologies, een divisie van Helix Medical, LLC / Blom-Singer é marca registrada da Hansa Medical Products, Inc., InHealth é marca registrada da International Healthcare Technologies, uma Divisão da Helix Medical, LLC.