Blom-Singer®

voice restoration systems

37-469-01 Rev. A

Manometer (with manometer adaptor)	4	ENG
Manometer (mit Manometeradapter)	8	DEU
Manómetro (con adaptador de manómetro)	13	ESP
Manomètre (avec adaptateur de manomètre)	18	
Manometro (con adattatore per manometro)	23	ITA
Manometer (met manometer adapter)	28	NED
Manómetro (com adaptador de manómetro)	33	POR


INHEALTH®
T E C H N O L O G I E S

User Information Data Sheet / Informationsblatt für Benutzer / Hoja de información del Usuario / Fiche Technique d'Information Utilisateur / Scheda informativa per il paziente / Informatieblad voor gebruikers / Folheto Informativo do Utente

1


Manometer

Manómetro / Manomètre / Manometro

2

Tracheostoma connector
Tracheostomaverbindung / Conector del traqueostoma / Connecteur de la trachéotomie / Connettore per stoma tracheale / Tracheostoma aansluiting / Conector de traqueostomia

INHEALTH® MANOMETER

(with manometer adaptor)

The InHealth manometer (with manometer adaptor) is designed as a training aid for clinicians to teach new and current patients how to identify and self-regulate their own respiratory pressure requirements during tracheoesophageal speech. Proper airflow control results in better voice quality, less fatigue, and longer duration of the tracheostoma adhesive housing seal.

INDICATIONS

The InHealth manometer (with manometer adaptor) is for use by qualified clinicians as a training aid to teach laryngectomized patients how to recognize and self-regulate their intratracheal pressure during tracheoesophageal speech.

CONTRAINDICATIONS

The InHealth manometer (with manometer adaptor) is only indicated to assess and monitor intratracheal speech airflow pressure produced by laryngectomized patients who speak using a tracheoesophageal voice prosthesis.

HOW SUPPLIED

The InHealth manometer (with manometer adaptor) is supplied non-sterile in a kit. Replacement manometer adaptors are available separately. The manometer is reusable, however, the manometer adaptor is disposable and intended for **single patient use only**. **The manometer adaptor must not be reused by more than one patient to prevent cross contamination.**

REF	DESCRIPTION
BE8000	Manometer Kit <i>includes 1 manometer, 1 manometer adaptor and 1 TruSeal adhesive housing</i>
BE8020	Replacement Manometer Adaptor Pack <i>includes 5 manometer adaptors and 5 TruSeal adhesive housings</i>

INSTRUCTIONS FOR USE

Manometer Adaptor Attachment

Please refer to the diagrams located at the front of this data sheet. The clinician's hands should always be thoroughly cleaned to help avoid introducing contaminants.

The manometer adaptor should be connected to the manometer before attaching the manometer adaptor to the tracheostoma adhesive housing.

1. Attach the manometer adaptor to the manometer by inserting the female luer lock connection of the manometer adaptor into the male luer lock connection of the manometer (diagram 1). Twist the two components together until securely connected.
2. Attach the manometer adaptor to the patient by inserting the tracheostoma connector end to the tracheostoma adhesive housing (diagram 2).
3. When the manometer adaptor is securely attached to the patient, and the manometer securely connected to the manometer adaptor, the clinician may proceed with the following assessment of intratracheal speech pressure.

Intratracheal Speech Pressure Assessment Procedure

The following procedure is provided by Eric D. Blom, Ph.D.

Failure to maintain the tracheostoma adhesive housing seal for a practical length of time is frequently the result of one or more of the following factors: (a) excessive intratracheal speech pressure; (b) excessive phlegm discharge; (c) ineffective valve seal application; or, (d) irregular tracheostomal architecture.

Excessive intratracheal speech pressure is a direct expression of resistance to airflow through a valved voice prosthesis, a tracheoesophageal vocal tract, or a combination of both. This "backpressure" pushes the tracheostoma adhesive housing seal loose from within, causing an air leak that can be durably repaired only by replacing the seal. Before initiation of therapy to teach a patient to attach and use a tracheostoma valve, a pressure measurement is recommended.

With the manometer adaptor inserted into the tracheostoma adhesive housing, which has been securely taped over the stoma, the patient takes a breath, lightly occludes the adaptor with a finger and counts to 15 at a normal conversational speech loudness level. Optimal levels are designated by color markings on the manometer: the green range indicates a pressure reading of 10 – 25 cm H₂O, which is excellent, and the blue range indicates a pressure reading of 26 – 40 cm H₂O, which is good. (See chart below.) Higher levels contribute to proportionally shorter duration of the tracheostoma adhesive housing seal.

PRESSURE RANGES (CM OF H₂O)

10 — 25 cm H ₂ O	Excellent – Green
26 — 40 cm H ₂ O	Good – Blue
41 — 120 cm H ₂ O	High

For some patients, simply monitoring the manometer, as they vary their loudness from softest (usually registering 10 – 15 cm H₂O) to loudest (often measuring 90cm H₂O or more), teaches the importance of controlling their voice loudness. Many tracheoesophageal speakers attempt to speak too loudly resulting in excessive intratracheal “backpressure,” deteriorating voice quality, and shortened duration of the tracheostoma adhesive housing seal.

WARNINGS AND PRECAUTIONS

If there are tears, cracks, or structural damage, do not continue to use the product. The manometer adaptor is disposable and intended for single patient use only. The manometer adaptor must not be reused by more than one patient to prevent cross contamination.

The InHealth manometer is pre-calibrated at the factory. Tampering with the instrument, e.g., unscrewing the housing, may result in malfunction, inaccurate pressure readings and/or leaks of the gauge. Tampering with the connection fittings on the manometer adaptor may result in inaccurate pressure readings.

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician or licensed speech pathologist.

SPECIAL ORDER PRODUCTS

If this data sheet accompanies a Special Order Product there will possibly be differences in the physical characteristics between the product enclosed and the product descriptions in this data sheet. These differences will not affect the safety or efficacy of the special order product. Special order products are nonreturnable.

ORDERING INFORMATION

USA & Canada

InHealth products can be ordered directly from International Healthcare Technologies at the following toll free number: (800) 477-5969, Monday—Thursday, 6:30 am—4:00 pm; and Friday, 6:30 am—2:00 pm, Pacific Time. (The toll free number can be used from the continental USA, Alaska, Hawaii, Puerto Rico, and the Virgin Islands. The toll free number for Canada is (800) 461-0991.) Products also may be ordered at the InHealth website: www.inhealth.com

If you have any questions or dissatisfaction with the product please contact customer service via telephone, fax, post, or E-mail. Tel: (800) 477-5969; Fax: (888)371-1530; Post: 1110 Mark Avenue, Carpinteria, CA 93013-2918 USA; E-mail: order@inhealth.com

International

For a list of international distributors, please contact International Healthcare Technologies, 1110 Mark Avenue, Carpinteria, CA 93013-2918, USA. Tel: +01(805) 684-9337, Fax: +01(805) 684-8594.

RETURNED GOODS POLICY

All returned merchandise must have a Return Merchandise Authorization number (RMA), and must be unopened and undamaged. RMA numbers are obtained from InHealth Customer Service. Returns without an RMA will not be accepted. Products must be returned in unopened packages, with manufacturer’s tamper evident seals intact to be accepted for replacement or credit. Products will not be accepted for replacement or credit if they have been in possession of the customer for more than 90 days. Special order products are nonreturnable. There is a 20% restocking fee on all returned merchandise.

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website da InHealth: www.inhealth.com

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BIBLIOGRAPHY

1. Blom, E.D.: "Tracheoesophageal Valve Fitting and Instruction." In E.D. Blom, M.I. Singer and R.C. Hamaker (Eds.), **Tracheoesophageal Voice Restoration Following Total Laryngectomy** (pp. 103-108). San Diego: Singular Publishing Group, 1998.

InHealth is a registered trademark of International Healthcare Technologies, A Division of Helix Medical, Inc., © 2003 International Healthcare Technologies / InHealth ist ein eingetragenes Warenzeichen von International Healthcare Technologies, ein Unternehmen der Helix Medical, Inc. / InHealth es una marca comercial registrada de International Healthcare Technologies, División de Helix Medical, Inc. / InHealth est une marque déposée de International Healthcare Technologies, filiale de Helix Medical, Inc. / InHealth è un marchio depositato di International Healthcare Technologies, una divisione di Helix Medical, Inc. / InHealth is een gedeponeerd handelsmerk van International Healthcare Technologies, een divisie van Helix Medical, Inc. / InHealth é uma marca registada da International Healthcare Technologies, uma Divisão da Helix Medical, Inc.

The CE Mark is only applicable to the Manometer Adaptor. The absence of the CE Mark for the Manometer is in conformance with PED 97/23/EC.

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A Division of Helix Medical Inc.