

INHEALTH
TECHNOLOGIES®

INTRODUCING

The Blom-Singer® Electrolarynx EL 2000

Innovation in Every Word

ORDER ONLINE: WWW.INHEALTH.COM



Features:

- Activator Button (red)
- Pitch/volume rocker signs (+) (-)
- Built-in battery
- Up to 30 hrs speaking time
- ON/OFF Switch
- Magnetic charging USB port
- Intraoral accessory (EL2002 Oral Adapter*)

*Sold separately



Quick Tips:

- Make sure that the electrolarynx head is in full contact with your skin to ensure adequate sound transmission. Otherwise, you may hear an unpleasant buzzing noise.
- Gain the attention of your listener before you speak. If possible, face them when speaking. Seeing your lips when you speak can help your listener understand you better.
- Speak slowly and overarticulate your speech by exaggerating your mouth movements.
- Find the best placement for the electrolarynx. Try the cheek as well as the neck, and try placing the device at different angles.
- Take natural pauses in conversation. Push and hold the “on” button while you are speaking, release during pauses.
- Practice! The more you practice, the easier this type of communication will be for you and your listeners.