INHEALTH

TECHNOLOGIES®

Blom-Singer® Tracheoesophageal Puncture Set

MEDICAL PROFESSIONAL Instructions For Use

Blom-Singer®

voice restoration systems

R1

37729-01G

Blom-Singer® and **InHealth Technologies®** are registered trademarks of **Freudenberg Medical, LLC**.

Patent(s): www.inhealthpatents.com

IFU: <https://inhealth.com/ifu>

Copyright © 2020 Freudenberg Medical, LLC. All rights reserved.

TABLE OF CONTENTS

Inhaltsverzeichnis / Índice / Sommaire / Sommario / Inhoud / Obsah / Spis treści

| | | |
|-----|--|-----------|
| ENG | Blom-Singer Tracheoesophageal Puncture Set | 5 |
| DEU | Blom-Singer tracheoösophageale Punktionsset | 10 |
| SPA | Equipo para punción traqueoesofágica Blom-Singer | 16 |
| FRA | Kit pour ponction trachéo-œsophagienne Blom-Singer | 22 |
| ITA | Kit per puntura tracheoesofagea Blom-Singer | 28 |
| NLD | Blom-Singer tracheo-oesofageale punctieset | 34 |
| POR | Kit para punção traqueoesofágica Blom-Singer | 40 |
| CES | Sada pro tracheoezofageální punkci | 46 |
| POL | Zestaw do nakłucia tchawiczo-przelykowego Blom-Singer | 51 |

Blom-Singer®

voice restoration systems

MEDICAL PROFESSIONAL Instructions For Use
Gebrauchsanweisung für MEDIZINISCHE FACHKRÄFTE
Instrucciones de uso para PROFESIONALES MÉDICOS
Instructions d'utilisation destinées aux PROFESSIONNELS DE SANTÉ
Istruzioni per l'uso destinate ALL'OPERATORE SANITARIO
Gebruiksaanwijzingen voor MEDISCHE PROFESSIONALS
Instruções de uso para PROFISSIONAIS MÉDICOS
Návod k použití pro ZDRAVOTNICKÝ PERSONÁL
Instrukcja użytkowania dla PERSONELU MEDYCZNEGO

BLOM-SINGER® TRACHEOESOPHAGEAL PUNCTURE SET

Introduction

The Blom-Singer Tracheoesophageal Puncture Set is designed to be used during surgery for placement of a tracheoesophageal puncture.

This Blom-Singer medical device may be used only once by one patient. It may not be reused. See Warnings, Precautions and Complications.

Do not use the device if it becomes contaminated; discard and replace it with a newly packaged product.

INDICATIONS

The Blom-Singer Tracheoesophageal Puncture Set is used to puncture and stent the tracheoesophageal septum in preparation for placement of a valved voice prosthesis.

CONTRAINDICATIONS

Tracheoesophageal puncture may be contraindicated in patients whose medical status increases the risk of uncontrolled dilation of the puncture and aspiration of esophageal contents around the catheter or voice prosthesis. These may include: irradiation to the tracheostoma exceeding 6500 rads; uncontrolled diabetes mellitus; concurrent chemotherapy.

Tracheoesophageal puncture may also be contraindicated in the small tracheostoma in which presence of the catheter or a voice prosthesis may compromise respiration.

HOW SUPPLIED

The Blom-Singer Tracheoesophageal Puncture Set is provided **sterile**, in a double-wrapped tray ready to be introduced into the surgical field. **Do not resterilize.**

PRODUCT DESCRIPTION

The Blom-Singer Tracheoesophageal Puncture Set is designed to be used during surgery for placement of a tracheoesophageal

puncture. This set includes: (1) **Introducer Needle**, a disposable 16 gauge curved needle fitted with a 14 gauge Teflon® sheath; (2) **Leader Cable and Catheter**, a 14 Fr., 26" silicone catheter with tapered tip attached to a 24" Teflon-coated stainless-steel leader cable; (3) **Drape**, a fenestrated surgical drape with adhesive around the aperture.

WARNINGS AND PRECAUTIONS

Hazards and risk factors associated with the Blom-Singer Tracheoesophageal Puncture Set include patient candidacy and puncture complications.

Recommended use of this product is for trained medical physicians only.

This procedure is only recommended with use of an esophagoscope to allow for adequate visualization of the tracheoesophageal wall and to provide protection of the posterior esophagus.

The trained medical doctor must puncture into the esophagoscope to prevent puncturing into the esophagus which could result in dysphagia, bleeding, and/or spinal cord complications.

If the catheter is not attached properly to the neck, the catheter may become dislodged from the puncture and no longer act as a stenting device.

If the stent is not replaced immediately, the patient's puncture tract may close and another puncture operation would be warranted.

Patient candidacy for use of the device and procedure should be determined by a trained medical physician.

Once the puncture is created using the device and the catheter is placed, the individual should be monitored closely and the individual should be fitted with the appropriate size voice prosthesis 72 hours post-operatively.

If there are tears, cracks, or structural damage to the device, discontinue use.

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician or a licensed speech pathologist.

COMPLICATIONS

Although rare, the following complications have been identified to occur with tracheoesophageal puncture. They include:

- aberrant perforation of the posterior esophagus with subsequent mediastinal infection or abscess;
- peristomal infection or cellulitis with catheter usage, which may be prevented by appropriate use of prophylactic antibiotics;
- aspiration around the catheter that is used to stent the tracheoesophageal puncture.

INSTRUCTIONS FOR USE

The following procedural instructions are provided by Daniel Deschler, M.D. Please refer to the diagrams located at the back of this instruction manual.

The Blom-Singer Tracheoesophageal Puncture Set should only be used with careful tissue handling techniques and endoscopy with a rigid esophagoscope as described below. Placement of an esophagoscope stabilizes the esophageal mucosa against the posterior tracheal wall and facilitates easier establishment of the puncture with the advantage of direct visual monitoring. Puncturing into the esophagoscope also shields adjacent structures and minimizes the risk of misdirected perforation.

- 1.** Under anesthesia, introduce a 7x30 cm, or appropriately sized, rigid esophagoscope to the level of the tracheostoma (diagram 1). Orient its tip to align the short aspect of the bevel (or the fenestration on the instrument if it has been modified) adjacent to the tracheal wall immediately inferior to the mucocutaneous junction of the tracheostoma. Palpate or transilluminate the membranous trachea to confirm endoscope positioning.
- 2.** The puncture should be placed midline 5 mm from the mucocutaneous border of the tracheostoma. Introduce the sheathed needle through the posterior tracheoesophageal wall into the esophagoscope with simultaneous monitoring by the endoscopist (diagram 2).
- 3.** Stabilize the sheath in the puncture and withdraw the needle (diagram 3).
- 4.** Thread the Teflon coated leader cable (attached to the tapered

white silicone catheter) into the sheath until it is retrieved in the proximal end of the esophagoscope (diagram 4).

5. Prior to advancing the catheter into the puncture, moisten with saline to facilitate atraumatic insertion. Withdraw the catheter while simultaneously trailing the esophagoscope until the tip of the catheter exists the mouth (diagram 5). Cut the catheter at the black dots with a pair of scissors.

6. Withdraw the catheter at the tracheostoma until the previously cut tip resides at the level of the base of the tongue. Reintroduce the esophagoscope or a laryngoscope and direct the tip of the catheter into the distal esophagus (diagram 6). Suture the catheter to the peristomal skin immediately lateral to the tracheostoma to prevent accidental dislodgment. The catheter is now directed down the esophagus.

SPECIAL ORDER PRODUCTS

If this instruction manual accompanies a Special Order Product, there may be differences in the physical characteristics between the product and the product descriptions described in this instruction manual. These differences will not affect the safety or desired beneficial effect of the special order product. Special Order products are nonreturnable.

ORDERING INFORMATION

USA

Blom-Singer products may be ordered directly from InHealth Technologies. TELEPHONE: Toll-Free (800)477-5969 or (805)684-9337, Monday–Friday, 9:30am–7:00pm, Eastern Standard Time. FAX: Toll-Free (888)371-1530 or (805)684-8594. EMAIL: order@inhealth.com ORDER ON-LINE: www.inhealth.com POST: InHealth Technologies, 1110 Mark Avenue, Carpinteria, CA 93013-2918, USA, Attention: Customer Service

International

Please contact our customer service department for a distributor referral.

Consumer Affairs

If you have any questions or dissatisfaction with a product please

contact our customer service department via telephone, fax, post, or Email: productcomplaints@inhealth.com

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




Additional references available upon request.

Weitere Referenzen aus der Literatur sind auf Anfrage erhältlich. / Se pueden solicitar más referencias. / Autres références disponibles sur demande. / Riferimenti aggiuntivi sono disponibili a richiesta. / Meer informatie is op aanvraag beschikbaar. / Referências adicionais mediante solicitação. / Další odkazy na literaturu jsou k dispozici na požádání. / Dodatkowe informacje dostępne na żądanie.

1. Blom, E.D.: "Some Comments on the Escalation of Tracheoesophageal Voice Prosthesis Dimensions," **Arch Otolaryngol**, 129: 500-502, 2003.
2. Eerenstein, S.E., Grolman, W., Schouwenburg, P.F.: "Downsizing of Voice Prosthesis Diameter in Patients with Laryngectomy," **Arch Otolaryngol**, 128: 838-841, 2003.
3. Hamaker, R.C., Blom, E.D.: "Botulinum Neurotoxin for Pharyngeal Constrictor Muscle Spasm in Tracheoesophageal Voice Restoration," **Laryngoscope**, 113:1479-1482, 2003.
4. Evans, P.H.R., Blom, E.D.: "Functional Restoration of Speech," **Principles and Practice of Head and Neck Oncology**, Evans, P.H.R., Montgomery, P.W., Gullane, P.J. (Eds.) Taylor and Francis Group, London, 2003.
5. Blom, E.D., Singer M.I., and Hamaker, R.C., **Tracheoesophageal Voice Restoration Following Total Laryngectomy**. San Diego: Singular Publishing Group, Inc.1998.
6. Blom, E.D., Hamaker, R.C.: "Tracheoesophageal Voice Restoration Following Total Laryngectomy," **Cancer of the Head and Neck**, E.N. and Suen, J. (Eds.) W.B. Saunders Publishers. 839-852, 1996.
7. Blom, E.D., "Tracheoesophageal Speech," **Seminars in Speech and Language**. New York: Thieme Medical Publishers. 16:3, 191-204, 1995.
8. Pauloski, B.R., Blom, E.D., Logemann, J.A., Hamaker, R.C. "Functional Outcome After Surgery for Prevention of Pharyngospasms in Tracheoesophageal Speakers. Part II: Swallow Characteristics," **Laryngoscope**, 105:10: 1104-1110, 1995.
9. Blom, E.D., Pauloski, B.R., Hamaker, R.C. "Functional Outcome After Surgery for Prevention of Pharyngospasms in Tracheoesophageal Speakers. Part I: Speech Characteristics," **Laryngoscope**, 105:10: 1093-1103, 1995.

-
10. Blom, E.D., Hamaker, R.C., Freeman, S.B.: "Postlaryngectomy Voice Restoration," Lucente, F.E. (Ed.): **Highlights of the Instructional Courses**, Vol. 7. St. Louis: Mosby-Year Book. pp 3-10, 1994.
 11. Blom, E.D., Singer M.I., and Hamaker, R.C. "A Prospective Study of Tracheoesophageal Speech," **Arch Otolaryngol**, 112: 440-447, 1986.
 12. Blom, E.D., Singer, M.I. "Disinfection of Silicone Voice Prostheses," **Arch Otolaryngol**, 112: 1986.
 13. Blom, E.D., Singer M.I., and Hamaker, R.C. "An Improved Esophageal Insufflation Test," **Arch Otolaryngol**, 111: 211-212, 1985.
 14. Hamaker, R.C., Singer M.I., Blom E.D. and Daniels, H.A. "Primary Voice Restoration at Laryngectomy," **Arch Otolaryngol**, 111: 182-186, 1985.
 15. Singer, M.I. and Blom, E.D. "A Selective Myotomy for Voice Restoration After Total Laryngectomy," **Arch Otolaryngol**, 107: 670-673, 1981.
 16. Singer, M.I. and Blom, E.D. "An Endoscopic Technique for Restoration of Voice After Laryngectomy," **Ann Otol Rhinol Laryngol**, 89:6: 529-533, 1980.

INHEALTH

TECHNOLOGIES®

UNITED STATES

1110 Mark Avenue
Carpinteria, CA 93013-2918 USA
Toll-Free: 800.477.5969
Toll-Free Fax: 888.371.1530

INTERNATIONAL

Tel: +1.805.684.9337
Fax: +1.805.684.8594
E-mail: order@inhealth.com

ORDER ONLINE

www.inhealth.com

Freudenberg Medical, LLC
1110 Mark Avenue
Carpinteria, CA 93013-2918 USA
www.inhealth.com

EMERGO EUROPE
Westervoortsedijk 60
6827 AT Arnhem
The Netherlands