
INHEALTH

TECHNOLOGIES®

Blom-Singer® Dilation/Sizing System

INSTRUCTIONS FOR USE

Blom-Singer®

voice restoration systems

Blom-Singer® and InHealth Technologies® are registered trademarks of **Freudenberg Medical, LLC**.

IFU: <https://inhealth.com/ifu>

Copyright © 2025 Freudenberg Medical, LLC. All rights reserved.

INSTRUCTIONS FOR USE

使用说明 / 使用說明書 / Udhëzime për përdorim / إرشادات الاستخدام / Инструкции за употреба /
Upute za korištenje / Návod k použití / Brugsanvisning / Gebruiksaanwijzing / Kasutusjuhend /
Käyttöohjeet / Instructions d'utilisation / Gebrauchsanweisung / Οδηγίες χρήσης / הוראות שימוש /
Használati utasítás / Notkunarleiðbeiningar / Istruzioni per l'uso / 사용 지침 / Lietošanas
norādījumi / Naudojimo instrukcijos / Bruksanvisning / Instrukcja użytkowania / Instruções de
uso / Instrucțiuni de utilizare / Инструкции по применению / Uputstvo za upotrebu / Návod na
použitie / Instrucciones de uso / Bruksanvisningar / Kullanım Talimatlar

TABLE OF CONTENTS


目录 / 目次 / Tabela e përmbajtjes / جدول المحتويات / Съдържание / Sadržaj / Obsah /
 Indholdsfortegnelse / Inhoud / Sisukord / Sisällysluettelo / Table des matières /
 Inhaltsverzeichnis / Πίνακας περιεχομένων / תוכן עניינים / Tartalomjegyzék / Efnisyfirlit /
 Sommario / 목차 / Satura rādītājs / Turinys / Innholdsfortegnelse / Spis treści / Índice /
 Cuprins / Содержание / Sadržaj / Obsah / Índice / Innehållsförteckning / İçindekiler

EN	Dilation / Sizing System	8	HU	Tágító/méretező rendszer	92
ZH	扩张/尺寸测量系统	13	IS	Vikkunar- /mælikerfi	97
JA	拡張/サイジングシステム	17	IT	Sistema di dilatazione/ misurazione	102
SQ	Sistemi i dilatimit/matjes	22	KO	확장/크기 측정 시스템	108
AR	نظام التوسيع/تحديد الحجم	28	LV	Dilatatoru/šablonu sistēma	113
BG	Система за дилатация/ оразмеряване	32	LT	Išplėtimo / dydžio matavimo sistema	118
HR	Sustav za dilataciju/mjerenje veličine	38	NO	System for dilatasjon/ størrelsesbestemmelse	123
CS	Systém pro dilataci / měření velikosti	43	PL	System rozszerzania/ wymiarowania	128
DA	Udvidelses-/ dimensioneringssystem	48	PT	Sistema de dilatação/ dimensionamento	134
NL	Dilatatie-/ maatbepalingssysteem	53	RO	Sistem de dilatare/măsurare	140
ET	Dilatatsiooni- / suuruse mõõtmise süsteem	59	RU	Система дилатации/измерения	145
FI	Laajennus-mitoitusjärjestelmä	64	SR	Sistem za dilataciju/merenje	151
FR	Système de dilatation/mesure	69	SK	Dilatačný systém/systém na určenie veľkosti	156
DE	Dilatations-/ Größeneinteilungssystem	75	ES	Sistema de dilatación/medición	162
EL	Σύστημα διαστολής / μέτρησης μεγέθους	81	SV	System för dilatations-/ storleksbestämning	167
HE	מערכת הרחבה/מידת גודל	87	TR	Dilatasyon / Boyutlandırma Sistemi	172

DIAGRAMS

插图 / 図 / Diagramet / الرسم التوضيحي / Схеми / Dijagrami / Schémata / Diagrammer / Afbeeldingen / Diagrammid / Kaaviot / Diagrammes / Abbildungen / Διαγράμματα / תרשימים / Ábrák / Skýringarmyndir / Diagrammi / 그림 / Diagrammas / Diagramos / Diagrammer / Schematy / Diagramas / Diagrame / Схемы / Dijagrami / Grafické znázornenia / Diagramas / Diagram / Diyagramlar


1



2



3







4

5

BLOM-SINGER® DILATION/SIZING SYSTEM**PRODUCT DESCRIPTION**

Please refer to the diagrams located at the front of this instruction manual.

The Blom-Singer® Dilation/Sizing System consists of three separately available, non-sterile products, as shown in the chart below. These products may be cleaned or sterilized and reused per instructions specific to each device in this manual.

Dilation/ Sizing System Product	Clinician Use	Patient Use	Size	Use With
Puncture Dilator	●	●	18 Fr	16 Fr Voice Prostheses
			22 Fr	20 Fr Voice Prostheses
Voice Prosthesis Sizer	●		16 Fr	16 Fr Voice Prostheses
			20 Fr	20 Fr Voice Prostheses
Dilator-Sizer	●		18 Fr	16 Fr Voice Prostheses
			22 Fr	20 Fr Voice Prostheses

Puncture Dilator

The Blom-Singer® Puncture Dilator (**diagram 1**) is available for use by the patient and/or clinician. It is a tapered, solid silicone stent that is inserted into the tracheoesophageal puncture (TEP) upon removal of a voice prosthesis. It is intended to prevent the TEP from closing; prevent leakage of esophageal contents through an open TEP into the trachea during swallowing; and dilate the TEP in preparation for voice prosthesis insertion or replacement. The device can also be used to temporarily stent (up to a maximum of 30 days) the TEP after voice prosthesis removal/dislodgement.

Voice Prosthesis Sizer

The Blom-Singer® Voice Prosthesis Sizer (**diagram 2**) is only available for use by a qualified clinician. It is used to measure the length of the TEP tract to determine the appropriate voice prosthesis length. Incremental 2-millimeter markings and numbers on the device correspond to the lengths of Blom-Singer® Voice Prostheses. The product is packaged with three sizer inserters. Replacement inserters are available separately from InHealth Technologies®.

Dilator-Sizer

The Blom-Singer® Dilator-Sizer (**diagram 3**) is only available for use by a qualified clinician. It is a combination Puncture Dilator and Voice Prosthesis Sizer. The Dilator-Sizer is used to measure the length of the TEP tract to determine the appropriate voice prosthesis length and to dilate the TEP for ease of voice prosthesis insertion or replacement. Incremental 2-millimeter markings and numbers on the device correspond to the lengths of Blom-Singer® Voice Prostheses.

INDICATIONS

The Blom-Singer® Dilation/Sizing System is indicated for dilation, measuring, and maintaining the TEP for tracheoesophageal voice restoration with a voice prosthesis following total laryngectomy. The Blom-Singer® Sizing System is only used by clinicians for measurement of the TEP tract to select the appropriate voice prosthesis length.

CONTRAINDICATIONS

The products of the Blom-Singer® Dilation/Sizing System are medical devices and should only be used by individuals with experience and training in their use and care.

These devices should not be placed in a small tracheostomal airway if presence obstructs breathing.

Placement of the 18 Fr Puncture Dilator or Dilator-Sizer is contraindicated in a tracheoesophageal puncture that has dilated substantially beyond the diameter of a conventional 16 Fr Blom-Singer® voice prosthesis. Likewise, the 22 Fr Puncture Dilator or Dilator-Sizer should not be used in a puncture that has dilated substantially beyond the diameter of a 20 Fr Blom-Singer® voice prosthesis.

WARNINGS AND PRECAUTIONS

The products of the Blom-Singer® Dilation/Sizing System are medically prescribed devices. Initial sizing and instructions in their use must be provided by a qualified clinician trained in the use of these particular devices.

The Blom-Singer® Voice Prosthesis Sizer and Dilator-Sizer are medical devices and should be used only by a clinician trained in their use. The Voice Prosthesis Sizer is not a valved voice prosthesis and cannot function as such.

If there are any tears, cracks, or structural damage, discontinue use. Care should be exercised to avoid inhaling prosthetic devices into the airway. In the unlikely event that this should occur, the patient should attempt to cough the device out of the windpipe and immediately seek medical attention if unsuccessful.

Never clean the device while it is in position in the TEP or tracheostoma. Only clean the device after it has been removed from the TEP and tracheostoma.

Do not use solvents or petroleum-based products for cleaning these devices. These materials may damage the devices.

Do not use the device in the event of the packaging being damaged or unintentionally opened prior to use; discard and replace it with a new device.

COMPLICATIONS

Although rare, the following complications may occur: stoma contamination or infection, which may require removal of the device and/or use of appropriate antibiotics; accidental inhalation of a device into the airway, which may require removal by a physician; inhalation of liquids through a TEP into the airway, which may cause the patient to cough or lead to aspiration; inflammatory reaction around the TEP site and formation of granulation tissue; TEP dilation resulting in leakage of fluids around the device; dislodgment of the device and subsequent closure of the TEP; intractable leakage around the device, requiring surgical revision or closure of the TEP; dysphagia, or difficulty swallowing.

INSTRUCTIONS FOR USE

Please refer to the diagrams located at the front of this instruction manual. The Blom-Singer® Dilation/Sizing device and user's hands should always be clean (gloves recommended for clinicians) to avoid introducing contaminants into the tracheostoma TEP, or esophagus.

Removal and insertion of devices in the TEP should only be done with a clear view of the stoma (e.g. in front of a mirror, with a bright light focused directly into the stoma). Avoid

touching the portion of the devices that will enter the TEP.

Puncture Dilator (Patient/Clinician)

Dilate the TEP with the Puncture Dilator of the appropriate diameter (see sizing chart).

1. Position the tapered tip of the dilator in the TEP.
2. Gently insert the device completely into the neck strap (**diagram 4**), following the angle of the TEP tract: if the tract angles straight back or slightly downward, insert the dilator so its curvature is also directed downward. If the puncture angles upward, insert the dilator so its curvature is also in an upward direction.
3. Secure the dilator near the edge of the stoma by holding in place with clean hands or with medical adhesive tape across the neck strap. Remove the dilator if gagging or discomfort is experienced; users should consult with their clinician if persistent gagging or discomfort is noted.
4. Leave the dilator in place for 5 minutes, or longer if necessary, in order to sufficiently dilate the TEP and allow for insertion of the voice prosthesis/size* (*clinician only).
5. If medical adhesive tape was used, remove the tape from the neck strap. Grasp the dilator securely, and carefully withdraw it from the TEP.
6. Immediately insert a voice prosthesis/size* per Instructions for Use supplied with the selected device (*clinician only).

Dilator-Sizer (Clinician Only)

To dilate and measure the TEP tract, use the Dilator-Sizer of the appropriate diameter (see sizing chart).

Always measure the TEP tract to confirm correct voice prosthesis length before selecting and inserting a voice prosthesis.

1. Gently insert the Dilator-Sizer of the appropriate diameter (see sizing chart) completely into the TEP (**diagram 5**), following the angle of the TEP tract: if the TEP tract angles straight back or slightly downward, insert the Dilator-Sizer so its curvature is also directed downward; if the puncture angles upward, insert the Dilator-Sizer so its curvature is also in an upward direction. Remove the Dilator-Sizer if the patient experiences gagging or discomfort.
2. Carefully withdraw the Dilator-Sizer until gentle resistance is detected, indicating that the esophageal wings are seated against the anterior wall of the esophagus (**diagram 6**). Pulling too firmly on the Dilator-Sizer may compress the TEP tract, resulting in a shorter, erroneous measurement.
3. Note the marker and size number that is closest to the entrance of the TEP and select a Blom-Singer® voice prosthesis of corresponding size. If the measurement is between two markers, select a voice prosthesis of the longer length.
4. Fully re-insert the Dilator-Sizer and secure the neck strap by holding it in place with clean hands or with medical adhesive tape to the patient's neck. Leave the Dilator-Sizer in the TEP to dilate 5 minutes, or longer if necessary, to allow for insertion of the selected voice prosthesis.
5. If medical adhesive tape was used, remove the tape from the neck strap. Grasp the Dilator-Sizer securely, and carefully withdraw it from the TEP.
6. Immediately insert the voice prosthesis per Instructions for Use supplied with the selected voice prosthesis.

Voice Prosthesis Sizer (Clinician Only)

To measure the TEP tract, use the Voice Prosthesis Sizer of the appropriate diameter (see

sizing chart). Use only the inserters packaged with the product (or in replacement packs). Always measure the TEP tract to confirm correct voice prosthesis length before selecting and inserting a voice prosthesis.

1. Place the tip of the inserter into the open end of the Voice Prosthesis Sizer (**diagram 7**). Avoid touching the portion that inserts into the TEP.
2. Attach the neck strap on the Voice Prosthesis Sizer securely over the safety peg on the inserter (**diagram 8**). This should prevent inadvertent dislodgment of the device during insertion into the TEP.
3. Hold the end of the inserter stick and insert the Voice Prosthesis Sizer completely into the TEP until the tracheal collar meets the posterior tracheal wall (**diagram 9**).
4. Carefully withdraw the Voice Prosthesis Sizer until very gentle resistance is detected indicating that the esophageal collar is seated against the anterior wall of the esophagus (**diagram 10a**). Pulling too firmly on the Voice Prosthesis Sizer may compress the TEP tract resulting in a shorter, erroneous measurement.
5. Note the marker and size number that is closest to the entrance of the TEP (**diagram 10b**) and select a Blom-Singer® Voice Prosthesis of corresponding size. If the measurement is between two markers, select a voice prosthesis of the longer length.
6. Carefully remove the Voice Prosthesis Sizer and replace with a Puncture Dilator.
7. Proceed with dilation and immediate voice prosthesis insertion as per Instructions for Use supplied with the selected voice prosthesis.

DEVICE CLEANING: Puncture Dilator (for Single Patient Reuse)

The Puncture Dilator may be cleaned and reused by a single patient. Clean the Puncture Dilator thoroughly with a mild detergent and rinse under a strong flush of tap water. Dry the dilator with sterile gauze. Once a dilator has been cleaned and dried completely, it should be stored in a clean, resealable plastic bag. Always handle with clean hands and avoid touching portions of the dilator that will enter the TEP. Inspect the device routinely for structural damage. If there are tears, cracks, or structural damage, discontinue use and replace device.

INSTITUTION STERILIZATION AND REUSE INSTRUCTIONS

The following devices of the Blom-Singer® Dilation/Sizing System may be cleaned and sterilized for reuse between patients to a limit of **60 uses**: Puncture Dilator, Dilator-Sizer, Voice Prosthesis Sizer (NOT the sizer inserter stick).

Prior to sterilization, clean the device using one of the following procedures:

Manual Cleaning

1. Rinse the device under cold running utility water for 90 seconds. While rinsing, use a soft bristled brush to brush the device for 5 seconds. Use an appropriately sized lumen or channel brush to brush the lumen of the device for 5 seconds, as applicable.
2. Immerse the device in a prepared enzymatic detergent solution for 1 minute (following manufacturer's instruction for dilution). While immersed, use a soft bristled brush to brush the device for 5 seconds. Use an appropriately sized lumen or channel brush to brush the lumen of the device for 5 seconds, as applicable.
3. Remove the device from the bath of detergent solution and rinse under cold running utility water for 1 minute.
4. Dry the device using pressurized air.

Automated Cleaning

1. If utilizing a washer disinfectant for cleaning use the following settings:

Stage	Recirculation Time (Minutes)	Temperature	Detergent Type (if applicable)
Pre-Wash 1	2:00	Cold Tap Water	N/A
Wash 1	1:00	43°C Tap Water (Set Point)	Enzymatic Detergent
Neutralization Wash	2:00	43°C Tap Water (Set Point)	pH neutral Detergent
Rinse 1	1:00	43°C Tap Water (Set Point)	N/A
Pure Water Rinse	1:00	43°C Tap Water (Set Point)	N/A
Dry Time	7:00	115°C	N/A

Note: use recommended detergent concentrations based on specific washer-disinfector and/or detergent manufacturer's recommendations.

To sterilize after cleaning:

Device to be sterilized must be placed in a single sterilization pouch.

Two separate parameters using pre-vacuum steam sterilization autoclave have been validated to a limit of **60 uses**:

- A. Sterilization at 132°C for 4 minutes after performing 4 preconditioning pulses
OR
- B. Sterilization at 134°C for 3 minutes after performing 4 preconditioning pulses

After sterilization cycle is complete:

Device requires a minimum drying time of 20 minutes.

Inspect the device routinely for damage. If there are tears, cracks, or unreadable markings, discontinue use and replace device.

SPECIAL STORAGE AND/OR HANDLING CONDITIONS

There are no special storage and/or handling conditions for the device.

DISPOSAL INSTRUCTIONS

The device is not biodegradable and can be contaminated when used. Carefully dispose of the device per local guidelines.

ORDERING INFORMATION

USA

Blom-Singer® products may be ordered directly from InHealth Technologies®.

Telephone: Toll-Free 800.477.5969 or +1.805.684.9337

Fax: Toll-Free 888.371.1530 or +1.805.684.8594

Email: order@inhealth.com Order Online: www.inhealth.com

Mail: InHealth Technologies, 1110 Mark Avenue, Carpinteria, CA 93013-2918, USA

International

Please contact our customer service department for a distributor referral.

PRODUCT COMPLAINTS/EU SERIOUS INCIDENTS

If you are dissatisfied with the device or have any questions, please contact
productcomplaints@inhealth.com.

Telephone: +1-800-477-5969

Fax: +1-888-371-1530

Any serious incident that has occurred in relation to the device should be reported to Freudenberg Medical, LLC as noted above and the competent authority of the EU Member State in which the user and/or patient is established.

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