Blom-Singer® FirstFit™ Surgical Kit

Instructions for Use

Table of Contents

目录 / 目次 / Tabela e përmbajtjes / جدول المحتويات / Съдържание / Sadržaj / Obsah / Indholdsfortegnelse / Inhoud / Sisukord / Sisällysluettelo / Table des matières / Inhaltsverzeichnis / Πίνακας περιεχομένων / תוכן עניינים / Tartalomjegyzék / Efnisyfirlit / Sommario / 목차 / Saturs / Turinys / Innholdsfortegnelse / Spis treści / Índice / Cuprins / Содержание / Sadržaj / Obsah / Índice / Innehållsförteckning / İçindekiler

EN	FirstFit Surgical Kit	FI	FirstFit-kirurgiasarja	NO	FirstFit kirurgisk sett
ZH	FirstFit 外科手术套件	FR	FirstFit Kit chirurgical	PL	Zestaw chirurgiczny FirstFit
JA	FirstFit 外科用キット	DE	FirstFit Chirurgie-Set	PT	Kit cirúrgico FirstFit
SQ	Kompleti kirurgjikal FirstFit	EL	Χειρουργικό κιτ FirstFit	RO	Trusă chirurgicală FirstFit
AR	FirstFit مجموعة الأدوات الجراحية	HE	ערכת ניתוח FirstFit	RU	Хирургический набор FirstFit
BG	Хирургичен комплект FirstFit	HU	FirstFit sebészeti készlet	SR	FirstFit Hirurški komplet
HR	Kirurški komplet FirstFit	IS	FirstFit skurðaðgerðarsett	SK	Chirurgická súprava FirstFit
CS	Chirurgická souprava FirstFit	IT	Kit chirurgico FirstFit	ES	Kit quirúrgico FirstFit
DA	FirstFit Kirurgisk sæt	KO	FirstFit 수술 키트	SV	FirstFit operationssats
NL	FirstFit chirurgische kit	LV	FirstFit Kirurgiskais komplekts	TR	FirstFit Cerrahi Kiti
ET	FirstFit Kirurgiline komplekt	LT	„FirstFit“ chirurginis rinkinys		

Blom-Singer® and InHealth Technologies® are registered trademarks of **Freudenberg Medical, LLC**.

Patent(s): www.inhealthpatents.com

IFU: <https://inhealth.com/ifu>


Copyright © 2023 Freudenberg Medical, LLC. All rights reserved.


Instructions for Use

使用说明 / 使用說明書 / Udhëzime për përdorim / إرشادات الاستخدام / Инструкции за употреба / Upute za korištenje / Návod k použití / Brugsanvisning / Gebruiksaanwijzing / Kasutusjuhend / Käyttöohjeet / Instructions d'utilisation / Gebrauchsanweisung / Οδηγίες χρήσης / שימוש / הוראות / Használati utasítás / Notkunarleiðbeiningar / Istruzioni per l'uso / 사용 지침 / Lietošanas norādījumi / Naudojimo instrukcijos / Bruksanvisning / Instrukcja użytkowania / Instruções de uso / Instrucțiuni de utilizare / Инструкции по применению / Uputstvo za upotrebu / Návod na použitie / Instrucciones de uso / Bruksanvisningar / KϫYürϫlϫk Tarihi

Diagrams

插图 / 図 / Diagramet / الرسم التوضيحي / Схеми / Dijagrami / Schémata / Diagrammer / Afbeeldingen / Diagrammid / Kaaviot / Diagrammes / Abbildungen / Διαγράμματα / תרשימים / Ábrák / Skýringarmyndir / Diagrammi / 그림 / Diagrammas / Diagramos / Diagrammer / Schematy / Diagramas / Diagrame / Схеми / Crteži / Grafické znázornenia / Diagramas / Diagram / Diyagramlar

ENGLISH

Blom-Singer® FirstFit™ Surgical Kit Instructions for Use (IFU)

Please refer to the diagrams located at the front of this instruction manual.

DESCRIPTION

The Blom-Singer® FirstFit™ Surgical Kit is designed for primary or secondary tracheoesophageal puncture (TEP) with retrograde voice prosthesis insertion. This set includes: one (1) 13-gauge curved cored puncture needle; one (1) 36-inch guidewire with rounded tips; one (1) pharynx protector tool with handle and notched cylinder; one (1) voice prosthesis inserter pre-loaded with a Blom-Singer Indwelling Voice Prosthesis of the diameter and length specified on the package label. The voice prosthesis is made of medical-grade silicone and consists of a one-way silicone flap valve, an esophageal flange, a body that holds the valve assembly, a tracheal flange, and a safety strap. The device and its components are provided sterile (by ethylene oxide) in a sealed tray, ready to be introduced into the surgical field.

Tracheoesophageal puncture diameter selection is typically based on surgeon's practice and preference. Tracheoesophageal puncture and subsequently voice prosthesis length can vary by an individual patient's anatomy and should be considered when selecting which surgical kit to use.

INDICATIONS

The device is indicated for use during surgical creation of primary or secondary tracheoesophageal puncture, dilation of the tracheoesophageal wall, and to guide the placement of the voice prosthesis for tracheoesophageal voice restoration following total laryngectomy.

CONTRAINDICATIONS

Tracheoesophageal puncture and therefore the device is contraindicated:

- in patients whose medical status or anatomical abnormalities increase the risk of uncontrolled dilation of the puncture and aspiration of esophageal contents around the voice prosthesis, including those with significant irradiation to the tracheostoma, uncontrolled diabetes mellitus, concurrent chemotherapy, bleeding disorders, and severe malnutrition;
- in a small tracheostoma where presence of a voice prosthesis may compromise respiration;
- in patients who have severe arthritis or blindness, are not able to care for a voice prosthesis, or are not interested in speaking with a voice prosthesis.

WARNINGS AND PRECAUTIONS

Use of this product is for trained physicians only.

This device and its components are for single patient use only. It may not be reused due to risk for infection. Use of a device that is not sterile may cause infection. Do not resterilize the product. Re-processing and/or resterilization of the device may alter its performance and device failure may occur. Do not use the device if the package is opened or in any way damaged, is expired or becomes contaminated; discard and replace it with an unopened product. If there are tears, cracks, or structural damage to any components of the kit, discontinue use.

A primary TEP procedure with this device is recommended only with use of the included pharynx protector to guard the posterior esophagus. The physician must puncture into the pharynx protector to prevent puncturing the posterior esophagus, which could result in dysphagia, bleeding, and/or spinal complications.

A secondary TEP procedure is recommended only with use of a rigid esophagoscope to allow for adequate visualization of the tracheoesophageal wall and to provide protection of the posterior esopha-

gus. Do not attempt to use the pharynx protector for secondary TEP. The physician must puncture into the esophagoscope to prevent puncturing directly into the esophagus, which could result in dysphagia, bleeding, and/or spinal complications.

In either primary or secondary TEP procedures, the guidewire must be placed and kept in the puncture tract before the needle is removed. If the guidewire is not kept within the puncture tract after the needle is removed, the puncture tract may be lost and another puncture procedure may be warranted. Additionally, failure to follow specific insertion and placement instructions may result in unsuccessful puncture or possible damage to tissue, which may require an additional TEP procedure. Once the puncture is completed using the device and the indwelling voice prosthesis is placed, the individual should be monitored for any complications.

This device and its components are indicated for retrograde voice prosthesis placement only. Anterograde placement would result in improper voice prosthesis function.

Voice Prosthesis

The clinician is to instruct patients on the use and care of the voice prosthesis as appropriate. Please refer to Instructions for Use for Indwelling Voice Prosthesis and Voice Prosthesis Cleaning system for further information.

Changes in the anatomy or medical status of the patient may lead to improper fitting and/or function of the voice prosthesis. Dislodgement or extrusion of the voice prosthesis from the TEP may occur and result in ingestion, aspiration (inhalation) or tissue damage. A foreign body such as a voice prosthesis in the airway may cause complications such as acute respiratory distress and/or respiratory arrest. Prosthesis size/length selection can vary by individual patient's anatomy and should be considered when selecting which surgical kit to use. If the voice prosthesis is improperly sized, it may cause tissue necrosis (tissue death) and/or device dislodgement.

The attached safety strap on an indwelling voice prosthesis should only be removed after the

esophageal flange has been verified to have deployed in the esophagus. Never attempt to insert or reinsert an indwelling voice prosthesis that has the safety strap removed or unsecured from the safety peg.

Voice Prosthesis Leakage

When the flap valve of the voice prosthesis fails to close completely, fluid may pass through the device from the esophagus to the trachea, which may cause coughing or aspiration. Use of a smaller diameter voice prosthesis than the tracheoesophageal puncture size may result in peripheral leakage (leakage around the device). Recurrent leakage of the voice prosthesis should be evaluated by a clinician as leakage could cause aspiration pneumonia. Selection of a different device size, model/option may be indicated. Gentle handling and pressure should always be used when cleaning the voice prosthesis to avoid device damage or dislodgement.

Voice Production

To prevent post-operative complications, the patient should not begin speaking with the voice prosthesis until the clinician has indicated that it is safe to do so. The lumen of the voice prosthesis must be kept clear of blockage for it to function properly to allow the patient to voice.

In some users, the inability to relax the muscles of the throat may account for an inability to speak fluently and with minimal effort. This problem requires professional assessment.

Patients requiring post-operative radiation may have transient interruption of voice in the third or fourth week of treatment. The device can remain in place as determined by the clinician.

Voice Prosthesis Dislodgement

Care must be exercised during device insertion, removal, or use of cleaning devices to avoid injury to the TEP or accidental displacement of the device, which could result in aspiration of the device into the trachea. Should aspiration occur, the patient should attempt to cough the device out of the trachea. Further medical attention may be necessary if coughing the device out is unsuccessful. If the voice prosthesis is dislodged from the TEP post-operatively, a Blom-Singer Puncture Dilator or suitable device of the appropriate diameter should be immediately placed in the puncture to keep it from closing and leaking fluids. A replacement device should be inserted within 24 hours. Foreign objects should not be inserted into the voice prosthesis. Inserting objects other than the Blom-Singer cleaning devices may cause dislodgment and subsequent aspiration or ingestion of the voice prosthesis or its components.

Voice Prosthesis Replacement

The indwelling voice prosthesis is not intended to be inserted or removed by the patient. The indwelling voice prosthesis is not a permanent device and requires replacement periodically. The device may be left in place in the TEP until it has persistent leakage or is not providing adequate voice for speech, requires resizing, or as otherwise indicated by the clinician. Removal of the indwelling voice prosthesis should only be done by grasping the tracheal flange of the device securely with a locking hemostat. Please refer to the Instructions for Use manual packaged with the selected replacement device for further information.

COMPLICATIONS

Although rare, the following complications may occur:

- Aberrant perforation of the posterior esophagus with subsequent mediastinal infection or abscess
- Peristomal infection or cellulitis, which may be prevented by appropriate use of prophylactic antibiotics
- Aspiration around the indwelling voice prosthesis positioned in the tracheoesophageal puncture
- Accidental aspiration of the voice prosthesis into the airway, which may require removal by a physician; puncture dilation resulting in leakage of fluids around the voice prosthesis; inflammatory reaction around the puncture site and formation of granulation tissue; dislodgment of the voice prosthesis and subsequent closure of the TEP; intractable (uncontrollable) leakage around the voice prosthesis, requiring surgical revision or closure of the puncture; dysphagia (difficulty swallowing); microbial growth deposits causing voice prosthesis leakage or valve incompetence; accidental ingestion of the voice prosthesis into the esophagus.

INSTRUCTIONS FOR USE

The clinician must carefully determine the device size and model to address the clinical needs of the individual patient.

PRIMARY TRACHEOESOPHAGEAL PUNCTURE WITH VOICE PROSTHESIS PLACEMENT

Primary tracheoesophageal puncture and voice prosthesis placement is undertaken following total laryngectomy.

1. Insert the pharynx protector deep enough into the esophagus to ensure protection of the posterior esophageal wall at the level of the puncture (**Diagram 1a**).
2. Remove and discard the blue needle protector from the tip of the puncture needle.
3. Use the puncture needle to carefully create the TEP in the tracheoesophageal wall at the midline 5-10mm below the mucocutaneous juncture into the lumen of the esophagus.
4. Penetration into the esophagus is carefully performed and immediately terminated when the leading edge of the puncture needle is observed.
5. While holding the puncture needle in place, pass one end of the guidewire through the puncture needle hub until it exits the pharynx protector (**Diagram 2**).
6. Withdraw the puncture needle while maintaining grasp on the esophageal end of the guidewire to ensure the guidewire does not pass completely through the TEP.
7. Remove pharynx protector.
8. While maintaining grasp on the tracheal end of the guidewire to secure it in the TEP, pass the esophageal end of the guidewire into the nose of the voice prosthesis inserter until it exits one of the holes of the device (**Diagram 3**).
9. To lock the guidewire in place on the voice prosthesis inserter, push the guidewire tip back into the device via the adjacent open hole and pull the guidewire from the nose of the inserter until flush to create a looped lock (**Diagram 4**).
10. Firmly grasp the tracheal end of the guidewire and with continuous force pull the guidewire back through the puncture until the tip of the voice prosthesis inserter is observed on the tracheal side of the TEP.
11. Maintain grasp on the guidewire with one hand while grasping the nose of the voice prosthesis inserter with the other.

12. Slowly but firmly pull the nose of the voice prosthesis inserter through the TEP until the tracheal flange is visualized in its entirety and deployed against the posterior tracheal wall. The inserter device should release from the prosthesis EXCEPT for at its safety strap (**Diagram 5**).
13. Visually confirm appropriate seating of the esophageal and tracheal flanges against the tracheoesophageal wall. Once confirmed, the voice prosthesis can be released from the device by cutting the safety strap (**Diagram 6**). If leaving the safety strap is desired, it can be carefully released from the inserter by lifting the safety strap from the peg of the inserter.

align the short aspect of the bevel adjacent to the tracheal wall 5mm between the mucocutaneous junction of the tracheostoma. Palpate or transilluminate the membranous trachea to confirm endoscope positioning.

2. Remove and discard the blue needle protector from the tip of the puncture needle.
3. Use the puncture needle to carefully create the TEP in the tracheoesophageal wall at the midline 5-10mm below the mucocutaneous juncture into the lumen of the esophagus.
4. Penetration into the esophagus is carefully performed and immediately terminated when the leading edge of the puncture needle is observed.
5. While holding the puncture needle in place, pass one end of the guidewire through the puncture needle hub until it exits the needle tip into the esophagoscope and out of the mouth (**Diagram 2**).
6. Withdraw the puncture needle while maintaining grasp on the esophageal end of the guidewire to ensure the guidewire does not pass completely through the TEP.
7. Withdraw the esophagoscope while maintaining grasp on the tracheal end of the guidewire to ensure the guidewire does not pass completely through the TEP.
8. While maintaining grasp on the tracheal end of the guidewire to secure it in the TEP, pass the esophageal end of the guidewire into the nose of the voice prosthesis inserter until it exits one of the holes in the body of the device (**Diagram 3**).
9. To lock the guidewire in place on the voice prosthesis inserter, push the guidewire tip back into the device via the adjacent open hole and pull the guidewire from the nose of the inserter until flush to create a looped lock (**Diagram 4**).
10. Firmly grasp the tracheal end of the guidewire and with continuous force pull the guidewire back through the puncture until the tip of the voice prosthesis inserter is observed on the tracheal side of the TEP.
11. Maintain grasp on the guidewire with one hand while grasping the nose of the voice prosthesis inserter with the other.
12. Slowly but firmly pull the nose of the voice prosthesis inserter through the TEP until the

POST-PRIMARY PROCEDURE INSTRUCTIONS

Closely monitor the voice prosthesis post-operatively and carefully clean as needed using only Blom-Singer cleaning devices (sold separately). Use of the placed voice prosthesis for speaking is strictly deferred until deemed appropriate by treating physician to avoid the potential of voicing related airflow causing a pharyngoesophageal fistula.

Once any oral intake is initiated, the voice prosthesis can be checked for leakage by having the patient drink while the clinician observes the device with a bright light.

SECONDARY TRACHEOESOPHAGEAL PUNCTURE WITH VOICE PROSTHESIS PLACEMENT

Prior to consideration of a secondary TEP, confirm a tracheostoma diameter of 1.5 cm or larger. A barium swallow may also be indicated to verify that pharyngoesophageal dimensions are adequate for rigid endoscopy. Insufflation testing (kit sold separately) can be completed to assess for tracheoesophageal speech candidacy.

NOTE: The pharynx protector is not used during secondary TEP procedures. A lighted rigid esophagoscope should be used to ensure protection of the posterior esophageal wall at the level of the puncture.

1. Under anesthesia, introduce a lighted rigid esophagoscope (at minimum: a size 6 mm diameter and 25 cm length) to the level of the tracheostoma (**Diagram 1b**). Orient its tip to

tracheal flange is visualized in its entirety and deployed against the posterior tracheal wall. The inserter device should release from the prosthesis EXCEPT for at its safety strap (**Diagram 5**).

13. Visually confirm appropriate seating of the esophageal and tracheal flanges against the tracheoesophageal wall (reinsertion of esophagoscope may be needed). Once confirmed, the voice prosthesis can be released from the device by cutting the safety strap (**Diagram 6**). If leaving the safety strap is desired, it can be carefully released from the inserter by lifting the safety strap from the inserter peg.

POST-SECONDARY PROCEDURE INSTRUCTIONS

Closely monitor the voice prosthesis post-operatively and carefully clean as needed in situ using only Blom-Singer cleaning devices (sold separately). Use of the placed voice prosthesis for speaking is strictly deferred until deemed appropriate by treating physician.

If oral intake is initiated, the voice prosthesis should be checked for leakage by having the patient drink while the clinician observes the device with a bright light.

RELEASING/RELOADING THE INSERTER

In the event the voice prosthesis insertion portion of the procedure is unsuccessful, the inserter can be released and reloaded onto the guidewire. The TEP tract and guidewire should first be secured:

1. Pass the available end of the guidewire through the TEP tract into the esophagus. If during a secondary procedure, reintroduce the esophagoscope as needed.
2. Firmly hold the voice prosthesis inserter while pushing the guidewire toward the voice prosthesis to release the guidewire from the looped lock (**Diagram 7**).
3. Remove the inserter from the guidewire. It can now be re-locked on the other end of the guidewire as per steps 8 and 9 in above instructions (**Diagram 8**).

To re-load the voice prosthesis onto the inserter:

1. Slide the central sleeve of the inserter towards the nose of the device to reveal the safety peg.

2. Pass the tip of the voice prosthesis safety strap through the blue ring on the inserter (**Diagram 9**) and advance it into the central sleeve of the inserter (**Diagram 10**).
3. Press the safety strap onto the safety peg, then release the central sleeve so that it covers the pegged portion of the inserter (**Diagram 11**).
4. Advance the strap through the inserter ring and gently fold the edges of the tracheal flange of the voice prosthesis to guide it through the inserter ring into place (**Diagram 12**). Do not apply pressure near the valved esophageal side of the prosthesis when reloading.

Once the voice prosthesis is re-loaded and the inserter is re-attached to the guidewire for voice prosthesis insertion, refer back to step 10 in either the Primary or Secondary TEP instructions above and resume the given procedure and post-procedure instructions.

SPECIAL STORAGE AND/OR HANDLING CONDITIONS

There are no special storage and/or handling conditions for this device.

DISPOSAL INSTRUCTIONS

Caution: Exercise care when handling the Blom-Singer FirstFit™ Surgical Kit, as it contains sharp instruments. The device is not biodegradable and can be considered contaminated when used. Carefully dispose of this device per local guidelines (including for contaminated devices/sharps).

PRODUCT COMPLAINTS/EU SERIOUS INCIDENTS

If you are dissatisfied with the device or have any questions, please contact productcomplaints@inhealth.com.

Telephone: +1-800-477-5969

Fax: +1-888-371-1530

Any serious incident that has occurred in relation to the device should be reported to Freudenberg Medical, LLC as noted above and the competent authority of the EU Member State in which the user and/or patient is established.

INTERNATIONAL CUSTOMERS

InHealth Technologies products are available world-wide through authorized InHealth Technologies distributors. For assistance or additional product information please contact your distributor.

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SYMBOLS GLOSSARY (USA only)

SYMBOL	STANDARD No.	SYMBOL TITLE & REF No.	MEANING OF SYMBOL
	EN ISO 15223-1:2016	Manufacturer; 5.1.1	Indicates the medical device manufacturer, as defined in EU Directives or Regulations
	EN ISO 15223-1:2016	Authorized representative in the European Community; 5.1.2	Indicates the authorized representative in the European Community
	EN ISO 15223-1:2016	Catalogue number; 5.1.6	Indicates the manufacturer's catalogue number so that the medical device can be identified
	EN ISO 15223-1:2016	Batch code; 5.1.5	Indicates the manufacturer's batch code so that the batch or lot can be identified
	EN ISO 15223-1:2016	Use-by date; 5.1.4	Indicates the date after which the medical device is not to be used
	EN ISO 15223-1:2016	Caution; 5.4.4	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself
	EN ISO 15223-1:2016	Do not use if package is damaged; 5.2.8	Indicates a medical device that should not be used if the package has been damaged or opened
	EN ISO 15223-1:2016	Do not re-use; 5.4.2	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure
	EN ISO 15223-1:2016	Do not resterilize; 5.2.6	Indicates a medical device that is not to be resterilized
	EN ISO 15223-1:2016	Consult electronic instructions for use; 5.4.3	Indicates the need for the user to consult the electronic instructions for use
	EN ISO 15223-1:2016	Sterilized using ethylene oxide; 5.2.3	Indicates a medical device that has been sterilized using ethylene oxide gas

INHEALTH

TECHNOLOGIES®

UNITED STATES

1110 Mark Avenue
Carpinteria, CA 93013-2918 USA
Toll-Free: 800.477.5969
Toll-Free Fax: 888.371.1530

INTERNATIONAL

Tel: +1.805.684.9337
Fax: +1.805.684.8594
E-mail: order@inhealth.com

ORDER ONLINE

www.inhealth.com

Freudenberg Medical, LLC

1110 Mark Avenue
Carpinteria, CA 93013-2918 USA
www.inhealth.com