INHEALTH

TECHNOLOGIES®

Nasal Septal Perforation Sizing Tool

INSTRUCTIONS FOR USE

Blom-Singer®

voice restoration systems

InHealth Technologies® are registered trademarks of Freudenberg Medical, LLC.

IFU: <https://inhealth.com/ifu>


Copyright © 2024 Freudenberg Medical, LLC. All rights reserved.







INSTRUCTIONS FOR USE

使用说明 / 使用説明書 / Udhëzime për përdorim / إرشادات الاستخدام / Инструкции за употреба / Upute za korištenje / Návod k použití / Brugsanvisning / Gebruiksaanwijzing / Kasutusjuhend / Käyttööhjeet / Instructions d'utilisation / Gebrauchsanweisung / Οδηγίες χρήσης / הוראות שימוש / Használati utasítás / Notkunarleiðbeiningar / Istruzioni per l'uso / 사용 지침 / Lietošanas norādījumi / Naudojimo instrukcijos / Bruksanvisning / Instrukcja użytkowania / Instruções de uso / Instrucțiuni de utilizare / Инструкции по применению / Uputstvo za upotrebu / Návod na použitie / Instrucciones de uso / Bruksanvisningar / Kullanım Talimatları





DIAGRAMS

插图 / 図 / Diagramet / الرسم التوضيحي / Схеми / Dijagrami / Schémata / Diagrammer / Afbeeldingen / Diagrammid / Kaaviot / Diagrammes / Abbildungen / Διαγράμματα / תרשימים / Ábrák / Skýringarmyndir / Diagrammi / 그림 / Diagrammas / Diagramos / Diagrammer / Schematy / Diagramas / Diagrame / Схемы / Crteži / Grafické znázornenia / Diagramas / Diagram / Diyagramlar

Nasal Septal Perforation Sizing Tool Instructions for Use



TABLE OF CONTENTS

目录 / 目次 / Tabela e përmbajtjes / جدول المحتويات / Съдържание / Sadržaj / Obsah /
 Indholdsfortegnelse / Inhoud / Sisukord / Sisällysluettelo / Table des matières /
 Inhaltsverzeichnis / Πίνακας περιεχομένων / תוכן עניינים / Tartalomjegyzék / Efnisyfirlit /
 Sommario / 목차 / Saturs / Turinys / Innholdsfortegnelse / Spis treści / Índice / Cuprins /
 Содержание / Sadržaj / Obsah / Índice / Innehållsförteckning / İçindekiler

EN	Nasal Septal Perforation Sizing Tool	7	HU	Eszköz orrsövény-perforáció méretének meghatározására	58
ZH	鼻中隔穿孔尺寸測量器	10	IS	Stærðarmælir fyrir rof á nefskip	61
JA	鼻中隔穿孔サイジングツール	12	IT	Strumento di dimensionamento della perforazione del setto nasale	64
SQ	Instrument matjeje për perforimin e septumit nazal	15	KO	비중격 천공 크기 측정 도구	67
19	أداة قياس حجم ثقب الحاجز الأنفي	AR	LV	Deguna starpsienas perforācijas mērīšanas rīks	70
BG	Инструмент за оразмеряване на перфорация на назалния септум	22	LT	Nosies pertvaros perforacijos dydžio nustatymo įrankis	73
HR	Alat za određivanje veličine perforacije nazalnog septuma	25	NO	Verktøy for størrelsesmåling av nasal septumperforasjon	76
CS	Nástroj ke stanovení velikosti perforace nosního septa	28	PL	Przemyiar perforacji przegrody nosa	79
DA	Redskab til måling af nasal septumperforation	31	PT	Instrumento para dimensionamento da perfuração no septo nasal	82
NL	Meetinstrument voor neuseptumperforatie	35	RO	Instrument de dimensionare pentru perforație de sept nazal	85
ET	Nina vaheseina perforatsiooni suuruse mõõtmisvahend	38	RU	Инструмент для определения размера перфорации носовой перегородки	88
FI	Nenän väliseinän perforaation mittaustyökalu	41	SR	Alatka za merenje nazalne septalne perforacije	92
FR	Outil de dimensionnement pour perforations du septum nasal	44	SK	Nástroj na určenie veľkosti nazálnej septálnej perforácie	95
DE	Vorrichtung zur Größenbestimmung einer Nasenseptumperforation	48	ES	Instrumento de medición para la perforación del tabique nasal	98
EL	Εργαλείο προσδιορισμού μεγέθους διάτρησης ρινικού διαφράγματος	51	SV	Storleksverktyg för perforation av nässeptum	102
55	כלי למדידת נקב במחיצת האף	HE	TR	Nazal Septal Perforasyon Boyutlandırma Aracı	105

ENGLISH

NASAL SEPTAL PERFORATION SIZING TOOL

Please refer to the diagrams at the front of this instruction manual. Keep the Instructions for Use (IFU), or a copy of it, until all devices have been utilized. IFU is also available online at www.inhealth.com/ifu.

DESCRIPTION

The Nasal Septal Perforation Sizing Tool is a graded handheld device designed to measure nasal septal perforations. The device is available in oval and round configurations with gradations that coincide with the offered sizes of Blom-Singer® Nasal Septal Perforation Prostheses. The device provides a method for measurement of size and thickness of nasal perforations to aid in selection of the appropriate Blom-Singer® Nasal Septal Perforation Prosthesis. The device is non-sterile, single-use, and disposable.

INDICATIONS

The device is indicated for use in an outpatient or clinical setting by qualified and trained medical professionals for the purpose of measuring nasal septal perforations.

CONTRAINDICATIONS

The device is contraindicated for use in patients with nasal abnormalities, beyond a nasal septal perforation, or for use in patients with an allergy to polypropylene plastic.

WARNINGS AND PRECAUTIONS

Use of this device is for qualified and trained professionals only. This device is for sizing reference only, aligned with the Blom-Singer® Nasal Septal Perforation Prostheses. Actual gradations measurements may be within 0.127 mm of claimed size. Do not insert more than one device into the same nostril. Do not insert the device and visualization tool into the same nostril.

This device is for single patient use only. It may not be reused due to risk of infection. Do not use the device if the packaging is opened or in any way damaged, is expired, or becomes contaminated upon arrival; discard and replace it with an unopened product. Do not use if device is broken in any way upon arrival; discontinue use and replace it with an unopened product.

COMPLICATIONS

Although rare, the following complications may occur:

- Enlargement of septal perforation
- Damage to the nasal septum causing fracture, laceration, or heightened problems with airflow
- Nasal tissue irritation

INSTRUCTIONS FOR USE

NOTE: A topical decongestant and/or anesthetic may be sprayed into the nose prior to sizing to aid in visualization and use of endoscope. The use of personal protective equipment is up to the discretion of the qualified and trained professional.


1. Insert a lighted endoscope in either nostril to visualize the perforation shape and select the appropriate sizing tool, round or oval (Diagram 1).
2. In the nostril opposite the lighted endoscope, insert the selected sizing tool with the gradations toward the septum wall (Diagram 1, 2A, 2B). Gradations are on either side of the handle to accommodate for inserting in either nostril.
3. Orient the tool into the perforation to align with the perforation shape. If the initially selected tool shape does not align with perforation, use the alternate shape tool (Diagram 3A, 3B).
4. To determine the size of the perforation and which Blom-Singer® Nasal Septal Perforation Prosthesis to select: Visualize which gradation is visible in the perforation to match shape and size. The entire sizing tool gradation should be visible inside the perforation to ensure an appropriate fit (Diagram 3A - 5B). Table 1 shows the correspondence of gradation to prosthesis size. To further aid in differentiating between gradations, the top and bottom gradations have a smooth surface while the middle gradation has a textured surface.

Table 1. Septal Perforation Sizing Guide (not to scale)

Round			Oval		
Sizing Tool	Size	Prosthesis	Sizing Tool	Size	Prosthesis
	5 mm			11 x 13 mm	
	7 mm			13 x 17 mm	
	10 mm			15 x 21 mm	

5. To determine the thickness of the septum wall and which Blom-Singer® Nasal Septal Perforation Prosthesis to select, use Table 2. The top gradation is equivalent to 2 mm thickness. The middle gradation is an additional 1 mm. The entire sizing tool gradation(s) should be visible inside the perforation to ensure an appropriate fit (Diagram 6A - 7B). If top gradation is the only visible, conclude the septum wall thickness is 2 mm. If the top and middle gradations are both visible, conclude the septum wall thickness is 3 mm. If between two gradations, choose the 3 mm to ensure a comfortable fit.

Table 2. Septum Wall Thickness Sizing Guide (not to scale)

Sizing Tool		Prosthesis Thickness
Round	Oval	Thickness
		 2 mm
		 3 mm

- Dispose of opened tools when complete with sizing.

For more information regarding the Blom-Singer® Nasal Septal Perforation Prosthesis, please refer to the Instructions for Use Manual packaged with the selected device.

SPECIAL STORAGE AND/OR HANDLING CONDITIONS

There are no special storage and/or handling conditions for this device.

DISPOSAL INSTRUCTIONS

The device is not biodegradable and can be considered contaminated when used. Dispose of device after use per local guidelines. Recycle packaging per local guidelines.

PRODUCT COMPLAINTS/EU SERIOUS INCIDENTS

If you are dissatisfied with the device or have any questions, please contact productcomplaints@inhealth.com

Telephone: +1-800-477-5969

Fax: +1-888-371-1530

Any serious incident that has occurred in relation to the device should be reported to Freudenberg Medical, LLC as noted above and the company authority of the EU Member State in which the user and/or patient is established.

INTERNATIONAL CUSTOMERS

InHealth Technologies products are available worldwide through authorized InHealth Technologies distributors. For assistance or additional product information, please contact your distributor.

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SYMBOLS GLOSSARY (USA only)

SYMBOL	STANDARD No.	SYMBOL TITLE & REF No.	MEANING OF SYMBOL
	EN ISO 15223-1:2016	Manufacturer; 5.1.1	Indicates the medical device manufacturer, as defined in EU Directives or Regulations
	93/42/EEC EU 2017/745	CE Marking; 93/42/EEC Annex XII, EU 2017/745 Annex V	CE marking of conformity
	EN ISO 15223-1:2016	Authorized representative in the European Community; 5.1.2	Indicates the authorized representative in the European Community
	EN ISO 15223-1:2016	Batch code; 5.1.6	Indicates the manufacturer's catalogue number so that the medical device can be identified
	EN ISO 15223-1:2016	Batch code; 5.1.5	Indicates the manufacturer's batch code so that the batch or lot can be identified
	EN ISO 15223-1:2016	Caution; 5.4.4	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself
	EN ISO 15223-1:2016	Do not use if package is damaged; 5.2.8	Indicates a medical device that should not be used if the package has been damaged or opened
	EN ISO 15223-1:2016	Consult instructions for use; 5.4.3	Indicates the need for the user to consult the instructions for use
	EN ISO 15223-1:2016	Do not re-use; 5.4.2	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure
	EN ISO 15223-1:2016	Non-sterile; 5.2.7	Indicates a medical device that has not been subjected to a sterilization process.

INHEALTH

TECHNOLOGIES®

UNITED STATES

1110 Mark Avenue
Carpinteria, CA 93013-2918 USA
Toll-Free: +1-800-477-5969
Toll-Free Fax: +1-888-371-1530

INTERNATIONAL

Tel: +1-805-684-9337
Fax: +1-805-684-8594
Email: order@inhealth.com

ORDER ONLINE

www.inhealth.com

Freudenberg Medical, LLC
1110 Mark Avenue
Carpinteria, CA 93013-2918 USA
www.inhealth.com

EMERGO EUROPE
Westervoortsedijk 60
6827 AT Arnhem
The Netherlands