INHEALTH

TECHNOLOGIES®

Blom-Singer® Patient Changeable Voice Prostheses

INSTRUCTIONS FOR USE

Blom-Singer®

voice restoration systems

Blom-Singer® and **InHealth Technologies®** are registered trademarks of **Freudenberg Medical, LLC**.

Patent(s): www.inhealthpatents.com

IFU: <https://inhealth.com/ifu>






Copyright © 2020 Freudenberg Medical, LLC. All rights reserved.

INSTRUCTIONS FOR USE

使用说明 / 使用說明書 / Udhëzime për përdorim / إرشادات الاستخدام / Инструкции за употреба / Upute za korištenje / Návod k použití / Brugsanvisning / Gebruiksaanwijzing / Kasutusjuhend / Käyttöohjeet / Instructions d'utilisation / Gebrauchsanweisung / Οδηγίες χρήσης / הוראות שימוש / Használati utasítás / Notkunarleiddbeiningar / Istruzioni per l'uso / 사용 지침 / Lietošanas norādījumi / Naudojimo instrukcijos / Bruksanvisning / Instrukcja użytkowania / Instruções de uso / Instrucțiuni de utilizare / Инструкции по применению / Uputstvo za upotrebu / Návod na použitie / Instrucciones de uso / Bruksanvisningar / Kullanım Talimatlar

DIAGRAMS

插图 / 図 / Diagramet / الرسم التوضيحي / Схеми / Dijagrami / Schémata / Diagrammer / Afbildingen / Diagrammid / Kaaviot / Diagrammes / Abbildungen / Διαγράμματα / תרשימים / Ábrák / Skýringarmyndir / Diagrammi / 그림 / Diagrammas / Diagramos / Diagrammer / Schematy / Diagramas / Diagreame / Схеми / Crteži / Grafické znázornenia / Diagramas / Diagram / Diyagramlar




TABLE OF CONTENTS

目録 / 目次 / Tabela e përmbajtjes / جدول المحتويات / Съдържание / Sadržaj / Obsah / Indholdsfortegnelse / Inhoud / Sisukord / Sisällysluettelo / Table des matières / Inhaltsverzeichnis / Πίνακας περιεχομένων / תוכן עניינים / Tartalomjegyzék / Efnisyfirlit / Sommario / 목차 / Saturs / Turinys / Innholdsfortegnelse / Spis treści / Índice / Cuprins / Содержание / Sadržaj / Obsah / Índice / Innehållsförteckning / İçindekiler

EN	Patient Changeable Voice Prosthesis	6	HU	Beteg által cserélhető hangprotézis	75
ZH	患者可更换发音假体	10	IS	Talventill sem sjúklingur getur sjálfur skipt um	79
JA	患者交換型ボイスプロテーゼ	14	IT	Protesi fonatoria sostituibile dal paziente	83
SQ	Proteza e zëvendësueshme e zërit e pacientit	18	KO	환자 교체가능 인공후두	88
22	جهاز الحنجرة الهوائية القابل للتغيير	AR	LV	Pacienta maināmā balsu protēze	92
BG	Сменяема гласова протеза	26	LT	Paciento keičiamas balso protezas	96
HR	Glasovna proteza koju može mijenjati pacijent	31	NO	Utskiftbar stemmeprotese	100
CS	Pacientem vyměnitelná hlasová protéza	35	PL	Wymienna proteza głosowa	104
DA	Stemmeprotese, som kan udskiftes af patienten	39	PT	Prótese vocal substituível pelo paciente	109
NL	Door de patiënt verwisselbare stemprothese	44	RO	Proteză fonatorie înlocuibilă de către pacient	113
ET	Patsiendi vahetatav hääleprotees	48	RU	Заменяемый пациентом голосовой протез	118
FI	Potilaan vaihdettavissa oleva puheproteesi	52	SR	Glasovna proteza koju menja pacijent	123
FR	Prothèse vocale changeable par le patient	56	SK	Hlasová protéza prispôsobiteľná pacientovi	127
DE	Vom Patienten auswechselbare Stimmprothese	61	ES	Prótesis fonatoria reemplazable por el paciente	132
EL	Φωνητική πρόθεση που αλλάζεται από τον ασθενή	66	SV	Röstprotes som kan bytas av patient	136
71	ותנדביבור הניתן להחליפה על ידי המטופל	HE	TR	Hasta Tarafından Değiştirilebilir Ses Protezi	140

ENGLISH

BLOM-SINGER® PATIENT CHANGEABLE VOICE PROSTHESES

PRODUCT DESCRIPTION

Please refer to the diagrams located at the front of this instruction manual. The Blom-Singer patient changeable voice prosthesis is available in a Low Pressure style (diagram 1A) or a Duckbill style (diagram 1B). Each device consists of a one-way valve (diagram 1A(a) and 1B(a)), an esophageal retention flange (diagram 1A(b) and 1B(b)), a body that holds the valve assembly (diagram 1A(c) and 1B(c)), a tracheal retention flange (diagram 1A(d) and 1B(d)), and a neck strap with safety hole (diagram 1A(e) and 1B(e)). The device is packaged with an inserter tool. The Blom-Singer voice prosthesis is made of silicone.

INDICATIONS (Reasons to prescribe the device or procedure)

The Blom-Singer voice prosthesis is indicated for tracheoesophageal voice restoration following total laryngectomy.

CONTRAINDICATIONS (Reasons that make it inadvisable to prescribe the particular device or procedure)

The Blom-Singer voice prosthesis is a medical device that should be used only by individuals with experience and training in its use and care.

WARNINGS AND PRECAUTIONS

This product is intended for use and reuse by a single patient only. A qualified, trained clinician is to instruct patients on the use and care of this device and to provide patients with the instructions for use supplied with the product. Initial sizing of the voice prosthesis must be performed by a qualified clinician trained in the use of this particular device. Any handling or management of the voice prosthesis should be performed with clean, dry hands. Handle the voice prosthesis carefully to avoid damaging the prosthesis. Only use a lint-free cloth or tissue to dry the device. Use of non lint-free materials may leave debris which can be aspirated into the airway. Do not use solvents or petroleum-based products for cleaning or lubricating the device. These materials may damage silicone or cause the device to not work properly. Insertion, removal, and cleaning of the prosthesis in the puncture should only be done with the user positioned in front of a mirror with a bright light focused directly on the stoma. Caution should be exercised to avoid accidental displacement of the voice prosthesis which could result in aspiration (inhalation) of the prosthesis. In the unlikely event that this should occur, the user should attempt to cough the voice prosthesis out of the windpipe and immediately seek medical attention if unsuccessful. If there are tears, cracks, or structural damage to the device, discontinue use. Should the device show signs of malfunction, contact your clinician as soon as possible for evaluation. Changes in your anatomy or medical status may lead to improper fitting and/or function of the device. It is recommended that your voice prosthesis and TEP (an opening made by a surgeon between the trachea and the esophagus) be evaluated by your clinician at regular intervals as determined by your clinician. Users with a voice prosthesis who undergo postoperative radiation may experience temporary interruption of voice during treatment. Consult with your clinician for further evaluation if this persists. The device can remain in place throughout radiation treatments. It should be noted that insertion or removal of the voice prosthesis may cause slight bleeding at the TE puncture. You should consult your healthcare provider if bleeding is persistent. If using the Blom-Singer gel caps for insertion, always wait at least three minutes following insertion of the Low Pressure voice prosthesis into the TE puncture before detaching the strap from the safety peg and withdrawing the inserter.

Voice Production

To prevent post-operative complications, do not begin speaking with the voice prosthesis until your doctor has informed you that it is safe for you to do so. The voice prosthesis lumen must be kept clear of blockage to enable patients to voice.

Voice Prosthesis Dislodgement

If the voice prosthesis is accidentally dislodged from the TE puncture, immediately place a Blom-Singer puncture dilator, or a suitable device as recommended by your clinician, of the appropriate diameter, in the puncture to keep it from closing and leaking fluids. The voice prosthesis should be reinserted as soon as possible. Contact your clinician if you are unable to insert the dilator, reinsert or replace your voice prosthesis, or if dislodgement recurs. Never insert foreign objects into the voice prosthesis. Inserting objects other than the Blom-Singer Flushing Device may cause dislodgment of the voice prosthesis or its components and may cause these objects to be swallowed or inhaled.

Voice Prosthesis Leakage

When the valve (diagram 1A(a)) fails to close completely, a few drops of fluid may pass back through the valve from the esophagus to the trachea, which may cause coughing or lead to aspiration. Recurrent leakage of fluids through the voice prosthesis could cause aspiration pneumonia. Gentle handling should always be used when cleaning the voice prosthesis to avoid damage, which could cause leakage.

Device Lifetime

The voice prosthesis may be used until it ceases to function properly, that is, until it leaks or is not providing adequate voice for speech. The voice prosthesis may be used for up to 6 months but this duration should be determined per evaluation and appropriate consultation with your clinician.

COMPLICATIONS

Although rare, the following complications may occur with silicone prostheses of the Blom-Singer type:

- stoma contamination or sepsis, which may require removal of the voice prosthesis and/or appropriate antibiotics; accidental aspiration of the voice prosthesis into the airway, which may require removal by a physician; occasional extrusion of the prosthesis, requiring replacement after dilation of the TE puncture and additional supervision of the stoma care regimen; puncture dilation resulting in leakage of fluids around the voice prosthesis - if leakage is unmanageable, surgical revision or puncture closure may be required; inflammatory reaction around the puncture site and formation of granulation tissue; dislodgment of the voice prosthesis and subsequent closure of the TE puncture; dysphagia (or difficulty swallowing); tearing or other damage to the voice prosthesis from improper use; microbial growth deposits causing leakage through the voice prosthesis or valve incompetence; accidental ingestion of the voice prosthesis into the esophagus and/or GI tract.

INSTRUCTIONS FOR USE

Please refer to the diagrams located at the front of this instruction manual.

Your clinician will instruct you on the use and routine care of the Blom-Singer voice prosthesis. Cleaning of the voice prosthesis should only be done in front of a mirror with a bright light focused directly on the stoma so that the open end of the prosthesis is clearly visible. The following procedural instructions are provided by Eric D. Blom, Ph.D. Once the user has been instructed by a qualified, trained clinician on how to remove, clean and reinsert the voice prosthesis, a maintenance schedule is selected. Some users prefer to remove and clean their voice prosthesis once a week,

while others may leave it in for a longer period of time. It is suggested that users alternate between two voice prostheses.

Voice Prosthesis Removal

Each time the voice prosthesis is removed, a clean Blom-Singer Puncture Dilator, of the appropriate diameter, should be used to prevent leakage of esophageal contents into the trachea. The dilator also serves to slightly expand the puncture, making prosthesis reinsertion easier. For instructions on using the puncture dilator, refer to the Blom-Singer Dilation Sizing System Instructions for Use.

1. Remove the tape from the voice prosthesis neck strap.
2. Grasp the voice prosthesis securely and pull it gently out of the puncture. Do not discard the voice prosthesis if you plan to reinsert it.
3. Insert a clean Blom-Singer Puncture Dilator, of the appropriate diameter, by positioning the tapered tip into the puncture, and gently pushing the dilator into the puncture (diagram 7).
4. Secure the dilator by fastening the handle with medical adhesive tape to the skin around the stoma so it does not slide out.

Voice Prosthesis Cleaning

Please refer to the Voice Prosthesis Cleaning System instructions for use (supplied separately with Blom-Singer Flushing Device) for complete instructions on how to safely clean the voice prosthesis while it is still in your TE puncture. The purpose of cleaning the voice prosthesis is to remove blockage of its lumen that impairs the patient's ability to voice.

Additional Voice Prosthesis Care Instruction

1. Remove the device from the TE puncture.
2. Rinse the device with water.
3. Soak the device in 3% hydrogen peroxide overnight.
4. Rinse thoroughly with water and allow to air dry or dry with a lint-free cloth before reinserting into TE puncture.

Insertion of the Voice Prosthesis

Please refer to the diagrams located at the front of this instruction manual.

Prior to use of the voice prosthesis, check the valve mechanism to ensure it is intact and working properly: the valve of the Low Pressure Voice Prosthesis should close flat inside the prosthesis. The slit valve on the Duckbill Voice Prosthesis should be closed without evidence of gapping.

1. Place the tip of the inserter into the open end of the voice prosthesis (see diagram 2). Avoid squeezing or damaging the valve of the voice prosthesis. Handle the voice prosthesis by the base of the neck strap to avoid touching the portion that inserts into the puncture.
2. Attach the neck strap of the voice prosthesis securely over the safety peg of the inserter (diagram 3). This should prevent inadvertent dislodgment of the voice prosthesis off the inserter during insertion into the puncture.
3. Hold the end of the inserter (diagram 4), and align the top of the voice prosthesis partially in the puncture with the neck strap oriented upwards.
4. Insert the voice prosthesis until the esophageal flange can be felt to “snap” into the esophagus. If the voice prosthesis does not insert easily on the first attempt, do not continue to try to insert it. Instead, reinsert a clean puncture dilator, then try again. If you experience repeated difficulty inserting the voice prosthesis, consult your clinician.

5. While firmly holding the inserter, lift the neck strap off the safety peg on the inserter (diagram 5). Gently withdraw the inserter with a twisting motion from the fully inserted voice prosthesis. Caution should be exercised when removing the inserter to avoid accidental displacement of the voice prosthesis.
6. Secure the neck strap to the neck with a strip of medical adhesive tape to prevent strap movement (diagram 6).

Alternate Method of Insertion of Low Pressure Voice Prosthesis

An alternate method of inserting the Low Pressure Voice Prosthesis is also available. The esophageal flange of the voice prosthesis can be inserted into a gel cap, using the Blom-Singer Gel Cap Insertion System, supplied separately. Note: The Gel Cap Insertion System should not be used with special order low pressure voice prostheses. Refer to the Instructions for Use packaged with the Blom-Singer Gel Cap Insertion System for instructions on loading the Low Pressure Voice Prosthesis into a gel cap. It is recommended that the user discuss the use of this insertion method with a qualified, trained clinician.

1. Handle the voice prosthesis by the base of the neck strap to avoid touching the portion that inserts into the puncture. Place the tip of the voice prosthesis inserter into the open end of the prosthesis (diagram 2).
2. Attach the strap of the prosthesis securely over the safety peg on the inserter (diagram 3). This should prevent inadvertent dislodgment of the prosthesis off the inserter during insertion into the puncture.
3. Dip the gel-cap tip of the voice prosthesis in a light coating of an oil-free, water-soluble lubricant (diagram 8).
4. Hold the end of the inserter (diagram 4), and align the gel-cap tip of the voice prosthesis partially in the puncture with the neck strap oriented upwards. Firmly insert the voice prosthesis until the tracheal flange is flush against the posterior tracheal wall. If the voice prosthesis does not insert easily on the first attempt, do not continue to try to insert it. Instead, reinsert a clean puncture dilator, then try again. This may require use of a new gel cap. If you experience repeated difficulty inserting the voice prosthesis, consult your clinician.
5. Hold the prosthesis in the fully inserted position for at least three minutes to allow the gel cap to dissolve completely. The esophageal retention flange, on an appropriately sized prosthesis that is positioned correctly, will unfold and the prosthesis will be in its proper position in the esophagus.
6. While firmly holding the inserter, lift the neck strap off the safety peg on the inserter (diagram 5). Gently withdraw the inserter with a twisting motion from the fully inserted voice prosthesis. Caution should be exercised when removing the inserter to avoid accidental displacement of the voice prosthesis.
7. Secure the neck strap to the neck with a small strip of medical adhesive tape to prevent strap movement (diagram 6).

SPECIAL STORAGE AND/OR HANDLING CONDITIONS

There are no special storage and/or handling conditions for this device.

DISPOSAL INSTRUCTIONS

This product is not biodegradable and can be considered contaminated when used. Carefully dispose of this device per local guidelines.

ORDERING INFORMATION

Online Ordering
www.inhealth.com

USA

InHealth Technologies
1110 Mark Avenue
Carpinteria, CA 93013-2918
Customer Service: +1-800-477-5969
Fax: +1-888-371-1530
Email: order@inhealth.com

INTERNATIONAL

Contact InHealth Technologies Customer Service at +1-805-684-9337 for international distributor contact.

PRODUCT COMPLAINTS/EU SERIOUS INCIDENTS

If you are dissatisfied with the device or have any questions, please contact

productcomplaints@inhealth.com.

Telephone: +1-800-477-5969

Fax: +1-888-371-1530

Any serious incident that has occurred in relation to the device should be reported to Freudenberg Medical, LLC as noted above and the competent authority of the EU Member State in which the user and/or patient is established.

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SYMBOLS GLOSSARY (USA only)

SYMBOL	STANDARD No.	SYMBOL TITLE & REF No.	MEANING OF SYMBOL
	EN ISO 15223-1:2016	Manufacturer; 5.1.1	Indicates the medical device manufacturer, as defined in EU Directives or Regulations
	93/42/EEC EU 2017/745	CE Marking; 93/42/EEC Annex XII, EU 2017/745 Annex V	CE marking of conformity
	EN ISO 15223-1:2016	Authorized representative in the European Community; 5.1.2	Indicates the authorized representative in the European Community
	EN ISO 15223-1:2016	Batch code; 5.1.6	Indicates the manufacturer's catalogue number so that the medical device can be identified
	EN ISO 15223-1:2016	Batch code; 5.1.5	Indicates the manufacturer's batch code so that the batch or lot can be identified
	EN ISO 15223-1:2016	Use-by date; 5.1.4	Indicates the date after which the medical device is not to be used
	EN ISO 15223-1:2016	Caution; 5.4.4	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself
	EN ISO 15223-1:2016	Do not use if package is damaged; 5.2.8	Indicates a medical device that should not be used if the packaged has been damaged or opened
	EN ISO 15223-1:2016	Consult instructions for use; 5.4.3	Indicates the need for the user to consult the instructions for use
	ISO 7000:2019	Single patient multiple use; 3706	To indicate that the medical device may be used multiple times (multiple procedures) on a single patient

INHEALTH

TECHNOLOGIES®

UNITED STATES

1110 Mark Avenue
Carpinteria, CA 93013-2918 USA
Toll-Free: 800.477.5969
Toll-Free Fax: 888.371.1530

INTERNATIONAL

Tel: +1.805.684.9337
Fax: +1.805.684.8594
E-mail: order@inhealth.com

ORDER ONLINE

www.inhealth.com

Freudenberg Medical, LLC
1110 Mark Avenue
Carpinteria, CA 93013-2918 USA
www.inhealth.com

CE 0344

EMERGO EUROPE
Westervoortsedijk 60
6827 AT Arnhem
The Netherlands