# INHEALTH

TECHNOLOGIES®

## Blom-Singer® Plug Insert

### INSTRUCTIONS FOR USE

**Blom-Singer®**

voice restoration systems

**Blom-Singer®** and **InHealth Technologies®** are registered trademarks of **Freudenberg Medical, LLC**.

Patent(s): [www.inhealthpatents.com](http://www.inhealthpatents.com)

IFU: <https://inhealth.com/ifu>


Copyright © 2020 Freudenberg Medical, LLC. All rights reserved.

## INSTRUCTIONS FOR USE

使用说明 / 使用說明書 / Udhëzime për përdorim / إرشادات الاستخدام / Инструкции за употреба / Upute za korištenje / Návod k použití / Brugsanvisning / Gebruiksaanwijzing / Kasutusjuhend / Käyttöohjeet / Instructions d'utilisation / Gebrauchsanweisung / Οδηγίες χρήσης / הוראות שימוש / Használati utasítás / Notkunarleiðbeiningar / Istruzioni per l'uso / 사용 지침 / Lietošanas norādījumi / Naudojimo instrukcijos / Bruksanvisning / Instrukcja użytkowania / Instruções de uso / Instrucțiuni de utilizare / Инструкции по применению / Uputstvo za upotrebu / Návod na použitie / Instrucciones de uso / Bruksanvisningar / Kullanım Talimatlar

## DIAGRAMS

插图 / 図 / Diagramet / الرسم التوضيحي / Схеми / Dijagrami / Schémata / Diagrammer / Afbeeldingen / Diagrammid / Kaaviot / Diagrammes / Abbildungen / Διαγράμματα / תרשימים / Ábrák / Skýringarmyndir / Diagrammi / 그림 / Diagrammas / Diagramos / Diagrammer / Schematy / Diagramas / Diagrame / Схемы / Crteži / Grafické znázornenia / Diagramas / Diagram / Diyagramlar

## TABLE OF CONTENTS

目録 / 目次 / Tabela e përmbajtjes / جدول المحتويات / Съдържание / Sadržaj / Obsah /  
 Indholdsfortegnelse / Inhoud / Sisukord / Sisällysluettelo / Table des matières /  
 Inhaltsverzeichnis / Πίνακας περιεχομένων / תוכן עניינים / Tartalomjegyzék / Efnisyfirlit /  
 Sommario / 목차 / Saturs / Turinys / Innholdsfortegnelse / Spis treści / Índice / Cuprins /  
 Содержание / Sadržaj / Obsah / Índice / Innehållsförteckning / İçindekiler

|           |                    |           |           |                        |    |
|-----------|--------------------|-----------|-----------|------------------------|----|
| <b>EN</b> | Plug Insert        | 05        | <b>HU</b> | Betétdugó              | 51 |
| <b>ZH</b> | 插塞                 | 08        | <b>IS</b> | Innskotstengi          | 54 |
| <b>JA</b> | プラグ・インサート          | 10        | <b>IT</b> | Tappo                  | 57 |
| <b>SQ</b> | Inserti bllokues   | 13        | <b>KO</b> | 플러그 삽입물                | 60 |
| 16        | Plug Insert        | <b>AR</b> | <b>LV</b> | Spraudņa ieliktnis     | 63 |
| <b>BG</b> | Втулка             | 19        | <b>LT</b> | Kaištis                | 65 |
| <b>HR</b> | Čep                | 22        | <b>NO</b> | Plugginnsats           | 68 |
| <b>CS</b> | Zátková vložka     | 24        | <b>PL</b> | Wtyk                   | 71 |
| <b>DA</b> | Stikindsats        | 27        | <b>PT</b> | Suplemento com plugues | 74 |
| <b>NL</b> | Insteekstop        | 30        | <b>RO</b> | Dop                    | 77 |
| <b>ET</b> | Pistikvahetükk     | 33        | <b>RU</b> | Зарлушка               | 80 |
| <b>FI</b> | Tulppasisäke       | 36        | <b>SR</b> | Čep                    | 83 |
| <b>FR</b> | Insert d'occlusion | 39        | <b>SK</b> | Upchávací zátká        | 86 |
| <b>DE</b> | Stopfeneinsatz     | 42        | <b>ES</b> | Obturador              | 89 |
| <b>EL</b> | Ένθετο πώμα        | 45        | <b>SV</b> | Plugginnsats           | 92 |
| 49        | אטם ההחדרה         | <b>HE</b> | <b>TR</b> | Geçme Tapa             | 95 |

## ENGLISH

**BLOM-SINGER® PLUG INSERT**

Please refer to the diagrams located at the front of this instruction manual.

**PRODUCT DESCRIPTION**

The Blom-Singer Plug Insert (diagram 1) is designed to prevent leakage through a voice prosthesis until the voice prosthesis and tracheoesophageal puncture (TEP) are evaluated by a clinician; it does not allow voicing. The device is non-sterile, made of silicone, and is available in two sizes: 16 Fr and 20 Fr. The device is for use with the Blom-Singer Classic and Advantage Indwelling Voice Prostheses. The Plug Insert consists of a plug tip (diagram 1a) and neck strap (diagram 1b). An inserter tool is included with each device. These products may be used multiple times by a single patient only. Not for use between patients.

**INDICATIONS**

The Blom-Singer Plug Insert may be used in conjunction with a Blom-Singer indwelling voice prosthesis to temporarily prevent leakage of fluids through the voice prosthesis into the trachea during swallowing.

**CONTRAINDICATIONS**

The Plug Insert is a medical product and should be used only by individuals with experience and training in its use and care. The Plug Insert is not intended to permanently replace a properly functioning voice prosthesis. It is not intended to be placed directly in a surgically created fistula between the esophagus and trachea. It is not intended to be used in combination with any other voice prosthesis or device other than the Blom-Singer ADVANTAGE or CLASSIC Indwelling Voice Prostheses. Use of all tracheostoma accessories is contraindicated (must not be used) when the Plug Insert is placed within the voice prosthesis. Use of tracheostoma accessories in conjunction with the Plug Insert may result in accidental displacement of the device or voice prosthesis, which could result in inhalation of the device or voice prosthesis. Do not use the Plug Insert with a 4 mm length voice prosthesis. The Plug Insert is not compatible with 4 mm length voice prostheses.

**WARNINGS AND PRECAUTIONS**

The qualified, trained clinician is to instruct patients on the use and care of this device and to provide patients with the instructions for use supplied with the product. Should the device show signs of malfunction, contact your clinician as soon as possible to have the device performance evaluated. If there are tears, cracks, or structural damage to the device, discontinue use and replace the device. Do not use the device in the event of the packaging being damaged or unintentionally opened prior to use; discard and replace it with a new device. Changes in your anatomy or medical status may lead to improper fitting and/or function of the device. It is recommended that your voice prosthesis and TEP be evaluated by your clinician at regular intervals. The Plug Insert is designed for temporary use only (maximum of 30 days) and is not intended to replace the indwelling voice prosthesis. After placement of the plug insert, contact your clinician as soon as possible. Do not detach or modify the neck strap on the Plug Insert. Do not use solvents or petroleum-based products for cleaning or lubricating the device; these materials may damage silicone or cause the device to not work properly. Removal of debris from the Plug Insert must only be done after the Plug Insert has been removed from the voice prosthesis. Only use a lint-free cloth to dry the device; use of non-lint free materials may leave debris which can be aspirated into the airway. Never force the Plug Insert into the voice prosthesis; this could result in dislodging the voice prosthesis into the esophagus. Should you experience repeated difficulty in managing the Plug Insert, contact your clinician.

### **Voice Prosthesis Dislodgement**

Care should be exercised when inserting or removing the Plug Insert or inserter tool to avoid accidental displacement of the voice prosthesis which could result in aspiration (inhalation) of the prosthesis. In the unlikely event that this should occur, you should attempt to cough the voice prosthesis out of the windpipe; immediately seek medical attention if unsuccessful. The neck strap should not be cut or removed from the Plug Insert. The neck strap is intended to help prevent accidental aspiration (inhalation) of the device should the device become dislodged.

If the voice prosthesis is accidentally dislodged from the TE puncture, immediately place a Blom-Singer Puncture Dilator, or a suitable device as recommended by your doctor, of the appropriate diameter, in the puncture to keep it from closing and leaking fluids. A replacement voice prosthesis should be reinserted within 24 hours. Inserting objects other than the plug insert may cause dislodgment of the voice prosthesis or its components and may cause these objects to be swallowed or inhaled. Do not insert any cleaning device into the voice prosthesis while the Plug Insert is in place.

### **Complications**

Although rare, the following complications may occur with Blom-Singer Plug Inserts:

- stoma or device contamination or bacterial infection, which may require removal of the prosthesis and or treatment with appropriate antibiotics;
- accidental aspiration of devices into the airway, which may require removal by a physician;
- aspiration of liquids through the prosthesis into the airway, which may cause the patient to cough;
- dislodgment of the voice prosthesis and subsequent closure of the TE puncture;
- tearing or other damage to the devices from improper use;
- accidental swallowing of devices into the esophagus and/or gastrointestinal (GI) tract.

### **INSTRUCTIONS FOR USE**

*Please refer to the diagrams located at the front of this instruction manual.*

Your clinician will instruct you on the use and care of this Blom-Singer device. Ensure that the size of the Plug Insert correctly matches the size of the voice prosthesis. Insertion and removal of the Plug Insert should be done in front of a mirror with a bright light focused directly into the stoma. Handle the device by the base of the strap to avoid touching the portion that inserts into the voice prosthesis. The device should always be handled with clean hands to prevent contaminants from entering into the tracheostoma, puncture, or esophagus. The following procedural instructions are provided by Eric D. Blom, Ph.D.

#### **Insertion of the Plug Insert**

1. Place the tip of the inserter tool into the open end of the Plug Insert and attach the hole in the neck strap on the Plug Insert securely over the safety peg on the inserter tool (diagram 2). This should prevent inadvertent dislodgment of the Plug Insert off the inserter tool during insertion into the voice prosthesis.
2. Hold the end of the inserter tool (diagram 3) and align the top of the Plug Insert tip partially inside the voice prosthesis with the neck strap oriented upwards (diagram 4).
3. Carefully and gently push the Plug Insert into the inner lumen of the voice prosthesis until the Plug Insert does not advance any further into the voice prosthesis (diagram 5). When used with most compatible voice prostheses, the Plug Insert will be fully inserted and positioned flush with the tracheal flange of the voice prosthesis (diagram 5a). Do not force the Plug Insert into the voice prosthesis. NOTE: The 20 Fr Plug Insert will only partially insert into 20 Fr voice prostheses that are 5mm or 6mm in length (diagram 5b).

4. If the Plug Insert is inserted and functioning correctly, there should no longer be leakage through the voice prosthesis. After confirming that the Plug Insert is securely in place and inserted, the inserter tool should be carefully removed.
5. To remove the inserter tool: hold the inserter tool and the neck strap together firmly in one hand. Carefully lift the neck strap off the safety peg of the inserter tool. Hold the neck strap against the neck while slowly and gently withdrawing the inserter tool with a twisting motion (diagram 6).
6. Apply a strip of medical adhesive tape over the neck strap to secure the device to the skin above the tracheostoma to prevent strap movement (diagram 7).

#### **Removal of the Plug Insert**

Grasp the neck strap of the Plug Insert and gently pull it out of the voice prosthesis. If the Plug Insert does not come out easily, stabilize/brace the voice prosthesis per your clinician's instructions. Then, re-attempt removal of the Plug Insert.

#### **Plug Insert Care**

After use, the Plug Insert should be inspected and any visible debris should be removed. Remove visible debris by holding the Plug Insert under running water. Air dry or use a lint-free cloth to dry. Once the Plug Insert has dried completely, it should be stored in a clean, resealable plastic bag.

#### **SPECIAL STORAGE AND/OR HANDLING CONDITIONS**

There are no special storage and/or handling conditions for the device.

#### **DISPOSAL INSTRUCTIONS**

The device is not biodegradable and can be contaminated when used. Carefully dispose of the device per local guidelines.

#### **ORDERING INFORMATION**

Online Ordering **www.inhealth.com**  
USA InHealth Technologies  
1110 Mark Avenue  
Carpinteria, CA 93013-2918

Customer Service: +1-800-477-5969  
Fax: +1-888-371-1530  
Email: **order@inhealth.com**

International Contact InHealth Technologies Customer Service at +1-805-684-9337 for international distributor contact.

#### **PRODUCT COMPLAINTS/EU SERIOUS INCIDENTS**

If you are dissatisfied with the device or have any questions, please contact **productcomplaints@inhealth.com**.

Telephone: +1-800-477-5969  
Fax: +1-888-371-1530

Any serious incident that has occurred in relation to the device should be reported to Freudenberg Medical, LLC as noted above and the competent authority of the EU Member State in which the user and/or patient is established.

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**SYMBOLS GLOSSARY (USA only)**

| SYMBOL                                                                             | STANDARD No.             | SYMBOL TITLE & REF No.                                     | MEANING OF SYMBOL                                                                                                                                                                                                          |
|------------------------------------------------------------------------------------|--------------------------|------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|    | EN ISO 15223-1:2016      | Manufacturer; 5.1.1                                        | Indicates the medical device manufacturer, as defined in EU Directives or Regulations                                                                                                                                      |
|    | 93/42/EEC<br>EU 2017/745 | CE Marking; 93/42/EEC Annex XII, EU 2017/745 Annex V       | CE marking of conformity                                                                                                                                                                                                   |
|    | EN ISO 15223-1:2016      | Authorized representative in the European Community; 5.1.2 | Indicates the authorized representative in the European Community                                                                                                                                                          |
|    | EN ISO 15223-1:2016      | Catalogue number; 5.1.6                                    | Indicates the manufacturer's catalogue number so that the medical device can be identified                                                                                                                                 |
|    | EN ISO 15223-1:2016      | Batch code; 5.1.5                                          | Indicates the manufacturer's batch code so that the batch or lot can be identified                                                                                                                                         |
|   | EN ISO 15223-1:2016      | Use-by date; 5.1.4                                         | Indicates the date after which the medical device is not to be used                                                                                                                                                        |
|    | EN ISO 15223-1:2016      | Caution; 5.4.4                                             | Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself |
|   | EN ISO 15223-1:2016      | Do not use if package is damaged; 5.2.8                    | Indicates a medical device that should not be used if the package has been damaged or opened                                                                                                                               |
|  | EN ISO 15223-1:2016      | Consult instructions for use; 5.4.3                        | Indicates the need for the user to consult the instructions for use                                                                                                                                                        |
|  | ISO 7000:2019            | Single patient multiple use; 3706                          | To indicate that the medical device may be used multiple times (multiple procedures) on a single patient                                                                                                                   |

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