
INHEALTH

TECHNOLOGIES®

Blom-Singer® Tracheostoma Valve Housings

INSTRUCTIONS FOR USE

Blom-Singer®

voice restoration systems

Blom-Singer® and **InHealth Technologies®** are registered trademarks of **Freudenberg Medical, LLC**.

Patent(s): www.inhealthpatents.com

IFU: <https://inhealth.com/ifu>




Copyright © 2020 Freudenberg Medical, LLC. All rights reserved.

INSTRUCTIONS FOR USE

使用说明 / 使用說明書 / Udhëzime për përdorim / إرشادات الاستخدام / Инструкции за употреба / Upute za korištenje / Návod k použití / Brugsanvisning / Gebruiksaanwijzing / Kasutusjuhend / Käyttöohjeet / Instructions d'utilisation / Gebrauchsanweisung / Οδηγίες χρήσης / הוראות שימוש / Használati utasítás / Notkunarleiðbeiningar / Istruzioni per l'uso / 사용 지침 / Lietošanas norādījumi / Naudojimo instrukcijos / Bruksanvisning / Instrukcja użytkowania / Instruções de uso / Instrucțiuni de utilizare / Инструкции по применению / Uputstvo za upotrebu / Návod na použitie / Instrucciones de uso / Bruksanvisningar / Kullanım Talimatlar

DIAGRAMS

插图 / 図 / Diagramet / الرسم التوضيحي / Схеми / Dijagrami / Schémata / Diagrammer / Afbeeldingen / Diagrammid / Kaaviot / Diagrammes / Abbildungen / Διαγράμματα / תרשימים / Ábrák / Skýringarmyndir / Diagrammi / 그림 / Diagrammas / Diagramos / Diagrammer / Schematy / Diagramas / Diagrame / Схеми / Crteži / Grafické znázornenia / Diagramas / Diagram / Diyagramlar



TABLE OF CONTENTS

目録 / 目次 / Tabela e përmbajtjes / جدول المحتويات / Съдържание / Sadržaj / Obsah /
 Indholdsfortegnelse / Inhoud / Sisukord / Sisällysluettelo / Table des matières / Inhaltsverzeichnis /
 Πίνακας περιεχομένων / תוכן עניינים / Tartalomjegyzék / Efnisyfirlit / Sommario / 목차 / Saturs /
 Turinys / Innholdsfortegnelse / Spis treści / Índice / Cuprins / Содержание / Sadržaj / Obsah / Índice /
 Innehållsförteckning / İçindekiler

EN	Tracheostoma Housings	6	HU	Tracheostoma szelepházak	38
ZH	气管造口外壳	8	IS	Hettur fyrir barkarauf	40
JA	気管孔ハウジング	10	IT	Alloggiamenti per tracheostomia	42
SQ	Mbajtëset e trakeostomës	12	KO	기관절개창 하우징	44
14	مبيئات الفغرة الرغامية	AR	LV	Traheostoma korpusi	46
BG	Трахеостома корпуси	16	LT	Tracheostomos korpusai	48
HR	Ležište traheostome	18	NO	Tracheostoma Housings	50
CS	Tracheostomické kryty	20	PL	Tracheostomabehuizingen	52
DA	Tracheostomické kryty	22	PT	Base para traqueostoma	54
NL	Tracheostomabehuizingen	24	RO	Carcase pentru traheostomă	56
ET	Trahheostoomi mansetid	26	RU	Carcase pentru traheostomă	58
FI	Trakeostoomakotelot	28	SR	Tracheostómické kryty	60
FR	Logements pour trachéostomie	30	SK	Tracheostómické kryty	62
DE	Tracheostomické kryty	32	ES	Alojamientos para traqueostoma	64
EL	Θήκες τραχειοστομίας	34	SV	Trakeostomihöljen	66
36	תושבות פיום קנה	HE	TR	Trakeostoma Muhafazaları	68

ENGLISH

BLOM-SINGER® TRACHEOSTOMA HOUSINGS

PRODUCT DESCRIPTION

The Blom-Singer Tracheostoma Housing (diagram 1) attaches to the skin surrounding the tracheostoma and acts as a holder for tracheostomal devices, such as: the Blom-Singer HME/ HumidiFilter System, the Blom-Singer Adjustable Tracheostoma Valve, or the Blom-Singer Shower Guard. The Blom-Singer Tracheostoma Housing is available in multiple styles and sizes. The device is reusable for a single patient. The valve housing requires user assembly with attachment of adhesive discs (supplied separately).

INDICATIONS (Reasons to prescribe the device or procedure)

The Blom-Singer Tracheostoma Housing is intended for use by laryngectomy patients with a tracheostoma. The Blom-Singer Tracheostoma Housing is designed to be applied to the skin area surrounding the tracheostoma and acts as a holder for tracheostomal devices, including but not limited to: the Blom-Singer HME/ HumidiFilter System, the Blom-Singer Adjustable Tracheostoma Valve, or the Blom-Singer Shower Guard.

WARNINGS AND PRECAUTIONS

Wash hands with soap and water prior to using to avoid contamination. Do not use solvents or petroleum-based products for cleaning the housings, as this could damage the devices. Voice prosthesis neck straps may obstruct the adhesive seal around the stoma. To reduce the possibility of air leaks, excess strap length may be trimmed as needed. Do not trim below the portion of the strap with the safety hole (diagram 2). Rare occurrence of skin irritation, rash or minor infection can occur. If you experience excessive skin irritation or rash, stop use immediately and contact your clinician. The Tracheostoma Valve Housing is for re-use by a single patient. Cross contamination leading to patient infection can occur if the device is shared between patients. Inspect the housing for physical damage, such as cracks or tears, and do not use if damaged. The device may not adhere securely to skin if: the skin is not completely dry after cleaning, if these instructions are not followed, if an appropriate model/size is not selected, if attempting to talk too loudly, or due to patient anatomy. Contact your clinician for assistance if needed.

INSTRUCTIONS FOR USE

Please refer to the diagrams located at the front of this instruction manual.

Removal and attachment of any devices worn over the tracheostoma should only be done with the user positioned in front of a mirror with a bright light focused directly on the stoma. Tracheostoma Valve Housing should be applied to clean, dry skin. Apply a protective skin barrier as needed (supplied separately). If a stronger seal is desired, apply a thin coat of medical grade liquid silicone adhesive (supplied separately) to the skin prior to attaching the housing. Allow the silicone adhesive to dry completely before attaching the assembled housing.

Blom-Singer Tracheostoma Valve Housing + Adhesive Disc Assembly

1. Peel the backing paper from one side of the disc only (diagram 3).
2. Place the Tracheostoma Valve Housing on top of the adhesive disc.
3. Flatten the Tracheostoma Valve Housing against the adhesive disc surface, pressing firmly around the entire housing to ensure complete contact with the adhesive disc (diagram 4).
4. Remove the paper backing on the reverse side of the disc to expose the second adhesive surface (diagram 5).

Attaching Assembled Housing to Skin

1. Position the assembled housing over the tracheostoma.
2. Press the assembled housing flush against the skin surface to avoid air pockets between the skin and the housing.
3. Press the assembled housing firmly against the skin around the entire housing to achieve maximum adhesive contact (diagram 6).

Removal

1. Carefully remove the assembled housing from the skin. Medical grade adhesive remover (supplied separately) may be used to assist with device removal. **Extreme care should be used to avoid getting adhesive remover or fumes in the tracheostoma when cleaning the surrounding skin area.**
2. Peel back the adhesive disc from the Tracheostoma Valve Housing. Retain the housing for cleaning and reuse. Discard the used adhesive disc.

Cleaning

The Tracheostoma Valve Housing may be cleaned and reused. Remove residual adhesive from the surface of the Tracheostoma Valve Housing with medical grade adhesive remover or by soaking the device overnight in water and then rubbing off the adhesive under warm tap water.

Device Replacement

Replace Tracheostoma Valve Housing when its flexibility decreases.

SPECIAL STORAGE AND/OR HANDLING CONDITIONS

There are no special storage and/or handling conditions for this device.

DISPOSAL INSTRUCTIONS

This device is not biodegradable and can be contaminated when used. Carefully dispose of the device per local guidelines.

ORDERING INFORMATION

Blom-Singer products may be ordered directly from InHealth Technologies. TELEPHONE: Toll-Free (800)477-5969 or (805)684-9337. FAX: Toll-Free (888)371-1530 or (805)684-8594. EMAIL: order@inhealth.com ORDER ONLINE: www.inhealth.com POST: InHealth Technologies, 1110 Mark Avenue, Carpinteria, CA 93013-2918, USA, Attention: Customer Service

INTERNATIONAL

Please contact our customer service department for a distributor referral.

PRODUCT COMPLAINTS/EU SERIOUS INCIDENTS

If you are dissatisfied with the device or have any questions, please contact **productcomplaints@inhealth.com**.

Telephone: +1-800-477-5969

Fax: +1-888-371-1530

Any serious incident that has occurred in relation to the device should be reported to Freudenberg Medical, LLC as noted above and the competent authority of the EU Member State in which the user and/or patient is established.

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SYMBOLS GLOSSARY (USA only)

SYMBOL	STANDARD No.	SYMBOL TITLE & REF No.	MEANING OF SYMBOL
	EN ISO 15223-1:2016	Manufacturer; 5.1.1	Indicates the medical device manufacturer, as defined in EU Directives or Regulations
	93/42/EEC EU 2017/745	CE Marking; 93/42/EEC Annex XII, EU 2017/745 Annex V	CE marking of conformity
	EN ISO 15223-1:2016	Authorized representative in the European Community; 5.1.2	Indicates the authorized representative in the European Community
	EN ISO 15223-1:2016	Batch code; 5.1.6	Indicates the manufacturer's catalogue number so that the medical device can be identified
	EN ISO 15223-1:2016	Batch code; 5.1.5	Indicates the manufacturer's batch code so that the batch or lot can be identified
	EN ISO 15223-1:2016	Use-by date; 5.1.4	Indicates the date after which the medical device is not to be used
	EN ISO 15223-1:2016	Caution; 5.4.4	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself
	EN ISO 15223-1:2016	Do not use if package is damaged; 5.2.8	Indicates a medical device that should not be used if the packaged has been damaged or opened
	EN ISO 15223-1:2016	Consult instructions for use; 5.4.3	Indicates the need for the user to consult the instructions for use
	ISO 7000:2019	Packaging unit; 2794	To indicate the number of pieces in the package
	ISO 7000:2019	Single patient multiple use; 3706	To indicate that the medical device may be used multiple times (multiple procedures) on a single patient
	EN ISO 15223-1:2016	Do not re-use; 5.4.2	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure

INHEALTH

TECHNOLOGIES®

UNITED STATES

1110 Mark Avenue
Carpinteria, CA 93013-2918 USA
Toll-Free: 800.477.5969
Toll-Free Fax: 888.371.1530

INTERNATIONAL

Tel: +1.805.684.9337
Fax: +1.805.684.8594
E-mail: order@inhealth.com

ORDER ONLINE

www.inhealth.com

Freudenberg Medical, LLC
1110 Mark Avenue
Carpinteria, CA 93013-2918 USA
www.inhealth.com

EMERGO EUROPE
Westervoortsedijk 60
6827 AT Arnhem
The Netherlands