# INHEALTH<sup>®</sup>

TECHNOLOGIES

## Blom-Singer<sup>®</sup> Indwelling Voice Prostheses

USER Instructions For Use

**Blom-Singer<sup>®</sup>**

voice restoration systems

**R1**

37639-01H

37639-01H Effective June 2018 / Gültig ab Juni 2018 / Efectivo a partir de junio de 2018 / Date d'entrée en vigueur : juin 2018 / Valido da giugno 2018 / Met ingang van juni 2018 / Vigente em junho de 2018 / V platnosti od června 2018 / Wchodzi w życie w czerwcu 2018 r.

**Blom-Singer** and **InHealth Technologies** are registered trademarks of **Freudenberg Medical, LLC**

Patent(s): [www.inhealthpatents.com](http://www.inhealthpatents.com)

# TABLE OF CONTENTS

Inhaltsverzeichnis / Índice / Sommaire / Sommario / Inhoud / Índice / Obsah / Spis treści

<b>ENG</b>	<b>Indwelling Voice Protheses</b>	<b>4</b>
<b>DEU</b>	<b>Verweil-Stimmprothesen</b>	<b>10</b>
<b>SPA</b>	<b>Prótesis fonatorias fijas</b>	<b>18</b>
<b>FRA</b>	<b>Prothèses vocales à demeure</b>	<b>26</b>
<b>ITA</b>	<b>Protesi fonatoria fissa</b>	<b>34</b>
<b>NLD</b>	<b>Verblijfstemprothese</b>	<b>42</b>
<b>POR</b>	<b>Próteses vocais fixas</b>	<b>48</b>
<b>CES</b>	<b>Vnitřní hlasové protézy</b>	<b>56</b>
<b>POL</b>	<b>Stałe protezy głosowe</b>	<b>62</b>

## **Blom-Singer®**

voice restoration systems

### **USER Instructions For Use**

Gebrauchsanweisung für BENUTZER  
Instrucciones para el USUARIO  
Instructions d'utilisation destinées à L'UTILISATEUR  
Istruzioni per l'uso destinate ALL'UTENTE  
Gebruiksaanwijzing voor de GEBRUIKER  
Instruções de uso para o USUÁRIO  
Návod k použití pro UŽIVATELE  
Instrukcja użytkowania dla UŻYTKOWNIKA

## **BLOM-SINGER® INDWELLING VOICE PROSTHESES**

### **INDICATIONS**

The Blom-Singer indwelling voice prosthesis is indicated for tracheoesophageal voice restoration following total laryngectomy, when placement, or replacement, of an indwelling voice prosthesis is performed by a qualified, trained medical professional.

### **CONTRAINDICATIONS**

The Blom-Singer indwelling voice prosthesis is a medical product that should be used only by a qualified, medical professional with experience and training in its use and care.

### **HOW SUPPLIED**

The Blom-Singer indwelling voice prosthesis is supplied in sterile or non-sterile options, as a complete voice restoration system for primary, secondary, or replacement procedures. Sterile voice prostheses are sterilized by ethylene oxide. Included in each package is a Voice Prosthesis Cleaning System for the patient, two (2) gel caps, flange introducer and inserter.

### **PRODUCT DESCRIPTION**

Please refer to the diagram located at the back of this instruction manual. Four distinct styles and a large selection of sizes are available to allow your clinician to choose a voice prosthesis best suited for your needs.

The Blom-Singer voice prosthesis is made of silicone, a biocompatible material used in many medical products. The gel cap used for insertion of the voice prosthesis is made of a vegetable base gelatin. Consult with your physician as individual reactions to these materials may occur.

The Blom-Singer voice prosthesis, (diagram 1), consists of a one-way silicone flap valve (diagram 1a), an esophageal retention flange (diagram 1b), a body that holds the valve assembly (diagram 1c), a tracheal retention flange (diagram 1d), and a neck strap (diagram 1e).

## WARNINGS AND PRECAUTIONS

The qualified, trained medical professional is to instruct patients on the use and care of this device and to provide patients with the Blom-Singer indwelling voice protheses instructions for use supplied with the product.

**This Blom-Singer medical device may be used only once by one patient. It may not be reused.** Reuse of this device is not possible after the initial placement because the strap is removed from the voice prosthesis at the time of initial placement.

Should the voice prosthesis show signs of malfunction, contact your physician/clinician as soon as possible to have the device performance evaluated.

If there are tears, cracks, or structural damage to the voice prosthesis or cleaning devices, discontinue use.

Changes in your anatomy or medical status may lead to improper fitting and/or function of the device. It is recommended that your voice prosthesis and TEP be evaluated by your clinician at regular intervals.

### **Voice Production**

To prevent post-operative complications, do not begin speaking with the voice prosthesis until your doctor has informed you that it is safe for you to do so.

The voice prosthesis lumen must be kept clear of blockage to enable patients to voice.

### **Voice Prosthesis Dislodgement**

If the voice prosthesis is accidentally dislodged from the TE puncture, immediately place a Blom-Singer Puncture Dilator, or a suitable device as recommended by your doctor, of the appropriate diameter, in the puncture to keep it from closing and leaking fluids. Contact your doctor immediately. Care should be exercised when inserting or removing cleaning devices to avoid accidental displacement of the voice prosthesis which could result in aspiration (inhalation) of the voice prosthesis. In the unlikely event that this should occur, you should attempt to cough the voice prosthesis out of the windpipe. Immediately seek medical attention if unsuccessful.

---

Never insert foreign objects into the voice prosthesis. The cleaning devices are designed to remove materials that cause blockage of the lumen of a voice prosthesis. Inserting objects other than the supplied cleaning devices may cause dislodgment of the voice prosthesis or its components and may cause these objects to be swallowed or aspirated..

### **Voice Prosthesis Leakage**

Always use gentle handling and pressure when cleaning the voice prosthesis to avoid deterioration of the device and damaging internal components which could cause leakage.

Recurrent leakage of fluids through the voice prosthesis should be evaluated by your clinician soon after it is noticed. When the flap valve (diagram 1a) fails to close completely a few drops of fluid may pass back through the valve from the esophagus to the trachea, which may cause coughing.

### **Microbial Growth**

Most users have microbial growth in their mouth, throat, or esophagus; this is a commonly occurring condition. Microbial growth deposits can cause valve deformation and failure, i.e., liquid leakage through the voice prosthesis and an increase in the pressure necessary to voice. Replacement of the voice prosthesis may be required.

### **Device Removal**

Never attempt to remove the device by yourself or allow anyone other than a qualified, trained medical professional to remove the voice prosthesis.

The indwelling voice prosthesis is not a permanent device and requires replacement periodically by your doctor.

The indwelling voice prosthesis may be left in place in the TE puncture until it ceases to function correctly, that is, until it leaks or is not providing adequate voice for speech. Many clinicians recommend routine follow-up evaluations not longer than 6 months from initial placement.

**Caution:** Federal (USA) law restricts this device to sale by or on the order of a physician.

## **COMPLICATIONS**

Although rare, the following complications have been identified

to occur with silicone prostheses of the Blom-Singer type. They include: stoma contamination or bacterial infection, which may require removal of the voice prosthesis and/or appropriate antibiotics; accidental aspiration of the voice prosthesis into the airway, which may require removal by a physician; inhalation of liquids through the prosthesis into the airway which may cause the patient to cough; loss of prosthesis retention requiring device replacement; puncture dilation resulting in leakage of fluids around the voice prosthesis; inflammatory reaction around the puncture site and formation of granulation tissue; dislodgment of the voice prosthesis and subsequent closure of the TE puncture; unmanageable leakage around the voice prosthesis, requiring surgical revision or closure of the puncture; difficulty speaking; dysphagia, tearing or other damage to the voice prosthesis from improper use; microbial growth deposits causing leakage through the voice prosthesis or valve incompetence; accidental ingestion of the voice prosthesis into the esophagus.

## **VOICE PROSTHESIS CLEANING AND CARE INSTRUCTIONS**

The following instructions should be made clear to the patient as part of the routine care and cleaning of the indwelling voice prosthesis.

The following information is a brief overview of the proper care and maintenance of the voice prosthesis. Please refer to the Voice Prosthesis Cleaning System instructions for use supplied with this product for complete instructions on how to use the cleaning devices.

The purpose of cleaning the voice prosthesis is to remove blockage of its lumen that impairs the patient's ability to voice.

The cleaning device provides a means for careful cleaning of the voice prosthesis while in the patient's TE puncture.

**Caution:** Use **only** the supplied cleaning devices. **Do not** insert foreign objects into the voice prosthesis. Inserting objects other than the supplied cleaning device may cause dislodgment of the voice prosthesis or its components.

The user's hands should always be washed with soap and water prior to using the cleaning system.

---

Cleaning of the indwelling voice prosthesis should only be done in front of a mirror with a bright light focused directly on the stoma so that the open end of the prosthesis is clearly visible.

Only use a lint free cloth or tissue to dry the device. Use of non-lint free materials may leave debris which can be aspirated into the airway.

Do not use solvents or petroleum-based products for cleaning or lubricating the device. These materials may damage silicone or cause the device to not work properly.

## **SPECIAL ORDER PRODUCTS**

If this instruction manual accompanies a Special Order Product, there may be differences in the physical characteristics between the product and the product descriptions described in this instruction manual. These differences will not affect the safety or desired beneficial effect of the special order product. Special Order products are nonreturnable.

## **ORDERING INFORMATION**

### **USA**

Blom-Singer products may be ordered directly from InHealth Technologies. TELEPHONE: Toll-Free (800)477-5969 or (805)684-9337

FAX: Toll-Free (888)371-1530 or (805)684-8594. EMAIL: [order@inhealth.com](mailto:order@inhealth.com)

ORDER ON-LINE: [www.inhealth.com](http://www.inhealth.com) POST: InHealth Technologies, 1110 Mark Avenue, Carpinteria, CA 93013-2918, USA, Attention: Customer Service

### **International**

Please contact our customer service department for a distributor referral.

### **Consumer Affairs**

If you have any questions or dissatisfaction with a product please contact our customer service department via telephone, fax, post, or Email: [productcomplaints@inhealth.com](mailto:productcomplaints@inhealth.com)

## GLOSSARY

### **Clinician**

A health professional, such as a physician, licensed speech-language pathologist, or nurse who is involved in clinical practice.

### **Contraindications**

A reason that makes it inadvisable to prescribe a particular procedure or treatment.

### **Esophagus**

The muscular membranous tube for the passage of food from the pharynx to the stomach; the “food pipe.”

### **Indications**

A reason to prescribe a drug or perform a procedure.

### **Microbial Growth**

Most users have microbial growth in their mouth, throat, or esophagus; this is a commonly occurring condition. Microbial growth deposits can cause valve deformation and failure, i.e., liquid leakage through the voice prosthesis and an increase in the pressure necessary to voice. Replacement of the voice prosthesis may be required.

### **Tracheoesophageal puncture (TEP)**

An opening made by a surgeon between the trachea and the esophagus. The opening is made at the time of initial surgery or later so that a voice prosthesis can be placed into this passage.

### **Trachea**

The airway through which respiratory air passes to and from the lungs (windpipe).

### **Tracheoesophageal (TE)**

Relating to the trachea and esophagus.

### **Voice prosthesis**

A one-way speech valve that is inserted into the TEP to allow air to pass from the trachea through the prosthesis and into the esophagus to produce sound.

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
Additional references available upon request.

Weitere Referenzen aus der Literatur sind auf Anfrage erhältlich. /  
Se pueden solicitar más referencias. / Autres références disponibles  
sur demande. / Riferimenti aggiuntivi sono disponibili a richiesta. /  
Meer informatie is op aanvraag beschikbaar. / Referências adicionais  
mediante solicitação. / Další odkazy na literaturu jsou k dispozici na  
požádání. / Dodatkowe informacje dostępne na żądanie.

1. Blom, E.D.: "Some Comments on the Escalation of Tracheoesophageal Voice Prosthesis Dimensions," **Arch Otolaryngol**, 129: 500-502, 2003.
2. Eerenstein, S.E., Grolman, W., Schouwenburg, P.F.: "Downsizing of Voice Prosthesis Diameter in Patients with Laryngectomy," **Arch Otolaryngol**, 128: 838-841, 2003.

- 
3. Hamaker, R.C., Blom, E.D.: "Botulinum Neurotoxin for Pharyngeal Constrictor Muscle Spasm in Tracheoesophageal Voice Restoration," **Laryngoscope**, 113:1479-1482, 2003.
  4. Evans, P.H.R., Blom, E.D.: "Functional Restoration of Speech," **Principles and Practice of Head and Neck Oncology**, Evans, P.H.R., Montgomery, P.W., Gullane, P.J. (Eds.) Taylor and Francis Group, London, 2003.
  5. Blom, E.D., Singer M.I., and Hamaker, R.C., **Tracheoesophageal Voice Restoration Following Total Laryngectomy**. San Diego: Singular Publishing Group, Inc.1998.
  6. Blom, E.D., Hamaker, R.C.: "Tracheoesophageal Voice Restoration Following Total Laryngectomy," **Cancer of the Head and Neck**, E.N. and Suen, J. (Eds). W.B. Saunders Publishers. 839-852, 1996.
  7. Blom, E.D., "Tracheoesophageal Speech," **Seminars in Speech and Language**. New York: Thieme Medical Publishers. 16:3, 191-204, 1995.
  8. Pauloski, B.R., Blom, E.D., Logemann, J.A., Hamaker, R.C. "Functional Outcome After Surgery for Prevention of Pharyngospasms in Tracheoesophageal Speakers. Part II: Swallow Characteristics," **Laryngoscope**, 105:10: 1104-1110, 1995.
  9. Blom, E.D., Pauloski, B.R., Hamaker, R.C. "Functional Outcome After Surgery for Prevention of Pharyngospasms in Tracheoesophageal Speakers. Part I: Speech Characteristics," **Laryngoscope**, 105:10: 1093-1103, 1995.
  10. Blom, E.D., Hamaker, R.C., Freeman, S.B.: "Postlaryngectomy Voice Restoration," Lucente, F.E. (Ed.): **Highlights of the Instructional Courses**, Vol. 7. St. Louis: Mosby-Year Book. pp 3-10, 1994.
  11. Blom, E.D., Singer M.I., and Hamaker, R.C. "A Prospective Study of Tracheoesophageal Speech," **Arch Otolaryngol**, 112: 440-447, 1986.
  12. Blom, E.D., Singer, M.I. "Disinfection of Silicone Voice Prostheses," **Arch Otolaryngol**, 112: 1986.
  13. Blom, E.D., Singer M.I., and Hamaker, R.C. "An Improved Esophageal Insufflation Test," **Arch Otolaryngol**, 111: 211-212, 1985.
  14. Hamaker, R.C., Singer M.I., Blom E.D. and Daniels, H.A. "Primary Voice Restoration at Laryngectomy," **Arch Otolaryngol**, 111: 182-186, 1985.
  15. Singer, M.I. and Blom, E.D. "A Selective Myotomy for Voice Restoration After Total Laryngectomy," **Arch Otolaryngol**, 107: 670-673, 1981.
  16. Singer, M.I. and Blom, E.D. "An Endoscopic Technique for Restoration of Voice After Laryngectomy," **Ann Otol Rhinol Laryngol**, 89:6: 529-533, 1980.

# INHEALTH®

TECHNOLOGIES

---

## UNITED STATES

1110 Mark Avenue  
Carpinteria, CA 93013-2918 USA  
Toll-Free: 800.477.5969  
Toll-Free Fax: 888.371.1530

## INTERNATIONAL

Tel: +1.805.684.9337  
Fax: +1.805.684.8594  
E-mail: [order@inhealth.com](mailto:order@inhealth.com)

## ORDER ONLINE

[www.inhealth.com](http://www.inhealth.com)

**Freudenberg Medical, LLC**  
1110 Mark Avenue  
Carpinteria, CA 93013-2918 USA  
[www.inhealth.com](http://www.inhealth.com)

CE 0344

**EMERGO EUROPE**  
Prinsessegracht 20  
2514 AP The Hague  
The Netherlands