# INHEALTH<sup>®</sup>

## TECHNOLOGIES

### Blom-Singer<sup>®</sup> Indwelling TEP Occluder

**MEDICAL PROFESSIONAL Instructions For Use**

**Blom-Singer<sup>®</sup>**

voice restoration systems

R1  
37720-01G

37720-01G Effective October 2017 / Gültig ab Oktober 2017 /  
Efectivo a partir de octubre de 2017 / Date d'entrée en vigueur :  
octobre 2017 / Valido da ottobre 2017 / Met ingang van oktober  
2017 / Vigente em outubro de 2017 / V platnosti od října 2017 /  
Wchodzi w życie w październiku 2017 r.

**Blom-Singer and InHealth Technologies** are registered  
trademarks of **Freudenberg Medical, LLC**

## TABLE OF CONTENTS

Inhaltsverzeichnis / Índice / Sommaire / Sommario / Inhoud / Obsah / Spis treści

ENG	<b>Blom-Singer Indwelling TEP Occluder</b>	5
DEU	<b>Blom-Singer TEP-Verweil-Okkluder</b>	17
SPA	<b>Dispositivo de oclusión fijo para TEP Blom-Singer</b>	31
FRA	<b>Dispositif d'occlusion PTO à demeure Blom-Singer</b>	44
ITA	<b>Occlusore TEP fisso Blom-Singer</b>	58
NLD	<b>Blom-Singer afsluitbare TEP-verblijfstemprothese</b>	72
POR	<b>Oclusor de PTE fixo Blom-Singer</b>	86
CES	<b>Nástroj pro uzavření vnitřní TEP</b>	99
POL	<b>Okluder stały NTP Blom-Singer</b>	111

**Blom-Singer®**

voice restoration systems

### MEDICAL PROFESSIONAL Instructions For Use

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**4 | 37720-01G**

## **BLOM-SINGER® INDWELLING TEP OCCLUDER**

The Blom-Singer Indwelling TEP (Tracheoesophageal Puncture) Occluder is a non-valved silicone device designed for placement in and maintenance of the TE (tracheoesophageal) puncture following total laryngectomy.

The Indwelling TEP Occluder is intended to:

1. Maintain the TEP fistula by keeping the puncture tract open.
2. Prevent leakage of esophageal contents into the trachea during swallowing.

**This Blom-Singer medical device may be used only once by one patient. It may not be reused.** See Warnings, Precautions and Complications.

**Do not use the device if it becomes contaminated;** discard and replace it with a newly packaged device.

This device is for anterograde placement only. Not intended for retrograde placement.

InHealth relies on the qualified, trained medical professional to instruct patients on the use and care of this device and to provide patients with the Blom-Singer instructions for use supplied with the product.

### **INDICATIONS**

The Blom-Singer Indwelling TEP Occluder is indicated for anterograde placement in and maintenance of the TE puncture following total laryngectomy when placement, or replacement, of the indwelling silicone prosthesis is performed by a qualified, trained medical professional. This device is for use with healed, intact tracheoesophageal puncture fistulas only.

### **CONTRAINDICATIONS**

The Blom-Singer Indwelling TEP Occluder is a medical product that should be used only by a qualified, medical professional with experience and training in its use and care.

---

This device should not be placed in a small tracheostomal airway if its presence obstructs breathing. This device is not a valved voice prosthesis and does not function like one. Tracheoesophageal speech is not possible with an Indwelling TEP Occluder.

Placement of the Indwelling prosthesis is contraindicated in a TE puncture that has dilated substantially beyond the diameter of a conventional 16 French or 20 French Blom-Singer Indwelling prosthesis.

## **HOW SUPPLIED**

The Indwelling TEP Occluder is available by special order in 16 French or 20 French diameters. It is supplied nonsterile; do not sterilize or place in an autoclave as damage to the device may occur. Provided in each package is the appropriately sized gel caps and inserters.

## **WARNINGS AND PRECAUTIONS**

The Blom-Singer Indwelling TEP Occluder is a medically prescribed device. Initial sizing and instructions for use must be provided by a qualified medical professional trained in the use of this particular device.

InHealth relies on the qualified, trained medical professional to instruct patients on the use and care of this device and to provide patients with the Blom-Singer instructions for use supplied with the product.

If the device is not functioning properly, the patient should have it evaluated by a qualified, trained medical professional soon as possible.

If there are tears, cracks, or structural damage to the device or accessories, discontinue use. Signs of a defective device include, but are not limited to,

- occluder degradation due to patient chemistry caused by microbial growth colonization on device;
- improper sizing caused by a compromised tracheoesophageal wall status;
- device pushed through TEP caused by improper use of inserter tool, improper sizing, or a compromised tracheoesophageal wall status;

- leakage around the device;
- tissue reaction or infection of puncture tract.

Handle the device and the gel cap loading tool carefully to avoid damaging the products.

Use only the specially sized gel caps provided by InHealth Technologies. The Gel Cap Insertion System is designed for use with Blom-Singer devices only.

Do not use petroleum-based products such as Vaseline® to lubricate the device as these products can damage silicone.

Changes in the anatomy or medical status of the patient may lead to improper fitting and/or function of the device. It is recommended that the device and TEP be evaluated by the clinician at regular intervals.

### **Device Dislodgement**

Care must be exercised during device insertion or removal to avoid injury to the TEP or accidental displacement of the device which could result in aspiration of the device into the trachea. Should aspiration occur, the patient should attempt to cough the device out of the trachea. Further medical attention may be necessary if coughing the device out is unsuccessful.

Confirm gel cap dissolution and deployment of the esophageal retention flange per procedures provided to ensure that the device is securely retained in the TEP.

If the device is accidentally dislodged from the puncture, the patient should immediately place a Blom-Singer puncture dilator, of the appropriate diameter, in the puncture to keep it from closing and leaking fluids. The device should be reinserted within 24 hours.

Foreign objects should never be inserted into the device by the patient.

### **Management of Indwelling TEP Occluder During Planned or Unplanned Airway Events**

During a planned (e.g., elective) or unplanned (e.g., emergency) airway event that involves stomal manipulation, the device should not be removed.

---

If the device becomes dislodged, it should be replaced as soon as possible with an all-purpose catheter of a compatible diameter to prevent aspiration of esophageal contents, and to maintain the opening between the trachea and the esophagus so that an Indwelling occluder can be inserted at a later time.

Standard resuscitation procedures for neck breathers should be utilized in a medical emergency.

Please instruct your patient to consider obtaining a medical alert bracelet or card to identify them as a neck breather. This should be kept in a clearly visible location.

### **Device Removal**

The attached neck strap on an Indwelling occluder should only be removed after the esophageal retention flange position within the esophagus has been verified. Never attempt to insert or reinsert an Indwelling occluder that has the neck strap removed.

Never remove one Indwelling occluder and insert another occluder without first dilating the TE puncture and remeasuring the tract to confirm correct device length.

The patient should never attempt to remove the device or allow anyone other than a qualified, trained medical professional to remove the device.

The Indwelling occluder is not a permanent device and requires replacement periodically. Many clinicians recommend routine follow-up evaluations not longer than 6 months. InHealth Technologies makes no warranty, either expressed or implied, as to the lifetime of the Indwelling TEP Occluder. The device lifetime may vary with individual use and biological conditions.

Removal of the Indwelling occluder should only be done by grasping the outer tracheal flange of the device securely with a locking hemostat. Pull gently and firmly until the device is fully withdrawn.

**Caution:** Federal (USA) law restricts this device to sale by or on the order of a physician or a licensed speech pathologist.

## **COMPLICATIONS**

Although rare, the following complications have been identified to occur with silicone devices of the Blom-Singer type. They include:

stoma contamination or sepsis, which may require removal of the device and/or appropriate antibiotics; accidental aspiration of the device into the airway, which may require removal by a physician; occasional extrusion of the device, requiring replacement after dilation of the TE puncture and additional supervision of the stoma care regimen; puncture dilation resulting in leakage of fluids around the device; inflammatory reaction around the puncture site and formation of granulation tissue; dislodgment of the device and subsequent closure of the TE puncture; intractable leakage around the device, requiring surgical revision or closure of the puncture; difficulty swallowing; tearing or other damage to the device from improper use; accidental ingestion of the device into the esophagus and/or GI tract.

## **PRODUCT DESCRIPTION**

The Blom-Singer Indwelling TEP Occluder is offered nonsterile. Nonsterile product is not intended for primary placement procedures.

The Blom-Singer Indwelling TEP Occluder is made of silicone, a biocompatible material used in many medical products. This device is latex-free.

The gel cap used for insertion of the device is made of a vegetable base gelatin. Individual patient reactions to these materials may occur.

The Blom-Singer Indwelling Occluder, (diagram 1) consists of an esophageal retention flange (diagram 1a), a tracheal retention flange (diagram 1b), and a neck strap (diagram 1c).

**Note:** Since this is a non-speaking device, cleaning of the device is not necessary. However, on the portion of the device that shows out of the puncture, the patient may need to remove dried mucus secretions using a tweezer.

In order to manage dryness and secretions of the airway, a heat and moisture exchange system is recommended. If dry secretions persist and require removal, it is recommended that the area of crusting be humidified to prevent trauma to the tracheal tissue.

## **INSTRUCTIONS FOR USE**

*The following procedure is provided by Eric D. Blom, Ph.D. Please refer to the diagrams located at the back of this instruction manual.*

---

Use of gloves, protective eye wear, and a respiratory face mask is recommended for the following procedure.

The Indwelling Occluder and its accessories should always be thoroughly clean and dry to help avoid introducing contaminants into the tracheostoma, TE puncture, or esophagus.

**Note:** Insertion or removal of the device in the TE puncture should only be done with a bright light focused directly on the stoma.

Have the following Blom-Singer items prepared for aseptic insertion of the device:

- Blom-Singer Indwelling Occluder
- Blom-Singer Gel Cap Insertion System
- Blom-Singer Dilator-Sizer
- Blom-Singer Voice Prosthesis Sizer (or Dilator-Sizer)
- Blom-Singer Puncture Dilator (or Dilator-Sizer)
- Blom-Singer Prosthesis Inserter

## Dilate and Measure

Please refer to the Blom-Singer Dilation/Sizing System Instructions For Use for complete product details.

**Optional:** *The Blom-Singer Sizer or Puncture Dilator may be used as an alternate procedure. Refer to the instructions for use supplied with the product.*

**1.** Dilate the TEP with a Blom-Singer Dilator-Sizer of the appropriate diameter. Gently insert the Dilator-Sizer into the puncture completely to the neck strap (diagram 2) and leave in place for five minutes.

**2.** Carefully withdraw the Dilator-Sizer until very gentle resistance is detected indicating that the esophageal wings are seated against the anterior wall of the esophagus (diagram 2a). Note: Pulling too firmly on the Dilator-Sizer may squeeze or compress the TE tract resulting in an erroneous measurement that indicates it is shorter than it actually is.

**3.** Note the marker that is closest to the entrance of the TE puncture. Read the size number on the Dilator-Sizer and select a Blom-Singer Indwelling Occluder of corresponding size. If the

measurement is between two markers, select the longer one.

- 4.** Fully insert the Dilator-Sizer and fasten the neck strap with adhesive tape on to the patient's neck. Leave the Dilator-Sizer in the puncture to dilate 5 minutes, or longer if necessary, in order to allow for easy insertion of the selected Indwelling Occluder. Prior to insertion of the Blom-Singer Indwelling Occluder, remove the tape from the neck strap, grasp the dilator securely, and carefully withdraw it from the puncture. Immediately insert the Indwelling Occluder.

## **Device Preparation**

- 1.** Prepare the device for easy,atraumatic insertion by utilizing the Blom-Singer Gel Cap Insertion System. Use only the specially sized gel caps provided with this product.

The Gel Cap Insertion System is designed for use with Blom-Singer Prostheses only; do not use with any other type devices. The following steps describe the hand folding method of gel cap insertion.

- a.** Remove the device from the package. Always handle the device with clean, gloved hands. Be certain that your gloves and the device are completely dry before loading the gel cap. Failure to do so may cause premature dissolving of the gel cap.
- b.** Remove the Gel Cap Insertion System items from the package: gel cap, flange introducer, and inserter tool.
- c.** Fold and load the device into the unprinted gel cap (discard the printed end) according to the following steps.
- d.** Fold the esophageal retention flange tightly in half toward the center of the device (diagram 4a).
- e.** Fold the two outer edges of the esophageal retention flange tightly against the folded flange (diagram 4b) and hold the device in this folded position.
- f.** Insert the folded portion of the esophageal retention flange into the gel cap (diagram 4c). Gently tuck and push the folded device as far as possible into the gel cap.
- g.** Use the flange introducer to gently tuck the remaining portion of the folded esophageal retention flange completely

---

into the gel cap (diagram 4d). Do not use any sharp or serrated instruments to avoid damaging the device.

- h.** The esophageal retention flange should be folded toward the center of the device and completely enclosed inside the gel cap (diagram 4e).

## Device Insertion

**Warning:** Caution should be exercised to avoid accidental displacement of the device which could result in aspiration of the device into the trachea. Should aspiration occur, the patient should attempt to cough the device out of the trachea. Further medical attention may be necessary if coughing the device out is unsuccessful.

- 1.** Place the device on the inserter and lock the neck strap onto the safety peg (diagram 3).
- 2.** Caution the patient not to swallow as you remove the previously inserted Dilator-Sizer to avoid aspiration of saliva into the trachea prior to inserting the device.
- 3.** Apply a light coating of water-soluble lubricant (oil-free) to the tip of the gel-capped end of the device and immediately insert the device completely into the TE puncture by aligning the tip of the device partially in the puncture with the neck strap oriented upwards (diagram 5). **Note:** Do not use petroleum-based products such as Vaseline® to lubricate the device as these products can damage silicone.
- 4.** Hold the device in this position of complete insertion for at least three minutes. This allows ample time for the gel cap to dissolve and release the esophageal retention flange (diagram 1a) within the esophagus. Patients should be instructed to swallow their saliva or warm water to facilitate dissolving the gel cap.

**Caution:** If the device does not insert easily on the first attempt, do not continue to try to insert it. Instead, insert a clean puncture dilator, of the appropriate diameter, for a few minutes to dilate the TE tract, and then try again.

- 5.** Detach the device neck strap (diagram 1c) from the safety peg on the inserter. Place a finger against the neck strap and carefully withdraw the inserter from the device with a twisting motion. **Warning:** Caution should be exercised when removing

the inserter to avoid accidental displacement of the device which could result in aspiration of the device. Secure the device with adhesive tape across the neck strap. It is recommended to orient the device with the neck strap up.

## NECK STRAP DETACHMENT

The Indwelling Occluder is designed to include the optional detachment of the neck strap following confirmation that the device esophageal retention flange is **securely** positioned within the esophageal lumen against the anterior wall of the esophagus. Removal of the neck strap is an **optional** procedure.

**The neck strap should not be detached from the Indwelling Occluder if:**

- use of a laryngectomy tube may potentially dislodge the device;
- a patient has previously been wearing a device that is larger in diameter than the Indwelling Occluder currently being inserted;
- the diameter of the TE puncture is dilated significantly larger than the device diameter;
- the health of the stoma tissue or physical condition of the patient is in question; or,
- direct flexible endoscopy is not utilized for confirmation of esophageal flange deployment.

**Note:** It is important **not** to detach the neck strap from the device in the instances cited above. Tape the neck strap to the peristomal skin and observe security of device retention for an extended period of time.

*The following procedure for confirming deployment of the esophageal retention flange is recommended by Eric D. Blom, Ph.D.*

After insertion of the device in the puncture tract and completion of the 3-minute waiting period for gel cap dissolution and esophageal retention flange deployment, repeatedly rotate the device on the inserter while it is in the puncture. If the esophageal retention flange has deployed and has properly seated against the mucosa of the anterior esophageal wall, the device on the inserter will rotate easily in the puncture tract. If the gel cap has

---

dissolved and incorrectly deployed the esophageal retention flange inside the puncture tract rather than inside the lumen of the esophagus, considerable resistance to device rotation will be felt.

A second method for confirming esophageal retention flange deployment is to utilize direct flexible endoscopy confirmation of deployment.

If the assessment methods described above fail to provide satisfactory verification of correct esophageal retention flange deployment, remove the device, dilate and remeasure the puncture tract, and repeat the device insertion and confirmation process.

Once deployment of the esophageal retention flange has been confirmed, detach the neck strap from the tracheal retention flange (diagram 1b) at the area of reduced neck strap width, where it meets the tracheal flange, by cutting carefully with scissors (diagram 6).

**Note:** if you are uncertain that the esophageal retention flange has properly unfolded, and the device is in place inside the esophageal lumen, after using the steps described above, do **not** detach the neck strap from the device. Securely tape the neck strap to the neck and observe the patient for a few days for retention ability.

## **Removal of the Device**

The Indwelling Occluder is not a permanent device and requires replacement periodically. Many clinicians recommend routine follow-up evaluations not longer than 6 months. InHealth Technologies makes no warranty, either expressed or implied, as to the lifetime of the Indwelling TEP Occluder. The device lifetime may vary with individual use and biological conditions.

**Removal of the Indwelling Occluder should only be done by grasping the outer tracheal flange of the device securely with a locking hemostat.** Pull gently and firmly until the device is fully withdrawn.

Insert a Dilator-Sizer of the appropriate diameter, and tape it in position for five minutes prior to inserting a new Indwelling Occluder.

**Warning:** Clearly inform the patient **never** to attempt to remove the device or allow anyone other than a qualified, trained medical professional to remove the Indwelling Occluder.

**Never** remove one device and insert another device without first dilating the TE puncture and remeasuring the tract to confirm device length.

## **SPECIAL ORDER PRODUCTS**

If this instruction manual accompanies a Special Order Product, there may be differences in the physical characteristics between the product and the product descriptions described in this instruction manual. These differences will not affect the safety or efficacy of the special order product. Special Order products are nonreturnable.

## **ORDERING INFORMATION**

### **USA**

Blom-Singer products may be ordered directly from InHealth Technologies. TELEPHONE: Toll-Free (800)477-5969 or (805)684-9337, Monday–Friday, 9:30am–7:00pm, Eastern Standard Time. FAX: Toll-Free (888)371-1530 or (805)684-8594. EMAIL: [order@inhealth.com](mailto:order@inhealth.com) ORDER ON-LINE: [www.inhealth.com](http://www.inhealth.com) POST: InHealth Technologies, 1110 Mark Avenue, Carpinteria, CA 93013-2918, USA, Attention: Customer Service

### **International**

Please contact our customer service department for a list of international distributors.

### **Consumer Affairs**

If you have any questions or dissatisfaction with a product please contact our customer service department via telephone, fax, post, or Email: [productcomplaints@inhealth.com](mailto:productcomplaints@inhealth.com)

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






Additional references available upon request.

Weitere Referenzen aus der Literatur sind auf Anfrage erhältlich. / Se pueden solicitar más referencias. / Autres références disponibles sur demande. / Riferimenti aggiuntivi sono disponibili a richiesta. / Meer informatie is op aanvraag beschikbaar. / Referências adicionais mediante solicitação. / Další odkazy na literaturu jsou k dispozici na požádání. / Dodatkowe informacje dostępne na żądanie.

1. Blom, E.D.: "Some Comments on the Escalation of Tracheoesophageal Voice Prosthesis Dimensions," **Arch Otolaryngol**, 129: 500-502, 2003.
2. Eerenstein, S.E., Grolman, W., Schouwenburg, P.F.: "Downsizing of Voice Prosthesis Diameter in Patients with Laryngectomy," **Arch Otolaryngol**, 128: 838-841, 2003.
3. Hamaker, R.C., Blom, E.D.: "Botulinum Neurotoxin for Pharyngeal Constrictor Muscle Spasm in Tracheoesophageal Voice Restoration," **Laryngoscope**, 113:1479-1482, 2003.
4. Evans, P.H.R., Blom, E.D.: "Functional Restoration of Speech," **Principles and Practice of Head and Neck Oncology**, Evans, P.H.R., Montgomery, P.W., Gullane, P.J. (Eds.) Taylor and Francis Group, London, 2003.
5. Blom, E.D., Singer M.I., and Hamaker, R.C., **Tracheoesophageal Voice Restoration Following Total Laryngectomy**. San Diego: Singular Publishing Group, Inc.1998.
6. Blom, E.D., Hamaker, R.C.: "Tracheoesophageal Voice Restoration Following Total Laryngectomy," **Cancer of the Head and Neck**, E.N. and Suen, J. (Eds). W.B. Saunders Publishers. 839-852, 1996.
7. Blom, E.D., "Tracheoesophageal Speech," **Seminars in Speech and Language**. New York: Thieme Medical Publishers. 16:3, 191-204, 1995.
8. Pauloski, B.R., Blom, E.D., Logemann, J.A., Hamaker, R.C. "Functional Outcome After Surgery for Prevention of Pharyngospasms in Tracheoesophageal Speakers. Part II: Swallow Characteristics," **Laryngoscope**, 105:10: 1104-1110, 1995.
9. Blom, E.D., Pauloski, B.R., Hamaker, R.C. "Functional Outcome After Surgery for Prevention of Pharyngospasms in Tracheoesophageal Speakers. Part I: Speech Characteristics," **Laryngoscope**, 105:10: 1093-1103, 1995.
10. Blom, E.D., Hamaker, R.C., Freeman, S.B.: "Postlaryngectomy Voice Restoration," Lucente, F.E. (Ed.): **Highlights of the Instructional Courses**, Vol. 7. St. Louis: Mosby-Year Book. pp 3-10, 1994.

11. Blom, E.D., Singer M.I., and Hamaker, R.C. "A Prospective Study of Tracheoesophageal Speech," **Arch Otolaryngol**, 112: 440-447, 1986.
12. Blom, E.D., Singer, M.I. "Disinfection of Silicone Voice Prostheses," **Arch Otolaryngol**, 112: 1986.
13. Blom, E.D., Singer M.I., and Hamaker, R.C. "An Improved Esophageal Insufflation Test," **Arch Otolaryngol**, 111: 211-212, 1985.
14. Hamaker, R.C., Singer M.I., Blom E.D. and Daniels, H.A. "Primary Voice Restoration at Laryngectomy," **Arch Otolaryngol**, 111: 182-186, 1985.
15. Singer, M.I. and Blom, E.D. "A Selective Myotomy for Voice Restoration After Total Laryngectomy," **Arch Otolaryngol**, 107: 670-673, 1981.
16. Singer, M.I. and Blom, E.D. "An Endoscopic Technique for Restoration of Voice After Laryngectomy," **Ann Otol Rhinol Laryngol**, 89:6: 529-533, 1980.

# **INHEALTH®**

## TECHNOLOGIES

---

### **UNITED STATES**

1110 Mark Avenue  
Carpinteria, CA 93013-2918 USA  
Toll-Free: 800.477.5969  
Toll-Free Fax: 888.371.1530

### **INTERNATIONAL**

Tel: +1.805.684.9337  
Fax: +1.805.684.8594  
E-mail: [order@inhealth.com](mailto:order@inhealth.com)

### **ORDER ONLINE**

[www.inhealth.com](http://www.inhealth.com)

**Freudenberg Medical, LLC**  
1110 Mark Avenue  
Carpinteria, CA 93013-2918 USA  
[www.inhealth.com](http://www.inhealth.com)

**CE 0344**

**EMERGO EUROPE**  
Prinsessegracht 20  
2514 AP The Hague  
The Netherlands