INHEALTH

TECHNOLOGIES®

Blom-Singer® Indwelling TEP Occluder

INSTRUCTIONS FOR USE

Blom-Singer®

voice restoration systems

37720-01H Effective May 2020 / 2020 年 5 月生效 / 2020年5月発効 / Në fuqi nga maji 2020 / 2020 اعتبارًا من مايو / В сила от май 2020 г. / Stupa na snagu u svibnju 2020. / V platnosti od května 2020 / Gældende fra maj 2020 / Met ingang van mei 2020 / Kehtiv alates maist 2020 / Voimassa toukokuusta 2020 alkaen / Date d'entrée en vigueur : mai 2020 / Gültig ab Mai 2020 / Σε ισχύ από τον Μάιο του 2020 / 2020 בתוקף החל ממאי / hatályba lépés 2020. május / Tekur gildi í maí 2020 / Valido da maggio 2020 / 2020년 5월부터 발효 / Spēkā no 2020. gada maija / Įsigalioja 2020 m. gegužės mėn / Gjelder fra mai 2020 / Obowiązuje od maja 2020 / Vigente em maio de 2020 / Cu intrare în vigoare în mai 2020 / Версия от мая 2020 г. / Na snazi od maja 2020. / Účinné k máju 2020 / Efectivo a partir de mayo de 2020 / Gäller från maj 2020 / Geçerlilik Tarihi Mayıs 2020

Blom-Singer® and InHealth Technologies® are registered trademarks of Freudenberg Medical, LLC.

Patent(s): www.inhealthpatents.com







Copyright © 2020 Freudenberg Medical, LLC. All rights reserved.

INSTRUCTIONS FOR USE

使用说明 / 使用說明書 / Udhëzime për përdorim / إرشادات الاستخدام / Инструкции за употреба / Upute za korištenje / Návod k použití / Brugsanvisning / Gebruiksaanwijzing / Kasutusjuhend / Käyttöohjeet / Instructions d'utilisation / Gebrauchsanweisung / Οδηγίες χρήσης / הוראות שימוש / Használati utasítás / Notkunarleiðbeiningar / Istruzioni per l'uso / 사용 지침 / Lietošanas norādījumi / Naudojimo instrukcijos / Bruksanvisning / Instrukcja użytkowania / Instruções de uso / Instrucțiuni de utilizare / Инструкции по применению / Uputstvo za upotrebu / Návod na použitie / Instrucciones de uso / Bruksanvisningar / Kullanım Talimatlar

DIAGRAMS

插图 / 図 / Diagramet / الرسم التوضيحي / Схеми / Dijagrami / Schémata / Diagrammer /
 Afbeeldingen / Diagrammid / Kaaviot / Diagrammes / Abbildungen / Διαγράμματα / תרשימים /
 Ábrák / Skýringarmyndir / Diagrammi / 그림 / Diagrammas / Diagramos / Diagrammer /
 Schematy / Diagramas / Diagramme / Схеми / Crteži / Grafické znázornenia / Diagramas /
 Diagram / Diyagramlar





TABLE OF CONTENTS

目录 / 目次 / Tabela e përmbajtjes / جدول المحتويات / Съдържание / Sadržaj / Obsah /
 Indholdsfortegnelse / Inhoud / Sisukord / Sisällysluttelo / Table des matières /
 Inhaltsverzeichnis / Πίνακας περιεχομένων / תוכן עניינים / Tartalomjegyzék / Efnisyfirlit /
 Sommario / 목차 / Saturs / Turinys / Innholdsfortegnelse / Spis treści / Índice / Cuprins /
 Содержание / Sadržaj / Obsah / Índice / Innehållsförteckning / İçindekiler

EN	Indwelling TEP Occluder	6	HU	Állandó TEP okkluder	74
ZH	留置式气管食管穿刺封堵器	10	IS	Eiginleg loka fyrir barka- og vélandaásthungu	78
JA	留置型TEPオクルーダー	13	IT	Occlusore TEP fisso	82
SQ	Tapa e TEP e fiksuar	18	KO	인체삽입형 TEP 폐색기	87
22	مطبّق البزل الرغامي المريني المستقر	AR	LV	levietojamais TEP aizsēgs	91
BG	Импантируем ТЕП оклудър	26	LT	[vidų implantuojamas gerklės-stemplės punkcijos vietos okliuzorius	95
HR	Trajni TEP okluder	30	NO	TEP-tillukning til implantasjon	99
CS	Nástroj pro uzavření vnitřní TEP	34	PL	Okluder stały NTP	103
DA	Indlagt TEP-okkluderingsanordning	39	PT	Oclusor de PTE Fixo	107
NL	Afsluitbare TEP-verblijfstemprothese	43	RO	Obturator TEP permanent	112
ET	Sisemine TEP-i sulgur	48	RU	Постоянный obturator ТПП	117
FI	Pysyvä TEP-suilin	52	SR	Trajni TEP okluder	122
FR	Dispositif d'occlusion PTO à demeure	56	SK	Zavedený oklúder TEP	126
DE	TEP-Verweil-Okkluder	60	ES	Dispositivo de oclusión fijo para TEP	130
EL	Μόνιμο εργαλείο απόφραξης ΤΟΠ	65	SV	Kvarliggande TEP-ockluderare	134
70	תפנית TEP	HE	TR	Kalıcı TEP Tikacı	138

ENGLISH

BLOM-SINGER® INDWELLING TEP OCCLUDER**PRODUCT DESCRIPTION**

Please refer to the diagrams located at the front of this instruction manual.

The Blom-Singer Indwelling TEP (Tracheoesophageal Puncture: opening made by a surgeon between trachea or “windpipe” and esophagus for placement of a voice prosthesis) Occluder is a nonsterile, non-valved silicone device designed for placement in and maintenance of the TEP following total laryngectomy. The device consists of a large esophageal flange, a tracheal flange, and a safety strap. The device is available in 16 Fr or 20 Fr diameters. Two (2) gel caps of the appropriate size, one (1) flange introducer, and one (1) inserter tool are also provided.

INDICATIONS (Reasons to prescribe the device or procedure)

The device is indicated for secondary placement in and maintenance of the TEP following total laryngectomy when placement, or replacement, of the indwelling silicone prosthesis is performed by a qualified, trained medical professional. This device is for use with healed, intact tracheoesophageal puncture fistulas only. This device is for anterograde placement only. The device helps to prevent leakage of esophageal contents into the trachea during swallowing.

CONTRAINDICATIONS (Reasons that make it inadvisable to prescribe the particular device or procedure)

The device is a medical product that should be used only by a qualified clinician with experience and training in its use and care. This device should not be placed in a small tracheostomal airway if its presence obstructs breathing. This device is not a valved voice prosthesis and does not function like one. Tracheoesophageal speech is not possible with an Indwelling TEP Occluder. Placement of the device is contraindicated in a TEP that has dilated substantially beyond the diameter of a conventional 16 Fr or 20 Fr Blom-Singer indwelling prosthesis.

WARNINGS AND PRECAUTIONS

Initial sizing and instructions for use must be provided by a qualified clinician trained in the use of this particular device. The clinician is to instruct patients on the use and care of this device and to provide patients with these instructions for use supplied with the product. This device is for single use only; reuse of this device could result in infection. If the device is not functioning properly, the patient should have it evaluated by a clinician as soon as possible. The device and accessories must be carefully removed from its packaging and handled in a manner to prevent contamination of the device. **Do not use the device if it becomes contaminated; discard and replace it with a newly packaged device.** Handle the device and accessories carefully to avoid damage. If there are tears, cracks, or structural damage to the device or accessories, discontinue use and contact InHealth Technologies Product Complaints. Do not use the device in the event of the packaging being damaged or unintentionally opened prior to use; discard and replace it with a new device. Do not sterilize or place in an autoclave as damage to the device may occur. Use only the appropriately-sized gel caps. Do not use petroleum-based products such as Vaseline® to lubricate the device as these can damage silicone. Changes in the anatomy or medical status of the patient may lead to improper fitting and/or function of the device. Dislodgement or extrusion of the device from the TEP may occur and result in ingestion, aspiration (inhalation) or tissue damage. A foreign body such as an Indwelling Occluder in the airway may cause complications such as acute respiratory (breathing) distress and/or respiratory arrest (stop breathing). Always measure the TEP tract when changing the device to select the proper device size. If the device is improperly sized, it may cause tissue necrosis (tissue death) and extrusion (device to protrude out of TEP). Individual patient reactions to the device materials may occur. Consult a clinician immediately if there are indications of tissue edema (swelling) and/or inflammation/infection. Evaluate patients with bleeding disorders or if they are undergoing anticoagulant (blood thinning) treatment for the risk of hemorrhage (bleeding) prior to placement or replacement of the device.

Device Dislodgement

Care must be exercised during device insertion or removal to avoid injury to the TEP or accidental displacement of the device, which could result in aspiration (inhalation) of the device into the trachea (windpipe). Should aspiration occur, the patient should attempt to cough the device out of the trachea. Further medical attention may be necessary if coughing the device out is unsuccessful. Confirm gel cap dissolution and deployment of the esophageal retention flange per procedures provided to ensure that the device is securely retained in the TEP. If the device is accidentally dislodged from the TEP, the patient or clinician should immediately place a Blom-Singer Puncture Dilator or suitable device of the appropriate diameter, in the TEP to keep it from closing and leaking fluids. A replacement device should be inserted within 24 hours. Foreign objects should never be inserted into the device by the patient. Inserting foreign objects may cause dislodgement of the device or its components and may cause these objects to be swallowed or aspirated.

Device Insertion/Removal

The attached safety strap on the device should only be removed after the esophageal flange has been verified to have deployed in the esophagus. Never attempt to insert or reinsert a device with the safety strap removed. Never remove one Indwelling Occluder and insert another one without first dilating the TEP and re-measuring the tract to confirm correct device length. The patient should never attempt to insert or remove the device or allow anyone other than a qualified, trained clinician to insert or remove the device. Removal of the device should only be done by grasping the tracheal flange of the device securely with a locking hemostat. The device is not for permanent use and requires replacement periodically.

COMPLICATIONS

Although rare, the following complications may occur with the device. They include:

- Contamination, infection or sepsis, which may require removal of the device and/or appropriate antibiotics.
- Accidental aspiration of the device into the airway, which may require removal by a physician.
- Occasional extrusion of the device, requiring replacement after dilation of the TEP and additional supervision of the stoma (opening through neck into trachea or wind pipe) care regimen.
- Puncture dilation resulting in leakage of fluids around the device.
- Inflammatory reaction around the TEP site and formation of granulation tissue.
- Dislodgment of the device and subsequent closure of the TEP.
- Intractable (uncontrollable) leakage around the device, possibly requiring surgical revision or closure of the TEP.
- Dysphagia (difficulty swallowing).
- Tearing or other damage to the device from improper use.
- Accidental ingestion of the device into the esophagus and/or GI tract.

INSTRUCTIONS FOR USE

NOTE: THIS DEVICE IS NOT INTENDED TO BE PLACED OR REMOVED BY THE PATIENT. All instructions related to device preparation, insertion, and removal are for the clinician only. If you are a patient, refer to the “Device Care” section.

The following procedure is provided by Eric D. Blom, Ph.D. Please also refer to the diagrams located at the front of this instruction manual. **NOTE:** Measure the TEP tract before inserting the device to confirm correct device length. **NOTE:** Insertion or removal of the device in the TEP must only be done using a bright light to illuminate the tracheostoma and TEP. **NOTE:** Dilate TEP before secondary placement of the device. Use of gloves, protective eye wear, and a respiratory face mask is recommended for the following procedure. The device and its accessories should always be thoroughly clean and dry to help avoid introducing contaminants into the tracheostoma, TEP, or esophagus.

Dilate and Measure (Diagrams 1 & 2)

Please refer to the Blom-Singer Dilator/Sizing System Instructions for Use for complete product details.

Device Preparation

Prepare the device for insertion by utilizing the hand folding method of gel cap insertion. Use only the appropriately sized gel caps provided with this product. Handle the device with clean, gloved hands. Remove the device and insertion accessories from the package: gel caps, flange introducer, and inserter tool. Be certain that your hands/gloves and the device are completely dry before loading the gel cap. Failure to do so may cause premature dissolving of the gel cap.

1. Use only the short end of the gel cap. Discard the gel cap body (longer end of gel cap). (Diagram 3)
2. Fold the esophageal flange tightly in half toward the center of the device (Diagram 4).
3. Fold the outer edges of the folded flange tightly against each other (Diagram 5) and hold the device in this folded position.
4. Insert the folded portion of the esophageal flange into the gel cap (Diagram 6). Gently tuck and push the folded esophageal flange as far as possible into the gel cap.
5. Use the curved side of the flange introducer and following the curve of the device body, gently tuck the remaining portion of the folded esophageal flange completely into the gel cap (Diagram 7). **Do not use any sharp or serrated instruments to avoid damaging the device/accessories.**
6. The esophageal flange should be folded in a forward position and completely enclosed inside the gel cap (Diagram 8).

Device Insertion

1. Place the device on the inserter stick and lock the safety strap onto the safety peg (Diagram 9).
2. Caution the patient not to swallow as you remove the Dilator-Sizer to avoid aspiration of saliva into the trachea prior to inserting the device.
3. Apply a light coating of water-soluble lubricant (not petroleum-based) to the tip of the gel-capped end of the device and insert the device completely into the TEP with the safety strap oriented upwards. Insert the device completely into the dilated TEP until the tracheal flange of the device is seated firmly against the posterior tracheal mucosa.
4. Hold the device in this position of complete insertion for at least three minutes. This allows ample time for the gel cap to dissolve and release the esophageal flange within the esophagus. Patients should be instructed to swallow their saliva to facilitate dissolving the gel cap.
Caution: If the device does not insert easily on the first attempt, do not continue to try to insert it. Instead, re-insert the puncture dilator of the appropriate diameter for a few minutes to dilate the TEP, and then try again.
5. Rotate the device on the inserter while it is in the TEP. The device will rotate easily on the inserter in the TEP if the esophageal retention flange has deployed and seated against the anterior esophageal wall. Considerable resistance to device rotation will be felt if the gel cap has not dissolved or incorrectly deployed inside the TEP rather than inside the lumen of the esophagus. If indicated, direct visualization of esophageal flange deployment can be obtained with flexible endoscopy. If the assessment methods described above fail to provide satisfactory verification of correct esophageal retention flange deployment, remove the device, dilate and re-measure the TEP, and repeat the device insertion and confirmation process.

Safety Strap Detachment

The device is designed to allow the optional detachment of the safety strap following confirmation that the esophageal flange is securely positioned within the esophageal lumen against the anterior wall of the esophagus. Removal of the safety strap is an **optional** procedure. Once deployment of the esophageal flange has been confirmed:

1. Detach the device neck strap from the safety peg on the inserter. Place a finger against the neck strap and carefully withdraw the inserter from the device with a twisting motion.
Warning: Caution should be exercised when removing the inserter to avoid accidental displacement of the device which could result in aspiration of the device.
2. Detach the safety strap from the tracheal flange at the area of reduced strap width, where it meets the tracheal flange, by first applying locking hemostats onto the strap and then cutting carefully with scissors (Diagram 10).

The safety strap should not be detached from the Indwelling Occluder if: Use of a laryngectomy tube or button may potentially dislodge the device; a patient has previously been wearing a device that is larger in diameter than the device currently being inserted; the diameter of the TEP is dilated significantly larger than the device diameter; the health of the stoma tissue or physical condition of the patient is in question. **Note: It is important not to detach the neck strap from the device in the instances cited above. Tape the safety strap to the peristomal skin and confirm security of the device retention.**

Removal of the Device

This device is intended to be used for 6 months but this duration should be determined per clinical evaluation and appropriate consultation with a skilled practitioner. Removal of the device should only be done by grasping the tracheal flange of the device securely with a locking hemostat. Pull firmly using caution until the device is fully withdrawn. Insert an appropriately sized dilator-sizer after the device is removed and tape it in position for five minutes prior to inserting a replacement device.

Device Care (for Patient)

Since this is a non-valved device, cleaning is not necessary. Careful removal of dried mucus from the tracheal flange may be completed as instructed by your clinician.

SPECIAL STORAGE AND/OR HANDLING CONDITIONS

Only the gel caps packaged with the device have special storage and/or handling conditions. They must be kept dry and out of direct sunlight.

DISPOSAL INSTRUCTIONS

This product is not biodegradable and can be considered contaminated when used. Carefully dispose of this device per local guidelines.

SPECIAL ORDER PRODUCTS

These devices are Special Order products and are nonreturnable.

ORDERING INFORMATION USA

Blom-Singer products may be ordered directly from InHealth Technologies.

ORDER ON-LINE: www.inhealth.com TELEPHONE: Toll-Free (800)477-5969 or (805)684-9337, Monday–Friday, 9:30am–7:00pm, Eastern Standard Time. FAX: Toll-Free (888)371-1530 or (805)684-8594. EMAIL: order@inhealth.com POST: InHealth Technologies, 1110 Mark Avenue, Carpinteria, CA 93013-2918, USA, Attention: Customer Service

International

Please contact our customer service department for a distributor referral.

PRODUCT COMPLAINTS/EU SERIOUS INCIDENTS

If you are dissatisfied with the device or have any questions, please contact productcomplaints@inhealth.com.

Telephone: +1-800-477-5969

Fax: +1-888-371-1530

Any serious incident that has occurred in relation to the device should be reported to Freudenberg Medical, LLC as noted above and the competent authority of the EU Member State in which the user and/or patient is established.

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SYMBOLS GLOSSARY (USA only)

SYMBOL	STANDARD No.	SYMBOL TITLE & REF No.	MEANING OF SYMBOL
	EN ISO 15223-1:2016	Manufacturer; 5.1.1	Indicates the medical device manufacturer, as defined in EU Directives or Regulations
	93/42/EEC EU 2017/745	CE Marking; 93/42/EEC Annex XII, EU 2017/745 Annex V	CE marking of conformity
	EN ISO 15223-1:2016	Authorized representative in the European Community; 5.1.2	Indicates the authorized representative in the European Community
	EN ISO 15223-1:2016	Catalogue number; 5.1.6	Indicates the manufacturer's catalogue number so that the medical device can be identified
	EN ISO 15223-1:2016	Batch code; 5.1.5	Indicates the manufacturer's batch code so that the batch or lot can be identified
	EN ISO 15223-1:2016	Use-by date; 5.1.4	Indicates the date after which the medical device is not to be used
	EN ISO 15223-1:2016	Caution; 5.4.4	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself
	EN ISO 15223-1:2016	Do not use if package is damaged; 5.2.8	Indicates a medical device that should not be used if the packaged has been damaged or opened
	EN ISO 15223-1:2016	Consult instructions for use; 5.4.3	Indicates the need for the user to consult the instructions for use
	EN ISO 15223-1:2016	Non-sterile; 5.2.7	Indicates a medical device that has not been subjected to a sterilization process
	EN ISO 15223-1:2016	Do not re-use; 5.4.2	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure

SYMBOLS GLOSSARY (USA only)

SYMBOL	STANDARD No.	SYMBOL TITLE & REF No.	MEANING OF SYMBOL
	EN ISO 15223-1:2016	Keep away from sunlight; 5.3.2	Indicates a medical device that needs protection from light sources
	EN ISO 15223-1:2016	Keep dry; 5.3.4	Indicates a medical device that needs to be protected from moisture

INHEALTH

TECHNOLOGIES®

UNITED STATES

1110 Mark Avenue
Carpinteria, CA 93013-2918 USA
Toll-Free: 800.477.5969
Toll-Free Fax: 888.371.1530

INTERNATIONAL

Tel: +1.805.684.9337
Fax: +1.805.684.8594
E-mail: order@inhealth.com

ORDER ONLINE

www.inhealth.com

Freudenberg Medical, LLC
1110 Mark Avenue
Carpinteria, CA 93013-2918 USA
www.inhealth.com

CE 0344

EMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands