INHEALTH

TECHNOLOGIES®

Blom-Singer® Shower Guard

INSTRUCTIONS FOR USE

Blom-Singer®

voice restoration systems

Blom-Singer® and InHealth Technologies® are registered trademarks of **Freudenberg Medical, LLC**.

Patent(s): www.inhealthpatents.com

IFU: <https://inhealth.com/ifu>

Copyright © 2020 Freudenberg Medical, LLC. All rights reserved.

INSTRUCTIONS FOR USE

使用说明 / 使用說明書 / Udhëzime për përdorim / إرشادات الاستخدام / Инструкции за употреба / Upute za korištenje / Návod k použití / Brugsanvisning / Gebruiksaanwijzing / Kasutusjuhend / Käyttöohjeet / Instructions d'utilisation / Gebrauchsanweisung / Οδηγίες χρήσης / הוראות שימוש / Használati utasítás / Notkunarleiðbeiningar / Istruzioni per l'uso / 사용 지침 / Lietošanas norādījumi / Naudojimo instrukcijos / Bruksanvisning / Instrukcja użytkowania / Instruções de uso / Instrucțiuni de utilizare / Инструкции по применению / Uputstvo za upotrebu / Návod na použitie / Instrucciones de uso / Bruksanvisningar / Kullanım Talimatlar

DIAGRAMS

插图 / 図 / Diagramet / الرسم التوضيحي / Схеми / Dijagrami / Schémata / Diagrammer /
 Afbeeldingen / Diagrammid / Kaaviot / Diagrammes / Abbildungen / Διαγράμματα / תרשימים /
 Ábrák / Skýringarmyndir / Diagrammi / 그림 / Diagrammas / Diagramos / Diagrammer /
 Schematy / Diagramas / Diagreame / Схеми / Crteži / Grafické znázornenia / Diagramas /
 Diagram / Diyagramlar



TABLE OF CONTENTS

目録 / 目次 / Tabela e përmbajtjes / جدول المحتويات / Съдържание / Sadržaj / Obsah /
 Indholdsfortegnelse / Inhoud / Sisukord / Sisällysluettelo / Table des matières /
 Inhaltsverzeichnis / Πίνακας περιεχομένων / תוכן עניינים / Tartalomjegyzék / Efnisyfirlit /
 Sommario / 목차 / Saturs / Turinys / Innholdsfortegnelse / Spis treści / Índice / Cuprins /
 Содержание / Sadržaj / Obsah / Índice / Innehållsförteckning / İçindekiler

EN	Shower Guard	5	HU	Zuhanyvédő	37
ZH	淋浴护罩	6	IS	Sturtuhlíf	39
JA	シャワーガード	8	IT	Protezione per doccia	41
SQ	Mbrojtësja për dush	11	KO	샤워 가드	43
13	واقى الاستحمام	AR	LV	Dušas aizsargs	45
BG	Shower Guard	15	LT	Dušo apsauga	47
HR	Štitnik pri tuširanju	17	NO	Dusjbeskyttelse	49
CS	Kryt pro sprchování	19	PL	Zabezpieczenie pod prysznic	51
DA	Badebeskytter	21	PT	Protetor para banho	53
NL	Douchebeveiliging	23	RO	Protecție pentru duș	55
ET	Dušikaitse	25	RU	Колпачок для защиты протеза в душе	57
FI	Suihkusuoja	27	SR	Štitnik pri tuširanju	59
FR	Protection pour la douche	29	SK	Sprchový chránič	61
DE	Duschschutz	31	ES	Protector para ducha	63
EL	Προστατευτικό για το μπάνιο	33	SV	Duschskydd	65
35	מגן מקלחת	HE	TR	Duş Koruyucu	67

ENGLISH

BLOM-SINGER® SHOWER GUARD**PRODUCT DESCRIPTION**

The Blom-Singer Shower Guard is designed to shield the tracheostoma when showering and washing your hair. When properly applied, it will permit you to breathe freely while preventing the shower spray and water that runs down the neck from entering the tracheostoma. It is supplied with one Blom-Singer Tracheostoma Valve Housing and 30 Blom-Singer Adhesive Tape Discs. The device is also compatible with the Blom-Singer StomaSoft Laryngectomy Tube.

INDICATIONS (Reasons to prescribe the device or procedure)

The device is designed for use by individuals who have a permanent opening in their neck (stoma) following total laryngectomy (removal of the voice box).

CONTRAINDICATIONS (Reasons that make it inadvisable to prescribe the device or procedure)

The device is designed only for use in protecting the airway opening (stoma) during showering and must be properly attached as described below. This device is not designed to protect the airway if the head and neck are submerged in water. It must not be used in any water sports.

WARNINGS AND PRECAUTIONS

The device is designed only for single-patient use in protecting the airway opening (tracheostoma) during showering and must be properly attached as described below. The device is not designed to protect the airway if the head and neck are submerged in water. The device must not be used in any water sports. If there are tears, cracks, or structural damage to the device, discontinue use.

COMPLICATIONS

The following complications may occur:

- The penetration of a small amount of shower spray or droplets of water into the airway (stoma) causing coughing, aspiration, or discomfort.
- Patient infection, airway or lung irritation, or tissue irritation.
- Drowning.
- Temporarily obstructed respiration.

INSTRUCTIONS FOR USE

The device should only be used by individuals who have been trained in its use and care.

Insertion, Removal and Cleaning

Before attaching the Shower Guard, a Blom-Singer tracheostoma housing or Blom-Singer StomaSoft Laryngectomy tube must first be in place. To place a housing or tube, follow the Instructions for Use provided with the selected product accessory.

Insertion

With clean hands, insert the Shower Guard partially into the lower rim of the previously applied tracheostoma housing or tube with the perforated breathing filter facing downward (diagram 1). Gently push the Shower Guard until it snaps into place. A properly inserted Shower Guard will fit securely inside the housing or tube. You should be able to breathe without difficulty with the Shower Guard in place.

Removal

The Shower Guard is designed for removal without also removing the tracheostoma housing or tube. Use one hand to hold the housing or tube in place. Grasp the Shower Guard by the outer rim with your other hand and pull to remove.

Cleaning

The tubular vent (diagram 2a) and the perforated breathing filter (diagram 2b) of the device can be disassembled for cleaning. Rinse with clean water and either air-dry or use a lint-free material to dry the device. Never use facial or toilet tissue to dry the device as particles of lint or fabric may be deposited onto the device that can be inhaled into the airway upon reinsertion. Once a device has been cleaned and dried completely, it should be reassembled and stored in a clean, resealable plastic bag until its next use.

SPECIAL STORAGE AND/OR HANDLING CONDITIONS

There are no special storage and/or handling conditions for this device.

DISPOSAL INSTRUCTIONS

This device is not biodegradable and can be contaminated when used. Carefully dispose of the device per local guidelines.

ORDERING INFORMATION

USA

Online Ordering www.inhealth.com

USA

InHealth Technologies

1110 Mark Avenue

Carpinteria, CA 93013-2918

Customer Service: +1-800-477-5969

Fax: +1-888-371-1530

Email: order@inhealth.com

INTERNATIONAL

Contact InHealth Technologies Customer Service at +1-805-684-9337 for international distributor contact.

PRODUCT COMPLAINTS/EU SERIOUS INCIDENTS

If you are dissatisfied with the device or have any questions, please contact

productcomplaints@inhealth.com.

Telephone: +1-800-477-5969

Fax: +1-888-371-1530

Any serious incident that has occurred in relation to the device should be reported to Freudenberg Medical, LLC as noted above and the competent authority of the EU Member State in which the user and/or patient is established.

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SYMBOLS GLOSSARY (USA only)

SYMBOL	STANDARD No.	SYMBOL TITLE & REF No.	MEANING OF SYMBOL
	EN ISO 15223-1:2016	Manufacturer; 5.1.1	Indicates the medical device manufacturer, as defined in EU Directives or Regulations
	93/42/EEC EU 2017/745	CE Marking; 93/42/EEC Annex XII, EU 2017/745 Annex V	CE marking of conformity
	EN ISO 15223-1:2016	Authorized representative in the European Community; 5.1.2	Indicates the authorized representative in the European Community
	EN ISO 15223-1:2016	Batch code; 5.1.6	Indicates the manufacturer's catalogue number so that the medical device can be identified
	EN ISO 15223-1:2016	Batch code; 5.1.5	Indicates the manufacturer's batch code so that the batch or lot can be identified
	EN ISO 15223-1:2016	Use-by date; 5.1.4	Indicates the date after which the medical device is not to be used
	EN ISO 15223-1:2016	Caution; 5.4.4	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself
	EN ISO 15223-1:2016	Do not use if package is damaged; 5.2.8	Indicates a medical device that should not be used if the packaged has been damaged or opened
	EN ISO 15223-1:2016	Consult instructions for use; 5.4.3	Indicates the need for the user to consult the instructions for use

INHEALTH

TECHNOLOGIES®

UNITED STATES

1110 Mark Avenue
Carpinteria, CA 93013-2918 USA
Toll-Free: 800.477.5969
Toll-Free Fax: 888.371.1530

INTERNATIONAL

Tel: +1.805.684.9337
Fax: +1.805.684.8594
E-mail: order@inhealth.com

ORDER ONLINE

www.inhealth.com

Freudenberg Medical, LLC
1110 Mark Avenue
Carpinteria, CA 93013-2918 USA
www.inhealth.com

EMERGO EUROPE
Westervoortsedijk 60
6827 AT Arnhem
The Netherlands