
INHEALTH

TECHNOLOGIES®

Blom-Singer® Insufflation Test Set

MEDICAL PROFESSIONAL INSTRUCTIONS FOR USE

Blom-Singer®

voice restoration systems

Blom-Singer® and InHealth Technologies® are registered trademarks of **Freudenberg Medical, LLC.**

Patent(s): www.inhealthpatents.com

IFU: <https://inhealth.com/ifu>





Copyright © 2020 Freudenberg Medical, LLC. All rights reserved.

MEDICAL PROFESSIONAL INSTRUCTIONS FOR USE

医学专业使用说明 / 医療専門家用取扱説明書 / Udhëzimet e përdorimit për mjekët profesionistë / إرشادات الاستخدام للمتخصصين في المجال الطبي / Инструкции за употреба за медицински специалисти / Upute za korištenje za zdravstveno osoblje / Návod k použití pro zdravotníky / Brugsanvisning til uddannede læger / Gebruiksaanwijzing voor de medische professional / Tervishoiutöötaja kasutusjuhend / Lääketieteen ammatillaisen käyttöohjeet / Instructions d'utilisation destinées aux professionnels de santé / Gebrauchsanweisung für medizinische Fachkräfte / Οδηγίες χρήσης για τον επαγγελματία υγείας / הוראות שימוש לצוות הרפואי / Egészségügyi szakembereknek szóló használati utasítás / Notkunarleiðbeiningar fyrir heilbrigðisstarfsmenn / Istruzioni per l'uso destinate al personale sanitario / 의료 전문가 사용 지침 / Medicīnas speciālistam paredzētas lietošanas instrukcijas / Medikams profesionalams skirta naudojimo instrukcija / Bruksanvisning for medisinske fagpersoner / Instrukcja obsługi dla personelu medycznego / Instruções de uso para profissionais médicos / Instrucțiuni de utilizare pentru cadre medicale specializate / Инструкции по применению для профессиональных медработников / Uputstvo za upotrebu za zdravstvenog radnika / Návod na použitie pre zdravotníckych pracovníkov / Instrucciones de uso para profesionales médicos / Sjukvårdspersonal, bruksanvisning / Sağlık Uzmanı Kullanım Talimatları

DIAGRAMS

插图 / 図 / Diagramet / الرسم التوضيحي / Схеми / Dijagrami / Schémata / Diagrammer /
 Afbeeldingen / Diagrammid / Kaaviot / Diagrammes / Abbildungen / Διαγράμματα / תרשימים /
 Ábrák / Skýringarmyndir / Diagrammi / 그림 / Diagrammas / Diagramos / Diagrammer /
 Schematy / Diagramas / Diagreame / Схеми / Crteži / Grafické znázornenia / Diagramas /
 Diagram / Diyagramlar



TABLE OF CONTENTS

目録 / 目次 / Tabela e përmbajtjes / جدول المحتويات / Съдържание / Sadržaj / Obsah / Indholdsfortegnelse / Inhoud / Sisukord / Sisällysluettelo / Table des matières / Inhaltsverzeichnis / Πίνακας περιεχομένων / תוכן עניינים / Tartalomjegyzék / Efnisyfirlit / Sommario / 목차 / Saturs / Turinys / Innholdsfortegnelse / Spis treści / Índice / Cuprins / Содержание / Sadržaj / Obsah / Índice / Innehållsförteckning / İçindekiler

EN	Insufflation Test Set	6	HU	Befúvás tesztkészlet	42
ZH	吹气测试仪	8	IS	Búnaður til prófunar á útblæstri	44
JA	通気検査セット	10	IT	Set di test per insufflazione	47
SQ	Kompleti i testit të insufacionit	12	KO	흡입제 테스트 세트	49
15	جهاز اختبار النفخ	AR	LV	Insufļācijas pārbaudes komplekts	51
BG	Комплект за тестване на инсуфлация	17	LT	Insufliacijos bandymo rinkinys	53
HR	Komplet za test insuflacije	19	NO	Testsett for insufflasjon	56
CS	Sada pro testování insuflace	21	PL	Zestaw do testu wdmuchiwania	58
DA	Insufflationstestsæt	24	PT	Conjunto de teste de insuflação	60
NL	Testset voor inblazing	26	RO	Set de testare prin insuflație	63
ET	Insufletsiooni testikomplekt	28	RU	Набор для проверки вдвувания	65
FI	Insufflaatiotestisarja	30	SR	Komplet za test insuflacije	68
FR	Ensemble de test d'insufflation	33	SK	Sada na testovanie insuflácie	70
DE	Insufflationstest-Set	35	ES	Kit de prueba para la insuflación	72
EL	Σετ δοκιμής εμφύσησης	38	SV	Set för insufflationstest	75
40	הוראות שימוש לצוות הרפואי	HE	TR	Üfleme Testi Seti	77

ENGLISH

BLOM-SINGER® INSUFFLATION TEST SET

PRODUCT DESCRIPTION

The Insufflation Test Set (diagram 1) is a disposable, one-time use system consisting of a 14 Fr. silicone catheter (diagram 1A) that is attached to a tracheostoma adaptor (diagram 1B). The catheter is 50 cm long and imprinted with a 25 cm length marker (diagram 1C). A Blom-Singer tracheostoma valve housing (diagram 1D) and adhesive tape discs (diagram 1E) are included to attach the tracheostoma adaptor to the skin surrounding the patient's tracheostoma. This product is for single patient use only. Not for use between patients as infection may result. The Insufflation Test Set is supplied Nonsterile; do not sterilize as damage to the device may occur. Included in each package is a silicone catheter with attached tracheostoma adaptor, a Blom-Singer Tracheostoma Valve Housing and Blom-Singer Adhesive Tape Discs.

INDICATIONS (Reasons to prescribe the device or procedure)

The Blom-Singer Insufflation Test Set is designed to aid the clinician in subjectively predicting the ability of the laryngectomized patient to produce tracheoesophageal voice.

WARNINGS AND PRECAUTIONS

The Insufflation Test Set should only be used by clinicians who have received proper training in its use. If there are tears, cracks, or structural damage to the device, discontinue use. Do not use the device if it becomes contaminated. Care must be exercised to ensure that the catheter does not coil in the pharynx as this may result in excessive coughing, gagging, and/or irritation. If significant spasm of the pharyngeal constrictor muscles with the opposition to airflow is suspected, the assessment should be abbreviated to avoid excessive gastric distension or difficulty in breathing.

INSTRUCTIONS FOR USE

The Insufflation Test Set and user's hands should be clean and dry to prevent infection or contaminants from entering into the patient's body. Use of gloves, protective eye wear, and a respiratory face mask by the administering clinician is recommended for the catheter insertion and insufflation test procedure. Removal and attachment of any devices worn over the tracheostoma should only be done with a bright light focused directly on the stoma. The included Blom-Singer tracheostoma housing must first be attached to the skin surrounding the patient's tracheostoma. The following procedural instructions are provided by Eric D. Blom, Ph.D. Please refer to the diagrams located at the front of this instruction manual. Tracheostoma Valve Housing should be applied to clean, dry skin. Apply a protective skin barrier as needed (supplied separately). If a stronger seal is desired, apply a thin coat of medical grade liquid silicone adhesive (supplied separately) to the skin prior to attaching the housing. Allow the silicone adhesive to dry completely before attaching the assembled housing.

Blom-Singer Tracheostoma Valve Housing + Adhesive Tape Disc Assembly

1. Peel the backing paper from one side of the disc only (diagram 2).
2. Place the Tracheostoma Valve Housing on top of the adhesive disc.
3. Flatten the Tracheostoma Valve Housing against the adhesive disc surface, pressing firmly around the entire housing to ensure complete contact with the adhesive disc (diagram 3).
4. Remove the paper backing on the reverse side of the disc to expose the second adhesive surface (diagram 4).

Attaching Assembled Housing to Skin

1. Position the assembled housing over the tracheostoma.
2. Press the assembled housing flush against the skin surface to avoid air pockets between the skin and the housing.
3. Press the assembled housing firmly against the skin around the entire housing to achieve maximum adhesive contact (diagram 5).

Catheter Placement

If indicated, topical anesthetic may be used. Apply a medical-grade, water soluble lubricant to the catheter tip to facilitate insertion. Introduce the catheter into one nostril and continue to pass until the 25 cm length marker reaches the nostril, thereby positioning the tip within the upper cervical esophagus (diagram 6A). Insert the tracheostoma adaptor (diagram 6B) into the tracheostoma valve housing (diagram 6C). Tape the catheter to the skin as shown to prevent movement during the test procedure (diagram 6D).

Patient Preparation

Prior to the assessment procedure, it may be useful to demonstrate to the patient, in front of a mirror, how to occlude the tracheostoma adaptor with a finger to attempt phonation on exhalation. Practice this until he/she is comfortable with the task.

Assessment Procedure

Instruct the patient to take a breath in, then occlude the opening in the adaptor with a finger to divert exhaled air to the esophagus as he or she attempts to produce voice (diagram 6D).

Recommended Test Stimuli

A minimum of five (5) separate attempts of sustained vowel production on a single exhalation; and five (5) attempts to count uninterrupted from 1 to 15.

Insufflation Test Interpretation

Positive Result

A positive insufflation test result is characterized by sustained phonation without interruption for eight (8) seconds or longer, and the ability to count fluently from 1 to 15. If a positive insufflation test result is achieved, the patient may be a suitable candidate for voice restoration via placement of a voice prosthesis.

Negative Result

A negative insufflation test result may vary from complete inability to produce voice, to brief 3 to 4 syllable utterances. Pharyngeal constrictor spasm is typically characterized by brief periods of phonation, opposition to air flow, and gastric filling. If the insufflation test result is negative, patients may require further assessment or intervention in order to be a suitable candidate for successful voice restoration via a voice prosthesis.

SPECIAL STORAGE AND/OR HANDLING CONDITIONS

There are no special storage and/or handling conditions for this device.

DISPOSAL INSTRUCTIONS

This device is not biodegradable and can be contaminated when used. Carefully dispose of the device per local guidelines.

ORDERING INFORMATION

Blom-Singer products may be ordered directly from InHealth Technologies. TELEPHONE: Toll-Free (800)477-5969 or (805)684-9337. FAX: Toll-Free (888)371-1530 or (805)684-8594. EMAIL: order@inhealth.com ORDER ONLINE: www.inhealth.com POST: InHealth Technologies, 1110 Mark Avenue, Carpinteria, CA 93013-2918, USA, Attention: Customer Service

INTERNATIONAL

Please contact our customer service department for a distributor referral.

PRODUCT COMPLAINTS/EU SERIOUS INCIDENTS

If you are dissatisfied with the device or have any questions, please contact productcomplaints@inhealth.com.

Telephone: +1-800-477-5969

Fax: +1-888-371-1530

Any serious incident that has occurred in relation to the device should be reported to Freudenberg Medical, LLC as noted above and the competent authority of the EU Member State in which the user and/or patient is established.

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SYMBOLS GLOSSARY (USA only)

SYMBOL	STANDARD No.	SYMBOL TITLE & REF No.	MEANING OF SYMBOL
	EN ISO 15223-1:2016	Manufacturer; 5.1.1	Indicates the medical device manufacturer, as defined in EU Directives or Regulations
	93/42/EEC EU 2017/745	CE Marking; 93/42/EEC Annex XII, EU 2017/745 Annex V	CE marking of conformity
	EN ISO 15223-1:2016	Authorized representative in the European Community; 5.1.2	Indicates the authorized representative in the European Community
	EN ISO 15223-1:2016	Batch code; 5.1.6	Indicates the manufacturer's catalogue number so that the medical device can be identified
	EN ISO 15223-1:2016	Batch code; 5.1.5	Indicates the manufacturer's batch code so that the batch or lot can be identified
	EN ISO 15223-1:2016	Use-by date; 5.1.4	Indicates the date after which the medical device is not to be used
	EN ISO 15223-1:2016	Caution; 5.4.4	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself
	EN ISO 15223-1:2016	Do not use if package is damaged; 5.2.8	Indicates a medical device that should not be used if the packaging has been damaged or opened
	EN ISO 15223-1:2016	Consult instructions for use; 5.4.3	Indicates the need for the user to consult the instructions for use
	ISO 7000:2019	Packaging unit; 2794	To indicate the number of pieces in the package
	EN ISO 15223-1:2016	Do not re-use; 5.4.2	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure

INHEALTH

TECHNOLOGIES®

UNITED STATES

1110 Mark Avenue
Carpinteria, CA 93013-2918 USA
Toll-Free: 800.477.5969
Toll-Free Fax: 888.371.1530

INTERNATIONAL

Tel: +1.805.684.9337
Fax: +1.805.684.8594
E-mail: order@inhealth.com

ORDER ONLINE

www.inhealth.com

Freudenberg Medical, LLC
1110 Mark Avenue
Carpinteria, CA 93013-2918 USA
www.inhealth.com

EMERGO EUROPE
Westervoortsedijk 60
6827 AT Arnhem
The Netherlands