
INHEALTH[®]

TECHNOLOGIES

Blom-Singer[®] Adjustable Bi-Flanged Fistula Prosthesis

MEDICAL PROFESSIONAL Instructions For Use

Blom-Singer[®]

voice restoration systems

R1
37742-01D

37742-01D Effective June 2017 / Gültig ab Juni 2017 / En vigencia a partir de junio del 2017 / Date d'entrée en vigueur : Juin 2017 / Valido da Giugno 2017 / Met ingang van juni 2017 / Vigente em junho de 2017 / V platnosti od červen 2017 / Wchodzi w życie w czerwcu 2017 r.

Blom-Singer and **InHealth Technologies** are registered trademarks of **Freudenberg Medical, LLC**

Patent(s): www.inhealthpatents.com

Blom-Singer®

ENG	Blom-Singer® Adjustable Bi-Flanged Fistula Prosthesis	4
DEU	Blom-Singer® Anpassbare Doppelflanschfistelprothese	18
SPA	Prótesis de dos bridas regulable para fístula Blom-Singer®	34
FRA	Prothèse ajustable pour fistule à deux collettes Blom-Singer®	49
ITA	Protesi per fistola a doppia flangia regolabile Blom-Singer®	65
NLD	Blom-Singer® verstelbare fistelprothese met twee kragen	80
POR	Prótese para fístula ajustável com dois flanges Blom-Singer®	95
CES	Upravitelná dvoupřírubová fistulová protéza Blom-Singer®	110
POL	Regulowana protezę do zamykania przetoki z dwoma kołnierzami Blom-Singer®	123
	Diagrams	139
	Bibliography	142

MEDICAL PROFESSIONAL Instructions For Use

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BLOM-SINGER® ADJUSTABLE BI-FLANGED FISTULA PROSTHESIS

Introduction

Please refer to the diagrams located at the back of this instruction manual. A glossary of terms appears at the end of this chapter.

The Blom-Singer® Adjustable Bi-Flanged Fistula Prosthesis may be used only once by one patient. It may not be reused. See Warnings and Precautions, and Complications.

An accessory product, the Blom-Singer® Insertion and Removal Catheter is available separately, and is intended for use ONLY with the Blom-Singer® Adjustable Bi-Flanged Fistula Prosthesis. It is available for use in certain insertion and removal options as determined by the physician. It may be used only once on one patient, and may not be reused. See Warnings and Precautions, and Complications.

The physician should discuss the risks and benefits of using a Blom-Singer® Adjustable Bi-Flanged Fistula Prosthesis with the patient.

INDICATIONS

The Blom-Singer® Adjustable Bi-Flanged Fistula Prosthesis is indicated for the management of hypopharyngeal fistulae.

A fistula is an abnormal passage or communication leading from an internal space out to the surface of the body. Patients having hypopharyngeal (such as pharyngocutaneous) fistulae need a means to reduce leakage of saliva and food/drink or esophageal contents into soft tissue, or external to the body.

Traditional treatment is typically to pack the site until it spontaneously heals or with the placement of a salivary bypass tube. Depending on fistula size and location, the patient may require a feeding tube and hospitalization for multiple days until the site has healed enough for normal eating.

Some of the research-proven predictors of tissue breakdown leading to formation of fistulae may be:

- Hypothyroidism
- Diabetes Mellitus
- Salivary Pooling
- Radiation therapy effects on the wound bed
- Gastric Reflux
- Residual Disease
- Infection

The objective of the Blom-Singer® Adjustable Bi-Flanged Fistula Prosthesis is to provide an alternative to the traditional treatment of fistulae described above. This prosthesis may maintain a dry tract without actual obturation. The device may also allow post-laryngectomy speech, eating by mouth, and swallowing. This is a short-term medical device, intended for use for not more than 29 days.

CONTRAINDICATIONS

The Blom-Singer® Adjustable Bi-Flanged Fistula Prosthesis should be used only by a physician with experience and training in its use and care.

The Blom-Singer® Adjustable Bi-Flanged Fistula Prosthesis is contraindicated for use in a fistula tract length that is longer than 25mm or a diameter that exceeds 25mm.

The Blom-Singer® Adjustable Bi-Flanged Fistula Prosthesis is contraindicated for use in a fistula tract that is not demonstratively patent along its entire length.

The Blom-Singer® Adjustable Bi-Flanged Fistula Prosthesis is contraindicated for use in a fistula where an existing pharyngoesophageal stricture may complicate or restrict prosthesis insertion or removal through the pharynx or cervical esophagus.

The Blom-Singer® Adjustable Bi-Flanged Fistula Prosthesis is

contraindicated in a fistula tract that is in close, unsafe proximity to the carotid artery such that mechanical erosion is potentially possible.

The Blom-Singer® Adjustable Bi-Flanged Fistula Prosthesis is contraindicated in a fistula if the potential for airway obstruction or device aspiration exists.

The Blom-Singer® Adjustable Bi-Flanged Fistula Prosthesis is contraindicated for pediatric patients.

HOW SUPPLIED

The Blom-Singer® Adjustable Bi-Flanged Fistula Prosthesis is provided sterile (sterilized by gamma radiation). It is packaged within a thermoform tray with a Tyvek® lid and placed in a box with the Instructions for Use. An accessory product, the Blom-Singer® Fistula Prosthesis Insertion and Removal Catheter is available separately. It is also provided sterile (sterilized by gamma radiation). It is packaged within a Tyvek® pouch and placed in a box with the Instructions for Use.

PRODUCT DESCRIPTION

The Blom-Singer® Adjustable Bi-Flanged Fistula Prosthesis is made of medical-grade silicone rubber. It is not made with natural rubber latex.

The Blom-Singer® Adjustable Bi-Flanged Fistula Prosthesis is designed to seal a fistula when the prosthesis is placed, adjusted and a snug flange-to-flange fit is achieved. The Fistula Prosthesis is available in three flange diameters: 25mm, 38mm, and 50mm.

The flange diameter options allow overlap of the outside diameter of the fistula by a recommended 8-10mm for the 25mm and 38mm sizes, and by a recommended 12.5mm overlap for the 50mm size. The three (3) sizes can seal a fistula that varies between approximately 6mm diameter and no larger than 25mm diameter, and tract length from approximately 1.5mm to 25mm.

The flanges, which are flexible and translucent, allow the device

to conform to the contours of the peri-fistula tissues to reduce leakage through the fistula tract. They may be individually trimmed to other diameters at the judgment of a physician. The 50mm Fistula Prosthesis has dual stems to facilitate a sufficient flange-to-flange fit at the perimeter of the flanges. Greater flange overlap is recommended for use with exceptionally large fistulæ.

The Fistula Prosthesis flanges are connected to an 8Fr diameter stretchable (elastomeric) beaded stem. Both flanges are the same diameter for a given size device. The inner flange is flexible enough to allow folding for trans-fistula or trans-oral insertion by a physician.

The stem has beads at regular intervals to accommodate a wide variety of fistula tract lengths, and to allow for a snug fit between the outer and inner flanges to create a surface seal against the tissues. The beaded stem is molded from silicone with barium sulfate and is therefore radiopaque.

The distance between the Fistula Prosthesis flanges can be manually adjusted tighter or looser depending on the sealing needs. The flexibility, or stretch, of the stem also allows self-adjustment to maintain the seal.

The design includes a locking loop for each of the beaded stems. After placement by a physician, and after achievement of a snug flange-to-flange surface seal against the tissue, feed the tapered end of the stem through the locking loop to help prevent loss of the proper fit. The excess stem can then be shortened at the physician's discretion.

The accessory product, the Blom-Singer® Insertion and Removal Catheter is made of medical-grade thermoplastic polyurethane. The Blom-Singer® Insertion and Removal Catheter is a 41cm (16in) long, 11Fr rod device that tapers down to 8Fr diameter and incorporates attachment eyelets for the Blom-Singer® Adjustable Bi-Flanged Fistula Prosthesis. The Blom-Singer® Insertion and Removal Catheter is designed to aid in the insertion or removal of the Blom-Singer® Adjustable Bi-Flanged Fistula Prosthesis.

WARNINGS AND PRECAUTIONS

The Blom-Singer® Adjustable Bi-Flanged Fistula Prosthesis is a medically prescribed device. Initial sizing, insertion, adjustment, and instructions for use must be provided by a physician trained in the use of this particular device.

Patients should be instructed on the use and care of the device by their trained physician or clinician.

If the prosthesis is not functioning properly, the patient should have it evaluated by a physician as soon as possible.

Do not use the prosthesis or the insertion and removal catheter if the integrity of the sterile package has been compromised and/or if the device looks contaminated. If this is the case, discard and use a new prosthesis/catheter.

Do not use the prosthesis or the insertion and removal catheter if the product has expired beyond the allowable shelf life. If this is the case, discard and use a new prosthesis/catheter.

Inspect the prosthesis before, during and after insertion. If there are tears, cracks, or structural damage to the prosthesis/catheter, discontinue use.

Do not use petroleum-based products such as Vaseline® to lubricate the prosthesis as these products can damage silicone.

Changes in the anatomy or medical status of the patient may lead to improper fitting and/or function of the device. It is recommended that the prosthesis and fistula be evaluated by the physician at regular intervals, and the fit of the prosthesis adjusted as clinically indicated.

The prosthesis and the insertion and removal catheter are for single-use only. Should they be reused, there is a risk of infection due to microbial contamination. Do not resterilize the products. Reprocessing and/or resterilization of the devices may alter their performance and device failure may occur.

Do not use sharp instruments such as hemostats on the beaded stem since they may cause damage which could result in a broken beaded stem.

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

COMPLICATIONS

Although rare, the following complications could potentially occur with the Blom-Singer® Adjustable Bi-Flanged Fistula Prosthesis:

- fistula contamination or sepsis, which may require removal of the prosthesis and/or appropriate medical management of sepsis
- accidental aspiration of the prosthesis into the airway, which may require removal by a physician
- accidental ingestion of the prosthesis into the esophagus and/or GI tract
- occasional extrusion of the prosthesis, requiring replacement
- fistula dilation resulting in leaking of fluids around the prosthesis
- inflammatory reaction around the fistula
- intractable leakage around the prosthesis, requiring adjustment or re-fitting of the prosthesis; or other surgical/prosthetic or medical management options
- tearing or other damage to the prosthesis from improper use

INSTRUCTIONS FOR USE

The following instructions are provided by Eric D. Blom, Ph.D. and Stephen B. Freeman, M.D., F.A.C.S. for insertion and removal procedures.

Fistula Assessment

Prior to insertion of a Blom-Singer® Adjustable Bi-Flanged Fistula Prosthesis, gently probe the fistula tract with a cotton tip applicator or 8Fr soft catheter to grossly estimate its length, diameter, and direction (Fig. 1). Fistula tract length and direction may also be assessed radiographically using an approved contrast substance.

Select a fistula prosthesis flange diameter that provides at least an 8 to 10mm overlap circumferentially around the outside edges of the fistula. The distal and proximal flanges should be of equal size. In some instances, the flanges may need to be carefully trimmed with scissors to provide the appropriate flange shape and/or size, at the physician's discretion.

Fistula Prosthesis Insertion Options

NOTE: The Blom-Singer® Insertion and Removal Catheter is available separately, and is intended for use ONLY with the Blom-Singer® Adjustable Bi-Flanged Fistula Prosthesis. It is not intended for use in Insertion Option A below, but is intended for use in Option B, Trans-oral retrograde insertion.

A) Trans-fistula insertion

Hand-fold the inner flange of the prosthesis into as small a configuration as possible and grasp it tightly with a hemostat. Also gently grasp the end of the beaded stem, keeping slight tension on the stem for safety (Fig 2).

Gently insert the folded flange through the fistula tract releasing the flange from its folded configuration at the inner aspect of the fistula tract. Withdraw the hemostat while securely holding the end of the beaded stem (Fig 3).

Lightly tug outwardly on the beaded stem to confirm that the flange has deployed and is securely seated against the inner tissue surface (Fig. 4).

Insert the tapered tip on the beaded stem of the prosthesis through the center hole in the outer flange. Move the outer flange down the beaded stem until a "snug" inner-flange-to-outer-flange relationship is achieved that seals sufficiently enough to reduce fistula leakage. Once this flange adjustment is complete, feed the tapered end of the beaded stem through the locking loop, leaving a small amount of slack in the stem to form a curved contour (Fig. 5). Detach the excess stem, leaving 4 to 5 beads beyond the face of the locking loop (Fig. 6). These

remaining beads may be needed later for adjustment and/or for prosthesis removal. Test for leakage as described below under "Post-Insertion Inspection".

B) Trans-oral retrograde insertion

This insertion option is the method of choice for fistula diameters of approximately 6mm or less to avoid trans-fistula insertion dilation. In trans-oral retrograde insertion, a dedicated Blom-Singer® Fistula Prosthesis Insertion and Removal Catheter is used (available separately, FM1010). Insert the Insertion and Removal Catheter through the fistula, out through the mouth (Fig. 8). To securely attach the beaded stem of the prosthesis (Figs. 7, 9) thread the beaded stem though the center eyelet of the catheter then carefully slide the beaded stem up (Fig. 7, A) into the smaller locking eyelet of the catheter.

Gently pull the catheter with prosthesis attached retrograde until the beaded stem exits the fistula and the inner flange is securely seated against the inner tissue surface (Figs. 9, 10).

Carefully detach the beaded stem from the Insertion and Removal Catheter by sliding the stem back to the center eyelet. Take care to retain the end of the beaded stem in one hand. The outer flange is assembled onto the beaded stem and is "snugly" adjusted as previously described above. Test for leakage as described below under "**Post-Insertion Inspection**".

Caution: Inspect the Prosthesis before, during and after insertion. If there are tears, cracks, or structural damage, such as stem breakage, to the prosthesis replace the device with another prosthesis.

Caution: Ensure the beaded stem with the inner flange attached is properly secured through the outer flange and fed through the locking loop to avoid unintentional dislodgement.

Post-Insertion Inspection

Visually inspect the placement and contact of the prosthesis flange(s) around and against the peri-fistula tissues. Circumferential

contact is necessary for adequate sealing. Visually confirm that the peri-fistula tissue has not been dilated or irritated from the insertion procedure. Confirm with the patient that the prosthesis fits comfortably. Discontinue use of the prosthesis if proper fit or level of comfort are medically concerning.

Ask the patient to drink 120 to 180ml (4 to 6oz) of water. Confirm that the leakage, if any, is minimal and acceptable for the patient's condition. Repeat the above a second time and confirm.

If the amount of leakage is not acceptable, the prosthesis may not be inserted or adjusted properly. Adjust the flange-to-flange tightness against the peri-fistula tissues and reconfirm with a swallow of water.

If the amount of leakage is still not acceptable for the patient's condition, consider replacing it with a different sized prosthesis.

In the event acceptable leakage cannot be achieved, the patient's tissue condition and/or anatomy may be the cause, and in this case use of the prosthesis is contraindicated.

In-Situ Prosthesis Cleaning

Instruct the patient on the routine care and cleaning of the Blom-Singer® Adjustable Bi-Flanged Fistula Prosthesis.

The fistula prosthesis flanges are manufactured of transparent silicone. This provides the opportunity to visually monitor the status of the fistula. Movement of the outer flange disturbs the seal to the peri-fistula tissues and is discouraged. However, if necessary, remove any observed debris from behind the outer flange with a sterile, 6-inch cotton-tipped plastic applicator. Ensure the fistula prosthesis is repositioned properly to achieve a snug fit that seals against the peri-fistula tissues.

During evaluation of the fistula by the physician at regular intervals, visually monitor the condition of the peri-fistula tissues to confirm there is no allergic or toxic reaction, significant irritation, or tissue infection. If there is a medical concern, discontinue use of the prosthesis.

Instruct the patient that they may carefully clean around the Fistula Prosthesis while taking care to not disturb the seal. Always use gentle handling and pressure when cleaning in or around the prosthesis to avoid deterioration and damage of the device which could cause leakage.

Caution: The user's hands and any cleaning materials should always be thoroughly clean to prevent contaminants from entering in or around the fistula. Cleaning should only be done in front of a mirror with a bright light focused directly on the fistula so it is clearly visible. Never use facial or toilet tissue as particles of lint or fabric may be deposited which can disturb the seal of the device. Occasional and light contact of a mild soap and water solution is permitted, followed by drying any excess fluid on the prosthesis. Drying the excess fluid by dabbing or blotting are the preferred methods, rather than wiping. Do not wipe in the same direction as the beaded stem passes through the locking loop as this may cause the fit of the beaded stem to loosen. If the prosthesis is not functioning properly, the patient should have it evaluated by a physician as soon as possible.

Caution: If the beaded stem disengages from the locking loop, the patient should have it evaluated by a physician as soon as possible.

Fistula Prosthesis Removal Options

NOTE: The Blom-Singer® Insertion and Removal Catheter is available separately, and is intended for use ONLY with the Blom-Singer® Adjustable Bi-Flanged Fistula Prosthesis. It is not intended for use in Removal Option A below, but is intended for use in Option B, Antegrade trans-fistula removal.

A) Sedated endoscopic removal

In the event the physician elects this approach, the inner flange of the fistula prosthesis can be removed in a sedated endoscopic removal procedure. Detach the outer flange from the beaded stem while holding on to the beaded stem. Grasp the stem of the

prosthesis with laryngeal forceps or alligator forceps as close as possible to where the beaded stem attaches to the inner flange and withdraw the prosthesis out through the mouth.

B) Antegrade trans-fistula removal

In this non-sedated removal procedure, the outer flange of the fistula prosthesis is removed and then the beaded stem (and its attached inner flange) is securely attached to a dedicated Blom-Singer® Fistula Prosthesis Insertion and Removal Catheter (available separately, FM1010). To securely attach the beaded stem of the prosthesis to the removal catheter thread the beaded stem though the center eyelet of the catheter then carefully slide the beaded stem down (Fig.7, B) into the smaller locking eyelet of the catheter. Then push the fistula prosthesis out through the mouth (Fig. 7, 11, 12).

Caution: Do not use aggressive force on the beaded stem during the removal procedure to avoid stem breakage, trauma to the peri-fistula tissues, and to avoid accidental ingestion or aspiration.

Large Fistula Prosthesis Removal for Prosthesis Downsizing

When a large fistula has decreased in size such that a smaller prosthesis is indicated, the large prosthesis can be removed from a still-relatively large fistula by simply pulling it out with a hemostat. An appropriate size smaller device is then inserted by one of the previously described procedures. After downsizing and proper insertion, again follow the above instructions “**Post-Insertion Inspection**”.

ORDERING INFORMATION

USA

Blom-Singer products may be ordered directly from InHealth Technologies. TELEPHONE: Toll-Free (800)477-5969 or

(805)684-9337, Monday — Friday, 9:30 am — 7:00 pm, Eastern Standard Time. FAX: Toll-Free (888)371-1530 or (805)684-8594. EMAIL: order@inhealth.com ORDER ON-LINE: www.inhealth.com POST: InHealth Technologies, 1110 Mark Avenue, Carpinteria, CA 93013-2918, USA, Attention: Customer Service.

Consumer Affairs

If you have any questions or dissatisfaction with a product please contact our customer service department via telephone, fax, post, or Email: productcomplaints@inhealth.com

INTERNATIONAL CUSTOMERS

Blom-Singer products are available worldwide through authorized InHealth Technologies distributors. For assistance or additional product information please contact the distributor.

GLOSSARY

Contraindications

A reason that makes it inadvisable to prescribe a particular procedure or treatment.

Esophageal

Relating to the esophagus.

Esophagus

The muscular membranous tube for the passage of food from the pharynx to the stomach; the "food pipe."

Fistula

An abnormal duct or passage that connects an abscess or cavity to the body surface.

Hypopharyngeal

Located beneath the pharyngeal apparatus.

Indications

A reason to prescribe a drug or perform a procedure.

Pharyngeal

Relating to the Pharynx.

Pharynx

The superior expanded portion of the alimentary tract, between the mouth and nasal cavities (superiorly and anteriorly) and the esophagus (inferiorly); consisting of nasopharynx, oropharynx, and laryngopharynx, the first two being shared with the respiratory tract; the throat.

Sepsis

The presence of various pathogenic organisms, or their toxins, in the blood or tissues.

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




1. McLean, J.N., Nicholas, C., Duggal, P., Chen, A., Grist, W.G., Losken, A., Carlson, G.W.: **"Surgical Management of Pharyngocutaneous Fistula After Total Laryngectomy"**, Annals of Plastic Surgery, July 2011. From the Division of Plastic Surgery, Department of Surgery, Emory University School of Medicine, Atlanta GA; and Department of Otolaryngology, Emory University School of Medicine, Atlanta GA.
2. Aarts, M.C., Rovers, M.M., Grau, C., Grolman, W., van der Heijden, G.J.M.: **"Salvage Laryngectomy After Primary Radiotherapy: What Are Prognostic Factors for the Development of Pharyngocutaneous Fistulae?"**, Otolaryngology – Head and Neck Surgery, 144 (1) 5-9 , 2011.
3. Bohannon, I.A., Carroll, W.R., Magnuson, J.S., Rosenthal, E.L.: **"Closure of post-laryngectomy pharyngocutaneous fistulae"**, Head and Neck Oncology, 3:29, 2011.
4. Lorenz, K.J., Griesler, L., Ehrhart, T., Maier, H.: **"Role of Reflux in Tracheoesophageal Fistula Problems After Laryngectomy"**, Annals of Otology, Rhinology & Laryngology, 119 (11): 719-728, 2010.
5. Saki, N., Nikakhlagh, S., Kazemi, M.: **"Pharyngocutaneous Fistula after Laryngectomy: Incidence, Predisposing Factors, and Outcome"**, Archives of Iranian Medicine, 11 (3): 314-317, 2008.
6. Gehrking, E., Raap, M., Sommer, K.D.: **"Classification and Management of Tracheoesophageal and Tracheopharyngeal Fistulas After Laryngectomy"**, Laryngoscope, 117: 1943-1951, 2007.
7. Dedivitis, R.A., Ribeiro, K.C.B., Castro, M.A.F., Nascimento, P.C.: **"Pharyngocutaneous fistula following total laryngectomy"**, Acta Otorhinolaryngologica Italica, v.27 (1): 2007.

-
8. Paydarfar, J.A., Birkmeyer, N.J.: "**Complications in Head and Neck Surgery, A Meta-analysis of Postlaryngectomy Pharyngocutaneous Fistula**", Arch Otolaryngology Head Neck Surgery, 132: 2006.
 9. Lörken, A, Krampert, J., Kau, R.J., Arnold, W.: "**Experiences with the Montgomery Salivary Bypass Tube (MSBT)**", Dysphagia, 12 (2): 79-83: Spring 1997.

INHEALTH®

TECHNOLOGIES

UNITED STATES

1110 Mark Avenue
Carpinteria, CA 93013-2918 USA
Toll-Free: 800.477.5969
Toll-Free Fax: 888.371.1530

INTERNATIONAL

Tel: +1.805.684.9337
Fax: +1.805.684.8594
E-mail: order@inhealth.com

ORDER ONLINE

www.inhealth.com

Freudenberg Medical, LLC
1110 Mark Avenue
Carpinteria, CA 93013-2918 USA
www.inhealth.com

CE 0344

EMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands