
# INHEALTH

TECHNOLOGIES®

## Blom-Singer® Adjustable Bi-Flanged Fistula Prosthesis

### INSTRUCTIONS FOR USE

**Blom-Singer®**

voice restoration systems

*Blom-Singer® Adjustable Bi-Flanged Fistula Prosthesis*

37742-01E Effective April 2021 / 2021 年 4 月生效 / 2021年4月發効 / Në fuqi nga prilli 2021 / 2021 نافذ المفعول بدءاً من أبريل / В сила от април 2021 г. / Stupa na snagu u travnju 2021 / V platnosti od dubna 2021 / Træder i kraft april 2021 / Met ingang van april 2021 / Jöustub 2021. aasta aprilli / Voimaantulo: huhtikuu 2021 / Date d'entrée en vigueur : avril 2021 / Gültig ab april 2021 / Με ισχύ από τον Απρίλιο του 2021 / 2021 החל באפריל / Érvényes: 2021 áprilisától / Gildir frá april 2021 / Valido da aprile 2021 / 2021년 4월 발효 / Spēkā no 2021. gada aprīlis / Įsigalioja 2021 m. balandžio mėn. / Gjelder fra april 2021 / Wchodzi w życie w kwietniu 2021 r. / Vigente em abril de 2021 / Cu intrare în vigoare din aprilie 2021 / Версия от апреля 2021 г. / Stupa na snagu u aprilu 2021. / Platný od apríla 2021 / Vigente desde abril de 2021 / Gäller från april 2021 / Nisan 2021'de Yürürlüğe Girdi

**Blom-Singer® and InHealth Technologies®** are registered trademarks of **Freudenberg Medical, LLC**.

Patent(s): [www.inhealthpatents.com](http://www.inhealthpatents.com)





Copyright © 2021 Freudenberg Medical, LLC. All rights reserved.

## INSTRUCTIONS FOR USE

使用说明 / 使用說明書 / Udhëzime për përdorim / إرشادات الاستخدام / Инструкции за употреба / Upute za korištenje / Návod k použití / Brugsanvisning / Gebruiksaanwijzing / Kasutusjuhend / Käyttöohjeet / Instructions d'utilisation / Gebrauchsanweisung / Οδηγίες χρήσης / הוראות שימוש / Használati utasítás / Notkunarleiðbeiningar / Istruzioni per l'uso / 사용 지침 / Lietošanas norādījumi / Naudojimo instrukcijos / Bruksanvisning / Instrukcja użytkowania / Instruções de uso / Instrucțiuni de utilizare / Инструкции по применению / Uputstvo za upotrebu / Návod na použitie / Instrucciones de uso / Bruksanvisningar / Kullanım Talimatları

## DIAGRAMS

插图 / 図 / Diagramet / الرسم التوضيحي / Схеми / Dijagrami / Schémata / Diagrammer /  
 Afbeeldingen / Diagrammid / Kaaviot / Diagrammes / Abbildungen / Διαγράμματα / תרשימים /  
 Ábrák / Skýringarmyndir / Diagrammi / 그림 / Diagrammas / Diagramos / Diagrammer /  
 Schematy / Diagramas / Diagrame / Схеми / Crteži / Grafické znázornenia / Diagramas /  
 Obrázky / Diyagramlar

## TABLE OF CONTENTS

目录 / 目次 / Tabela e përmbajtjes / جدول المحتويات / Съдържание / Sadržaj / Obsah /  
 Indholdsfortegnelse / Inhoud / Sisukord / Sisällysluettelo / Table des matières /  
 Inhaltsverzeichnis / Πίνακας περιεχομένων / תוכן עניינים / Tartalomjegyzék / Efnisyfirlit /  
 Sommario / 목차 / Saturs / Turinys / Innholdsfortegnelse / Spis treści / Índice / Cuprins /  
 Содержание / Sadržaj / Obsah / Índice / Innehållsförteckning / İçindekiler

EN	Adjustable Bi-Flanged Fistula Prosthesis	5	HU	Állítható, kettős peremű fistulaprotézis	67
ZH	可调节式双凸缘瘘管假体	8	IS	Stíllanlegt tvíflansa stöðtaki fyrir fistil	71
JA	調節可能バイフランジ型瘻孔プロテーゼ	12	IT	Protesi per fistola a doppia flangia regolabile	75
SQ	Protezë e rregullueshme e fistulës, me dy flanxha	15	KO	조절 가능 이중 플랜지 누공 인공후두	79
19	بدلة ناسور ثنائية الشفة قابلة للتعديل	AR	LV	Pielāgojama divu atloku fistulas protēze	83
BG	Регулируема двуфлангова фистулна протеза	23	LT	Reguliuojamas dviejų jungių fistulės protezas	86
HR	Prilagodljiva proteza za fistulu s dvama obodima	27	NO	Justerbar fistelprotese med to flensskiver	90
CS	Upravitelná dvouprůrbová fistulová protéza	31	PL	Regulowana proteza do zamykania przetoki z dwoma kolierzami	94
DA	Justerbar fistelprotese med to flanger	35	PT	Prótese ajustável para fistula com dois flanges	98
NL	Verstelbare fistelprothese met twee kragen	38	RO	Proteză pentru fistule, reglabilă, cu două manșoane	102
ET	Reguleeritav kahe äärikuga fistuli protees	43	RU	Регулируемый двухфланцевый протез фистулы	106
FI	Säädettävä kaksilaipainen fisteliproteesi	47	SR	Podesiva proteza za fistulu sa dvostrukom prirubnicom	111
FR	Prothèse ajustable pour fistule à deux collerettes	50	SK	Nastavitelná dvojprůrbová fistulová protéza	115
DE	Anpassbare Doppelflanschfistelprothese	55	ES	Prótesis de dos bridas regulable para fistula	118
EL	Ρυθμιζόμενη πρόθεση σιργγίου διπλού παρεμβύσματος	59	SV	Justerbar biflansad fistelprotes	123
64	תותבת פיסטולה דו-אוגנית מתכווננת	HE	TR	Ayarlanabilir Çift Flanşlı Fistül Protezi	126

## ENGLISH

**BLOM-SINGER® ADJUSTABLE BI-FLANGED FISTULA PROSTHESIS  
PRODUCT DESCRIPTION**

The Blom-Singer® Adjustable Bi-Flanged Fistula Prosthesis is a sterile (by gamma radiation), medical-grade silicone rubber device, designed to seal a fistula when properly placed, adjusted and a snug flange to flange fit is achieved. The device is not made with natural rubber latex. The Fistula Prosthesis is available in three flange diameters: 25 mm, 38 mm, and 50 mm. The three (3) sizes can seal a fistula that varies between approximately 6 mm diameter and no larger than 25 mm diameter. The 8 Fr diameter stretchable beaded stem(s) attached to the inner flange can be manually adjusted to accommodate fistula tract lengths from approximately 1.5 mm to 25 mm. The 50 mm Fistula Prosthesis has dual stems to facilitate flange-to-flange fit at the perimeter of the flanges. The beaded stem is molded from silicone with barium sulfate. Barium sulfate is a radiopaque material.

The flexible and translucent flanges allow the Fistula Prosthesis to conform to the contours of the peri-fistula tissue. The flanges may be individually trimmed to other shapes/sizes at the judgement of the physician. Both flanges are the same diameter for a given device size. The inner flange is flexible, to allow folding for trans-fistula insertion by a physician. The external flange includes a locking loop(s) for the beaded stem(s).

**INDICATIONS**

The Blom-Singer® Adjustable Bi-Flanged Fistula Prosthesis is indicated as an option for the management of hypopharyngeal (such as pharyngocutaneous) fistulae. The device may provide a means to reduce leakage of saliva, food/drink, or esophageal contents external to the body. It also may maintain a dry tract without traditional obturation. This is a short-term medical device, intended for use for not more than 29 days.

**CONTRAINDICATIONS**

The Blom-Singer® Adjustable Bi-Flanged Fistula Prosthesis should be used only by a physician with experience and training in its use and care. It is contraindicated for use:

- in a fistula tract length that is longer than 25 mm or a diameter that exceeds 25 mm
- in a fistula tract that is not demonstratively patent along its entire length
- in a fistula tract that is in close proximity to the carotid artery such that mechanical erosion is possible
- if the potential for airway obstruction or device aspiration exists
- with a patient where an existing pharyngoesophageal stricture may complicate or restrict prosthesis insertion or removal through the pharynx or cervical esophagus
- with pediatric patients

**WARNINGS AND PRECAUTIONS**

The Blom-Singer® Adjustable Bi-Flanged Fistula Prosthesis is a medically prescribed device. The physician should discuss the risks and benefits of using a Blom-Singer® Adjustable Bi-Flanged Fistula Prosthesis with the patient. Sizing, insertion, adjustment, and instructions for use must be provided by a physician trained in the use of this particular device. Patients should be instructed on the use and care of the device by their trained physician or clinician. If the prosthesis is not functioning properly, the patient should have the device evaluated by a physician as soon as possible. Do not use the prosthesis in the event of the packaging or sterile packaging being damaged or unintentionally opened prior to use; discard the prosthesis and replace it with a new prosthesis. Do not use the prosthesis if the product has expired; discard the expired product and use a new prosthesis. Inspect the prosthesis before, during, and after insertion. If there are tears, cracks, or structural damage to the prosthesis, discontinue use.

### *Blom-Singer® Adjustable Bi-Flanged Fistula Prosthesis*

Do not use petroleum-based products such as Vaseline® to lubricate the prosthesis as these products can damage silicone. Changes in the anatomy or medical status of the patient may lead to improper fitting and/or function of the device. The prosthesis and fistula should be evaluated by the physician at regular intervals, and the fit of the prosthesis should be adjusted as clinically indicated. The prosthesis is for single-use only; do not reuse the devices, reuse could result in infection due to microbial contamination. Do not resterilize the products. Reprocessing and/or reesterilization of the devices may alter their performance and device failure may occur. Do not use sharp instruments (such as hemostats) on the beaded stem, as damage or breakage could occur. Do not use aggressive force on the beaded stem during the insertion and removal procedures to avoid stem breakage, trauma to the peri-fistula tissues, and to avoid accidental ingestion or aspiration. Always use gentle handling and pressure to avoid deterioration and damage of the device which could cause leakage. Ensure the beaded stem with the inner flange attached is properly secured through the outer flange and fed through the locking loop to avoid unintentional dislodgement. During evaluations of the fistula by the physician, visually monitor the condition of the peri-fistula tissues to confirm there is no allergic or toxic reaction, significant irritation, or tissue infection. If there is a medical concern, discontinue use of the prosthesis. If the beaded stem disengages from the locking loop, the patient should be evaluated by a physician as soon as possible.

### **COMPLICATIONS**

Although rare, the following complications may occur:

- fistula contamination or sepsis, which may require removal of the prosthesis and/or appropriate medical management of sepsis
- accidental aspiration of the prosthesis into the airway, which may require removal by a physician
- accidental ingestion of the prosthesis into the esophagus and/or GI tract
- occasional extrusion of the prosthesis, requiring replacement
- fistula dilation resulting in leaking of fluids around the prosthesis
- inflammatory reaction of or around the fistula
- intractable leakage around the prosthesis, requiring adjustment, re-fitting, or replacement of the prosthesis or other surgical/medical management options
- tearing or other damage to the prosthesis from improper use

### **INSTRUCTIONS FOR USE**

*The following instructions are provided by Eric D. Blom, Ph.D. and Stephen B. Freeman, M.D., F.A.C.S.*

*Note: This device is not intended to be placed or removed by the patient. All instructions for use are for the physician only.*

Any care, handling, or adjustment of the device should only be done with a direct view of the fistula so it is clearly visible (e.g. with a bright light focused directly on the fistula). The user's hands should be thoroughly clean to prevent contaminants from entering in or around the fistula. Gloves are recommended.

#### **Fistula Assessment**

Prior to insertion of a Blom-Singer® Adjustable Bi-Flanged Fistula Prosthesis, examine the fistula tract to grossly estimate its length, diameter, and direction. Fistula tract length and direction may also be assessed radiographically using an approved/appropriate contrast substance, as necessary. Select a fistula prosthesis flange diameter that provides at least an 8 to 10 mm overlap circumferentially around the outside edges of the widest aspect of the fistula. A 12.5 mm overlap is recommended when using the 50 mm (large) size fistula prosthesis.

#### **Fistula Prosthesis Insertion**

**Note:** For fistula diameters of 6 mm or less, user should exercise caution during insertion to

avoid fistula dilation.

1. Hand-fold the inner flange of the prosthesis into as small a configuration as possible and grasp it tightly with a hemostat (**Diagram 1**).
2. Gently grasp the end of the beaded stem with the other hand, keeping slight tension on the stem for safety.
3. Gently insert the folded flange through the fistula tract, releasing the flange from its folded configuration at the inner aspect of the fistula tract.
4. Withdraw the hemostat while securely holding the end of the beaded stem (**Diagram 2**).
5. Gently tug the beaded stem outward to confirm that the flange has deployed and is securely seated against the inner tissue surface (**Diagram 3**).
6. Insert the tapered tip of the beaded stem through the center hole in the outer flange. The outer flange should be positioned such that the locking loop is facing out.
7. Advance the outer flange along the beaded stem until a snug inner-flange-to-outer-flange relationship is achieved.
8. Feed the tapered end of the beaded stem through the locking loop, leaving a small amount of slack in the stem to form a curved contour (**Diagram 4**).
9. Inspect/adjust the device and test for leakage as described below under “**Post-Insertion Inspection.**”

#### **Post-Insertion Inspection**

Inspect the placement and contact of the prosthesis flanges around and against the peri-fistula tissue. Circumferential contact is necessary for adequate sealing. Examine and confirm that the peri-fistula tissue has not been dilated or irritated from the insertion procedure. Confirm with the patient that the prosthesis fits comfortably; discontinue use of the prosthesis if proper fit or level of comfort are medically concerning. Instruct the patient to swallow their saliva, or if oral intake is indicated, take small sips of water, advancing up to 120 to 180 ml (4 to 6 oz). Confirm that the leakage, if any, is minimal and acceptable for the patient’s condition. If the amount of leakage is not acceptable, the prosthesis may not be inserted or adjusted properly. Confirm that the flanges are flush against the peri-fistula tissues, adjusting as needed. Re-test for leakage as needed. If the amount of leakage is still not acceptable for the patient’s condition, a larger-sized prosthesis may be indicated. In the event an acceptable amount of leakage cannot be achieved, the patient’s tissue condition and/or anatomy may not be compatible with the device, and use of the prosthesis is contraindicated.

#### **In-Situ Prosthesis Care**

The translucent silicone flanges of the device provide the ability to visually monitor the status of the fistula. Movement of the outer flange is discouraged, as it may disturb the seal to the peri-fistula tissues; if necessary, carefully remove any observed debris from behind the outer flange with a sterile, 6-inch cotton-tipped applicator. After any debris removal or care is completed, ensure the fistula prosthesis seal has not been disrupted and that the fit of the device is maintained. At the physician’s discretion, instruct the patient on any recommended care (e.g. the patient may carefully clean the skin area around the Fistula Prosthesis while taking care to not disturb the seal).

**Caution:** Never use facial or toilet tissue as particles of lint or fabric may be deposited which can disturb the seal of the device. If wiping excess fluid from around the device, do not wipe in the same direction as the beaded stem passes through the locking loop as this may cause the fit of the beaded stem to loosen. If the prosthesis is not functioning properly (e.g. leakage, beaded stem loosens) after performing device care, the patient should have the device evaluated by a physician as soon as possible.

## *Blom-Singer® Adjustable Bi-Flanged Fistula Prosthesis*

### **Fistula Prosthesis Removal**

The device may be removed endoscopically.

1. Detach the outer flange from the beaded stem while maintaining hold of the beaded stem.
2. With forceps, grasp the beaded stem as close as possible to where it attaches to the inner flange of the prosthesis, and withdraw the device out through the mouth.

### **Prosthesis Removal for Prosthesis Downsizing**

In the event that a large fistula has decreased in size such that a smaller prosthesis is indicated or desired, the prosthesis in use may be removed from a still-relatively large fistula by gently pulling it out with a hemostat. A smaller device can then be inserted by the previously described insertion procedure. After insertion, refer to “**Post-Insertion Inspection**” instructions.

### **SPECIAL STORAGE AND/OR HANDLING CONDITIONS**

There are no special storage and/or handling conditions for these devices.

### **DISPOSAL INSTRUCTIONS**

These devices are not biodegradable and can be contaminated when used. Carefully dispose of the devices per local guidelines.

### **ORDERING INFORMATION**

#### **USA**

Blom-Singer products may be ordered directly from InHealth Technologies.

ORDER ON-LINE: [www.inhealth.com](http://www.inhealth.com) TELEPHONE: Toll-Free (800)477-5969 or (805)684- 9337, FAX: Toll-Free (888)371-1530 or (805)684-8594. EMAIL: [order@inhealth.com](mailto:order@inhealth.com) POST: InHealth Technologies, 1110 Mark Avenue, Carpinteria, CA 93013-2918, USA, Attention: Customer Service

#### **INTERNATIONAL CUSTOMERS**

Blom-Singer products are available worldwide through authorized InHealth Technologies distributors. For assistance or additional product information please contact the distributor.

#### **PRODUCT COMPLAINTS/EU SERIOUS INCIDENTS**

If you are dissatisfied with the device or have any questions, please contact

**[productcomplaints@inhealth.com](mailto:productcomplaints@inhealth.com)**.

Telephone: +1-800-477-5969

Fax: +1-888-371-1530

Any serious incident that has occurred in relation to the device should be reported to Freudenberg Medical, LLC as noted above and the competent authority of the EU Member State in which the user and/or patient is established.

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POSTA: InHealth Technologies, 1110 Mark Avenue, Carpinteria, CA 93013-2918, USA,  
Attention: Customer Service

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**SYMBOLS GLOSSARY (USA only)**

SYMBOL	STANDARD No.	SYMBOL TITLE & REF No.	MEANING OF SYMBOL
	EN ISO 15223-1:2016	Manufacturer; 5.1.1	Indicates the medical device manufacturer, as defined in EU Directives or Regulations
	93/42/EEC EU 2017/745	CE Marking; 93/42/EEC Annex XII, EU 2017/745 Annex V	CE marking of conformity
	EN ISO 15223-1:2016	Authorized representative in the European Community; 5.1.2	Indicates the authorized representative in the European Community
	EN ISO 15223-1:2016	Catalogue number; 5.1.6	Indicates the manufacturer's catalogue number so that the medical device can be identified
	EN ISO 15223-1:2016	Batch code; 5.1.5	Indicates the manufacturer's batch code so that the batch or lot can be identified
	EN ISO 15223-1:2016	Use-by date; 5.1.4	Indicates the date after which the medical device is not to be used
	EN ISO 15223-1:2016	Date of manufacture; 5.1.3	Indicates the date when the medical device was manufactured
	EN ISO 15223-1:2016	Caution; 5.4.4	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself
	EN ISO 15223-1:2016	Do not use if package is damaged; 5.2.8	Indicates a medical device that should not be used if the package has been damaged or opened
	EN ISO 15223-1:2016	Consult instructions for use; 5.4.3	Indicates the need for the user to consult the instructions for use
	EN ISO 15223-1:2016	Do not re-use; 5.4.2	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure

**SYMBOLS GLOSSARY (USA only)**

	EN ISO 15223-1:2016	Sterilized using irradiation; 5.2.4	Indicates a medical device that has been sterilized using irradiation
	EN ISO 15223-1:2016	Do not re-sterilize; 5.2.6	Indicates a medical device that is not to be re-sterilized
	ISO 7000:2019	Single sterile barrier system; 3707	To indicate that there is a single sterile barrier system

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