INHEALTH

TECHNOLOGIES®

Blom-Singer® Laryngectomy Tubes

(including **StomaSoft®**)

INSTRUCTIONS FOR USE

Blom-Singer®

voice restoration systems

Blom-Singer® and **InHealth Technologies®** are registered trademarks of **Freudenberg Medical, LLC**.

Patent(s): www.inhealthpatents.com

IFU: <https://inhealth.com/ifu>

Copyright © 2020 Freudenberg Medical, LLC. All rights reserved.


INSTRUCTIONS FOR USE

使用说明 / 使用説明書 / Udhëzime për përdorim / إرشادات الاستخدام / Инструкции за употреба / Upute za korištenje / Návod k použití / Brugsanvisning / Gebruiksaanwijzing / Kasutusjuhend / Käyttöohjeet / Instructions d'utilisation / Gebrauchsanweisung / Οδηγίες χρήσης / օրօտ / Käyttöohjeet / שימוש / Használati utasítás / Notkunarleiðbeiningar / Istruzioni per l'uso / 사용 지침 / Lietošanas norādījumi / Naudojimo instrukcijos / Bruksanvisning / Instrukcja użytkowania / Instruções de uso / Instrucțiuni de utilizare / Инструкции по применению / Uputstvo za upotrebu / Návod na použitie / Instrucciones de uso / Bruksanvisningar / Kullanım Talimatları

Note: All diagrams show the StomaSoft model. / 注意: 所有插图显示的都是 StomaSoft 型号。 / 注: すべての図は StomaSoft モデルを示しています。 / Shënim: Të gjitha diagramet tregojnë modelin "StomaSoft". / .StomaSoft / ملحوظة: تظهر كل الرسوم التوضيحية طراز StomaSoft. / Забелешка: Всички диаграми показват модела StomaSoft. / Napomena: svi dijagrami prikazuju model StomaSoft. / Poznámka: na všech schématech je znázorněn model StomaSoft. / Bemærk: Alle diagrammer viser StomaSoft-modellen. / Let op: alle diagrammen tonen het StomaSoft-model. / Märkus: Kõikidel diagrammidel on näha mudel StomaSoft. / Huomaa: kaikissa kaavioissa esitetään StomaSoft-malli. / Remarque: Tous les diagrammes représentent le modèle StomaSoft. / Hinweis: Alle Diagramme zeigen das StomaSoft-Modell. / Σημείωση: Όλα τα διαγράμματα παρουσιάζουν το μοντέλο StomaSoft. / .StomaSoft / 注意: 所有插图显示的都是 StomaSoft 型号。 / Athugið: Allar skýringarmyndir sýna StomaSoft gerðina. / Nota: tutte le figure mostrano il modello StomaSoft. / 참고: 모든 그림은 StomaSoft 모델을 보여줍니다. / Piezime: visās shēmās ir attēlots StomaSoft modelis. / Pastaba: Visose diagramose rodomas „StomaSoft“ modelis. / Merk: Alle diagrammene viser StomaSoft-modellen. / Uwaga: Wszystkie schematy przedstawiają model StomaSoft. / Observação: todos os diagramas mostram o modelo StomaSoft. / Notā: Toate diagramele indică modelul StomaSoft. / Примечание: на всех схемах показана модель StomaSoft. / Napomena: Svi dijagrami prikazuju StomaSoft model. / Poznámka: Všetky diagramy zobrazujú model StomaSoft. / Nota: Todos los diagramas muestran el modelo StomaSoft. / OBS! Alla diagram visar modellen StomaSoft. / Not: Tüm diyagramlar StomaSoft modelini göstermektedir.

DIAGRAMS

插图 / 図 / Diagramet / الرسم التوضيحي / Диаграми / Dijagrami / Schémata / Diagrammer / Afbeeldingen / Diagrammid / Kaaviot / Diagrammes / Abbildungen / Διαγράμματα / תרשימים / Ábrák / Skýringarmyndir / Diagrammi / 그림 / Diagrammas / Diagramos / Diagrammer / Schematy / Diagramas / Diagramme / Схемы / Slike / Grafické znázornenia / Diagramas / Diagram / Diyagramlar


TABLE OF CONTENTS

目録 / 目次 / Tabela e përmbajtjes / جدول المحتويات / Съдържание / Sadržaj / Obsah /
 Indholdsfortegnelse / Inhoud / Sisukord / Sisällysluettelo / Table des matières /
 Inhaltsverzeichnis / Πίνακας περιεχομένων / תוכן עניינים / Tartalomjegyzék / Efnisyfirlit /
 Sommario / 목차 / Saturs / Turinys / Innhaldsfortegnelse / Spis treści / Índice / Cuprins /
 Содержание / Sadržaj / Obsah / Índice / Innehållsförteckning / İçindekiler

EN	Laryngectomy Tubes	6	HU	Laringektómiás tubus	48
ZH	喉切除管	8	IS	Barkakýlisnámsslöngur	51
JA	喉頭切除チューブ	11	IT	Tubi per laringectomia	54
SQ	Tuba laringektomie	14	KO	Laryngectomy Tubes	56
16	الأنابيب المخصصة لحالات استئصال الحنجرة	AR	LV	Laringektomijas caurulītes	59
BG	Тръби за ларингектомия	19	LT	Laringektomijos vamzdeliai	61
HR	Cijevi za laringektomiju	22	NO	Laryngektomislange	64
CS	Kanyly pro laryngektomii	24	PL	Rurki laryngektomijne	67
DA	Laryngektomituber	27	PT	Tubos de laringectomia	69
NL	Laryngectomiebuisjes	29	RO	Tuburi pentru laringectomie	72
ET	Larüngotoomiatorud	32	RU	Ларингэктомические канюли	75
FI	Laryngektomiaputket	35	SR	Tubusi za laringektomiju	78
FR	Tubes de laryngectomie	37	SK	Laryngektomické trubice	81
DE	Laryngektomietubus	40	ES	Tubos de laringectomía	83
EL	Σωλήνες λαρυγγεκτομής	43	SV	Laryngektomikanyler	86
46	צינוריות לרינגקטומיה	HE	TR	Larenjektomi Tüpleri	89

ENGLISH

DEVICE DESCRIPTION

The Blom-Singer® and Blom-Singer® StomaSoft® Laryngectomy Tubes (tubes) are made of flexible silicone and are available in many sizes. The standard Blom-Singer® Laryngectomy Tube is available in a sterile (by gamma irradiation), non-fenestrated option, and does not allow the attachment of any HME or other stoma accessories.

The StomaSoft® is available non-sterile in both fenestrated and non-fenestrated models, and allows for attachment of HME and other stoma accessories.

A fenestrated (having small openings or holes) tube may be used simultaneously with a voice prosthesis to facilitate unrestricted airflow for tracheoesophageal voice.

INDICATIONS FOR USE


The tube is designed for use by a single patient to maintain tracheostoma patency after total laryngectomy. The StomaSoft® model also allows attachment of HME and other stoma accessories.

CONTRAINDICATIONS

The tube should only be used by clinicians and patients that have been trained in its use. Carefully monitor the tracheostoma tissue during radiation therapy. The tube should not be used when experiencing inflammation, infection, bleeding, discharge, or excessive coughing. The tube is not meant for use with mechanical ventilators. The tube should not be used if its presence obstructs breathing.

WARNINGS AND PRECAUTIONS

The clinician is to instruct the patient on the usage of the device and to provide the patient these Instructions For Use. Wash hands with soap and water prior to handling the tube. Only use a lint-free cloth to wipe the tube to avoid getting lint into the airway. Do not sterilize or re-sterilize the device. Do not use solvent or petroleum-based products for cleaning or lubricating the tube; these materials may damage the tube. Only Blom-Singer® stoma accessories can be used with the StomaSoft® Tube; use of other accessories may result in complications. Never insert any cleaning device into the tube while it is positioned in the neck; only clean the device after it has been removed from the tracheostoma. The tube is for single patient use only and must not be shared between patients; infection can occur with shared use, which could lead to illness or serious injury. Inspect the tube before each use for physical damage, such as cracks or tears, and do not use if damaged. If physical damage is found, contact InHealth Technologies Product Complaints. Do not make modifications to the tube. Do not use the tube in the event of the packaging or sterile packaging being damaged or unintentionally opened prior to use;

discard and replace it with a new tube. Patients should consult a clinician if the tube no longer fits properly (e.g. difficulty inserting the tube), if they have difficulty breathing with the tube in place, or if they experience inflammation, infection, bleeding, discharge, or excessive coughing while using the tube.

POTENTIAL COMPLICATIONS

Complications that require potential removal of the tube, cleaning, and/or medical treatment include but are not limited to: stoma contamination/infection, accidental inhalation of an inappropriately sized tube into the airway, structural failure from excessive usage and/or cleaning, and obstruction of the tube by mucus.

INSTRUCTIONS FOR USE

Please refer to the diagrams at the front of this instruction manual.

Wash hands with soap and water prior to inserting or removing the tube. The user should be positioned in front of a mirror with a bright light shining directly on the stoma (**Diagram 1**).

Insertion of the Laryngectomy Tube

1. If desired, apply a small amount of medical grade, water-soluble lubricant to the rim of the stoma to help with insertion (**Diagram 2**).
2. If needed, slightly fold or compress the insertion tip of the tube to conform to the contour of the tracheostoma (**Diagram 3**).
3. Insert the insertion tip of the tube to the lower rim of the stoma and insert completely into the stoma until its larger funnel shaped end rests against the peristomal skin (**Diagram 4**).

Attaching an Optional Laryngectomy or Tracheostomy Tube Holder (sold separately)

1. Attach the straps by passing them through the openings in the side of the tube such that the tabs on the straps point away from the tube (**Diagram 5**).
2. Pass the straps around the back of the neck and velcro or tie to secure. The straps should be loose enough to permit one finger width between the backside of the straps and the surface of the neck.

Removal of the Laryngectomy Tube

1. If using a tube holder, untie or release the holder straps (**Diagram 6**).
2. Remove the tube from the stoma (**Diagram 7**).

Attachment or Removal of Accessory Devices (StomaSoft® models only)

To **insert** an accessory device: hold the outer rim of the accessory device and press into the top of the tube until fully seated (**Diagram 8**).

To **remove** an accessory device: hold the top of the tube and grasp the outer rim of the accessory device. Pull until free from the tube (**Diagram 9**).

Cleaning and Disinfection of the Laryngectomy Tube

1. Remove the tube from the tracheostoma and remove any attached devices.
2. Rinse the tube under warm running tap water.
3. Soak the tube in hot water (122-140°F / 50-60°C) containing at least two drops of dish soap per 8 oz. (235mL) of water for at least 15 minutes.
4. Remove tube from the soap and water solution and rinse under warm running tap water. Gently clean the inside and outside of the tube with a cleaning foam swab (sold separately from InHealth Technologies). Do not use cotton swabs to clean the tube as fibers may be left behind (**Diagram 10**).
5. If the tube is fenestrated (StomaSoft® only), clean the fenestrations in a similar manner using the Blom-Singer Cleaning Brush (sold separately from InHealth Technologies) (**Diagram 11**).

6. Disinfect the tube by soaking in a 3% hydrogen peroxide solution overnight.
7. Rinse thoroughly under running tap water.
8. Allow tube to air dry.

SPECIAL STORAGE AND/OR HANDLING CONDITIONS

There are no special storage and/or handling conditions for the tube.

DISPOSAL INSTRUCTIONS

The tube is not biodegradable and can be contaminated when used. Carefully dispose of the tube per local guidelines.

ORDERING INFORMATION

Online Ordering www.inhealth.com
USA
InHealth Technologies
1110 Mark Avenue
Carpinteria, CA 93013-2918
Customer Service: +1-800-477-5969
Fax: +1-888-371-1530
Email: order@inhealth.com

INTERNATIONAL Contact InHealth Technologies Customer Service at +1-805-684-9337 for international distributor contact.

PRODUCT COMPLAINTS/EU SERIOUS INCIDENTS

If you are dissatisfied with the device or have any questions, please contact productcomplaints@inhealth.com.

Telephone: +1-800-477-5969

Fax: +1-888-371-1530

Any serious incident that has occurred in relation to the device should be reported to Freudenberg Medical, LLC as noted above and the competent authority of the EU Member State in which the user and/or patient is established.

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SYMBOLS GLOSSARY (USA only)

SYMBOL	STANDARD No.	SYMBOL TITLE & REF No.	MEANING OF SYMBOL
	EN ISO 15223-1:2016	Manufacturer; 5.1.1	Indicates the medical device manufacturer, as defined in EU Directives or Regulations
	93/42/EEC EU 2017/745	CE Marking; 93/42/EEC Annex XII, EU 2017/745 Annex V	CE marking of conformity
	EN ISO 15223-1:2016	Authorized representative in the European Community; 5.1.2	Indicates the authorized representative in the European Community
	EN ISO 15223-1:2016	Catalogue number; 5.1.6	Indicates the manufacturer's catalogue number so that the medical device can be identified
	EN ISO 15223-1:2016	Batch code; 5.1.5	Indicates the manufacturer's batch code so that the batch or lot can be identified
	EN ISO 15223-1:2016	Use-by date; 5.1.4	Indicates the date after which the medical device is not to be used
	EN ISO 15223-1:2016	Date of manufacture; 5.1.3	Indicates the date when the medical device was manufactured
	EN ISO 15223-1:2016	Caution; 5.4.4	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself
	EN ISO 15223-1:2016	Do not use if package is damaged; 5.2.8	Indicates a medical device that should not be used if the package has been damaged or opened
	EN ISO 15223-1:2016	Consult instructions for use; 5.4.3	Indicates the need for the user to consult the instructions for use
	EN ISO 15223-1:2016	Non-sterile; 5.2.7	Indicates a medical device that has not been subjected to a sterilization process

SYMBOLS GLOSSARY (USA only)

	EN ISO 15223-1:2016	Sterilized using irradiation; 5.2.4	Indicates a medical device that has been sterilized using irradiation
	EN ISO 15223-1:2016	Do not resterilize; 5.2.6	Indicates a medical device that is not to be resterilized
	ISO 7000:2019	Single patient multiple use; 3706	To indicate that the medical device may be used multiple times (multiple procedures) on a single patient
	ISO 7000:2019	Single sterile barrier system; 3707	To indicate that there is a single sterile barrier system

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