
INHEALTH

TECHNOLOGIES®

Blom-Singer® Adhesive Housing

INSTRUCTIONS FOR USE

Blom-Singer®

voice restoration systems

Blom-Singer® and **InHealth Technologies®** are registered trademarks of **Freudenberg Medical, LLC**.

Patent(s): www.inhealthpatents.com

IFU: <https://inhealth.com/ifu>

Copyright © 2020 Freudenberg Medical, LLC. All rights reserved.

INSTRUCTIONS FOR USE

使用说明 / 使用説明書 / Udhëzime për përdorim / إرشادات الاستخدام / Инструкции за употреба / Upute za korištenje / Návod k použití / Brugsanvisning / Gebruiksaanwijzing / Kasutusjuhend / Käyttöohjeet / Instructions d'utilisation / Gebrauchsanweisung / Οδηγίες χρήσης / הוראות שימוש / Használati utasítás / Notkunarleiðbeiningar / Istruzioni per l'uso / 사용 지침 / Lietošanas norādījumi / Naudojimo instrukcijos / Bruksanvisning / Instrukcja użytkowania / Instruções de uso / Instrucțiuni de utilizare / Инструкции по применению / Uputstvo za upotrebu / Návod na použitie / Instrucciones de uso / Bruksanvisningar / Kullanım Talimatlar

TABLE OF CONTENTS

目録/目次 / Tabela e përmbajtjes / جدول المحتويات / Съдържание / Sadržaj / Obsah / Indholdsfortegnelse / Inhoud / Sisukord / Sisällysluettelo / Table des matières / Inhaltsverzeichnis / Πίνακας περιεχομένων / תוכן עניינים / Tartalomjegyzék / Efnisyfirlit / Sommario / 목차 / Saturs / Turinys / Innehållsförtegnelse / Spis treści / Índice / Cuprins / Содержание / Sadržaj / Obsah / Índice / Innehållsförteckning / İçindekiler

EN	Adhesive Housing	4	HU	Öntapadós szelepház	36
ZH	胶粘剂外壳	6	IS	Sjálflímandi hlíf	38
JA	粘着ハウジング	8	IT	Alloggiamento adesivo	40
SQ	Mbajtëse ngjitëse	10	KO	접착 하우징	42
12	مبيت لاصق	AR	LV	Stiprinājuma plāksne ar limenti	44
BG	Адхезивна рамка	14	LT	Lipnus korpusas	46
HR	Ljepljivo ležište	16	NO	Heftende innfatning	48
CS	Nalepovací kryt	18	PL	Ośłona samoprzylepna	50
DA	Selvklæbende ventilhus	20	PT	Base adesiva	52
NL	Zelfklevende behuizing	21	RO	Carcasă adezivă	54
ET	Liimitav mansett	24	RU	Адгезивный корпус	56
FI	Liimakotelo	26	SR	Lepljivo kućište	58
FR	Logement adhésif	27	SK	Lepiaci kryt	60
DE	Haftgehäuse	30	ES	Alojamiento adhesivo	62
EL	Αυτοκόλλητη θήκη	32	SV	Självhäftande hölje	64
34	תושבת נדבקת	HE	TR	Yapışkan Muhafaza	66

ENGLISH

BLOM-SINGER® ADHESIVE HOUSINGS**DEVICE DESCRIPTION**

- Single-use and preassembled
- Available in several shapes and sizes
- For use by laryngectomy patients
- Acts as a holder for Blom-Singer tracheostoma accessories (e.g. Heat and Moisture Exchange (HME) System, Adjustable Tracheostoma Valve II (ATSV II), Shower Guard, etc.)

INDICATIONS (Reasons to prescribe the device or procedure)

The housings are designed to be applied to the skin surrounding the tracheostoma or stoma.

CONTRAINDICATIONS (Reasons that make it inadvisable to prescribe the particular device or procedure)

The housings should only be used by people who have been trained on their use. The housings are for single patient use.

PRECAUTIONS

Wash hands with soap and water prior to using to avoid contamination.

Do not use solvents or petroleum-based products for cleaning the housings, as this could damage the products.

Voice prosthesis neck straps may obstruct the seal around the stoma. To reduce the possibility of air leaks, cut any excess straps, as needed.

Do not trim below the portion of the strap with the small safety hole (Figure 1).


Figure 1

WARNINGS

Rare occurrence of skin irritation, rash or minor infection can occur if skin cleaning is not adequate. If you experience excessive skin irritation or rash, stop use immediately and contact your clinician. The housings must not be reused. A reused housing will not provide a proper seal resulting in the failure of the device. The device may not adhere properly if the skin is not completely dry after cleaning or if the use instructions are not followed. Cross contamination leading to patient infection can occur if the device is shared between different patients. Inspect the housing for physical damage, such as cracks or tears, and do not use if damaged. If physical damage is found, contact InHealth Technologies Product Complaints. Keep silicone adhesives, skin barrier, alcohol fumes and adhesive remover or any fumes out of the entryway of the tracheostoma as they may cause coughing, skin or lung irritation and/or skin rash.

INSTRUCTIONS FOR USE

1. Stand in front of a mirror with a bright light shining directly on the application site. Clean the skin area around the stoma with mild soap and warm water. Allow the skin to dry completely.
2. Gently wipe the skin area with an alcohol pad.
3. If a protective skin barrier is needed, apply 1-2 layers.
4. If a stronger seal is necessary, apply a thin coat of medical grade liquid silicone adhesive to the skin prior to attaching the housing. Allow the silicone adhesive to dry for at least 3 to 4 minutes before attaching the housing.
5. Remove the housing from the plastic bag and place on a clean, dry surface.
6. Using the pull-tab, slowly peel the release liner from the housing to expose the adhesive. Do not allow the adhesive film to fold onto itself.
7. Stretch the skin gently to position the housing directly over the stoma. Properly position the housing with the pull-tab down unless otherwise instructed by your clinician (Figure 2).


Figure 2

8. Press the housing against the skin surface and gently rub the surface to avoid air pockets.

How to remove the housing

1. Wash hands with soap and water.
2. Stand in front of a mirror with a bright light shining directly on the application site.
3. To avoid skin irritation, carefully and slowly remove the housing from the skin.
4. Remove any adhesive from the skin with medical-grade adhesive remover.
5. Clean the skin with soap and water.

COMMON PROBLEMS

Follow the instructions below for common problems. If not resolved, contact your clinician for assistance.

Problem

Adhesive housing will not adhere securely to skin during use.

Causes and remedies

1. Excessive resistance to airflow through the voice prosthesis and/or throat. See your clinician for assessment.
2. Insufficient adhesive application and/or use. Review instructions for use.
3. Failure to allow the liquid silicone adhesive to dry thoroughly, at least 3-4 minutes, prior to attaching the housing.
4. Failure to carefully remove the hands-free speech device or HME/ HumidiFilter System prior to coughing.
5. Excess internal pressure against the adhesive seal caused by voicing with excess effort.

Problem

Skin irritation

Causes and remedies

1. Avoid repeated application of the housing on the same day.
2. Remove adhesive slowly and carefully to avoid irritation.
3. Apply protective skin barrier to skin surrounding tracheostoma prior to applying silicone adhesive or housing.
4. If difficulty with your adhesive cannot be easily remedied, consult your clinician for assistance.

SPECIAL STORAGE AND/OR HANDLING CONDITIONS

There are no special storage and/or handling conditions for this device.

DISPOSAL INSTRUCTIONS

This device is not biodegradable and can be contaminated when used. Carefully dispose of the device per local guidelines.

ORDERING INFORMATION

Online Ordering

www.inhealth.com

USA

InHealth Technologies

1110 Mark Avenue

Carpinteria, CA 93013-2918

Customer Service: +1-800-477-5969

Fax: +1-888-371-1530

Email: order@inhealth.com

INTERNATIONAL

Contact InHealth Technologies Customer Service at +1-805-684-9337

for international distributor contact.

PRODUCT COMPLAINTS/EU SERIOUS INCIDENTS

If you are dissatisfied with the device or have any questions, please contact productcomplaints@inhealth.com.

Telephone: +1-800-477-5969

Fax: +1-888-371-1530

Any serious incident that has occurred in relation to the device should be reported to Freudenberg Medical, LLC as noted above and the competent authority of the EU Member State in which the user and/or patient is established.

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SYMBOLS GLOSSARY (USA only)

SYMBOL	STANDARD No.	SYMBOL TITLE & REF No.	MEANING OF SYMBOL
	EN ISO 15223-1:2016	Manufacturer; 5.1.1	Indicates the medical device manufacturer, as defined in EU Directives or Regulations
	93/42/EEC EU 2017/745	CE Marking; 93/42/EEC Annex XII, EU 2017/745 Annex V	CE marking of conformity
	EN ISO 15223-1:2016	Authorized representative in the European Community; 5.1.2	Indicates the authorized representative in the European Community
	EN ISO 15223-1:2016	Batch code; 5.1.6	Indicates the manufacturer's catalogue number so that the medical device can be identified
	EN ISO 15223-1:2016	Batch code; 5.1.5	Indicates the manufacturer's batch code so that the batch or lot can be identified
	EN ISO 15223-1:2016	Use-by date; 5.1.4	Indicates the date after which the medical device is not to be used
	EN ISO 15223-1:2016	Caution; 5.4.4	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself
	EN ISO 15223-1:2016	Do not use if package is damaged; 5.2.8	Indicates a medical device that should not be used if the packaged has been damaged or opened
	EN ISO 15223-1:2016	Consult instructions for use; 5.4.3	Indicates the need for the user to consult the instructions for use
	EN ISO 15223-1:2016	Do not re-use; 5.4.2	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure
	ISO 7000:2019	Packaging unit; 2794	To indicate the number of pieces in the package

INHEALTH

TECHNOLOGIES®

UNITED STATES

1110 Mark Avenue
Carpinteria, CA 93013-2918 USA
Toll-Free: 800.477.5969
Toll-Free Fax: 888.371.1530

INTERNATIONAL

Tel: +1.805.684.9337
Fax: +1.805.684.8594
E-mail: order@inhealth.com

ORDER ONLINE

www.inhealth.com

Freudenberg Medical, LLC
1110 Mark Avenue
Carpinteria, CA 93013-2918 USA
www.inhealth.com

EMERGO EUROPE
Westervoortsedijk 60
6827 AT Arnhem
The Netherlands