
INHEALTH

TECHNOLOGIES®

Blom-Singer® Adhesive Housing

INSTRUCTIONS FOR USE

Blom-Singer®

voice restoration systems

Blom-Singer® and InHealth Technologies® are registered trademarks of **Freudenberg Medical, LLC**.

IFU: <https://inhealth.com/ifu>








Copyright © 2024 Freudenberg Medical, LLC. All rights reserved.

INSTRUCTIONS FOR USE

使用说明 / 使用説明書 / Udhëzime për përdorim / إرشادات الاستخدام / Инструкции за употреба / Upute za korištenje / Návod k použití / Brugsanvisning / Gebruiksaanwijzing / Kasutusjuhend / Käyttööhjeet / Instructions d'utilisation / Gebrauchsanweisung / Οδηγίες χρήσης / הוראות שימוש / Használati utasítás / Notkunarleiðbeiningar / Istruzioni per l'uso / 사용 지침 / Lietošanas norādījumi / Naudojimo instrukcijos / Bruksanvisning / Instrukcja użytkowania / Instruções de uso / Instrucțiuni de utilizare / Инструкции по применению / Uputstvo za upotrebu / Návod na použitie / Instrucciones de uso / Bruksanvisningar / Kullanım Talimatlar

DIAGRAMS

插图 / 図 / Diagramet / الرسم التوضيحي / Схеми / Dijagrami / Schémata / Diagrammer / Afbeeldingen / Diagrammid / Kaaviot / Diagrammes / Abbildungen / Διαγράμματα / תרשימים / Ábrák / Skýringarmyndir / Diagrammi / 그림 / Diagrammas / Diagramos / Diagrammer / Schematy / Diagramas / Diagrame / Схемы / Crteži / Grafické znázornenia / Diagramas / Diagram / Diyagramlar




TABLE OF CONTENTS

目录/目次 / Tabela e përmbajtjes / جدول المحتويات / Съдържание / Sadržaj / Obsah / Indholdsfortegnelse / Inhoud / Sisukord / Sisällysluettelo / Table des matières / Inhaltsverzeichnis / Πίνακας περιεχομένων / תוכן עניינים / Tartalomjegyzék / Efnisyfirlit / Sommario / 목차 / Saturs / Turinys / Innholdsfortegnelse / Spis treści / Índice / Cuprins / Содержание / Sadržaj / Obsah / Índice / Innehållsförteckning / İçindekiler

EN	Adhesive Housing	6	HU	Öntapadós szelepház	41
ZH	胶粘剂外壳	8	IS	Hásing með límeftni	43
JA	粘着ハウジング	10	IT	Alloggiamento adesivo	45
SQ	Mbajtëse ngjyhtëse	12	KO	접착 하우징	47
14	غلاف لاصق	AR	LV	Stiprinājuma plāksne ar līmēti	49
BG	Адхезивна рамка	16	LT	Lipnusis korpusas	51
HR	Ljepljivo ležište	18	NO	Heftende innfatning	54
CS	Nalepovací kryt	20	PL	Ośłona samoprzylepna	56
DA	Selvklebende ventilhus	22	PT	Carcaça adesiva	58
NL	Zelfklevende behuizing	25	RO	Carcasă adezivă	60
ET	Liimitav mansett	27	RU	Адгезивный корпус	63
FI	Liimakotelo	29	SR	Lepljivo kućište	65
FR	Support adhésif	31	SK	Lepiaci kryt	67
DE	Haftgehäuse	34	ES	Alojamiento adhesivos	69
EL	Αυτοκόλλητη θήκη	36	SV	Självhäftande hölje	72
39	תושבת נדבקת	HE	TR	Yapışkan Muhafaza	74

ENGLISH

BLOM-SINGER® ADHESIVE HOUSINGS

DEVICE DESCRIPTION

- Single-use and preassembled
- Available in several shapes, sizes and styles
- For use by laryngectomy patients
- Acts as an attachment point for Blom-Singer tracheostoma accessories (e.g., Heat and Moisture Exchange (HME), Shower Guard)

INDICATIONS (Reasons to prescribe the device or procedure)

The housings are designed to be applied to the skin surrounding the tracheostoma or stoma.

CONTRAINDICATIONS (Reasons that make it inadvisable to prescribe the particular device or procedure)

The housings should only be used by people who have been trained on their use. The housings are for single patient use.

PRECAUTIONS

- Wash hands with soap and water prior to using to avoid contamination.
- Do not use solvents or petroleum-based products for cleaning the housings, as this could damage the products.
- Voice prosthesis neck straps may obstruct the seal around the stoma. To reduce the possibility of air leaks, cut any excess straps, as needed.
- Do not trim below the portion of the strap with the small safety hole (Diagram 1).

WARNINGS

Rare occurrence of skin irritation, rash or minor infection can occur if skin cleaning is not adequate. If you experience excessive skin irritation or rash, stop use immediately and contact your clinician. The housings must not be reused. A reused housing will not provide a proper seal resulting in the failure of the device. The device may not adhere properly if the skin is not completely dry after cleaning or if the use instructions are not followed. Cross contamination leading to patient infection can occur if the device is shared between different patients. Inspect the housing for physical damage, such as cracks or tears, and do not use if damaged. If physical damage is found, contact InHealth Technologies Product Complaints. Keep silicone adhesives, skin barrier, alcohol fumes and adhesive remover or any fumes out of the entryway of the tracheostoma as they may cause coughing, skin or lung irritation and/or skin rash.

INSTRUCTIONS FOR USE

Prior to using product, ensure area around the stoma is clean and dry.

1. Stand in front of a mirror with a bright light shining directly on the application site. Clean the skin area around the stoma with mild soap and warm water. Allow the skin to dry completely.
2. Gently wipe the skin area with an alcohol pad.
3. If a protective skin barrier is needed, apply 1-2 layers.
4. If a stronger seal is necessary, apply a thin coat of medical grade liquid silicone adhesive to the skin prior to attaching the housing. Allow the silicone adhesive to dry for at least 3 to 4 minutes before attaching the housing.
5. Remove the housing from the plastic bag and place on a clean, dry surface.
 - a. For AccuFit and HydroFit: Using the pull-tab, slowly peel the release liner from the housing to expose the adhesive (Diagram 2A). Be careful not to allow the adhesive housing to fold onto itself.
 - b. For TruFit: Tear away the center third of the release liner to expose this area of adhesive (Diagram 2B). Do not allow the adhesive housing to fold onto itself.
6. Stretch the skin gently to position the housing directly over the stoma. Properly position the housing with the pull-tab down unless otherwise instructed by your clinician (Diagram 3A, 3B).
 - a. For AccuFit and HydroFit: Press the housing against the skin surface and gently rub the surface to avoid air pockets (Diagram 4A).

- b. For TruFit: Press the housing against the skin surface and gently rub the center surface to avoid air pockets (Diagram 4B).
7. For TruFit only: One at a time, slowly peel the remaining 2 release liners from the housing to expose the adhesive (Diagram 5B, 6B).
8. For TruFit only: Press the newly exposed adhesive against the skin surface and gently rub the surfaces to avoid air pockets (Diagram 7B).

How to remove the housing

1. Wash hands with soap and water.
2. Stand in front of a mirror with a bright light shining directly on the application site.
3. To avoid skin irritation, carefully and slowly remove the housing from the skin.
4. Remove any adhesive from the skin with medical-grade adhesive remover.
5. Clean the skin with soap and water.

COMMON PROBLEMS

Follow the instructions below for common problems. If not resolved, contact your clinician for assistance.

Problem: Adhesive housing will not adhere securely to skin during use.

Causes and remedies

- Excessive resistance to airflow through the voice prosthesis and/or throat. See your clinician for assessment.
- Insufficient adhesive application and/or use. Review instructions for use.
- Failure to allow the liquid silicone adhesive to dry thoroughly, at least 3-4 minutes, prior to attaching the housing.
- Failure to carefully remove the hands-free speech device or HME/ HumidiFilter System prior to coughing.
- Excess internal pressure against the adhesive seal caused by voicing with excess effort.

Problem: Skin irritation

Causes and remedies

- Avoid repeated application of the housing on the same day.
- Remove adhesive slowly and carefully to avoid irritation.
- Apply protective skin barrier to skin surrounding tracheostoma prior to applying silicone adhesive or housing.
- If difficulty with your adhesive cannot be easily remedied, consult your clinician for assistance.

SPECIAL STORAGE AND/OR HANDLING CONDITIONS

There are no special storage and/or handling conditions for this device.

DISPOSAL INSTRUCTIONS

This device is not biodegradable and can be contaminated when used. Carefully dispose of the device per local guidelines.

ORDERING INFORMATION

USA

InHealth Technologies® products may be ordered directly from Freudenberg Medical, LLC, d/b/a "InHealth Technologies."

Telephone: Toll-Free 800.477.5969 or +1.805.684.9337

Fax: Toll-Free 888.371.1530 or +1.805.684.8594

Email: order@inhealth.com Order Online: www.inhealth.com

Mail: InHealth Technologies 1110 Mark Avenue Carpinteria, CA 93013-2918 USA

INTERNATIONAL

InHealth Technologies® products are available through authorized distributors worldwide. For international distributor information, please contact our customer service department.

PRODUCT COMPLAINTS/EU SERIOUS INCIDENTS

If you are dissatisfied with the device or have any questions, please contact
productcomplaints@inhealth.com

Telephone: +1.800.477.5969

Fax: +1.888.371.1530

Any serious incident that has occurred in relation to the device should be reported to Freudenberg Medical, LLC as noted above and the competent authority of the EU Member State in which the user and/or patient is established.

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SYMBOLS GLOSSARY (USA only)

SYMBOL	STANDARD No.	SYMBOL TITLE & REF No.	MEANING OF SYMBOL
	EN ISO 15223-1:2016	Manufacturer; 5.1.1	Indicates the medical device manufacturer, as defined in EU Directives or Regulations
	93/42/EEC EU 2017/745	CE Marking; 93/42/EEC Annex XII, EU 2017/745 Annex V	CE marking of conformity
	EN ISO 15223-1:2016	Authorized representative in the European Community; 5.1.2	Indicates the authorized representative in the European Community
	EN ISO 15223-1:2016	Batch code; 5.1.6	Indicates the manufacturer's catalogue number so that the medical device can be identified
	EN ISO 15223-1:2016	Batch code; 5.1.5	Indicates the manufacturer's batch code so that the batch or lot can be identified
	EN ISO 15223-1:2016	Use-by date; 5.1.4	Indicates the date after which the medical device is not to be used
	EN ISO 15223-1:2016	Caution; 5.4.4	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself
	EN ISO 15223-1:2016	Do not use if package is damaged; 5.2.8	Indicates a medical device that should not be used if the packaged has been damaged or opened
	EN ISO 15223-1:2016	Consult instructions for use; 5.4.3	Indicates the need for the user to consult the instructions for use
	EN ISO 15223-1:2016	Do not re-use; 5.4.2	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure
	ISO 7000:2019	Packaging unit; 2794	To indicate the number of pieces in the package

INHEALTH

TECHNOLOGIES®

UNITED STATES

1110 Mark Avenue
Carpinteria, CA 93013-2918 USA
Toll-Free: 800.477.5969
Toll-Free Fax: 888.371.1530

INTERNATIONAL

Tel: +1.805.684.9337
Fax: +1.805.684.8594
E-mail: order@inhealth.com

ORDER ONLINE

www.inhealth.com

Freudenberg Medical, LLC

1110 Mark Avenue
Carpinteria, CA 93013-2918 USA
www.inhealth.com

EMERGO EUROPE

Westervoortsedijk 60
6827 AT Arnhem
The Netherlands