
INHEALTH

TECHNOLOGIES®

Blom-Singer® SpeakFree® HME

Hands-Free Valve with *EasyFlow*® HME
-or-
Hands-Free Valve with *ClassicFlow*® HME

INSTRUCTIONS FOR USE

Blom-Singer®

voice restoration systems

Blom-Singer® and InHealth Technologies® are registered trademarks of **Freudenberg Medical, LLC**.

Patent(s): www.inhealthpatents.com

IFU: <https://inhealth.com/ifu>

Copyright © 2020 Freudenberg Medical, LLC. All rights reserved.

INSTRUCTIONS FOR USE

使用说明 / 使用說明書 / Udhëzime për përdorim / إرشادات الاستخدام / Инструкции за употреба / Upute za korištenje / Návod k použití / Brugsanvisning / Gebruiksaanwijzing / Kasutusjuhend / Käyttöohjeet / Instructions d'utilisation / Gebrauchsanweisung / Οδηγίες χρήσης / הוראות שימוש / Használati utasítás / Notkunarleiðbeiningar / Istruzioni per l'uso / 사용 지침 / Lietošanas norādījumi / Naudojimo instrukcijos / Bruksanvisning / Instrukcja użytkowania / Instruções de uso / Instrucțiuni de utilizare / Инструкции по применению / Uputstvo za upotrebu / Návod na použitie / Instrucciones de uso / Bruksanvisningar / Kullanım Talimatlar

DIAGRAMS

插图 / 図 / Diagramet / الرسم التوضيحي / Схеми / Dijagrami / Schémata / Diagrammer /
 Afbeeldingen / Diagrammid / Kaaviot / Diagrammes / Abbildungen / Διαγράμματα / תרשימים /
 Ábrák / Skýringarmyndir / Diagrammi / 그림 / Diagrammas / Diagramos / Diagrammer /
 Schematy / Diagramas / Diagrame / Схеми / Crteži / Grafické znázornenia / Diagramas /
 Diagram / Diyagramlar


TABLE OF CONTENTS

目录 / 目次 / Tabela e përmbajtjes / جدول المحتويات / Съдържание / Sadržaj / Obsah /
 Indholdsfortegnelse / Inhoud / Sisukord / Sisällysluettelo / Table des matières /
 Inhaltsverzeichnis / Πίνακας περιεχομένων / תוכן עניינים / Tartalomjegyzék / Efnisyfirlit /
 Sommario / 목차 / Saturs / Turinys / Innholdsfortegnelse / Spis treści / Índice / Cuprins /
 Содержание / Sadržaj / Obsah / Índice / Innehållsförteckning / İçindekiler

EN	SpeakFree® HME	06	HU	SpeakFree® HME	46
ZH	SpeakFree® HME	08	IS	SpeakFree® HME	48
JA	SpeakFree® HME	10	IT	SpeakFree® HME	50
SQ	SpeakFree® HME	13	KO	SpeakFree® HME	53
17	SpeakFree® HME	AR	LV	SpeakFree® HME	55
BG	SpeakFree® HME	19	LT	SpeakFree® HME	58
HR	SpeakFree® HME	22	NO	SpeakFree® HME	60
CS	SpeakFree® HME	24	PL	SpeakFree® HME	62
DA	SpeakFree® HME	26	PT	SpeakFree® HME	65
NL	SpeakFree® HME	28	RO	SpeakFree® HME	67
ET	SpeakFree® HME	30	RU	SpeakFree® HME	70
FI	SpeakFree® HME	33	SR	SpeakFree® HME	73
FR	SpeakFree® HME	35	SK	SpeakFree® HME	75
DE	SpeakFree® HME	38	ES	SpeakFree® HME	77
EL	SpeakFree® HME	41	SV	SpeakFree® HME	80
43	SpeakFree® HME	HE	TR	SpeakFree® HME	82

ENGLISH

BLOM-SINGER® SpeakFree® HME

DEVICE DESCRIPTION

The Blom-Singer SpeakFree® HME is a single use device that offers the benefits of a heat and moisture exchange system as well as a hands-free valve to achieve tracheostoma (stoma) closure for use of a voice prosthesis.

The device (**diagram 1**) is intended to be attached to either a Blom-Singer tracheostoma adhesive housing or a Blom-Singer StomaSoft® Laryngectomy Tube (**diagram 2**). The device contains a foam filter for heat and moisture exchange and a detachable crossbar. The purpose of the crossbar is to prevent unintended airway obstruction.

INDICATIONS (Reasons to prescribe the device or procedure)

The device is for individuals who have a permanent tracheostoma following total laryngectomy surgery. The device provides humidified and warmed air. The device allows for hands-free speech in patients who have undergone tracheoesophageal puncture (opening made by a surgeon between trachea and esophagus for placement of a voice prosthesis) with voice prosthesis insertion.

CONTRAINDICATIONS (Reasons that make it inadvisable to prescribe the device or procedure)

The device should not be placed over a tracheostoma if its presence obstructs respiration.

WARNINGS AND PRECAUTIONS

Wash hands with soap and water before and after handling the device. Only use a lint-free cloth or tissue on the device to avoid getting lint into the airway. This device is not reusable and may only be used by one patient. There is a risk of infection if this device is shared between patients or reused. Inspect the device before use for physical damage, such as cracks or tears, and do not use if damaged. If physical damage is found, contact InHealth Technologies Product Complaints.

Breathing

Breathing through the device should be comfortable and easy. If breathing is not easy, remove the device and consult your medical professional for assistance. When you begin using the device, a slight increase in airflow resistance may be felt. This resistance is a normal function of the device that may require some adjustment time. The device may need to be removed during instances of extreme physical activity and heavy breathing.

Coughing

The device must be removed before coughing or clearing the airway.

Sleeping

The device must be removed before sleeping to avoid the potential for airway obstruction.

POTENTIAL COMPLICATIONS

The following may occur:

- Labored breathing, which may occur if the device is not changed regularly
- Stoma contamination or infection that may require the use of appropriate antibiotics
- Irritation, including skin irritation around the stoma

INSTRUCTIONS FOR USE

Insertion, Removal and Care

Insertion and removal of any devices worn over the tracheostoma should only be done with the user positioned in front of a mirror with a bright light focused directly on the stoma.

A Blom-Singer tracheostoma housing or StomaSoft® Laryngectomy Tube must first be placed on/around the tracheostoma. Before inserting or removing the device, adjust it to its nearest-closed position by rotating the outer rim clockwise (**diagram 3**).

Insertion

To insert the device, grasp its outer rim and insert it partially into the lower rim of the previously applied housing or laryngectomy tube. Gently push the device into the attachment housing or tube until it snaps securely into place (**diagram 4**). Do not use the crossbar as a handle to insert the device.

Care

The device is designed for daytime/awake use only and is not designed to be disassembled, reassembled, or reused. If the device should become heavily soiled, remove it and replace with a new device. If mucus is removed from the device, never use facial or toilet tissue as particles of lint or fabric may be deposited onto the device that can be inhaled into the airway upon reinsertion. Do not wash the device as this will result in loss of function.

Removal

To remove the device, hold the tracheostoma housing or laryngectomy tube in place with one hand while pulling the device with the other hand until it releases (**diagram 5**). Do not use the crossbar as a handle to remove the device.

Attaching/Removing Crossbar

To attach the optional crossbar, hold the device and press the crossbar onto the top of the device until it snaps into place (**diagram 6**). To remove the crossbar, hold onto the device and pull the crossbar away from the device until it releases (**diagram 7**). Do not attempt to place or remove the crossbar while the device is attached to the tracheostoma.

Device Function

When properly adjusted, the valve remains in an open position during normal, quiet breathing and routine physical activity. As the user exhales to speak, the airflow pushes the valve into the closed position to provide a seal for speech. When the user stops speaking or inhales, the valve returns to its resting open position. This device can also be manually occluded to produce speech by placing a finger or thumb over the opening at the top of the device when the crossbar is detached. With the extremely high pressure of a cough, the valve may be forced partially through the opening in the device. If this happens, remove the device and reposition the valve with a clean fingertip.

Device Sensitivity Adjustment

The device is designed to accommodate the speech and breathing requirements of each individual user by rotating the outer ring of the device clockwise or counterclockwise. Clockwise rotation adjusts the valve to close more easily (**diagram 3**). Counterclockwise rotation adjusts the valve to require greater effort (respiratory airflow) to close (**diagram 8**). Set the device to a more closed position during speech and adjust it to a more open position during increased respiratory exertion (e.g. during physical exercise).

To adjust the device while seated in the tracheostoma housing, the user should first stabilize the base of the housing with one hand while simultaneously rotating the outer ring of the device with the other hand (**diagram 9**).

SPECIAL STORAGE AND/OR HANDLING CONDITIONS

These devices have special storage and/or handling conditions:

- Wash hands after handling the devices, as they contain substances that may be irritating to skin and eyes.
- The package containing the devices should be resealed after opening to preserve the effectiveness of the devices' heat and moisture exchange capabilities.

DISPOSAL INSTRUCTIONS

This product is not biodegradable and can be considered contaminated when used. Carefully dispose of this device per local guidelines.

ORDERING INFORMATION USA

Blom-Singer products may be ordered directly from InHealth Technologies.

ORDER ON-LINE: www.inhealth.com TELEPHONE: Toll-Free (800)477-5969 or (805)684-9337.
FAX: Toll-Free (888)371-1530 or (805)684-8594. EMAIL: order@inhealth.com POST: InHealth Technologies, 1110 Mark Avenue, Carpinteria, CA 93013-2918, USA, Attention: Customer Service

International

Please contact our customer service department for a distributor referral.

PRODUCT COMPLAINTS/EU SERIOUS INCIDENTS

If you are dissatisfied with the device or have any questions, please contact productcomplaints@inhealth.com.

Telephone: +1-800-477-5969

Fax: +1-888-371-1530

Any serious incident that has occurred in relation to the device should be reported to Freudenberg Medical, LLC as noted above and the competent authority of the EU Member State in which the user and/or patient is established.

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SYMBOLS GLOSSARY (USA only)

SYMBOL	STANDARD No.	SYMBOL TITLE & REF No.	MEANING OF SYMBOL
	EN ISO 15223-1:2016	Manufacturer; 5.1.1	Indicates the medical device manufacturer, as defined in EU Directives or Regulations
	93/42/EEC EU 2017/745	CE Marking; 93/42/EEC Annex XII, EU 2017/745 Annex V	CE marking of conformity
	EN ISO 15223-1:2016	Authorized representative in the European Community; 5.1.2	Indicates the authorized representative in the European Community
	EN ISO 15223-1:2016	Catalogue number; 5.1.6	Indicates the manufacturer's catalogue number so that the medical device can be identified
	EN ISO 15223-1:2016	Batch code; 5.1.5	Indicates the manufacturer's batch code so that the batch or lot can be identified
	EN ISO 15223-1:2016	Use-by date; 5.1.4	Indicates the date after which the medical device is not to be used
	EN ISO 15223-1:2016	Caution; 5.4.4	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself
	EN ISO 15223-1:2016	Do not use if package is damaged; 5.2.8	Indicates a medical device that should not be used if the package has been damaged or opened
	EN ISO 15223-1:2016	Consult instructions for use; 5.4.3	Indicates the need for the user to consult the instructions for use
	EN ISO 15223-1:2016	Do not re-use; 5.4.2	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure
	ISO 7000-2019	Packaging unit; 2794	To indicate the number of pieces in the package

INHEALTH

TECHNOLOGIES®

UNITED STATES

1110 Mark Avenue
Carpinteria, CA 93013-2918 USA
Toll-Free: 800.477.5969
Toll-Free Fax: 888.371.1530

INTERNATIONAL

Tel: +1.805.684.9337
Fax: +1.805.684.8594
E-mail: order@inhealth.com

ORDER ONLINE

www.inhealth.com

Freudenberg Medical, LLC
1110 Mark Avenue
Carpinteria, CA 93013-2918 USA
www.inhealth.com

EMERGO EUROPE
Westervoortsedijk 60
6827 AT Arnhem
The Netherlands