# INHEALTH

TECHNOLOGIES®

## Blom-Singer® Indwelling Voice Prostheses

for Primary, Secondary,  
Replacement TEP Procedures

**INSTRUCTIONS FOR USE**

**Blom-Singer®**

voice restoration systems

**Blom-Singer® and InHealth Technologies®** are registered trademarks of **Freudenberg Medical, LLC.**

Patent(s): [www.inhealthpatents.com](http://www.inhealthpatents.com)

IFU: <https://inhealth.com/ifu>







Copyright © 2020 Freudenberg Medical, LLC. All rights reserved.







## INSTRUCTIONS FOR USE

使用说明 / 使用說明書 / Udhëzime për përdorim / إرشادات الاستخدام / Инструкции за употреба / Upute za korištenje / Návod k použití / Brugsanvisning / Gebruiksaanwijzing / Kasutusjuhend / Käyttöohjeet / Instructions d'utilisation / Gebrauchsanweisung / Οδηγίες χρήσης / הוראות שימוש / Használati utasítás / Notkunarleiðbeiningar / Istruzioni per l'uso / 사용 지침 / Lietošanas norādījumi / Naudojimo instrukcijos / Bruksanvisning / Instrukcja użytkowania / Instruções de uso / Instrucțiuni de utilizare / Инструкции по применению / Uputstvo za upotrebu / Návod na použitie / Instrucciones de uso / Bruksanvisningar / Kullanım Talimatlar

## DIAGRAMS

插图 / 図 / Diagramet / الرسم التوضيحي / Схеми / Dijagrami / Schémata / Diagrammer / Afbeeldingen / Diagrammid / Kaaviot / Diagrammes / Abbildungen / Διαγράμματα / תרשימים / Ábrák / Skýringarmyndir / Diagrammi / 그림 / Diagrammas / Diagramos / Diagrammer / Schematy / Diagramas / Diagrame / Схеми / Crteži / Grafické znázornenia / Diagramas / Diagram / Diyagramlar

## TABLE OF CONTENTS

目録 / 目次 / Tabela e përmbajtjes / جدول المحتويات / Съдържание / Sadržaj / Obsah / Indholdsfortegnelse / Inhoud / Sisukord / Sisällysluettelo / Table des matières / Inhaltsverzeichnis / Πίνακας περιεχομένων / תוכן עניינים / Tartalomjegyzék / Efnisyfirlit / Sommarío / 목차 / Saturs / Turinys / Innholdsfortegnelse / Spis treści / Índice / Cuprins / Содержание / Sadržaj / Obsah / Índice / Innehållsförteckning / İçindekiler

<b>EN</b>	Indwelling Voice Protheses	6	<b>HU</b>	Beépített hangprotézisek	93
<b>ZH</b>	留置发音假体	11	<b>IS</b>	Innliggjandi raddventlar	99
<b>JA</b>	留置ボイスプロテーゼ	16	<b>IT</b>	Protesi fonatorie fisse	104
<b>SQ</b>	Proteza permanente zëri	21	<b>KO</b>	유치형 인공후두	110
27	جهاز الحنجرة الهوائية الدائم	<b>AR</b>	<b>LV</b>	levadāmās balss protēzes	115
<b>BG</b>	Имплантируеми гласови протези	31	<b>LT</b>	vidų implantuojamas balso protezas	121
<b>HR</b>	Trajne glasovne proteze	38	<b>NO</b>	Innoperert taleprotese	126
<b>CS</b>	Vnitřní hlasové protézy	43	<b>PL</b>	Stale protezy głosowe	131
<b>DA</b>	Indlagte stemmeproteses	48	<b>PT</b>	Próteses vocais fixas	137
<b>NL</b>	Verblijfstemprotheses	54	<b>RO</b>	Proteze fonatorii permanente	143
<b>ET</b>	Sisemine hääleprotees	60	<b>RU</b>	Постоянные голосовые протезы	149
<b>FI</b>	Pysyvät puheproteesit	65	<b>SR</b>	Trajna glasovna proteza	155
<b>FR</b>	Prothèses vocales à demeure	70	<b>SK</b>	Zavádzacie hlasové protézy	160
<b>DE</b>	Verweil-Stimmprothesen	76	<b>ES</b>	Prótesis fonatorias fijas	166
<b>EL</b>	Μόνιμες φωνητικές προθέσεις	82	<b>SV</b>	Kvarliggande röstproteser	172
89	תותבי דיבור פנימיים	<b>HE</b>	<b>TR</b>	Kalıcı Ses Protezleri	177

## ENGLISH

### **BLOM-SINGER® INDWELLING VOICE PROSTHESES**

for Primary, Secondary, Replacement Tracheoesophageal Puncture (TEP) Procedures

#### **PRODUCT DESCRIPTION**

*Please refer to the diagrams located at the front of this instruction manual.*

The Blom-Singer Indwelling Voice Prosthesis is provided in sterile or non-sterile options for primary, secondary, or replacement tracheoesophageal puncture (TEP: opening made by a surgeon between trachea or "windpipe" and esophagus for placement of a voice prosthesis) procedures. The device is packaged with one (1) flange introducer, one (1) inserter stick, two (2) appropriately-sized gel caps (made of vegetable-based gelatin), and a Voice Prosthesis Cleaning System. The device is made of medical-grade silicone and consists of a one-way silicone flap valve, an esophageal flange, a body that holds the valve assembly, a tracheal flange, and a safety strap. Devices are available in 16 Fr or 20 Fr diameter. Some devices may include silver oxide treatment to valve (Advantage and Dual Valve); a second one-way flap valve (Dual Valve); a titanium body ring (Advantage Hard Valve Assembly); large tracheal and/or esophageal flanges; or a valve with higher resistance (Increased Resistance). Sterile option voice prostheses are sterilized by ethylene oxide.

#### **INDICATIONS** (Reasons to prescribe the device or procedure)

The Blom-Singer Indwelling Voice Prosthesis is indicated for tracheoesophageal voice restoration following total laryngectomy, when placement or replacement of an indwelling voice prosthesis is performed by a qualified, trained clinician.

#### **CONTRAINDICATIONS** (Reasons that make it inadvisable to prescribe the particular device or procedure)

The Blom-Singer indwelling voice prosthesis is a medical product that should be used only by a qualified clinician with experience and training in its use and care.

#### **WARNINGS AND PRECAUTIONS**

Sizing of length and diameter for voice prosthesis selection must be determined by a qualified clinician trained in the use of this particular device. The clinician must carefully determine the device size and model to address the clinical needs of the individual patient and their TEP. The clinician is to instruct patients on the use and care of this device and to provide patients with these instructions for use supplied with the product. This medical device is for single use only. It may not be reused. Reuse of this device could result in infection. If the device is not functioning properly, the patient should have it evaluated by a clinician as soon as possible. The voice prosthesis and accessories must be carefully removed from its packaging and handled in a manner to prevent contamination of the device. Handle the device and accessories carefully to avoid damage. If there are tears, cracks, or structural damage to the prosthesis, cleaning device, or accessories, discontinue use and contact InHealth Technologies Product Complaints. Do not use the device in the event of the packaging or sterile packaging being damaged or unintentionally opened prior to use; discard and replace it with a new prosthesis. When using a sterile device, aseptic handling technique must be used. Use only the appropriately sized gel caps. Do not use petroleum-based products such as Vaseline® to lubricate the voice prosthesis as these products can damage the device. Changes in the anatomy or medical status of the patient may lead to improper fitting and/or function of the device. Dislodgement or extrusion of the device from the TEP may occur and result in ingestion, aspiration (inhalation) or tissue damage. A foreign body such as a voice prosthesis in the airway may cause complications such as acute respiratory (breathing) distress and/or respiratory arrest (stop breathing). Always measure the TEP when changing the device to select the proper device size. If the device is improperly sized, it may cause tissue necrosis (tissue death) and/or device dislodgement. If converting a patient from a large diameter voice prosthesis to a smaller diameter voice

prosthesis, the standard practice is to use a voice prosthesis with an enlarged esophageal or enlarged esophageal and tracheal flanges; use of a voice prosthesis with enlarged flanges may help reduce the risk of peripheral leakage/aspiration and voice prosthesis dislodgement. Individual patient reactions to the device materials may occur. Consult a clinician immediately if there are indications of tissue edema (swelling) and/or inflammation/infection. Evaluate patients with bleeding disorders or if they are undergoing anticoagulant (blood thinning) treatment for the risk of hemorrhage (bleeding) prior to placement or replacement of the device.

#### **Voice Production**

To prevent post-operative complications, the patient should not begin speaking with the voice prosthesis until the clinician has indicated that it is safe to do so. The lumen of the voice prosthesis must be kept clear of blockage for it to function properly to allow the patient to voice. In some users, the inability to relax the muscles of the throat may account for an inability to speak fluently and with minimal effort. This problem requires professional assessment. Patients requiring post-operative radiation may have transient interruption of voice in the third or fourth week of treatment. The device can remain in place as determined by the clinician.

#### **Voice Prosthesis Dislodgement**

Care must be exercised during device insertion, removal, or use of cleaning devices to avoid injury to the TEP or accidental displacement of the device, which could result in aspiration (inhalation) of the device into the trachea (windpipe). Should aspiration occur, the patient should attempt to cough the device out of the trachea. Further medical attention may be necessary if coughing the device out is unsuccessful. Confirm gel cap dissolution and deployment of the esophageal flange to ensure device is securely retained in TEP. If the voice prosthesis is dislodged from the TEP, a Blom-Singer Puncture Dilator or suitable device of the appropriate diameter should be immediately placed in the puncture to keep it from closing and leaking fluids. A replacement device should be inserted within 24 hours. Foreign objects should not be inserted into the device. Inserting objects other than the Blom-Singer cleaning devices may cause dislodgment and subsequent aspiration or ingestion of the voice prosthesis or its components.

#### **Voice Prosthesis Leakage**

When the flap valve fails to close completely, a few drops of fluid may pass through the device from the esophagus (food pipe) to the trachea, which may cause coughing or aspiration. Use of a smaller diameter voice prosthesis than the existing tracheoesophageal puncture size may result in peripheral leakage (leakage around the device). Recurrent leakage of the voice prosthesis should be evaluated by a clinician as leakage could cause aspiration pneumonia. Selection of a different device model/option may be indicated. Gentle handling and pressure should always be used when cleaning the voice prosthesis to avoid device damage, which could cause leakage.

#### **Microbial (Microscopic Organisms) Growth**

Microbial growth deposits on the device may cause valve deformation and failure, i.e., fluid leakage through or around the device and/or an increase in the pressure necessary to voice. Replacement of the device may be required.

#### **Device Insertion/Removal**

The attached safety strap on an indwelling voice prosthesis should only be removed after the esophageal flange has been verified to have deployed in the esophagus. Never attempt to insert or reinsert an indwelling voice prosthesis that has the safety strap removed.

Never remove one indwelling voice prosthesis and insert another without first dilating the TEP and re-measuring the tract to confirm correct voice prosthesis length. The patient should never attempt to insert or remove the device or allow anyone other than a qualified clinician to insert or remove the device. The indwelling voice prosthesis is not a permanent device and requires replacement periodically.

The device may be left in place in the TEP until it has persistent leakage or is not providing adequate voice for speech, or requires resizing. Removal of the indwelling voice prosthesis should only be done by grasping the tracheal flange of the device securely with a locking hemostat.

## COMPLICATIONS

Although rare, the following complications have been identified to occur with silicone protheses of the Blom-Singer type. They include: stoma (opening through neck into trachea or wind pipe) contamination or sepsis, which may require removal of the voice prosthesis and/or appropriate antibiotics; accidental aspiration of the voice prosthesis into the airway, which may require removal by a physician; occasional extrusion of the prosthesis, requiring replacement after dilation of the TEP and additional supervision of the stoma care regimen; puncture dilation resulting in leakage of fluids around the voice prosthesis; inflammatory reaction around the puncture site and formation of granulation tissue; dislodgment of the voice prosthesis and subsequent closure of the TEP; intractable (uncontrollable) leakage around the voice prosthesis, requiring surgical revision or closure of the puncture; dysphagia (difficulty swallowing); tearing or other damage to the voice prosthesis from improper use; microbial growth deposits causing voice prosthesis leakage or valve incompetence; accidental ingestion of the voice prosthesis into the esophagus.

## INSTRUCTIONS FOR USE

THIS DEVICE IS NOT INTENDED TO BE INSERTED OR REMOVED BY THE PATIENT.

*All instructions related to device preparation, insertion, and removal are for the clinician only. If you are a patient, refer to the section on "Device Cleaning and Care".* Prior to device preparation for placement, check the valve mechanism to be sure it is intact and working properly. The flap valve should close flat against the seating surface inside the voice prosthesis.

### Primary Voice Prosthesis Placement

The sterile option voice prosthesis may be placed by the physician at the time of total laryngectomy (primary puncture, primary placement). Nonsterile product is not intended for primary placement procedures.

### Secondary (or Replacement) Device Placement

*The following procedure is provided by Eric D. Blom, Ph.D. as recommended instructions for secondary or replacement procedures.*

**Note:** All secondary or replacement procedures must only be done with a bright light focused directly on the stoma and TEP.

Use of gloves and protective eyewear is recommended. Measure the TEP tract before inserting a voice prosthesis to confirm correct voice prosthesis sizing.

### Dilate and Measure (Diagrams 1 and 2)

Please refer to the Blom-Singer Dilation/Sizing System Instructions for Use for complete product details.

### Voice Prosthesis Preparation

Prepare the voice prosthesis for insertion by utilizing the hand folding method of gel cap insertion. Use only the appropriately sized gel caps provided with this product. Handle the device with clean, gloved hands.

Remove the voice prosthesis and insertion components (gel cap, flange introducer, and inserter tool) from the package. Be certain that your hands and the voice prosthesis are completely dry before loading the gel cap. Failure to do so may cause premature dissolution of the gel cap.

1. Use only the short end of the gel cap. Discard the longer end of the gel cap (diagram 3).
2. Fold the esophageal flange tightly in half toward the center of the device (diagram 4).
3. Fold the two outer edges of the esophageal flange tightly against the folded flange (diagram 5) and hold the prosthesis in this folded position.

4. Insert the folded portion of the esophageal flange into the gel cap (diagram 6). Gently tuck and push the folded esophageal flange as far as possible into the gel cap.
5. Use the curved side of the flange introducer and following the curve of the device body, gently tuck the remaining portion of the folded esophageal flange completely into the gel cap (diagram 7). Do not use any sharp or serrated instruments to avoid damaging the prosthesis or gel cap.
6. The esophageal flange should be folded in a forward position and completely enclosed inside the gel cap (diagram 8).

#### **Device Insertion**

1. Place the device on the inserter and lock the safety strap onto the safety peg (diagram 9).
2. Caution the patient not to swallow as you remove the dilator or sizer to avoid aspiration of saliva into the trachea prior to inserting the voice prosthesis.
3. Apply a light coating of water-soluble lubricant (not petroleum based) to the tip of the gel-capped end of the device and place the tip of the voice prosthesis in the TEP with the safety strap oriented upwards. Insert the device completely into the TEP until the tracheal flange of the device is seated firmly against the posterior tracheal mucosa.
4. Hold the device in this position of complete insertion for at least three minutes to allow the gel cap to dissolve and release the esophageal flange within the esophagus. Patients should be instructed to swallow their saliva to facilitate dissolving the gel cap.
5. Confirm the esophageal flange has deployed: Rotate the device on the inserter while it is in the TEP. The device will rotate easily on the inserter in the TEP if the esophageal flange has deployed and is seated against the anterior esophageal wall. Considerable resistance to device rotation will be felt if the gel cap has not dissolved or incorrectly deployed in the TEP instead of inside the lumen of the esophagus.
6. Detach the safety strap from the safety peg of the inserter. Place a finger against the safety strap and carefully withdraw the inserter from the device with a twisting motion. Secure the safety strap of the device to the neck with medical adhesive tape across the safety strap.  
**Warning:** Caution should be exercised when removing the inserter to avoid accidental displacement of the voice prosthesis, which could result in aspiration of the device.
7. Assess the capacity to produce tracheoesophageal voice. Digital occlusion of the tracheostoma will divert airflow from the trachea. A device that is correctly positioned with the gel cap fully dissolved should allow that airflow to open the flap valve and pass into the esophagus for tracheoesophageal voice production (diagram 10).
8. Ensure there is no leakage through or around the device by having the patient drink water (diagram 11).

#### **Confirmation of Valve Function and Esophageal Flange Deployment (Non-Dual Valve)**

If the proper deployment and positioning of the esophageal flange cannot be confirmed by methods described above, the clinician may need to determine gel cap dissolution and esophageal retention flange deployment by manually opening, with gentle contact, the one-way flap valve in the esophageal end of the voice prosthesis. Direct a bright light into the prosthesis to visualize the flap valve at the proximal end. **Slowly and carefully** advance the blunt end (non-cotton tip) of a cotton-tipped applicator into the voice prosthesis until **gentle contact** opens the flap valve. A flap valve that does not open with gentle contact may indicate an undeployed flange or inverted flap valve is impairing valve function. If indicated, direct visualization of esophageal flange deployment can be obtained with flexible endoscopy. If the assessment methods described above fail to provide satisfactory verification of correct esophageal retention flange deployment, remove the device, dilate and re-measure the puncture tract, and repeat the device insertion and confirmation process.

#### **Confirmation of Valve Function and Esophageal Flange Deployment (Dual Valve)**

To accommodate esophageal flange deployment of the Dual Valve Indwelling Voice Prosthesis, the inserter is designed so that it does not fully insert into the voice prosthesis. Within three

minutes, the inserter should visibly move forward into the voice prosthesis; the clinician should see this forward movement indicating that the gel cap has dissolved and that the esophageal flange has deployed. To confirm valve function, use a small wax curette or equivalent gauge instrument to gently open the tracheal valve while the patient swallows their saliva or water. The clinician should see the properly positioned esophageal valve at the other end and no leakage.

### **Safety Strap Detachment**

The device is designed to include the optional detachment of the safety strap following confirmation that the voice prosthesis esophageal flange is securely positioned within the esophageal lumen against the anterior wall of the esophagus. Removal of the safety strap is an **optional** procedure. Once deployment of the esophageal flange has been confirmed, detach the safety strap from the tracheal flange at the area of reduced safety strap width, where it meets the tracheal flange, by cutting carefully with scissors (diagram 12).

### **The safety strap should not be detached from the device if:**

- Use of a laryngectomy tube or button may potentially dislodge the device.
- A patient has previously had a device that is larger in diameter than the device being inserted
- The diameter of the TEP is dilated significantly larger than the device diameter.
- The health of the stoma tissue or physical condition of the patient is in question.

Do not detach the safety strap from the device in the instances cited above. Tape the safety strap to the peristomal skin and confirm security of the device.

### **Device Removal**

The device is not permanent and requires periodic replacement: it is intended for six (6) months of use, but this duration should be determined per clinical evaluation and appropriate consultation with the clinician.

Removal of the device **should only** be completed by grasping the tracheal flange of the device securely with a locking hemostat. Pull gently and firmly until the voice prosthesis is fully withdrawn. Insert a dilator or dilator-size of the appropriate diameter and tape it in position for five minutes prior to inserting a new device.

### **DEVICE CLEANING AND CARE**

The following information should be made clear to the patient regarding the routine care and cleaning of the indwelling voice prosthesis. The purpose of cleaning the prosthesis is to remove blockage of its lumen that impairs the patient's ability to voice while the voice prosthesis is in the patient's TEP. Cleaning should be done by using the Blom-Singer Cleaning Brush and the Blom-Singer Flushing Device. Please refer to the Voice Prosthesis Cleaning System instructions for use for complete instructions on how to use the cleaning devices.

**Caution:** Use **only** the Blom-singer cleaning devices. **Do not** insert objects other than the Blom-Singer cleaning devices into the voice prosthesis as this may cause damage or dislodgement of the voice prosthesis or its components. Cleaning of the indwelling voice prosthesis should only be done in front of a mirror with a bright light focused directly on the stoma so that the open end of the voice prosthesis is clearly visible. Only use a lint free cloth or tissue to dry the device. Use of non-lint free materials may leave debris, which can be aspirated into the airway. Do not use solvents or petroleum-based products for cleaning or lubricating the device. These materials may damage the silicone or cause the device to not work properly.

### **SPECIAL STORAGE AND/OR HANDLING CONDITIONS**

Only the **gel caps** packaged with the device have special storage and/or handling conditions. They must be **kept dry** and **out of direct sunlight**.

### **DISPOSAL INSTRUCTIONS**

This product is not biodegradable and can be considered contaminated when used. Carefully dispose of this device per local guidelines.

## SPECIAL ORDER PRODUCTS

Special Order products are nonreturnable.

## ORDERING INFORMATION

### USA

Blom-Singer products may be ordered directly from InHealth Technologies.

**ORDER ON-LINE:** www.inhealth.com **TELEPHONE:** Toll-Free (800)477-5969 or (805)684-9337.

**FAX:** Toll-Free (888)371-1530 or (805)684-8594. **EMAIL:** order@inhealth.com

**POST:** InHealth Technologies, 1110 Mark Avenue, Carpinteria, CA 93013-2918, USA,

Attention: Customer Service

### INTERNATIONAL

Please contact our customer service department for a distributor referral.

## PRODUCT COMPLAINTS/EU SERIOUS INCIDENTS

If you are dissatisfied with the device or have any questions, please contact

**productcomplaints@inhealth.com.**

Telephone: +1-800-477-5969

Fax: +1-888-371-1530

Any serious incident that has occurred in relation to the device should be reported to Freudenberg Medical, LLC as noted above and the competent authority of the EU Member State in which the user and/or patient is established.

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**POSTA:** InHealth Technologies, 1110 Mark Avenue, Carpinteria, CA 93013-2918, USA, Attention: Customer Service

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

**SYMBOLS GLOSSARY (USA only)**

<b>SYMBOL</b>	<b>STANDARD No.</b>	<b>SYMBOL TITLE &amp; REF No.</b>	<b>MEANING OF SYMBOL</b>
	EN ISO 15223-1:2016	Manufacturer; 5.1.1	Indicates the medical device manufacturer, as defined in EU Directives or Regulations
	93/42/EEC EU 2017/745	CE Marking; 93/42/EEC Annex XII, EU 2017/745 Annex V	CE marking of conformity
	EN ISO 15223-1:2016	Authorized representative in the European Community; 5.1.2	Indicates the authorized representative in the European Community
	EN ISO 15223-1:2016	Catalogue number; 5.1.6	Indicates the manufacturer's catalogue number so that the medical device can be identified
	EN ISO 15223-1:2016	Batch code; 5.1.5	Indicates the manufacturer's batch code so that the batch or lot can be identified
	EN ISO 15223-1:2016	Use-by date; 5.1.4	Indicates the date after which the medical device is not to be used
	EN ISO 15223-1:2016	Date of manufacture; 5.1.3	Indicates the date when the medical device was manufactured
	EN ISO 15223-1:2016	Caution; 5.4.4	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself
	EN ISO 15223-1:2016	Do not use if package is damaged; 5.2.8	Indicates a medical device that should not be used if the package has been damaged or opened
	EN ISO 15223-1:2016	Consult instructions for use; 5.4.3	Indicates the need for the user to consult the instructions for use
	EN ISO 15223-1:2016	Sterilized using ethylene oxide; 5.2.3	Indicates a medical device that has been sterilized using ethylene oxide

**SYMBOLS GLOSSARY (USA only)**

<b>SYMBOL</b>	<b>STANDARD No.</b>	<b>SYMBOL TITLE &amp; REF No.</b>	<b>MEANING OF SYMBOL</b>
	EN ISO 15223-1:2016	Non-sterile; 5.2.7	Indicates a medical device that has not been subjected to a sterilization process
	EN ISO 15223-1:2016	Do not re-use; 5.4.2	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure
	EN ISO 15223-1:2016	Keep away from sunlight; 5.3.2	Indicates a medical device that needs protection from light sources
	EN ISO 15223-1:2016	Keep dry; 5.3.4	Indicates a medical device that needs to be protected from moisture
	EN ISO 15223-1:2016	Do not re-sterilize; 5.2.6	Indicates a medical device that is not to be re-sterilized
	ISO 7000:2019	Single sterile barrier system; 3707	To indicate that there is a single sterile barrier system.

# INHEALTH

TECHNOLOGIES®

---

## UNITED STATES

1110 Mark Avenue  
Carpinteria, CA 93013-2918 USA  
Toll-Free: 800.477.5969  
Toll-Free Fax: 888.371.1530

## INTERNATIONAL

Tel: +1.805.684.9337  
Fax: +1.805.684.8594  
E-mail: [order@inhealth.com](mailto:order@inhealth.com)

## ORDER ONLINE

[www.inhealth.com](http://www.inhealth.com)

**Freudenberg Medical, LLC**  
1110 Mark Avenue  
Carpinteria, CA 93013-2918 USA  
[www.inhealth.com](http://www.inhealth.com)

CE 0344

**EMERGO EUROPE**  
Westervoortsedijk 60  
6827 AT Arnhem  
The Netherlands