
INHEALTH

TECHNOLOGIES®

Blom-Singer® Day&Night® HME

INSTRUCTIONS FOR USE

Blom-Singer®

voice restoration systems

Blom-Singer® and InHealth Technologies® are registered trademarks of **Freudenberg Medical, LLC**.

Patent(s): www.inhealthpatents.com

IFU: <https://inhealth.com/ifu>

Copyright © 2020 Freudenberg Medical, LLC. All rights reserved.

INSTRUCTIONS FOR USE

使用说明 / 使用說明書 / Udhëzime për përdorim / إرشادات الاستخدام / Инструкции за употреба / Upute za korištenje / Návod k použití / Brugsanvisning / Gebruiksaanwijzing / Kasutusjuhend / Käyttöohjeet / Instructions d'utilisation / Gebrauchsanweisung / Οδηγίες χρήσης / הוראות שימוש / Használati utasítás / Notkunarleiðbeiningar / Istruzioni per l'uso / 사용 지침 / Lietošanas norādījumi / Naudojimo instrukcijos / Bruksanvisning / Instrukcja użytkowania / Instruções de uso / Instrucțiuni de utilizare / Инструкции по применению / Uputstvo za upotrebu / Návod na použitie / Instrucciones de uso / Bruksanvisningar / Kullanım Talimatlar

TABLE OF CONTENTS

目录 / 目次 / Tabela e përmbajtjes / جدول المحتويات / Съдържание / Sadržaj / Obsah /
 Indholdsfortegnelse / Inhoud / Sisukord / Sisällysluettelo / Table des matières /
 Inhaltsverzeichnis / Πίνακας περιεχομένων / תוכן עניינים / Tartalomjegyzék / Efnisyfirlit /
 Sommario / 목차 / Saturs / Turinys / Innhaldsfortegnelse / Spis treści / Índice / Cuprins /
 Содержание / Sadržaj / Obsah / Índice / Innehållsförteckning / İçindekiler


EN	Day&Night® HME	06	HU	Day&Night® HME	51
ZH	Day&Night® HME	08	IS	Day&Night® HME	53
JA	Day&Night® HME	11	IT	Day&Night® HME	56
SQ	Day&Night® HME	13	KO	Day&Night® HME	59
16	Day&Night® HME	AR	LV	Day&Night® HME	62
BG	Day&Night® HME	19	LT	Day&Night® HME	64
HR	Day&Night® HME	22	NO	Day&Night® HME	67
CS	Day&Night® HME	25	PL	Day&Night® HME	70
DA	Day&Night® HME	28	PT	Day&Night® HME	73
NL	Day&Night® HME	30	RO	Day&Night® HME	76
ET	Day&Night® HME	34	RU	Day&Night® HME	79
FI	Day&Night® HME	36	SR	Day&Night® HME	82
FR	Day&Night® HME	39	SK	Day&Night® HME	84
DE	Day&Night® HME	42	ES	Day&Night® HME	87
EL	Day&Night® HME	45	SV	Day&Night® HME	90
48	Day&Night® HME	HE	TR	Day&Night® HME	93

SYMBOLS GLOSSARY

SYMBOL	TITLE
	Manufacturer
	CE Marking
	Authorized representative in the European Community
	Catalog number
	Batch code
	Use-by date
	Caution
	Do not use if package is damaged
	Consult instructions for use
	Do not re-use
	Packaging unit
	Medical device

DIAGRAMS

插图 / 図 / Diagramet / الرسم التوضيحي / Схеми / Dijagrami / Schémata / Diagrammer / Afbeeldingen / Diagrammid / Kaaviot / Diagrammes / Abbildungen / Διαγράμματα / תרשימים / Ábrák / Skýringarmyndir / Diagrammi / 그림 / Diagrammas / Diagramos / Diagrammer / Schematy / Diagramas / Diagramme / Схемы / Crteži / Grafické znázornenia / Diagramas / Diagram / Diyagramlar

BLOM-SINGER® Day&Night® HME**PRODUCT DESCRIPTION**

The Blom-Singer® Day&Night® Heat and Moisture Exchange (HME) cartridge is a single, 24-hour use device that offers the benefits of a heat and moisture exchange (HME) system.

Attachment of the HME to the tracheostoma can be achieved through the use of Blom-Singer® tracheostoma adhesive housing or a Blom-Singer® StomaSoft® Laryngectomy tube, sold separately. The Day&Night® HME cartridge contains a foam filter for heat and moisture exchange, side airflow vents to maximize airflow, and a locking mechanism to secure the device in the open position to prevent occlusion as desired.

INDICATIONS (Reasons to prescribe the device or procedure)

The Blom-Singer® Day&Night® HME is indicated for individuals who have a permanent tracheostoma (opening through neck into trachea or “windpipe”) following total laryngectomy surgery. It is designed to provide laryngectomees with humidified, warmed, and filtered air. It may be used for heat and moisture exchange by esophageal speakers, electrolarynx users, and tracheoesophageal voice prosthesis users.

CONTRAINDICATIONS (Reasons that make it inadvisable to prescribe the particular device or procedure)

The Blom-Singer® Day&Night® HME should not be placed over a tracheostoma airway if its presence obstructs respiration.

WARNINGS AND PRECAUTIONS

Wash hands with soap and water before and after handling the device, as the foam contains substances that may be irritating to skin and eyes. Only use a lint-free cloth if wiping the device to avoid getting lint into the airway. This device is not reusable and may only be used by one patient. There is a risk of infection if this device is shared between patients or reused. Inspect the device before use for physical damage (such as cracks or tears); do not use if damaged. If physical damage is found, contact InHealth Technologies Product Complaints.

Care should be exercised when inserting or removing HME devices to avoid aspiration (inhalation) of the voice prosthesis or HME. In the unlikely event that this should occur, you should attempt to cough the device out of the windpipe: immediately seek medical attention if unsuccessful.

Breathing

Breathing through the device should be comfortable and easy. If breathing is not easy, remove the device and consult your medical professional for assistance. When you begin using the device, a slight increase in airflow resistance may be felt. This resistance is a normal function of the device that may require some adjustment time. The device may need to be removed or engaged in the locked-open position during instances of increased physical activity/heavy breathing.

Coughing

It is recommended that the device be removed prior to coughing.

Sleeping

The optional locking mechanism should be engaged on the device when sleeping to prevent accidental occlusion.

POTENTIAL COMPLICATIONS

The following complications may occur with devices of this type, but are infrequent. They include:

- Stoma contamination or infection which may require use of appropriate antibiotics.
- Irritation, including skin irritation around the stoma.
- Labored breathing, which may occur if the device is heavily soiled, clogged with mucus, not changed regularly, or accidentally occluded.

INSTRUCTIONS FOR USE

Please refer to the diagrams located at the front of this instruction manual.

Insertion and removal of any devices worn over the tracheostoma should only be done with the user positioned in front of a mirror with a bright light focused directly on the stoma.

A Blom-Singer® tracheostoma housing or StomaSoft® Laryngectomy Tube must first be placed on/around the tracheostoma. See instructions for use manual for the specific instructions on insertion or placement of the chosen tracheostoma attachment.

Insertion

To insert the device, hold its outer rim and insert it partially into the lower rim of the previously applied tracheostoma housing or tube (**diagram 1**). Gently push the device into the applied attachment until it snaps into place. A properly inserted HME cartridge will sit flush and fit snugly inside the housing or tube (**diagram 2**).

Device Adjustment/Locking Mechanism

The Day&Night® HME has an optional twist-lock feature to keep it in the open position if desired.

To lock: hold the tracheostoma housing or tube in place with one hand while gently turning the upper lip of the device **clockwise** with the other hand (**diagram 3a**).

NOTE: When in the open/locked position, the user will be unable to occlude the device for tracheoesophageal speech with a voice prosthesis.

To unlock: hold the tracheostoma housing or tube in place with one hand while gently turning the upper lip of the device **counterclockwise** with the other hand (**diagram 3b**).

Tracheoesophageal Speech with a Voice Prosthesis

Tracheoesophageal speech with a voice prosthesis and Day&Night® HME is produced by pressing against the face of the device (**diagram 4a**). This will compress the foam and seal the device for voicing. Upon manual release, the Day&Night® HME face returns to its open, breathing position (**diagram 4b**).

Removal

The device can be removed in either the locked or unlocked position. To remove the device, hold the tracheostoma housing or tube in place with one hand while gently pulling the device with the other hand until it releases (**diagram 5**).

Device Care

The Day&Night® HME cartridge is designed for 24-hour use and is not designed to be cleaned, disassembled/reassembled, or reused. Mucus may be wiped away with a lint-free cloth as able. If the HME cartridge should become heavily soiled, remove it and replace it with a new HME cartridge.

SPECIAL STORAGE AND/OR HANDLING CONDITIONS

Only the foam filter portions of these devices have special storage and/or handling conditions:

- Wash hands after handling the cartridge, as it contains substances that may be irritating to skin and eyes.
- The packages containing the HME cartridges should be resealed after opening to preserve the effectiveness of the foam.

DISPOSAL INSTRUCTIONS

This device is not biodegradable and can be contaminated when used. Carefully dispose of the device per local guidelines.

ORDERING INFORMATION

USA InHealth Technologies® products may be ordered directly from Freudenberg Medical, LLC, d/b/a "InHealth Technologies."

Telephone: Toll-Free 800.477.5969 or +1.805.684.9337

Fax: Toll-Free 888.371.1530 or +1.805.684.8594

Email: order@inhealth.com Order Online: www.inhealth.com

Mail: InHealth Technologies 1110 Mark Avenue Carpinteria, CA 93013-2918 USA

INTERNATIONAL

InHealth Technologies products are available through authorized distributors worldwide. For international distributor information, please contact our customer service department.

PRODUCT COMPLAINTS/EU SERIOUS INCIDENTS

If you are dissatisfied with the device or have any questions, please contact productcomplaints@inhealth.com

Telephone: +1-800-477-5969

Fax: +1-888-371-1530

Any serious incident that has occurred in relation to the device should be reported to Freudenberg Medical, LLC as noted above and the competent authority of the EU Member State in which the user and/or patient is established.

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ÜRÜN ŞİKAYETLERİ/AB CİDDİ OLAYLAR