# INHEALTH

TECHNOLOGIES®

## Blom-Singer® Nasal Septal Perforation Prosthesis

MEDICAL PROFESSIONAL INSTRUCTIONS FOR USE

**Blom-Singer®**

voice restoration systems

37805-01D Effective January 2020 / 2020 年 1 月生效 / 2020年1月発効 / Në fuqi nga janari 2020 / 2020 اعتبارًا من شهر يناير / В сила от януари 2020 г. / Stupa na snagu u siječnju 2020. / V platnosti od ledna 2020 / Gældende fra januar 2020 / Met ingang van januari 2020 / Kehtib alates jaanuarist 2020 / Voimassa tammikuusta 2020 alkaen / Date d'entrée en vigueur : janvier 2020 / Gültig ab Januar 2020 / Σε ισχύ από τον Ιανουάριο του 2020 / 2020 בתוקף החל מינואר / Hatályos 2020. januártól / Gildir frá janúar 2020 / Valido da gennajo 2020 / 2020년 1월 발효 / Spēkā no 2020. gada janvāra / Įsigalioja 2020 m. sausio mėn. / Gjelder fra januar 2020 / Wchodzi w życie w styczniu 2020 r. / Vingente em janeiro de 2020 / Cu intrare în vigoare din ianuarie 2020 / Версия от января 2020 г. / Važi od januara 2020. / Platné od januára 2020 / Efectivo a partir de enero de 2020 / Gäller från januari 2020 / Geçerlilik Tarihi: Ocak 2020

Blom-Singer® and InHealth Technologies® are registered trademarks of Freudenberg Medical, LLC.

Patent(s): [www.inhealthpatents.com](http://www.inhealthpatents.com)

Copyright © 2020 Freudenberg Medical, LLC. All rights reserved.

## MEDICAL PROFESSIONAL INSTRUCTIONS FOR USE


医学专业使用说明 / 医療専門家用取扱説明書 / Udhëzimet e përdorimit për mjekët profesionistë / إرشادات الاستخدام للمتخصصين في المجال الطبي / Инструкции за употреба за медицински специалисти / Upute za korištenje za zdravstveno osoblje / Návod k použití pro zdravotníky / Brugsanvisning til uddannede læger / Gebruiksaanwijzing voor de medische professional / Tervishoitutöötaja kasutusjuhend / Lääketieteen ammattilaisen käyttöohjeet / Instructions d'utilisation destinées aux professionnels de santé / Gebrauchsanweisung für medizinische Fachkräfte / Οδηγίες χρήσης για τον επαγγελματία υγείας / הוראות שימוש לצוות הרפואי / Egészségügyi szakembereknek szóló használati utasítás / Notkunarleiðbeiningar fyrir heilbrigðisstarfsmenn / Istruzioni per l'uso destinate al personale sanitario / 의료 전문가 사용 지침 / Medicīnas speciālistam paredzētas lietošanas instrukcijas / Medikams profesionalams skirta naudojimo instrukcija / Bruksanvisning for medisinske fagpersoner / Instrukcja obsługi dla personelu medycznego / Instruções de uso para profissionais médicos / Instrucțiuni de utilizare pentru cadre medicale specializate / Инструкции по применению для профессиональных медработников / Uputstvo za upotrebu za zdravstvenog radnika / Návod na použitie pre zdravotníckych pracovníkov / Instrucciones de uso para profesionales médicos / Sjukvårdspersonal, bruksanvisning / Sağlık Uzmanı Kullanım Talimatları

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
### **BLOM-SINGER® NASAL SEPTAL PERFORATION PROSTHESIS**

#### **DEVICE DESCRIPTION**

The Blom-Singer® Nasal Septal Perforation Prosthesis is a bi-flanged elastomeric device. Each prosthesis is two (2) pieces which are attached to one another using an engineered magnetic circuit. For a given catalog item, both flanges are the same diameter and shape, as are the diameter and shape of the connecting shaft which most closely fits the septal perforation. An accurately fitting shaft shape and diameter minimize prosthesis movement and collection of debris. A qualified, trained medical professional will assess the diameter and shape of the septal perforation and select the proper device to ensure an exact fit is achieved to close the septal perforation.

#### **INDICATIONS FOR USE**

The Nasal Septal Perforation Prosthesis is indicated for the non-surgical closure of nasal septal perforations. Environment of use: Hospital, sub-acute, physician office settings.

#### **CONTRAINDICATIONS**

- Septal deformities are relative contraindications for prosthesis placement, as the flange on the convex side of the deformity may increase nasal obstruction.
- Women who are pregnant should consult their doctor before using a medical device containing a magnet.
- People who have magnetic or electronic medical equipment, such as pacemakers, defibrillators, and diabetic pumps, should not use the Nasal Septal Perforation Prosthesis.
- Should not be in situ during MRI examination (Magnetic resonance imaging).
- Should not be in situ during radiation therapy.

#### **WARNINGS AND PRECAUTIONS**

- The Nasal Septal Perforation Prosthesis is for single patient use only.
- Duration of use is up to 29 days.
- Initial sizing and instructions for use must be provided by a qualified medical professional knowledgeable in the use of this device.
- If the prosthesis is not functioning or sized properly (causing discomfort, swallowed, aspirated, etc.), the patient should have it evaluated by a trained medical professional as soon as possible.
- If suggested maintenance is not followed, it could result in irritation or infection.
- If the prosthesis is removed, it should only be reinserted by a qualified clinician. The patient should not attempt to reinsert it themselves.
- If there are tears, cracks, or structural damage to the prosthesis or packaging, do not use the device.
- Device is provided Non Sterile.
- Do not sterilize in radiation or any apparatus that would exceed 150° F (66° C).
- Do not use petroleum-based products such as Vaseline® to lubricate the prosthesis as these products can damage silicone.

#### **ADDITIONAL INFORMATION**

**Air Travel:** This prosthesis contains magnets. This prosthesis may set off metal detectors at airports. To avoid incidents, it is suggested to advise security in advance about your prosthesis.

## **INSTRUCTIONS FOR USE**

- The physician must evaluate the condition of the patient and this procedure based on his/her own medical training and experience and assess if this device meets the needs of the patient. The following procedure is for informational purposes only.
- The Nasal Septal Perforation Prosthesis should be sized so that there is no nasal tissue captured between the magnets.

## **EVALUATION OF SEPTAL PERFORATION** (common office supplies)

1. Visually estimate the septal thickness.
2. Gently insert a tongue blade or instrument that can be easily marked on in one nostril covering the septal perforation.
3. Use a marking instrument in opposite nostril to trace out the perforation.
4. If available, compare the outline of the perforation to the Nasal Septal Perforation Prosthesis reference guide.
5. Assess the shape to be oval or round.
6. Measure the perforation outline with a ruler to determine the required septal prosthesis.

## **PLACEMENT**

1. Securely grasp the edge of one (1) of the two (2) flange components gently with forceps or a locking hemostat.
2. Place one flange inside of nose so the magnet shaft inserts into the perforation.
3. Securely grasp the edge of the second flange component with forceps or locking hemostat and place in opposite nostril until the magnetic shafts of both flange components align and attract to each other in the septal perforation.
  - a. If using an oval Nasal Septal Perforation Prosthesis, minor alignment may be necessary once the Nasal Septal Perforation Prosthesis is magnetically coupled.
  - b. Ensure that there is no tissue between the magnets.
4. Detach and remove forceps or locking hemostat
5. Assess the security of the prosthesis placement by having the patient forcefully blow his/her nose. The Nasal Septal Perforation Prosthesis should not move or uncouple with this action.

## **REMOVAL**

1. Grasp both flange component edges simultaneously with forceps
2. Gently pull to remove the Nasal Septal Perforation Prosthesis. (DO NOT USE EXCESSIVE FORCE.)

## **SUGGESTED MAINTENANCE** (to be performed by a qualified clinician)

1. Visible debris removal with prosthesis in place:
  - a. Use a cotton tip applicator that is dampened with 0.9% saline solution to gently move around and under the flanges.
2. Visible debris removal, if prosthesis is removed:
  - a. If the prosthesis is removed, the physician may try to remove visible debris with mild soap and water. Rinse with water and air dry prior to reinsertion.

## SPECIAL STORAGE AND/OR HANDLING CONDITIONS

The device shall not be exposed to temperatures exceeding 66° C.

## DISPOSAL INSTRUCTIONS

This device is not biodegradable and can be contaminated when used. Carefully dispose of the device per local guidelines.

## ORDERING INFORMATION

### USA

Blom-Singer products may be ordered directly from InHealth Technologies. TELEPHONE: Toll-Free (800)477-5969 or (805)684-9337, Monday — Friday, 9:30 am — 7:00 pm, Eastern Standard Time. FAX: Toll-Free (888)371-1530 or (805)684-8594. EMAIL: [order@inhealth.com](mailto:order@inhealth.com) ORDER ON-LINE: [www.inhealth.com](http://www.inhealth.com) POST: InHealth Technologies, 1110 Mark Avenue, Carpinteria, CA 93013-2918, USA, Attention: Customer Service.

### INTERNATIONAL

Blom-Singer products are available worldwide through authorized InHealth Technologies distributors. For assistance or additional product information please contact the distributor.

## PRODUCT COMPLAINTS/EU SERIOUS INCIDENTS

If you are dissatisfied with the device or have any questions, please contact

**[productcomplaints@inhealth.com](mailto:productcomplaints@inhealth.com).**

Telephone: +1-800-477-5969

Fax: +1-888-371-1530

Any serious incident that has occurred in relation to the device should be reported to Freudenberg Medical, LLC as noted above and the competent authority of the EU Member State in which the user and/or patient is established.

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## SYMBOLS GLOSSARY (USA only)

SYMBOL	STANDARD No.	SYMBOL TITLE & REF No.	MEANING OF SYMBOL
	EN ISO 15223-1:2016	Manufacturer; 5.1.1	Indicates the medical device manufacturer, as defined in EU Directives or Regulations
	93/42/EEC EU 2017/745	CE Marking; 93/42/EEC Annex XII, EU 2017/745 Annex V	CE marking of conformity
	EN ISO 15223-1:2016	Authorized representative in the European Community; 5.1.2	Indicates the authorized representative in the European Community
	EN ISO 15223-1:2016	Catalogue number; 5.1.6	Indicates the manufacturer's catalogue number so that the medical device can be identified
	EN ISO 15223-1:2016	Batch code; 5.1.5	Indicates the manufacturer's batch code so that the batch or lot can be identified
	EN ISO 15223-1:2016	Use-by date; 5.1.4	Indicates the date after which the medical device is not to be used
	EN ISO 15223-1:2016	Caution; 5.4.4	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself
	EN ISO 15223-1:2016	Do not use if package is damaged; 5.2.8	Indicates a medical device that should not be used if the packaging has been damaged or opened
	EN ISO 15223-1:2016	Consult instructions for use; 5.4.3	Indicates the need for the user to consult the instructions for use
	ISO 7000:2019	Packaging unit; 2794	To indicate the number of pieces in the package

**SYMBOLS GLOSSARY (USA only)**

<b>SYMBOL</b>	<b>STANDARD No.</b>	<b>SYMBOL TITLE &amp; REF No.</b>	<b>MEANING OF SYMBOL</b>
	EN ISO 15223-1:2016	Upper limit of temperature; 5.3.6	Indicates the upper limit of temperature to which the medical device can be safely exposed.
	ISO 7000:2019	Single patient multiple use; 3706	To indicate that the medical device may be used multiple times (multiple procedures) on a single patient

# INHEALTH

TECHNOLOGIES®

---

## UNITED STATES

1110 Mark Avenue  
Carpinteria, CA 93013-2918 USA  
Toll-Free: 800.477.5969  
Toll-Free Fax: 888.371.1530

## INTERNATIONAL

Tel: +1.805.684.9337  
Fax: +1.805.684.8594  
E-mail: [order@inhealth.com](mailto:order@inhealth.com)

## ORDER ONLINE

[www.inhealth.com](http://www.inhealth.com)

### Freudenberg Medical, LLC

1110 Mark Avenue  
Carpinteria, CA 93013-2918 USA  
[www.inhealth.com](http://www.inhealth.com)

### EMERGO EUROPE

Prinsessegracht 20  
2514 AP The Hague  
The Netherlands