


# Introducing the New Blom-Singer® Dual Valve Voice Prosthesis




Our new dual valve Indwelling voice prosthesis uses leading-edge design to deliver on a very simple, but powerful benefit: two valves are better than one!

Should the primary valve fail due to Candida (yeast) colonization, the secondary valve is intended to remain viable to help prevent leakage of fluids.

As an added advantage, each valve is imbedded with silver oxide, an antimicrobial agent that provides significant protection to the device for laryngectomees who are susceptible to Candida growth. Blom-Singer voice prostheses that contain silver oxide have been proven to eliminate or reduce the need for antifungal agents, such as Nystatin™ or Mycelex®.

Included with the Dual Valve Voice Prosthesis is a unique pop-through insertion system, designed to assist the clinician in the deployment of the esophageal flange during placement.



The Dual Valve voice prostheses are offered in non-sterile 20 Fr. indwelling packages, in five lengths:

LENGTH	20 FR DIAMETER
6mm	DV2006-NS
8mm	DV2008-NS
10mm	DV2010-NS
12mm	DV2012-NS
14mm	DV2014-NS

Call us today to place an order at **(805) 684-9337** or visit our website at [www.inhealth.com](http://www.inhealth.com)

**INHEALTH**  
TECHNOLOGIES

**Blom-Singer®**  
voice restoration systems  
*The Voice of Experience  
Since 1978*