

A photograph of two women, one older with white curly hair and one younger with blonde hair, both smiling and petting a small black and white dog in a grassy outdoor setting. The older woman is wearing a grey zip-up top and large hoop earrings. The younger woman is wearing a black top and has a tattoo on her arm. The background shows greenery and a building.

**INHEALTH**  
TECHNOLOGIES®

**Blom-Singer®**

**2025/26 PRODUCT CATALOG**

## CONTENTS

<b>Pre-Operative Education</b>	<b>2</b>
“What to Expect” Brochure	2
Educational Videos	2
<b>Surgical Solutions</b>	<b>3</b>
FirstFit™ Surgical Kit (Sterile)	3
Indwelling Voice Protheses (Sterile)	3
<b>Post-Operative Care</b>	<b>4</b>
Laryngectomy Tubes	4
Marpac Tracheostomy Tube Holder/Collar™	4
StomaSoft® Cleaning Swabs	4
Electrolarynx	6
Post-Operative Care Kits (For Total Laryngectomy)	7
<b>Next Steps</b>	<b>8</b>
Voice Protheses (Clinician Placed)	8
Special Order Voice Protheses (Clinician Placed)	10
Voice Protheses (Patient Changeable)	12
Voice Prothesis Accessories	13
Heat & Moisture Exchange (HME)	15
Adhesive Housings	16
Tracheostoma Attachments	17
<b>And Beyond</b>	<b>18</b>
Hands-Free Speech Valves	18
Skin Care	18
Stoma Protection and Accessories	19
<b>Professional Products</b>	<b>21</b>
Katz Extractor® Oto-Rhino Foreign Body Remover	25
Nasal Septal Perforation Sizing Tool	26
Nasal Septal Perforation Protheses	26
<b>Ordering Information</b>	<b>27</b>
<b>New Patient Enrollment Instructions</b>	<b>27</b>
<b>Free Educational Resources</b>	<b>27</b>

# InHealth Technologies® offers a TOTAL Solution for your total laryngectomy from pre- to post-op and beyond

TOTAL Solution serves patients at all phases of voice restoration, from products that deliver solutions for pre- and post-operative care and beyond. This continuum of care is enhanced by the InHealth Technologies® customer-focused team, including leading clinical specialists, experienced patient consultants, and dedicated customer service members.

Cover Image:

Left: Joanne  
Year of TL: 2015  
Wearing AccuFit® Adhesive Housing  
and Day&Night® HME

Right: Meagon  
Year of TL: 2016  
Wearing HydroFit® Adhesive Housing  
and Day&Night® HME

# TOTAL Solution

PHASE

1

## Pre-Operative Education

- "What to Expect" Brochure
- Educational Videos

PHASE

2

## Surgical Solutions

- FirstFit™ Surgical Kit (Sterile)
- Sterile Option Voice Prosthesis

PHASE

3

## Post-Operative Care

- Post-Operative Care Kits
- Laryngectomy Tubes
- Electrolarynx

PHASE

4

## Next Steps

- Next Steps Kit
- Voice Prosthesis
- HMEs
- Adhesive Housings
- Educational Programs

PHASE

5

## And Beyond

- Hands-Free Speech Valve
- Skin Care





## ● PRE-OPERATIVE EDUCATION

### "What to Expect" Brochure

Available online and free to download at:

<https://inhealth.com/educational-materials>



### Educational Videos

Available at InHealth Technologies®

<https://inhealth.com/videos>



## ● SURGICAL SOLUTIONS

### FirstFit™ Surgical Kit (Sterile) (Available Only in the US)

CAUTION: FEDERAL (USA) LAW RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN

Designed for primary or secondary tracheoesophageal puncture (TEP) with retrograde voice prosthesis insertion. This set includes: one (1) 13-gauge curved cored puncture needle; one (1) 36-inch guidewire with rounded tips; one (1) pharynx protector tool with handle and notched cylinder; one (1) voice prosthesis inserter pre-loaded with an indwelling voice prosthesis of the diameter and length specified on the package label.

Ref	Description
TP 1608	FirstFit™ Surgical Kit + Classic™ 16 Fr 8 mm Voice Prosthesis
TP 1610	FirstFit™ Surgical Kit + Classic™ 16 Fr 10 mm Voice Prosthesis
TP 1612	FirstFit™ Surgical Kit + Classic™ 16 Fr 12 mm Voice Prosthesis
TP 1614	FirstFit™ Surgical Kit + Classic™ 16 Fr 14 mm Voice Prosthesis
TP 2008	FirstFit™ Surgical Kit + Classic™ 20 Fr 8 mm Voice Prosthesis
TP 2010	FirstFit™ Surgical Kit + Classic™ 20 Fr 10 mm Voice Prosthesis
TP 2012	FirstFit™ Surgical Kit + Classic™ 20 Fr 12 mm Voice Prosthesis
TP 2014	FirstFit™ Surgical Kit + Classic™ 20 Fr 14 mm Voice Prosthesis



### Blom-Singer® Classic™ Indwelling Voice Protheses (Sterile)

CAUTION: FEDERAL (USA) LAW RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN

Designed for primary TE puncture and placement procedures. Includes insertion and cleaning accessories. Clinician placed.

Length	16 Fr	20 Fr	Length	16 Fr	20 Fr
8mm	IN 1608-SO	IN 2008-SO	12mm	IN 1612-SO	IN 2012-SO
10mm	IN 1610-SO	IN 2010-SO	14mm	IN 1614-SO	IN 2014-SO



## ● POST-OPERATIVE CARE

### Blom-Singer® StomaSoft® Laryngectomy Tubes

Blom-Singer® StomaSoft® Laryngectomy Tubes are flexible, soft, and compatible with InHealth Technologies® HMEs and other accessories from InHealth Technologies®. They feature reinforced strap holders and are available with or without fenestrations.

Size/Length	I.D./O.D.	Non-Fenestrated	Fenestrated	Rx
8/36	9.5mm/12.0mm	BE 6698	BE 6698F	■
8/55	9.5mm/12.0mm	BE 6699	BE 6699F	■
9/36	10.5mm/13.5mm	BE 6700	BE 6700F	■
9/55	10.5mm/13.5mm	BE 6701	BE 6701F	■
10/36	12.0mm/15.0mm	BE 6702	BE 6702F	■
10/55	12.0mm/15.0mm	BE 6703	BE 6703F	■
12/36	13.5mm/17.0mm	BE 6704	BE 6704F	■
12/55	13.5mm/17.0mm	BE 6705	BE 6705F	■

### Blom-Singer® Laryngectomy Tubes (Sterile)

Blom-Singer® Laryngectomy Tubes (Sterile) are durable, lightweight silicone tubes that are used to keep the tracheostoma open after laryngectomy. The tubes do not hold an HME.

Size/Length	I.D./O.D.	Non-Fenestrated (Sterile)	Rx
8/36	9.5mm/12.0mm	BE 6398	■
8/55	9.5mm/12.0mm	BE 6399	■
9/36	10.5mm/13.5mm	BE 6400	■
9/55	10.5mm/13.5mm	BE 6401	■
10/36	12.0mm/15.0mm	BE 6402	■
10/55	12.0mm/15.0mm	BE 6403	■
12/36	13.5mm/17.0mm	BE 6404	■
12/55	13.5mm/17.0mm	BE 6405	■

### Marpac Tracheostomy Tube Holder/Collar™

Marpac Tracheostomy Tube Holder/Collar™ (for laryngectomy patients) is a single use, disposable Comfortable Collar™ used to secure a Blom-Singer® Laryngectomy Tube or a Blom-Singer® Laryngectomy Button. Latex free. Adult—one size fits all.

Ref	Description
204	Tracheostomy Tube Holder/Collar™

### StomaSoft® Cleaning Swabs

The StomaSoft® Cleaning Swabs are intended for cleaning laryngectomy tubes or buttons after removal from the patient. Available in packs of 10.

Ref	Description
TSWB-10	Foam Swab, pack of 10



Mike  
Year of TL: 2015  
Wearing StomaSoft®  
Laryngectomy Tube  
and Day&Night® HME



## Blom-Singer® Electrolarynx EL 2000

An electrolarynx is a speaking aid that produces sound vibrations that can be conducted through the throat and mouth to speak after a laryngectomy.

Ref	Description
EL 2000	Blom-Singer® Electrolarynx Digital Speech Aid
EL 2002	Oral Adapter

Oral Adaptor EL 2002



Electrolarynx EL 2000





## Blom-Singer® Post-Operative Care Kits (For Total Laryngectomy)

CAUTION: FEDERAL (USA) LAW RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN

The Blom-Singer® Post-Operative Care Kits were designed with valued input from clinicians and patients. These kits contain supplies to assist patients in their transition from the hospital to home after a total laryngectomy. They are available in several configurations to meet select needs. All configurations include your choice of StomaSoft® Laryngectomy Tube and HMEs as listed below.

The Post-Operative Care Kit is housed in a hard, clear storage case, which also includes information on how to obtain our Next Steps Kit. The Next Steps Kit includes additional sample products and educational items.

All kit configurations contain the following patient care supplies:

- 2 StomaSoft® Laryngectomy Tubes
- HMEs, pack of 30
- 7 Marpac Tracheostomy Tube Holder/Collars™
- BE 6048 Blom-Singer® Shower Guard
- Drawstring Bag
- HydroFit® Adhesive Housings, pack of 6
- Mirror with Stand
- Halcyon® LED Flashlight
- Next Steps Kit - Request Flyer



REF	Description
<b>PK 2100</b>	
BE 6701	Blom-Singer® StomaSoft® Laryngectomy Tube 9/55 Non-Fenestrated, 2 ea
BE 1091	Day&Night® HME Cartridge ClassicFlow® Heat and Moisture Exchange, 1 ea
<b>PK 2101</b>	
BE 6703	Blom-Singer® StomaSoft® Laryngectomy Tube 10/55 Non-Fenestrated, 2 ea
BE 1091	Day&Night® HME Cartridge ClassicFlow® Heat and Moisture Exchange, 1 ea
<b>PK 2102</b>	
BE 6701	Blom-Singer® StomaSoft® Laryngectomy Tube 9/55 Non-Fenestrated, 2 ea
BE 1092	Day&Night® HME Cartridge EasyFlow® Heat and Moisture Exchange, 1 ea
<b>PK 2103</b>	
BE 6703	Blom-Singer® StomaSoft® Laryngectomy Tube 10/55 Non-Fenestrated, 2 ea
BE 1092	Day&Night® HME Cartridge EasyFlow® Heat and Moisture Exchange, 1 ea
<b>PK 2104</b>	
EL 2000	Blom-Singer® Electrolarynx, 1 ea
EL 2002	Oral Adapter, 1 ea
BE 6703	Blom-Singer® StomaSoft® Laryngectomy Tube 10/55 Non-Fenestrated, 2 ea
BE 1091	Day&Night® HME Cartridge ClassicFlow® Heat and Moisture Exchange, 1 ea
<b>PK 2106</b>	
BE 6705	Blom-Singer® StomaSoft® Laryngectomy Tube 12/55 Non-Fenestrated, 2 ea
BE 1091	Day&Night® HME Cartridge ClassicFlow® Heat and Moisture Exchange, 1 ea
<b>PK 2107</b>	
BE 6705	Blom-Singer® StomaSoft® Laryngectomy Tube 12/55 Non-Fenestrated, 2 ea
BE 1092	Day&Night® HME Cartridge EasyFlow® Heat and Moisture Exchange, 1 ea

## ● NEXT STEPS

### VOICE PROSTHESES (CLINICIAN PLACED)

#### Blom-Singer® Classic™ Indwelling Voice Prostheses (Sterile)

Includes sterile prostheses for primary TE puncture and placement procedures. Includes insertion and cleaning accessories. Clinician placed.

Length	16 Fr	20 Fr	Rx	Length	16 Fr	20 Fr	Rx
8mm	IN 1608-SO	IN 2008-SO	■	12mm	IN 1612-SO	IN 2012-SO	■
10mm	IN 1610-SO	IN 2010-SO	■	14mm	IN 1614-SO	IN 2014-SO	■

#### Blom-Singer® Classic™ Indwelling Voice Prostheses

Includes nonsterile voice prosthesis for secondary or replacement TEP procedures, with insertion and cleaning accessories. Clinician placed.

Length	16 Fr	20 Fr	Rx	Length	16 Fr	20 Fr	Rx
4mm	IN 1604-NS	IN 2004-NS	■	14mm	IN 1614-NS	IN 2014-NS	■
6mm	IN 1606-NS	IN 2006-NS	■	16mm	IN 1616-NS	IN 2016-NS	■
8mm	IN 1608-NS	IN 2008-NS	■	18mm	IN 1618-NS	IN 2018-NS	■
10mm	IN 1610-NS	IN 2010-NS	■	20mm	IN 1620-NS	IN 2020-NS	■
12mm	IN 1612-NS	IN 2012-NS	■				

#### Blom-Singer® Dual Valve™ Indwelling Voice Prostheses

The Blom-Singer® Dual Valve™ Indwelling Voice Prosthesis has two valves embedded with silver oxide as a material preservative.\* When the primary valve fails, the secondary speech valve is intended to remain viable to help prevent leakage of fluids. The Blom-Singer® Dual Valve™ Voice Prosthesis is available with standard or large tracheal and esophageal flanges. Includes nonsterile voice prosthesis for secondary or replacement TEP procedures, with insertion/cleaning accessories. Clinician placed.

Length	20 Fr	20 Fr Large Flange	Rx	Length	20 Fr	20 Fr Large Flange	Rx
6mm	DV 2006-NS	DV 2006-LF	■	12mm	DV 2012-NS	DV 2012-LF	■
8mm	DV 2008-NS	DV 2008-LF	■	14mm	DV 2014-NS	DV 2014-LF	■
10mm	DV 2010-NS	DV 2010-LF	■				

#### Blom-Singer® Advantage® Indwelling Voice Prostheses

The Blom-Singer® Advantage® Indwelling Voice Prosthesis is embedded with silver oxide as a material preservative.\* Includes nonsterile voice prosthesis for secondary or replacement TEP procedures, with insertion and cleaning accessories.

##### HARD VALVE ASSEMBLY

Length	20 Fr	Rx
4mm	AD 2004-NS	■
6mm	AD 2006-NS	■
8mm	AD 2008-NS	■
10mm	AD 2010-NS	■
12mm	AD 2012-NS	■
14mm	AD 2014-NS	■

##### SOFT VALVE ASSEMBLY

Length	16 Fr	20 Fr	Rx
4mm	AS 1604-NS	AS 2004-NS	■
6mm	AS 1606-NS	AS 2006-NS	■
8mm	AS 1608-NS	AS 2008-NS	■
10mm	AS 1610-NS	AS 2010-NS	■
12mm	AS 1612-NS	AS 2012-NS	■
14mm	AS 1614-NS	AS 2014-NS	■



Classic™ Indwelling (Sterile)



Classic™ Indwelling



Dual Valve™



Dual Valve™ Large Flange



Advantage® Hard Valve Assembly



Advantage® Soft Valve Assembly

\*Silver oxide, a known anti-microbial substance, has been incorporated into the silicone flap valve of the Blom-Singer® Advantage® and Dual Valve™ Indwelling Voice Prostheses as a material preservative.



Ksenia  
Year of TL: 2019  
Wearing AccuFit®  
Adhesive Housing and  
Day&Night® HME



## SPECIAL ORDER VOICE PROSTHESES (CLINICIAN PLACED)


—DEVELOPED BY ERIC D. BLOM, PH.D

Unique problems call for unique solutions. InHealth Technologies® offers several voice prostheses to address those needs.

### Special Length

Available in one millimeter incremental lengths for a more precise TE puncture fit.


	5mm	7mm	9mm	Rx
16 Fr	IN 1605-SL	IN 1607-SL	IN 1609-SL	■
20 Fr	IN 2005-SL	IN 2007-SL	IN 2009-SL	■



### Increased Resistance

Designed for patients who experience intermittent valve opening with fluctuations in esophageal pressures, resulting in air ingestion and/or leakage.


	6mm	8mm	10mm	12mm	Rx
16 Fr	IN 1606-IR	IN 1608-IR	IN 1610-IR	IN 1612-IR	■
20 Fr	IN 2006-IR	IN 2008-IR	IN 2010-IR	IN 2012-IR	■



### Large Esophageal Flange


Provides a seal to help manage leakage around a voice prosthesis and enables downsizing from a wider diameter prosthesis to a smaller diameter. Effectiveness depends on an optimal prosthesis length fit that achieves a snug flange-to-flange seal against the wall surface(s).

	4mm	6mm	8mm	10mm	12mm	14mm	Rx
16 Fr	IN 1604-LEF	IN 1606-LEF	IN 1608-LEF	IN 1610-LEF	IN 1612-LEF	IN 1614-LEF	■
20 Fr	IN 2004-LEF	IN 2006-LEF	IN 2008-LEF	IN 2010-LEF	IN 2012-LEF	IN 2014-LEF	■



### LARGE ESOPHAGEAL FLANGE (SPECIAL LENGTH)

	5mm	7mm	Rx
16 Fr	IN 1605-LESL	IN 1607-LESL	■
20 Fr	IN 2005-LESL	IN 2007-LESL	■





## Large Esophageal and Tracheal Flanges

Provides a double seal in patients with enlarged puncture tracts and tear-drop or irregularly shaped punctures. Effectiveness depends on an optimal prosthesis length fit that achieves a snug flange-to-flange seal against the wall surface(s).

	4mm	6mm	8mm	10mm	12mm	14mm	Rx
16 Fr	IN 1604-LF	IN 1606-LF	IN 1608-LF	IN 1610-LF	IN 1612-LF	IN 1614-LF	■
20 Fr	IN 2004-LF	IN 2006-LF	IN 2008-LF	IN 2010-LF	IN 2012-LF	IN 2014-LF	■




### LARGE ESOPHAGEAL FLANGE (INCREASED RESISTANCE)

	4mm	6mm	8mm	Rx
16 Fr	IN 1604-LEIR	IN 1606-LEIR	IN 1608-LEIR	■
20 Fr	IN 2004-LEIR	IN 2006-LEIR	IN 2008-LEIR	■

## TEP Occluders

The Blom-Singer® TEP (Tracheoesophageal Puncture) Occluder is a non-valved silicone device designed for placement in, and maintenance of, the TEP following total laryngectomy. Not a voice prosthesis.

	4mm	6mm	8mm	10mm	12mm	14mm	Rx
16 Fr	IN 1604-LETO	IN 1606-LETO	IN 1608-LETO	IN 1610-LETO	IN 1612-LETO	IN 1614-LETO	■
20 Fr	IN 2004-LETO	IN 2006-LETO	IN 2008-LETO	IN 2010-LETO	IN 2012-LETO	IN 2014-LETO	■



## VOICE PROSTHESES (PATIENT CHANGEABLE)

### Blom-Singer® Duckbill Voice Prostheses

The Blom-Singer® Duckbill is a patient changeable, slit valve prosthesis.

Length	16 Fr	Rx	Length	16 Fr	Rx
6mm	DB 16-006	■	12mm	DB 16-012	■
8mm	DB 16-008	■	14mm	BE 6001	■
10mm	BE 6000	■	18mm	BE 6002	■

### Blom-Singer® Low Pressure Voice Prostheses

The Blom-Singer® Low Pressure is a patient changeable, one-way valve prosthesis.

Length	16 Fr	20 Fr	Rx	Length	16 Fr	20 Fr	Rx
6mm	BE 6009	BE 2014	■	18mm	BE 6012	BE 2026	■
8mm	LP 16-008	LP 20-008	■	22mm	BE 6013	BE 2030	■
10mm	BE 6010	BE 2018	■	25mm	BE 6014	BE 2033	■
12mm	LP 16-012	LP 20-012	■	28mm	BE 6015	—	■
14mm	BE 6011	BE 2022	■				

### Blom-Singer® Special Order Low Pressure Voice Prostheses

The Blom-Singer® Special Order Low Pressure is a patient changeable, one-way valve prosthesis.

Length	16 Fr Special Increased Resistance	16 Fr Special Length	20 Fr Special Length	20 Fr Special Increased Resistance	Rx
4mm	—	SPECIAL 11-902-06	SPECIAL 11-251-06	—	■
6mm	LP 16-006-IR	—	—	LP 20-006-IR	■
8mm	LP 16-008-IR	—	—	LP 20-008-IR	■
10mm	LP 16-010-IR	—	—	LP 20-010-IR	■
12mm	LP 16-012-IR	—	—	LP 20-012-IR	■
14mm	LP 16-014-IR	—	—	LP 20-014-IR	■
16mm	—	SPECIAL 11-902-04	SPECIAL 11-251-04	—	■
18mm	LP 16-018-IR	—	—	—	■
20mm	—	SPECIAL 11-902-05	SPECIAL 11-251-05	—	■
22mm	LP 16-022-IR	—	—	—	■
25mm	LP 16-025-IR	—	—	—	■
28mm	LP 16-028-IR	—	—	—	■

### Blom-Singer® Voice Prosthesis Sizers

CAUTION: FEDERAL (USA) LAW RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN

The Blom-Singer® Voice Prosthesis Sizer is a device used by clinicians to measure the length of the tracheoesophageal puncture to assist in selection of the appropriate voice prosthesis.

Ref	Description
BE 6036	Voice Prosthesis Sizer with 3 inserters (for 16 Fr voice prostheses)
BE 2036	Voice Prosthesis Sizer with 3 inserters (for 20 Fr voice prostheses)
BE 6096	Replacement Inserters (for use with both BE 6036 and BE 2036 Sizers), 10 ea



## Blom-Singer® Dilator-Sizers

CAUTION: FEDERAL (USA) LAW RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN

The Blom-Singer® Dilator-Sizer combines two tasks into one innovative product: sizing and dilation. In one step, the clinician can dilate and measure the tracheoesophageal puncture during voice prosthesis placement procedures.

Ref	Description
BE 2051	18 Fr Dilator-Sizer (for 16 Fr prostheses)
BE 2052	22 Fr Dilator-Sizer (for 20 Fr prostheses)

## Blom-Singer® Tracheoesophageal Puncture Dilators

The Blom-Singer® Puncture Dilator is available for use by the patient and/or clinician. It is a tapered, solid silicone stent that is inserted into the tracheoesophageal puncture (TEP) upon removal of a voice prosthesis. It is intended to prevent the TEP from closing; prevent leakage of esophageal contents through an open TEP into the trachea during swallowing; and dilate the TEP in preparation for voice prosthesis insertion or replacement.

Ref	Description	Rx
BE 6050	18 Fr Puncture Dilator-Sizer (for 16 Fr prostheses)	■
BE 2050	22 Fr Puncture Dilator-Sizer (for 20 Fr prostheses)	■

## VOICE PROSTHESIS ACCESSORIES

### Blom-Singer® Replacement Gel Caps

Use of the Blom-Singer® Replacement Gel Caps provides a smooth, rounded shape to the esophageal flange end of the voice prosthesis, enabling easy insertion into the tracheoesophageal puncture.

Ref	Description	Rx
BE 3190	Replacement Gel Caps, 16 Fr, 90 ea	■
BE 3195	Replacement Gel Caps, 18 Fr, 90 ea	■
BE 3290	Replacement Gel Caps, 20 Fr, 90 ea	■
BE 3295	Replacement Gel Caps, 22 Fr, 60 ea	■

### Blom-Singer® Flushing Device and Cleaning Brushes

For routine cleaning of a voice prosthesis as recommended by your clinician.

Ref	Description
IN 4050	Flushing Device, pack of 3
IN 4062*	Cleaning Brush, 16/20 Fr, Standard, 4mm–10mm, 3 ea
IN 4063*	Cleaning Brush, 16/20 Fr, Long, 12mm–20mm, 3 ea

\*Fits indwelling style voice prostheses only.

### Blom-Singer® Plug Inserts

The Blom-Singer® Plug Insert is designed for temporary management of leakage through a voice prosthesis. Not compatible with patient changeable voice prostheses or Dual Valve prostheses.


Ref	Description	Rx
IN 4069	16 Fr Plug Insert (Advantage® and Classic™, 6mm to 20mm)	■
IN 4070	20 Fr Plug Insert (Advantage® and Classic™, 6mm to 20mm)	■



Flushing Device



Cleaning Brush, Standard



Cleaning Brush, Long



NEW

# Blom-Singer® HMEs

## Simple choices to fit your lifestyle

### Day&Night® HME



- Wear all day and all night, including bedtime
- Attaches to the entire line of InHealth Technologies® baseplates, StomaSoft® tubes, and Blom-Singer® Laryngectomy Buttons



Day&Night®  
EasyFlow® HME and  
ClassicFlow® HME

### SpeakFree® HME

- Your HME and Hands-Free speaking valve, all in one
- Attaches to the entire line of InHealth Technologies® baseplates, StomaSoft® tubes, and Blom-Singer® Laryngectomy Buttons



SpeakFree®  
EasyFlow® HME and  
ClassicFlow® HME

Left: Joanne  
Year of TL: 2015  
Wearing AccuFit®  
Adhesive Housing and  
Day&Night® HME

Right: Meagon  
Year of TL: 2016  
Wearing HydroFit®  
Adhesive Housing and  
SpeakFree®  
Hands-Free Valve





## HEAT & MOISTURE EXCHANGE (HME)

Research shows that long-term HME use results in less mucus<sup>1</sup>, less coughing<sup>2</sup>, and decreased shortness of breath, fatigue, sleep problems, anxiety, and depression<sup>1,2</sup>.

### Blom-Singer® SpeakFree® Hands-Free Valve with ClassicFlow® HME

Blom-Singer® SpeakFree® HME Hands-Free Valve is the first heat and moisture exchange (HME) cartridge with a single-use fully integrated adjustable hands-free valve. The SpeakFree® HME is translucent and non-reflective for a natural look. The SpeakFree® ClassicFlow® HME is designed for increased moisture retention and airflow resistance.

Ref	Description
BE 1090EZ	SpeakFree® Hands-Free Valve with ClassicFlow® HME, 30 per pack



SpeakFree® Hands-Free  
Valve with  
ClassicFlow® HME

### Blom-Singer® SpeakFree® Hands-Free Valve with EasyFlow® HME

Blom-Singer® SpeakFree® HME Hands-Free Valve is the first heat and moisture exchange (HME) cartridge with a single-use fully integrated adjustable hands-free valve. The SpeakFree® HME is translucent and non-reflective for a natural look. With SpeakFree® EasyFlow®, the foam filter is more porous and allows you to breathe more freely to match your activity level and pulmonary needs.

Ref	Description
BE 1090EF	SpeakFree® Hands-Free Valve with EasyFlow® HME, 30 per pack



SpeakFree® Hands-Free  
Valve with  
EasyFlow® HME

### Blom-Singer® Day&Night® ClassicFlow® HME

The Blom-Singer® Day&Night® ClassicFlow® can be worn all day and night, including bedtime. The Day&Night® HME is translucent and non-reflective for a natural look. The Day&Night® ClassicFlow® HME is designed for increased moisture retention and airflow resistance.

Ref	Description
BE 1091	Day&Night® ClassicFlow® HME Cartridge, pack of 30 (24-hour use)



Day&Night®  
ClassicFlow® HME

### Blom-Singer® Day&Night® EasyFlow® HME

The Blom-Singer® Day&Night® EasyFlow® can be worn all day and night, including bedtime. The Day&Night® HME is translucent and non-reflective for a natural look. The Day&Night® EasyFlow® foam filter is more porous and allows you to breathe more freely to match your activity level and pulmonary needs.

Ref	Description
BE 1092	Day&Night® EasyFlow® HME Cartridge, pack of 30 (24-hour use)



Day&Night®  
EasyFlow® HME

1. Akerstaff AH, Hilgers FJ, et al. Improvements in respiratory and psychosocial functioning following total laryngectomy by the use of a heat and moisture exchanger. *Annals of Otolary, Rhinology & Laryngology*. 1993;102(11):878–883.
2. Akerstaff AH, Hilgers FJ, et al. Heat and moisture exchangers as a treatment option in the post-operative rehabilitation of laryngectomized patients. *Clinical Otolaryngology*. 1995;20:504–509.

## ADHESIVE HOUSINGS

InHealth Technologies® presents the FitSeries—adhesive housings designed to meet your skin and stoma needs after a total laryngectomy. Compatible with Day&Night® and SpeakFree® Hands-Free HME cartridges and Shower Guard.

### Blom-Singer® AccuFit® Adhesive Housings

The Blom-Singer® AccuFit® Adhesive Housing is designed to be our strongest adhesive option while still maintaining flexibility. It has increased elasticity and is less prone to tearing. Available in three shapes: round, oval, and oval extra.


Ref	Description
BE 6082	AccuFit® Adhesive Housing, Round, pack of 30
BE 6083	AccuFit® Adhesive Housing, Oval, pack of 30
BE 6084	AccuFit® Adhesive Housing, Oval Extra, pack of 30



### Blom-Singer® HydroFit® Adhesive Housings

The Blom-Singer® HydroFit® Adhesive Housing is a hydrocolloid adhesive baseplate made with a protective gel-like layer and a moisture-absorbent adhesive. The HydroFit® is gentle on the skin and flexible, but maintains strong adhesion. Available in three shapes: round, oval, and oval extra.

Ref	Description
BE 6085	HydroFit® Adhesive Housing, Round, pack of 30
BE 6086	HydroFit® Adhesive Housing, Oval, pack of 30
BE 6087	HydroFit® Adhesive Housing, Oval Extra, pack of 30



### Blom-Singer® TruFit™ Adhesive Housing (Available Only in the US and Canada)

The Blom-Singer® TruFit™ Adhesive Housing offers a robust baseplate compatible with deeper tracheostomas. Engineered with improved support to facilitate hands-free speech when used with a Blom-Singer® voice prosthesis and the SpeakFree® HME.

Ref	Description
BE 6088	TruFit™ Adhesive Housing, Oval Extra, pack of 30



## TRACHEOSTOMA ATTACHMENTS

### Blom-Singer® Tracheostoma Valve Housings

The Blom-Singer® Tracheostoma Valve Housing is a reusable, flexible housing that attaches to the Day&Night® HME and SpeakFree® Hands-Free HME cartridge or Shower Guard over the stoma. It is applied with double sided tape discs. Available in two sizes, standard and large.

Ref	Description
BE 6038	Tracheostoma Valve Housing, standard, 1 ea
BE 6039	Tracheostoma Valve Housing, large, 1 ea



### Blom-Singer® Adhesive Tape Discs

A variety of double-sided tape discs are available to provide attachment of the tracheostoma valve housing.

Ref	Description	Ref	Description
BE 6041	Regular Tape Discs, standard, pack of 30	BE 6044	Thin Tape Discs, large, pack of 30
BE 6042	Regular Tape Discs, large, pack of 30	BE 6034	Heavy Duty Tape Discs, standard, pack of 30
BE 6043	Thin Tape Discs, standard, pack of 30	BE 6035	Heavy Duty Tape Discs, large, pack of 30



### Blom-Singer® Laryngectomy Buttons (Available Only in the US)

The Blom-Singer® Laryngectomy Button is a self-retaining holder for devices intended for vocal and pulmonary rehabilitation after total laryngectomy. It can provide stomal support and help to keep the tracheostoma open.

Ref	Description	Size	Rx
E3131208	Laryngectomy Button, 12/8, 1 ea	OD 12, L 8	■
E3131408	Laryngectomy Button, 14/8, 1 ea	OD 14, L 8	■
E3131608	Laryngectomy Button, 16/8, 1 ea	OD 16, L 8	■
E3131808	Laryngectomy Button, 18/8, 1 ea	OD 18, L 8	■



NEW

## ● AND BEYOND

### HANDS-FREE SPEECH VALVE

#### Blom-Singer® SpeakFree® Hands-Free Valve with ClassicFlow® HME

Blom-Singer® SpeakFree® HME Hands-Free Valve is the first heat and moisture exchange (HME) cartridge with a single-use fully integrated adjustable hands-free valve. The SpeakFree® HME is translucent and non-reflective for a natural look. The SpeakFree® ClassicFlow® HME is designed for increased moisture retention and airflow resistance.

Ref	Description
BE 1090EZ	SpeakFree® Hands-Free Valve with ClassicFlow® HME, 30 per pack



SpeakFree® Hands-Free Valve with ClassicFlow® HME

#### Blom-Singer® SpeakFree® Hands-Free Valve with EasyFlow® HME

Blom-Singer® SpeakFree® HME Hands-Free Valve is the first heat and moisture exchange (HME) cartridge with a single-use fully integrated adjustable hands-free valve. The SpeakFree® HME is translucent and non-reflective for a natural look. With SpeakFree® EasyFlow®, the foam filter is more porous and allows you to breathe more freely to match your activity level and pulmonary needs.

Ref	Description
BE 1090EF	SpeakFree® Hands-Free Valve with EasyFlow® HME, 30 per pack



SpeakFree® Hands-Free Valve with EasyFlow® HME

### SKIN CARE

#### Skin-Prep™ Liquid Film-Forming Protective Barrier Wipes

Skin-Prep™ is a liquid film-forming dressing that, upon application to intact skin, forms a protective film to help reduce friction during removal of tapes and films.

Ref	Description
4204-00	Skin-Prep™ Liquid Film-Forming Protective Barrier Wipes, box of 50



#### Skin Tac™ Adhesive Barrier Wipes

Skin Tac™ liquid adhesive barrier wipes are latex-free, hypo-allergenic, clear, non-rubber adhesive barrier wipes.

Ref	Description
MS407W*	Skin Tac™ Adhesive Barrier Wipes, box of 50

\*May not be available in all international markets.



#### NFS Adhesive (Silicone Adhesive)

Pressure Sensitive Adhesive (PSA) is a silicone-based adhesive for use with all silicone prosthetic devices.

Ref	Description
B-530	NFS Adhesive, 1 oz





## Remove™ Adhesive Remover

Gently cleans hydrocolloid, acrylic and rubber-based adhesive residues from the skin without irritation, discomfort or trauma.

Ref	Description
59403125	Remove™ Adhesive Remover, box of 50



## STOMA PROTECTION AND ACCESSORIES

### ADDvox® Stoma Cover Mock Turtlenecks (Only available in the US)

The ADDvox® Stoma Cover Mock Turtleneck has inner layers of material that help keep debris from entering the stoma. They slip under a shirt like a fashion dickey and have an outer layer of soft 100% cotton knit. Each cover is 9-1/2" long x 8-3/4" wide. Velcro fastener adjusts to fit necks of 12-1/2" to 22". Each color sold separately.

Ref	Description	Color	Quantity
BE 6227-R2	ADDvox® Stoma Cover Mock Turtleneck	White	1 ea
BE 6228-R2	ADDvox® Stoma Cover Mock Turtleneck	Black	1 ea
BE 6229-R2	ADDvox® Stoma Cover Mock Turtleneck	Navy	1 ea



### ADDvox® Stoma Scarves (Only available in the US)

The ADDvox® Stoma Cover Scarf has inside layers to help filter out dust, dirt and pollen, is completely washable, and needs little care. This attractive neckline scarf looks just like a regular scarf or ascot and is available in a wide range of colors for men or women. All scarves have Velcro closures and adjust to fit necks from 17-3/4" to 22" (can be made smaller with tucks or stitching). Each color sold separately.

Ref	Description	Color	Quantity
BE 6232-R2	ADDvox® Stoma Cover Scarf	Black	1 ea
BE 6233-R2	ADDvox® Stoma Cover Scarf	White	1 ea
BE 6234-R2	ADDvox® Stoma Cover Scarf	Navy	1 ea



### ADDvox® Cloth Bib-Style Stoma Covers (Only available in the US)

The ADDvox® soft, long-lasting bib protects the tracheostoma and lets you breathe more easily. This "bib-style" Stoma Cover will look crisp and clean after repeated laundering. Five soft layers keep dust and dirt from reaching the trachea. Fits necks up to 22". Velcro closure. Each color sold separately.

Ref	Description	Color	Quantity
BE 6224-R2	ADDvox® Cloth Bib-Style Stoma Cover	White	1 ea
BE 6225-R2	ADDvox® Cloth Bib-Style Stoma Cover	Beige	1 ea
BE 6226-R2	ADDvox® Cloth Bib-Style Stoma Cover	Navy	1 ea



## ADDvox® Silicone Shower Shield (Only available in the US)

The one-piece ADDvox® Silicone Shower Shield is worn during showering to deflect the water away from your stoma. Extremely soft and pliable, it creates a soft, comfortable and conforming seal along the irregular contours of the neck, and vents on the underside allowing for airflow. Made of medical grade silicone, this latex-free Silicone Shower Shield is comfortable, durable, and not affected by temperature, ozone, sunlight, soap or detergents.

Ref	Description
BE 6079-R2	ADDvox® Silicone Shower Shield, 1 ea

## Blom-Singer® Shower Guard

The Blom-Singer® Shower Guard is a one-piece, right-angle plastic shower guard that comes with a Tracheostoma Standard Adhesive Housing and 30 double-sided Adhesive Tape Discs. Also compatible with the Blom-Singer® StomaSoft® Laryngectomy Tube and the Blom-Singer® Laryngectomy Button.

Ref	Description
BE 6048	Blom-Singer® Shower Guard—Includes Tracheostoma Valve Housing (BE 6038) and Regular Tape Discs (BE 6041), pack of 30

## Blom-Singer® Foam Stoma Protector

A simple and temporary (disposable) solution to help prevent debris from entering the stoma. Adhesive strip on one edge allows for easy attachment, access to the stoma to occlude for speech with a voice prosthesis, or for coughing. Measurement listed is approximate.

Ref	Description
BE 6220	Blom-Singer® Foam Stoma Protector, Tan, Large 2.7" L x 2.7" W (7cm x 7cm), pack of 30

## ADDvox® Stoma Filters with Microporous Adhesive

ADDvox® Stoma Filters offer protection against dust, dirt and pollen. Designed for single use with a hypoallergenic, non-irritating, self-adhesive microporous strip that is comfortable and easy to apply. Not recommended for cannula users. 2-1/4" wide x 2-1/2" high x 1/8" thick. Available in off-white and skin tone.

Ref	Description
BE 6222-R2	ADDvox® Stoma Filters, Skin tone, pack of 60
BE 6223-R2	ADDvox® Stoma Filters, Off-white, pack of 60





**Professional Products**

## Blom-Singer® Voice Prosthesis Sizers

CAUTION: FEDERAL (USA) LAW RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN

The Blom-Singer® Voice Prosthesis Sizer is a device used by clinicians to measure the length of the tracheoesophageal puncture to assist in selection of the appropriate voice prosthesis.

Ref	Description
BE 6036	Blom-Singer® Voice Prosthesis Sizer with 3 inserters (for 16 Fr voice prostheses)
BE 2036	Blom-Singer® Voice Prosthesis Sizer with 3 inserters (for 20 Fr voice prostheses)
BE 6096	Replacement Inserters (for use with both BE 6036 and BE 2036 Sizers), 10 ea

## Blom-Singer® Dilator-Sizers

CAUTION: FEDERAL (USA) LAW RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN

The Blom-Singer® Dilator-Sizer combines two tasks into one innovative product: sizing and dilation. In one step, the clinician can dilate and measure the tracheoesophageal puncture during voice prosthesis placement procedures.

Ref	Description
BE 2051	Blom-Singer® 18 Fr Dilator-Sizer (for 16 Fr prostheses)
BE 2052	Blom-Singer® 22 Fr Dilator-Sizer (for 20 Fr prostheses)

## Blom-Singer® Adjustable Bi-Flanged Fistula Prostheses (Sterile)

CAUTION: FEDERAL (USA) LAW RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN

Designed by Eric D. Blom, Ph.D to control leakage of saliva, food, and drink around a hypopharyngeal fistula. The prosthesis is adjustable in length and comes in two flange sizes, 25mm and 38mm, ensuring a proper fit.

Ref	Description
FM 1003	25mm (Small) Adjustable Bi-Flanged Fistula Prosthesis, 1 ea
FM 1004	38mm (Medium) Adjustable Bi-Flanged Fistula Prosthesis, 1 ea

## Manometer with Manometer Adaptor

A tool to teach patients how to self-regulate intratracheal air pressure. Proper airflow control may allow for better voice quality, less fatigue, and longer duration of the tracheostoma adhesive housing seal.

Ref	Description
BE 8000	Manometer Kit includes 1 manometer, 1 manometer adaptor and 1 adhesive housing
BE 8020	Replacement Manometer Adaptor includes 5 manometer adaptors and 5 adhesive housings

## Blom-Singer® Insufflation Test Set

The Blom-Singer® Insufflation Test Set is designed to aid the clinician in subjectively predicting the ability of the laryngectomized patient to produce tracheoesophageal voice.

Ref	Description
BE 6040	Blom-Singer® Insufflation Test Set, box of 6






Kevin  
Year of TL: 2001  
Wearing SpeakFree®  
Hands-Free HME



# Katz Extractor<sup>®</sup>

OTO-RHINO FOREIGN BODY REMOVER

The easy way out



## Katz Extractor® Oto-Rhino Foreign Body Remover

Designed in conjunction with David L. Katz, MD, MPH, the Katz Extractor® Oto-Rhino Foreign Body Remover is intended for extraction of impacted foreign bodies in the nasal passage. Single use only. This device may also be used in the external ear canal, based upon clinical judgment. Covered by one or more of the following patents and other patents pending: US#5,454,817. Katz Extractor® is a registered trademark of Freudenberg Medical, LLC.


Ref	Description
FB 5000	Katz Extractor® Oto-Rhino Foreign Body Remover, box of 5

### TO LEARN MORE

Scan QR code  
to visit the  
Katz Extractor® website







# Nasal Septal Perforation Prosthesis

A two-piece, magnetically coupled solution for non-surgical closure of nasal septal perforations.

## Nasal Septal Perforation Sizing Tool

With the Nasal Septal Perforation Sizing Tool, medical professionals can deliver better results with less complexity, all while prioritizing patient safety and comfort. Available in both oval and round configurations to accommodate diverse patient needs.

Ref	Description
SPT 002	Nasal Septal Perforation Sizing Tool, 1 box of 5 pairs, oval and round

## Blom-Singer® Nasal Septal Perforation Prostheses













The Blom-Singer® Nasal Septal Perforation Prosthesis is a two-piece, magnetically coupled solution for non-surgical closure of nasal septal perforations.

### OVAL

Ref	Description
SP 2112	Blom-Singer® Nasal Septal Perforation Prosthesis Oval 11 x 13 x 2mm
SP 2113	Blom-Singer® Nasal Septal Perforation Prosthesis Oval 11 x 13 x 3mm
SP 2132	Blom-Singer® Nasal Septal Perforation Prosthesis Oval 13 x 17 x 2mm
SP 2133	Blom-Singer® Nasal Septal Perforation Prosthesis Oval 13 x 17 x 3mm
SP 2152	Blom-Singer® Nasal Septal Perforation Prosthesis Oval 15 x 21 x 2mm
SP 2153	Blom-Singer® Nasal Septal Perforation Prosthesis Oval 15 x 21 x 3mm

### ROUND

Ref	Description
SP 1052	Blom-Singer® Nasal Septal Perforation Prosthesis Round 5 x 2mm
SP 1053	Blom-Singer® Nasal Septal Perforation Prosthesis Round 5 x 3mm
SP 1072	Blom-Singer® Nasal Septal Perforation Prosthesis Round 7 x 2mm
SP 1073	Blom-Singer® Nasal Septal Perforation Prosthesis Round 7 x 3mm
SP 1102	Blom-Singer® Nasal Septal Perforation Prosthesis Round 10 x 2mm
SP 1103	Blom-Singer® Nasal Septal Perforation Prosthesis Round 10 x 3mm



## ORDERING INFORMATION

### UNITED STATES

**Phone:** Toll-Free: 800.477.5969

**Fax:** Toll-Free: 888.371.1530

**E-mail:** [order@inhealth.com](mailto:order@inhealth.com)

**Prescription E-mail:** [ihtrx@inhealth.com](mailto:ihtrx@inhealth.com)

**Order Online:** [www.inhealth.com](http://www.inhealth.com)

**Mail:**

InHealth Technologies®  
1110 Mark Ave.  
Carpinteria, CA 93013 • USA

### CANADA

**Phone:** Toll-Free: 1.888.437.4123

**Fax:** Toll-Free: 1.226.641.5285

**E-mail:** [canadaorder@inhealth.com](mailto:canadaorder@inhealth.com)

**Website:** [www.inhealth.com/canada](http://www.inhealth.com/canada)

**Mail:**

InHealth Canada  
65 Spruce Street  
Tillsonburg ON N4G 5C4

**Rx: Products  
with the Rx  
notation require  
a prescription  
from a physician  
(US patients).**

## NEW PATIENT ENROLLMENT INSTRUCTIONS

### Enrolling as a new patient? It's easy! Here is what to do:

1. Obtain a **prescription** from your provider.
2. Fill out the **Patient Enrollment Form**. All forms are available under the **Privacy Policy and Forms** tab on the InHealth Technologies® website: [www.inhealth.com](http://www.inhealth.com)
3. Send back to InHealth Technologies® with a copy of your insurance card (front and back). You can send via email: [insurance@inhealth.com](mailto:insurance@inhealth.com) or fax: **888.371.1530**
4. Place your order by calling **1.800.477.5969** or emailing [order@inhealth.com](mailto:order@inhealth.com).

### Need insurance or Medicare help? Contact us!

Phone: **800.477.5969** option 3

Email: [medicare@inhealth.com](mailto:medicare@inhealth.com)

While you don't have to fill anything out on the **Patient Bill of Rights** form, please make sure you review it. We want you to know your rights.


## FREE EDUCATIONAL RESOURCES

- Available on [www.inhealth.com](http://www.inhealth.com) under the **Education** and **Resources** tabs.
- Visit InHealth Technologies® video channel at <https://inhealth.com/videos> to view our videos.



In 1978, Drs. Singer and Blom pioneered the methods of tracheoesophageal puncture and valved silicone voice prostheses that have been the international standard for voice restoration for over 45 years, allowing thousands of patients worldwide to regain their ability to speak. As a leading supplier of advanced voice restoration systems, our goal at InHealth Technologies® is to help laryngectomees regain and preserve a vital part of their lives—their voice. We are proud to continue the Blom-Singer® legacy today with clinical support, input, and leadership.

**Blom-Singer®**







Janice  
Year of TL: 2014  
Wearing AccuFit®  
Adhesive Housing  
with SpeakFree® HME  
Hands-Free Valve

## STAY CONNECTED...



PATIENTS  
SCAN QR CODE

### Speakers Club™

**Patients: Scan the QR code to sign up today!**

We created this digital E-News for you! **Speakers Club™ E-News** will provide you with information on upcoming patient events, new products and resources.



PROFESSIONALS  
SCAN QR CODE

### ProNET™

**Professionals: Scan the QR code to sign up today!**

We created this digital E-News specifically for professionals. **ProNET™ E-News** will deliver useful laryngectomy and voice prosthesis information for you and your practice, including material you may want to share with your patients.

**INHEALTH**  
TECHNOLOGIES®

United States  
Toll-Free: 800.477.5969  
Fax: 888.371.1530  
Email: [info@inhealth.com](mailto:info@inhealth.com)

Canada  
Toll-Free: +1.888.437.4123  
Fax: +1.226.641.5285  
Email: [canadaorder@inhealth.com](mailto:canadaorder@inhealth.com)

International  
Tel: +1.805.684.9337  
Fax: +1.805.684.8594  
E-mail: [international@inhealth.com](mailto:international@inhealth.com)

Terms & Conditions: <https://inhealth.com/terms-conditions>  
©2025 InHealth Technologies®—Manufactured by Freudenberg Medical, LLC (241104.00)

**Blom-Singer®**

The voice of experience since 1978



[www.inhealth.com](http://www.inhealth.com)