

# InHealth Technologies® Medicare Assignment Program Products

## Medicare Insurance Information

The InHealth Technologies® Medicare Assignment Program is available for patients with **Traditional Medicare Part B** as their primary insurance. If you qualify, you can enroll in the program; please allow up to 48 hours for enrollment before placing your order. With Medicare Assignment, **you are typically responsible only for your 20% co-payment** and any remaining deductible. If you have secondary insurance, your out-of-pocket cost may be \$0.

## Blom-Singer® StomaSoft® Laryngectomy Tubes

Blom-Singer® StomaSoft® Laryngectomy Tubes are flexible, soft, and compatible with InHealth Technologies® HMEs and other accessories from InHealth Technologies®. They feature reinforced strap holders and are available with or without fenestrations.

REF	Size	ID/OD	Fenestration	Rx
BE 6698	8/36	9.5/12.0	Non-fenestrated	■
BE 6699	8/55	9.5/12.0	Non-fenestrated	■
BE 6700	9/36	10.5/13.5	Non-fenestrated	■
BE 6701	9/55	10.5/13.5	Non-fenestrated	■
BE 6702	10/36	12.0/15.0	Non-fenestrated	■
BE 6703	10/55	12.0/15.0	Non-fenestrated	■
BE 6704	12/36	13.5/17.0	Non-fenestrated	■
BE 6705	12/55	13.5/17.0	Non-fenestrated	■
BE 6698F	8/36	9.5/12.0	Fenestrated	■
BE 6699F	8/55	9.5/12.0	Fenestrated	■
BE 6700F	9/36	10.5/13.5	Fenestrated	■
BE 6701F	9/55	10.5/13.5	Fenestrated	■
BE 6702F	10/36	12.0/15.0	Fenestrated	■
BE 6703F	10/55	12.0/15.0	Fenestrated	■
BE 6704F	12/36	13.5/17.0	Fenestrated	■
BE 6705F	12/55	13.5/17.0	Fenestrated	■

## Blom-Singer® Laryngectomy Tubes - Sterile

Blom-Singer® Laryngectomy Tubes (Sterile) are durable, lightweight silicone tubes that are used to keep the tracheostoma open after laryngectomy. The tubes do not hold an HME.

REF	Size	ID/OD	Fenestration	Rx
BE 6398	8/36	9.5/12.0	Non-fenestrated	■
BE 6399	8/55	9.5/12.0	Non-fenestrated	■
BE 6400	9/36	10.5/13.5	Non-fenestrated	■
BE 6401	9/55	10.5/13.5	Non-fenestrated	■
BE 6402	10/36	12.0/15.0	Non-fenestrated	■
BE 6403	10/55	12.0/15.0	Non-fenestrated	■
BE 6404	12/36	13.5/17.0	Non-fenestrated	■
BE 6405	12/55	13.5/17.0	Non-fenestrated	■

## Blom-Singer® Laryngectomy Buttons

The Blom-Singer® Laryngectomy Button is a self-retaining holder for devices intended for vocal and pulmonary rehabilitation after total laryngectomy. It can provide stomal support and help to keep the tracheostoma open.

REF	Description	Rx
E3131208	Laryngectomy Button 12/8	■
E3131408	Laryngectomy Button 14/8	■
E3131608	Laryngectomy Button 16/8	■
E3131808	Laryngectomy Button 18/8	■

## StomaSoft® Cleaning Swabs

The StomaSoft® Cleaning Swabs are intended for cleaning laryngectomy tubes or buttons after removal from the patient. Available in packs of 10.

REF	Description
TSWB-10	StomaSoft® Foam Swab, pack of 10

## Marpac Tracheostomy Tube Holder/Collar™

Marpac Tracheostomy Tube Holder/Collar™ (for laryngectomy patients) is a single use, disposable Comfortable Collar™ used to secure a Blom-Singer® Laryngectomy Tube or a Blom-Singer® Laryngectomy Button.

REF	Description
204	Marpac Tracheostomy Tube Holder/Collar™

## Blom-Singer® Duckbill Voice Prosthesis

The Blom-Singer® Duckbill is a patient changeable, slit valve prosthesis.

16 Fr.					
REF	Length	Rx	REF	Length	Rx
DB16-006	6mm	■	DB16-012	12mm	■
DB16-008	8mm	■	BE 6001	14mm	■
BE 6000	10mm	■	BE 6002	18mm	■

## Blom-Singer® Low Pressure Voice Prosthesis

The Blom-Singer® Low Pressure is a patient changeable, one-way valve prosthesis.

16 Fr.			20 Fr.		
REF	Length	Rx	REF	Length	Rx
BE 6009	6mm	■	BE 2014	6mm	■
LP16-008	8mm	■	LP20-008	8mm	■
BE 6010	10mm	■	BE 2018	10mm	■
LP16-012	12mm	■	LP20-012	12mm	■
BE 6011	14mm	■	BE 2022	14mm	■
BE 6012	18mm	■	BE 2026	18mm	■
BE 6013	22mm	■	BE 2030	22mm	■
BE 6014	25mm	■	BE 2033	25mm	■
BE 6015	28mm	■	—	—	—

## Blom-Singer® Electrolarynx

An electrolarynx is a speaking aid that produces sound vibrations that can be conducted through the throat and mouth to speak after a laryngectomy.

REF	Description
EL 2000	Electrolarynx

## Blom-Singer® SpeakFree® Hands-Free Valves

Blom-Singer® SpeakFree® Hands-Free Valve is the first heat and moisture exchange (HME) cartridge with a single-use fully integrated adjustable hands-free valve. The SpeakFree® HME is translucent and non-reflective for a natural look. The SpeakFree® ClassicFlow® HME is designed for increased moisture retention and airflow resistance. With SpeakFree® EasyFlow®, the foam filter is more porous and allows you to breathe more freely to match your activity level and pulmonary needs.

REF	Description
BE 1090EF	SpeakFree® EasyFlow® HME, pack of 30
BE 1090EZ	SpeakFree® ClassicFlow® HME, pack of 30

## Blom-Singer® HMEs (for 24-hour use)

The Blom-Singer® Day&Night® ClassicFlow® can be worn all day and night, including bedtime. The Day&Night® HME is translucent and non-reflective for a natural look. The Day&Night® ClassicFlow® HME is designed for increased moisture retention and airflow resistance. The Day&Night® EasyFlow® foam filter is more porous and allows you to breathe more freely to match your activity level and pulmonary needs.

REF	Description
BE 1091	Day&Night® ClassicFlow® HME, pack of 30
BE 1092	Day&Night® EasyFlow® HME, pack of 30

## Blom-Singer® AccuFit® Adhesive Housings

The Blom-Singer® AccuFit® Adhesive Housing is designed to be our strongest adhesive option while still maintaining flexibility. It has increased elasticity and is less prone to tearing. Available in three shapes: round, oval, and oval extra.

REF	Description
BE 6082	AccuFit® Adhesive Housing, Round, pack of 30
BE 6083	AccuFit® Adhesive Housing, Oval, pack of 30
BE 6084	AccuFit® Adhesive Housing, Oval Extra, pack of 30

## Blom-Singer® Hydrofit® Adhesive Housings

The Blom-Singer® Hydrofit® Adhesive Housing is a hydrocolloid adhesive baseplate made with a protective gel-like layer and a moisture absorbent adhesive. The Hydrofit® is gentle on the skin and flexible, but maintains strong adhesion. Available in three shapes: round, oval, and oval extra.

REF	Description
BE 6085	Hydrofit® Adhesive Housing, Round, pack of 30
BE 6086	Hydrofit® Adhesive Housing, Oval, pack of 30
BE 6087	Hydrofit® Adhesive Housing, Oval Extra, pack of 30

## Blom-Singer® TruFit™ Adhesive Housings

The Blom-Singer® TruFit™ Adhesive Housing offers a robust baseplate compatible with deeper tracheostomas. Engineered with improved support to facilitate hands-free speech when used with a Blom-Singer® voice prosthesis and the SpeakFree® HME.

REF	Description
BE 6088	TruFit™ Adhesive Housing, Oval Extra, pack of 30

## Skin-Prep™ Liquid Film-Forming Protective Barrier Wipes

Skin-Prep™ is a liquid film-forming dressing that, upon application to intact skin, forms a protective film to help reduce friction during removal of tapes and films.

REF	Description
4204-00	Skin-Prep™ Liquid Film-Forming Protective Barrier Wipes, box of 50

## Skin Tac™ Adhesive Barrier Wipes

Skin Tac™ liquid adhesive barrier wipes are latex-free, hypo-allergenic, clear, non-rubber adhesive barrier wipes.

REF	Description
MS407W	Skin Tac™ Adhesive Barrier Wipes, box of 50

## Remove™ Adhesive Remover

Gently cleans hydrocolloid, acrylic and rubber-based adhesive residues from the skin without irritation, discomfort or trauma.

REF	Description
59403125	Remove™ Adhesive Remover, box of 50

## Blom-Singer® Tracheoesophageal Puncture Dilators

The Blom-Singer® Puncture Dilator is available for use by the patient and/or clinician. It is a tapered, solid silicone stent that is inserted into the tracheoesophageal puncture (TEP) upon removal of a voice prosthesis. It is intended to prevent the TEP from closing; prevent leakage of esophageal contents through an open TEP into the trachea during swallowing; and dilate the TEP in preparation for voice prosthesis insertion or replacement.

REF	Description	Rx
BE 6050	18 Fr Puncture Dilator-Sizer (for 16 Fr prostheses)	■
BE 2050	22 Fr Puncture Dilator-Sizer (for 20 Fr prostheses)	■

## Blom-Singer® Plug Inserts

The Blom-Singer® Plug Insert is designed for temporary management of leakage through a voice prosthesis. Not compatible with patient changeable voice prostheses or Dual Valve prostheses.

REF	Description	Rx
IN 4069	16 Fr Plug Insert (Advantage® and Classic™, 6mm to 20mm)	■
IN 4070	20 Fr Plug Insert (Advantage® and Classic™, 6mm to 20mm)	■

## Blom-Singer® Cleaning Brushes

For routine cleaning of a voice prosthesis as recommended by your clinician. Fits indwelling style voice prosthesis only.

REF	Description
IN 4062	Cleaning Brush, 16/20 Fr, Standard, 4mm-10mm, 3 ea
IN 4063	Cleaning Brush, 16/20 Fr, Long, 12mm-20mm, 3 ea

**INHEALTH**  
TECHNOLOGIES®  
www.inhealth.com

1110 Mark Avenue  
Carpinteria, CA 93013  
Tel: 800.477.5969  
Fax: 734.354.5757  
insurance@inhealth.com