

HME STYLES & FEATURES

GUIDE

INHEALTH
TECHNOLOGIES®

| STYLES |  |  |  |  |
|---|---|---|---|---|
| FEATURES | Blom-Singer® SpeakFree® HME Hands-Free Valve with ClassicFlow® HME | Blom-Singer® SpeakFree® HME Hands-Free Valve with EasyFlow® HME | Blom-Singer® Day&Night® ClassicFlow® HME | Blom-Singer® Day&Night® EasyFlow® HME |
| Adjustable Hands Free Speech Valve | ✓ | ✓ | | |
| Single-Use Device | ✓ | ✓ | ✓ | ✓ |
| Easy Manual Occlusion for Voicing | ✓ | ✓ | ✓ | ✓ |
| Breathability | GOOD | GOOD | BEST | BEST |
| Moisture Exchange | BETTER | GOOD | BETTER | GOOD |
| Hygroscopic Filter | ✓ | ✓ | ✓ | ✓ |
| Day and Night Wear | DAYTIME | DAYTIME | 24 HR USE | 24 HR USE |
| Compatible with all Blom-Singer® Housings and StomaSoft® Tube | ✓ | ✓ | ✓ | ✓ |
| Optional Safety Feature | CROSSBAR | CROSSBAR | TWIST TO LOCK | TWIST TO LOCK |
| Made with Latex-free Materials | ✓ | ✓ | ✓ | ✓ |
| Non-reflective Surface Texture | ✓ | ✓ | ✓ | ✓ |
| Reference Number | BE 1090EZ 30 pk | BE 1090EF 30 pk | BE1091 30 pk | BE1092 30 pk |

Order online at www.inhealth.com

Toll-free: 800.477.5969

Email: order@inhealth.com

TOTAL
Solution