

Introducing the **InHealth 22.5 Fr Indwelling Voice Prosthesis**

*Advantage Hard Valve Assembly with
Large Esophageal Flange*



InHealth Technologies® is committed to providing you with voice prosthesis solutions to meet your patients' needs. Introducing the newest Blom-Singer® voice prosthesis, the InHealth 22.5 Fr Indwelling Voice Prosthesis, designed to aid in troubleshooting both central and peripheral leakage.

Features:

- **Advantage Hard Valve Assembly**

Valve is seated with a titanium ring and embedded with silver oxide, a known material preservative, **to help promote longer device life**

- **Large Esophageal Flange**

Larger esophageal flange to help **manage potential leakage around the voice prosthesis**

Sizes Available:

Ref	Description
VP 2206-LEF	22.5 Fr Advantage Hard Valve Assembly LEF Indwelling Voice Prosthesis, Sterile, 6mm
VP 2208-LEF	22.5 Fr Advantage Hard Valve Assembly LEF Indwelling Voice Prosthesis, Sterile, 8mm
VP 2210-LEF	22.5 Fr Advantage Hard Valve Assembly LEF Indwelling Voice Prosthesis, Sterile, 10mm
VP 2212-LEF	22.5 Fr Advantage Hard Valve Assembly LEF Indwelling Voice Prosthesis, Sterile, 12mm