



Katz Extractor[®]

OTO-RHINO FOREIGN BODY REMOVER

The easy way out

The fast, flexible way to remove foreign bodies from noses and ears

- **Efficient and Safe**
Great for restless patients.
- **Reliable**
Trusted by healthcare professionals worldwide.
- **Convenient**
Lightweight, fully assembled and ready to use.

**Order the
Katz Extractor[®]
Five-pack today!**

Katz Extractor®

OTO-RHINO FOREIGN BODY REMOVER

A quick and simple solution for foreign body removal from the nose or ear.




The Katz Extractor® Oto-Rhino Foreign Body Remover is a small, lightweight, single-use device for quick and easy extraction of foreign bodies from the nose or external ear canal. It comes ready to use in a five-pack dispenser box. There are no expensive supplies, such as forceps, hooks or balloon tip catheters/syringes to locate and assemble. No additional suction equipment is required.

To see a demonstration video, visit:
www.katzextractor.com




Please Scan
to access website
and IFU

It takes just three simple steps to complete an extraction with the Katz Extractor!



Insert past foreign body



Depress the syringe plunger to inflate balloon



Extract catheter out of canal with balloon inflated

ORDER INFORMATION

When ordering, please reference

Product Number:
FB 5000

Tel: **1.888.437.4123**

Email: Canadaorder@inhealth.com

CUSTOMER SERVICE

Our knowledgeable representatives are available to assist you
Monday through Thursday 9:00 a.m. to 4:00 p.m.
and Friday 9:00 a.m. to 2:00 p.m., EST

Designed in conjunction with David L. Katz, M.D., MPH, the Katz Extractor is intended principally for extraction of impacted foreign bodies in the nasal passages. This device may also be used in the external ear canal, based upon clinical judgement.

INHEALTH
CANADA

Stay Connected. Like, Share, Follow.



©2025 InHealth Technologies

Manufactured by Freudenberg Medical, LLC (210412.03)