

**INHEALTH**  
TECHNOLOGIES®

# Electrolarynx EL 1000



Order Online at  
[www.InHealth.com](http://www.InHealth.com)

**New charcoal gray color**

**Electronic volume and pitch push buttons have replaced mechanical wheels**

*– Easier operation and improved reliability*

**Fixed Device Head has replaced adjustable head**

*– Improved reliability, no adjustment needed*

**High contrast bubble labels**

*– Easier to read volume and pitch buttons*

**Raised On/Off button**

*– Easier to find and press*

**Anti-Roll Feature**

**Oral Adapter, Lanyard, and One 9-volt batteries included**

**3-Year Warranty**