

HME STYLES & FEATURES

GUIDE

INHEALTH
TECHNOLOGIES®

STYLES				
FEATURES	Blom-Singer® Day&Night® ClassicFlow® HME	Blom-Singer® Day&Night® EasyFlow® HME	Blom-Singer® SpeakFree® HME Hands-Free Valve with ClassicFlow® HME	Blom-Singer® SpeakFree® HME Hands-Free Valve with EasyFlow® HME
Wear Time	24 HR USE	24 HR USE	DAYTIME	DAYTIME
Single-Use Device	✓	✓	✓	✓
Easy Manual Occlusion for Voicing	✓	✓	✓	✓
Adjustable Hands-Free Speech Valve			✓	✓
Airflow				
Moisture Retention				
Hygroscopic Filter	✓	✓	✓	✓
International Standard 22mm Hub (Compatible with all InHealth products)	✓	✓	✓	✓
Optional Safety Feature	TWIST TO LOCK	TWIST TO LOCK	CROSSBAR	CROSSBAR
Not Made with Natural Rubber Latex	✓	✓	✓	✓
Nonreflective Surface Texture	✓	✓	✓	✓
Reference Number	BE 1091 30 pk	BE 1092 30 pk	BE 1090EZ 30 pk	BE 1090EF 30 pk

www.inhealth.com
Toll-free: 800.477.5969

TOTAL

Solution

Blom-Singer®
The voice of experience since 1978