

HME STYLES & FEATURES

GUIDE

INHEALTH
TECHNOLOGIES®

STYLES				
FEATURES	Blom-Singer® SpeakFree® HME Hands-Free Valve with ClassicFlow® HME	Blom-Singer® SpeakFree® HME Hands-Free Valve with EasyFlow® HME	Blom-Singer® Day&Night® ClassicFlow® HME	Blom-Singer® Day&Night® EasyFlow® HME
Adjustable Hands Free Speech Valve	✓	✓		
Single-Use Device	✓	✓	✓	✓
Easy Manual Occlusion for Voicing	✓	✓	✓	✓
Breathability	GOOD	GOOD	BEST	BEST
Moisture Exchange	BETTER	GOOD	BETTER	GOOD
Hygroscopic Filter	✓	✓	✓	✓
Day and Night Wear	DAYTIME	DAYTIME	24 HR USE	24 HR USE
Compatible with all Blom-Singer® Housings and StomaSoft® Tube	✓	✓	✓	✓
Optional Safety Feature	CROSSBAR	CROSSBAR	TWIST TO LOCK	TWIST TO LOCK
Made with Latex-free Materials	✓	✓	✓	✓
Non-reflective Surface Texture	✓	✓	✓	✓
Reference Number	BE 1090EZ 30 pk	BE 1090EF 30 pk	BE1091 30 pk	BE1092 30 pk

Order online at www.inhealth.com

Toll-free: 800.477.5969

Email: order@inhealth.com

TOTAL

Solution