

# The Blom Singer® Nasal Septal Perforation Prosthesis

A two-piece, magnetically coupled solution for non-surgical closure of nasal septal perforations.

A nasal septal perforation can vary in size and location, and it may be asymptomatic or cause a variety of aggravating and uncomfortable symptoms. The Blom-Singer® Nasal Septal Perforation Prosthesis comes in a variety of shapes, sizes, and thicknesses for a comfortable fit.



2-Piece Magnetically  
Coupled Prosthesis



Multiple Options  
to help select the  
right product fit



Reduce or eliminate  
symptoms<sup>1</sup>



Non-Surgical  
Closure

## MAGNETIC COUPLING

Neodymium magnets encased in medical grade silicone coupled to provide minimally invasive closure of the nasal septal perforation.<sup>2</sup>






ROUND






OVAL

## PRODUCT INFORMATION

### ROUND



### OVAL



**A better fit means less movement and a reduction in bleeding and crusting.**

## BILLING CODES

2026 Blom-Singer® Nasal Septal Perforation Prosthesis		Ambulatory Surgery Center	Outpatient Surgery at Hospital	Physician Fee Schedule
CPT Code	Procedure Description	Facility Payment (National Medicare Average)	APC (Ambulatory Payment Classifications)	Payment (National Medicare Average)
30220	Insertion Nasal Septal Prosthesis	\$659.17	\$686.92	\$261.85 - \$407.70
31231	Nasal Endoscopy, Diagnostic, Unilateral or Bilateral	\$104.73	\$139.29	\$207.40-\$300.00
HCPCS Code	HCPCS Short Description	Description		Fee Schedule (National Medicare Average)
L8047	Nasal Septal Prosthesis	A removable superficial prosthesis which occludes a hole in the nasal septum but which does not include superficial nasal tissue.		80% of product retail price
HCPCS Code	HCPCS Short Description	Description		Fee Schedule (National Medicare Average)
L8699	Prosthetic Implant NOS	Prosthesis implant, not otherwise specified		80% of product retail price

### Blom-Singer® Nasal Septal Perforation Prosthesis

According to PDAC and CMS, HCPCS L8047 is allowed once every five years.

If the physician changes the device sooner, they are able to submit the claim, along with a Medical Necessity explanation.

The physician is able to submit a claim for their service, using CPT 30220, as needed (no limitation).

Reimbursement information is applicable solely to the U.S. market and is subject to change.

\* the symptoms are, loss of airflow, crusting, whistling, epistaxis, and infections

1. Mullace, M, et al. *Acta Otorhinolaryngol Ital* 26, 2006; 216-218

2. Teschner M, et al. *Eur Arch Otorhinolaryngol*. 2012. Mar;269(3):861-5

**INHEALTH**  
TECHNOLOGIES®

www.inhealth.com

United States  
Toll-Free: 800.477.5969  
Fax: 888.371.1530  
email: info@inhealth.com

International  
Tel: +1.805.684.9337  
Fax: +1.805.684.8594  
email: international@inhealth.com

Canada  
Toll-Free: +1.888.437.4123  
Fax: +1.226.641.5285  
email: canadaorder@inhealth.com