

Introducing: The **NEW TruFit™** Adhesive Housing

The TruFit™ Adhesive Housing offers a robust baseplate compatible with deeper tracheostomas.



NEW
TruFit™
One size & one shape

NEW **TruFit™**

OVAL EXTRA
BE 6088

Key Benefits of TruFit™:

- Designed for patients with deep set stomas.
- Engineered with improved support to facilitate hands-free speech when used with a Blom-Singer® voice prosthesis and the SpeakFree® HME.
- Same reliable, strong adhesion as AccuFit®.
- Easy application with three piece peel-off backing.

AccuFit®



Key Benefits of AccuFit®:

- Designed with strong adhesion, while maintaining flexibility and stretch.
- Increased elasticity and higher resistance to tearing.

Three Shapes:

ROUND BE 6082	OVAL BE 6083	OVAL EXTRA BE 6084
-------------------------	------------------------	------------------------------

HydroFit®



Key Benefits of HydroFit®:

- Hydrocolloid adhesive baseplate made with a protective gel-like layer for sensitive skin.
- Absorbs sweat while still maintaining adhesion.

Three Shapes:

ROUND BE 6085	OVAL BE 6086	OVAL EXTRA BE 6087
-------------------------	------------------------	------------------------------