

Blom-Singer® StomaSoft® Laryngectomy Tube



StomaSoft® Laryngectomy Tube:
HME-compatible with and
without fenestrations

Features include:

- Sizes 8, 9, 10, and 12
- Available in two lengths (36mm and 55mm)
- Designed with a reinforced strap holder
- Compatible with InHealth Technologies® HME systems and accessories
- Available with and without fenestrations

For ordering information, go to www.inhealth.com

Size/Length	I.D./O.D.	Non-Fenestrated	Fenestrated	Rx
8/36	9.5mm/12.0mm	BE 6698	BE 6698F	■
8/55	9.5mm/12.0mm	BE 6699	BE 6699F	■
9/36	10.5mm/13.5mm	BE 6700	BE 6700F	■
9/55	10.5mm/13.5mm	BE 6701	BE 6701F	■
10/36	12.0mm/15.0mm	BE 6702	BE 6702F	■
10/55	12.0mm/15.0mm	BE 6703	BE 6703F	■
12/36	13.5mm/17.0mm	BE 6704	BE 6704F	■
12/55	13.5mm/17.0mm	BE 6705	BE 6705F	■