

The Blom-Singer® Laryngectomy Button

**Breathe with confidence.
Speak with ease.**

The Blom-Singer® Laryngectomy Button is a self-retaining holder for devices intended for vocal and pulmonary rehabilitation after total laryngectomy. It can provide stomal support and help to keep the tracheostoma open.



**The Blom-Singer®
Laryngectomy Button
enhances communication
and overall quality of life.**

Its ergonomic design offers a comfortable, airtight fit, helping patients maintain an open tracheostoma for breathing and communication.

**Compatible with Blom-Singer®
SpeakFree® and Day&Night® Heat
and Moisture Exchange (HME)**
systems to support lung health.



Day&Night® HME



SpeakFree® HME



Shower Guard

BENEFITS:

- **Stoma Support** - Designed to maintain an open tracheostoma and prevent stomal stenosis
- **Pulmonary Rehabilitation** - Compatible with all Blom-Singer® HMEs and attachments
- **Custom Fit** - Available in multiple sizes for a secure fit, with wings for attachment points, if needed
- **Safety Features** - Sterile and individually packaged, made of biocompatible silicone, and MRI-safe

PRODUCT ORDERING INFORMATION

REF NUMBER	DESCRIPTION	SIZE	Rx
E3131208	Laryngectomy Button 12/8	OD 12, L 8	■
E3131408	Laryngectomy Button 14/8	OD 14, L 8	■
E3131608	Laryngectomy Button 16/8	OD 16, L 8	■
E3131808	Laryngectomy Button 18/8	OD 18, L 8	■

INHEALTH
TECHNOLOGIES®

United States
1110 Mark Ave.,
Carpinteria, CA 93013
Toll-Free: 800.477.5969
Toll-Free Fax: 888.371.1530
www.inhealth.com