

HME STYLES & FEATURES

GUIDE

INHEALTH
TECHNOLOGIES®

STYLES				
FEATURES	Blom-Singer® SpeakFree™ HME Hands-Free Valve with ClassicFlow ® HME	Blom-Singer® SpeakFree™ HME Hands-Free Valve with EasyFlow ® HME	Blom-Singer® Day&Night ClassicFlow ® HME	Blom-Singer® EasyFlow ® HME
Adjustable Hands Free Speech Valve	✓	✓		
Single-Use Device	✓	✓	✓	✓
Easy Manual Occlusion for Voicing	✓	✓	✓	✓
Breathability	GOOD	GOOD	BEST	BEST
Moisture Exchange	BETTER	GOOD	BETTER	GOOD
Hygroscopic Filter	✓	✓	✓	✓
Day and Night Wear	DAYTIME	DAYTIME	✓	✓
Compatible with all Blom-Singer Housings and StomaSoft Tube	✓	✓	✓	✓
Optional Safety Feature	CROSSBAR	CROSSBAR	TWIST TO LOCK	TWIST TO LOCK
MucusShield™ Protection				
Made with Latex-free Materials	✓	✓	✓	✓
Non-reflective Surface Texture	✓	✓	✓	✓
Reference Number	BE 1090EZ 30 pk	BE 1090EF 30 pk	BE1091 30 pk	BE1092 30 pk