

# The Blom-Singer® Laryngectomy Button\*

**Breathe with confidence.  
Speak with ease.**

The Blom-Singer® Laryngectomy Button is a self-retaining holder for devices intended for vocal and pulmonary rehabilitation after total laryngectomy. It can provide stomal support and help to keep the tracheostoma open.



## The Blom-Singer® Laryngectomy Button enhances communication and overall quality of life.

Its ergonomic design offers a comfortable, airtight fit, helping patients maintain an open tracheostoma for breathing and communication.

**Compatible with Blom-Singer®  
SpeakFree® and Day&Night® Heat  
and Moisture Exchange (HME)**  
systems to support lung health.



Day&Night® HME



SpeakFree® HME



Shower Guard

## BENEFITS:

- **Stoma Support** – Designed to maintain an open tracheostoma and prevent stomal stenosis
- **Pulmonary Rehabilitation** – Compatible with all Blom-Singer® HMEs and attachments
- **Custom Fit** – Available in multiple sizes for a secure fit, with wings for attachment points, if needed
- **Safety Features** – Sterile and individually packaged, made of biocompatible silicone, and MRI-safe

## PRODUCT ORDERING INFORMATION

REF NUMBER	DESCRIPTION	SIZE	Rx
E3131208	Laryngectomy Button 12/8	OD 12, L 8	■
E3131408	Laryngectomy Button 14/8	OD 14, L 8	■
E3131608	Laryngectomy Button 16/8	OD 16, L 8	■
E3131808	Laryngectomy Button 18/8	OD 18, L 8	■

\*Available Only in the U.S.