

# The Advantages of the Blom Singer® Nasal Septal Perforation Prosthesis



**Options in Fit:** The Blom-Singer® Nasal Septal Perforation Prosthesis is available in round and oval shapes, in several sizes, and thicknesses to close nasal septal perforations. A better fit means less movement and a reduction in bleeding and crusting.



**Magnetic coupling:** Neodymium magnets encased in medical grade silicone couple to provide minimally invasive closure of the nasal septal perforation.<sup>1</sup>




**Ease of insertion:** Studies have documented that patients using a two (2) piece magnetically coupled nasal septal prosthesis experienced no particular discomfort on insertion, removal or replacement while all cases reported a precise, tight fit. A magnetic-based technique ensures a more practical debris removal procedure for the attending ENT specialist.<sup>1</sup>



INHEALTH  
TECHNOLOGIES®

Blom-Singer®



## Want to learn more about the Blom-Singer® Nasal Septal Perforation Prosthesis?

Visit us at [www.nasalseptal.com](http://www.nasalseptal.com)  
Call us at 800-477-5969

<sup>1</sup> Teschner M, et al. Eur Arch Otorhinolaryngol. 2012 Mar;269(3):861-5

INHEALTH  
TECHNOLOGIES®  
[www.inhealth.com](http://www.inhealth.com)

United States  
Toll-Free: 800.477.5969  
Fax: 888.371.1530  
email: [info@inhealth.com](mailto:info@inhealth.com)

International  
Tel: +1.805.684.9337  
Fax: +1.805.684.8594  
email: [international@inhealth.com](mailto:international@inhealth.com)

Canada  
Toll-Free: +1.888.437.4123  
Fax: +1.226.641.5285  
email: [canadaorder@inhealth.com](mailto:canadaorder@inhealth.com)

## The Blom Singer® Nasal Septal Perforation Prosthesis

A two-piece, magnetically coupled solution for non-surgical closure of nasal septal perforations.

# The Blom Singer® Nasal Septal Perforation Prosthesis

A two-piece, magnetically coupled solution for non-surgical closure of nasal septal perforations.

**Simplified solution for nasal perforations:** A nasal septal perforation can vary in size and location, and it may be asymptomatic or cause a variety of symptoms.

For those with a nasal septal perforation, and those who live with people with a nasal septal perforation, the impact can be aggravating and uncomfortable.

## COMMON CAUSES OF NASAL SEPTAL PERFORATIONS:

- Complications from nasal surgery
- Chronic epistaxis
- Aggressive digital nasal cleaning (nose picking)
- Nasal piercings
- Chemical Inhalants
- Cocaine use


## POSSIBLE SYMPTOMS:

- Small anterior lesions may cause a whistling sound
- Larger perforations may cause crusting, bloody discharge, difficulty breathing, nasal pressure, and discomfort
- Rhinorrhea
- Epistaxis
- Infection

## PRODUCT INFORMATION

### ROUND C)))


5 mm



2 mm  
SP 1052

3 mm  
SP 1053


7 mm



2 mm  
SP 1072


3 mm  
SP 1073

10 mm



2 mm  
SP 1102


3 mm  
SP 1103



Flange width  
is 5mm for  
all prostheses

### OVAL C)))


11 x 13 mm



2 mm  
SP 2112

3 mm  
SP 2113


13 x 17 mm



2 mm  
SP 2132

3 mm  
SP 2133

15 x 21 mm



2 mm  
SP 2152

3 mm  
SP 2153

**A better fit means less movement and a reduction in bleeding and crusting.**