STYLES & FEATURES VOICE PROSTHESES



As a leading supplier of advanced voice restoration systems, our goal at InHealth Technologies® is to help laryngectomees regain and preserve a vital part of their lives—their voice.

We know that every person may have different voice prostheses needs. We understand those needs and offer a wide range of solutions with different prostheses types, lengths, large flanges, increased resistance, and combinations of these.

STYLES						9
			CLINICIAN - PLACE	:D		PATIENT - PLACED
	CLASSIC* AVAILABLE IN 16FR & 20FR	SPECIAL ORDER AVAILABLE IN 16FR & 20FR	DUAL VALVE AVAILABLE IN 20FR	ADVANTAGE HARD-VALVE AVAILABLE IN 20FR	ADVANTAGE SOFT-VALVE AVAILABLE IN 16FR & 20FR	LOW PRESSURE AVAILABLE IN 16FR & 20FR
	The Blom-Singer® Classic Indwelling: "our best seller" voice prosthesis. *Sterile option available.	Special order voice prostheses offer increased valve weight, specialty lengths, or large flanges to address specific clinical needs.	The Blom-Singer® Dual Valve is the first ever voice prosthesis with two valves in the shaft embedded with silver oxide material. The second valve provide an additional layer of protection against central leakage.	The Blom-Singer® Advantage Hard Valve Assembly voice prosthesis has a silver oxide material embedded in the valve of the prosthesis with the addition of a titanium ring to reinforce the biofilm resistance.	The Blom-Singer® Advantage Soft Valve Assembly voice prosthesis has a silver oxide material embedded in the valve for biofilm resistance.	The low-pressure voice prosthesis can be managed by the patient. This prosthesis is a soft material that can be placed with or without a gel cap.
FEATURES						
Anti-fungal equipment (silver oxide) - designed to prevent biofilm from embedding into the material of the prosthesis			•	•	•	
Special lengths - designed to address specific sizing needs		•				•



STYLES & FEATURES VOICE PROSTHESES



STYLES













CLINICIAN - PLACED

PAH	ENI	- P	LACED

FEATURES	CLASSIC AVAILABLE IN 16FR & 20FR	SPECIAL ORDER AVAILABLE IN 16FR & 20FR	DUAL VALVE AVAILABLE IN 20FR	ADVANTAGE HARD VALVE AVAILABLE IN 20FR	ADVANTAGE SOFT VALVE AVAILABLE IN 16FR & 20FR
Increased valve resistance - designed to prevent unintentional valve opening		•			
Large esophageal flange - designed to address peripheral leakage or prosthesis dislodgement		•			
Large esophageal and tracheal flanges - designed to address peripheral leakage or dislodgement		•	•		
Increased resistance with large esophageal flange - designed to address central and peripheral leakage		•			
Blom-Singer Gel Cap Insertion Method - designed to provide a comfortable and atraumatic voice prosthesis placement.	•	•	•	•	•

LOW PRESSURE AVAILABLE IN 16FR & 20FR
•
OPTIONAL

ACCESSORIES	CLASSIC AVAILABLE IN 16FR & 20FR	SPECIAL ORDER AVAILABLE IN 16FR & 20FR	DUAL VALUE AVAILABLE IN 20FR	ADVANTAGE HARD VALVE AVAILABLE IN 20FR	ADVANTAGE SOFT VALVE AVAILABLE IN 16FR & 20FR
Flushing Device - for routine cleaning of voice prostheses	•	•	•	•	•
Brush -for routine cleaning of voice prostheses	•	•	•	•	•
Plug Insert - for temporary management of leakage through voice prostheses	•	•		•	•

LOW PRESSURE AVAILABLE IN 16FR & 20FR
•