

# Katz Extractor<sup>®</sup>

OTO-RHINO FOREIGN BODY REMOVER

The easy way out!

The fast, flexible way to remove foreign bodies from noses and ears.

Order the  
Katz Extractor<sup>®</sup>  
Five-pack today!

Ready-to-use  
Atraumatic  
Quick extraction  
Economical



# Katz Extractor<sup>®</sup>

OTO-RHINO FOREIGN BODY REMOVER


A quick and simple solution for foreign body removal from the nose or ear.




The Katz Extractor<sup>®</sup> Oto-Rhino Foreign Body Remover is a small, lightweight, single-use device for quick and easy extraction of foreign bodies from the nose or external ear canal. It comes ready to use in a five-pack dispenser box. There are no expensive supplies, such as forceps, hooks or balloon tip catheters/syringes to locate and assemble. No additional suction equipment is required.

To see a demonstration video, Visit: [www.katzextractor.com](http://www.katzextractor.com)

*It takes just three simple steps to complete an extraction with the Katz Extractor!*



Insert past foreign body



Depress syringe to inflate balloon



Extract catheter out of canal with balloon inflated

## ORDER INFORMATION

Product Number:  
FB 5000

To Order Online Visit:  
[www.katzextractor.com](http://www.katzextractor.com)

Tel: 800.477.5969

Email: [order@inhealth.com](mailto:order@inhealth.com)

### CUSTOMER SERVICE

Our knowledgeable representatives are available to assist you Monday through Thursday 6:30 a.m. to 4:00 p.m. and Friday 6:30 a.m. to 2:00 p.m., PT.

### SHIPPING / PAYMENT

All orders are shipped standard five business-day service. Overnight shipping is available upon request for orders placed before 2:00 p.m., PT.

Designed in conjunction with David L. Katz, M.D., MPH, the Katz Extractor is intended principally for extraction of impacted foreign bodies in the nasal passages. This device may also be used in the external ear canal, based upon clinical judgement.

**INHEALTH**  
TECHNOLOGIES<sup>®</sup>



© 2022 InHealth Technologies Manufactured by Freudenberg Medical, LLC (210412.01)